Light Induced Fluorescence Evaluator

Diagnostic mode

Treatment mode






Healthy dentin





Healthy and thick enamel appears in blue at the cusps. Healthy dentin is green





End of the excavation




The floor of the cavity is green: the active caries is completely removed





Lack of fluorescence shows the presence of an enamodentinal carious lesion (demineralization)





Presence of green-orange hard tissues at the bottom of the cavity: slow evolution of decay completely removed while keeping the healthy tertiary dentin



CARIOUS LESIONS



A red signal surrounded by a lack of fluorescence (demineralization) indicates the presence of an enamodentinal carious lesion





Presence of green-orange hard tissues at the bottom of the cavity: slow evolution of decay completely removed while keeping the healthy tertiary dentin



A black zone surrounded by a demineralization (lack of fluorescence) shows the presence of an enamo-dentinal carious lesion






From the article of Pr. Tassery, Head of the Preventive and Restorative Department, Marseilles Dental School, University AMU, LBN Montpellier I. Published in Minim. Interv. Dent. 2011 and in Clin. Oral. Invest. 2013

The medical device for dental care SOPROLIFE is of class IIa and manufactured by SOPRO notified body LNE/GMED. This medical device is not reimbursed by health insurance organizations. Read carefully the instructions on the labelling before use.



SoproLife Caries Classification



From the article of Pr. Rechmann, Director of the Clinical Sciences Research Group, Department of Preventive and Restorative Dental Sciences, School of Dentistry, University of California, San Francisco. Published in the JBO in April 2012 From the article of Pr. Tassery, Head of the Preventive and Restorative Department, Marseilles Dental School, University AMU, LBN Montpellier I. Published in Minim. Interv. Dent. 2011 and in Clin. Oral. Invest. 2013

