

**3M** Science.  
Applied to Life.™

Elipar™ DeepCure-S  
LED Curing Light

**Have  
confidence  
in a deep,  
uniform cure.**



# Be sure of your cure.

## Elipar™ DeepCure-S LED Curing Light


Do you frequently use posterior or bulk fill materials? Research shows that 69% of dentists who place bulk fill restorations are not confident of polymerization deep in the cavity.<sup>1</sup> We took this as a challenge and created a new high-performing light.



### A deep, uniform cure.

The Elipar™ DeepCure-S LED Curing Light holds true to its name. Due to optimized optics, you can be confident that your restorations will have a deep, uniform cure—from center to rim, from surface to cavity bottom and at clinically relevant distances. Laboratory test results below prove why.

#### Comparison of light penetration



More homogeneous energy distribution throughout the restoration. Images comparing the light penetration of various light curing devices show that the Elipar™ DeepCure-S LED Curing Light produces a more collimated and uniform beam profile—even in deeper areas.


Source: 3M internal data

<sup>1</sup> Clinicians Report, October 2014, Volume 7, Issue 10

# What does this mean for your clinical results?


It means a high depth of cure, especially if the light positioning is difficult.

The Elipar™ DeepCure-S LED Curing Light helps to compensate for slight movements during curing, delivering a high depth of cure, as shown below.



## Designed for comfort.

The light guide of the Elipar DeepCure-S LED curing light is designed with an optimum angle and height that makes it more comfortable for patients—even those with limited mouth-opening abilities. Clinicians also find it comfortable to use, with easier access to all tooth surfaces, even in the hard-to-reach posterior area.



## Dentist satisfaction ratings\*

98% were satisfied or very satisfied with the light.

90% said it increases confidence of a full cure to the bottom of the proximal box.

\*Source: Application test conducted with dentists in the U.S.A., Germany, Turkey and Denmark.

## For all clinical procedures where light curing matters.

The Elipar DeepCure-S LED curing light can be used with all 3M materials to produce reliable results.\*\*


- Filtek™ Composites
- Scotchbond™ Universal Adhesive
- RelyX™ Cements
- Clinpro™ Sealant



\*\*Elipar™ DeepCure-S LED Curing Light is compatible with all campherquinone- and Ivocerin-based photoinitiators.

# Elipar™ DeepCure-S LED Curing Light

For dentists who enjoy the look, feel and durability of high-quality stainless steel



## Technical Performance Data

<b>Wavelength</b>	430–480 nm
<b>Light intensity</b>	1,470 mW/cm <sup>2</sup> (-10%/+20%)
<b>Power supply</b>	Lithium-ion battery Approx. 120 min. battery runtime (~720 10-sec. cures) with constant light output regardless of battery charge
<b>Operation</b>	Intuitive two-button and single-mode operation Pre-set cure times: 5, 10, 15 and 20 seconds, continuous mode (120 sec.) and tack-cure mode
<b>Curing time</b>	Refer to material instructions; 10 sec. for many composites
<b>Light guide</b>	10 mm; black coated; autoclavable; optimal intraoral reach due to user- and patient-friendly geometry

## Ordering Information—Stainless Steel

Item #	Product Information
76976	<b>Elipar™ DeepCure-S LED Curing Light</b> Contains: Handpiece (Cordless), Charging Base; Li-ion Battery; 10mm Light Guide; Eye Shield
76981	<b>Elipar™ DeepCure-S Light Guide, 10mm</b>
76984	<b>Elipar™ DeepCure Eye Shield</b>
76985	<b>Elipar™ DeepCure-S Rechargeable Li-ion Battery</b>

[www.3M.com/CuringLights](http://www.3M.com/CuringLights)



3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA  
Phone 1-800-634-2249  
Web 3M.com/dental

3M Canada  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
Phone 1-888-363-3685

3M, ESPE, Clinpro, Elipar, Filtek, RelyX and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2016. All rights reserved. All other trademarks are not trademarks of 3M. Please recycle. Printed in U.S.A.  
70-2013-0701-7