

PLANMECA

3D imaging



ENGLISH

Passion to innovate

An introduction from our President



“Welcome to the future of digital imaging. It gives me great pleasure to introduce you to our world-leading **Planmeca ProMax® 3D** X-ray units and **Planmeca Romexis®** imaging software – with a pioneering combination of 3D images that takes you closer for an even greater understanding of what your patients need.

I’m extremely proud of our product innovations, and for over 40 years we’ve worked closely with dental professionals to set new standards in our field. What makes us a bit different is that all core product development and manufacturing takes place in Finland – ensuring exceptional quality and unmatched attention to detail at every stage of the process.

This brings us to our **Planmeca ProMax®** product family, taking care of all your 2D and 3D imaging needs in a single unit. Each product is a true all-in-one unit, offering easy-to-use controls and incredible patient comfort. We have a dedicated team of in-house R&D professionals behind the scenes, all determined to make the best possible products for you and your patients. Therefore I am thrilled to invite you to discover our complete selection of advanced 3D solutions.”

*Heikki Kyösti
President and founder
Planmeca Group*

Fantastic five.....	4
Unique 3D combination – an industry first	6
CBCT	8
3D face photo	10
3D model scanning.....	12
Real-time jaw movement – in 3D	14
Planmeca ProMax® 3D family Key features	16
Ease of operation	18
Advanced technology.....	20
Pioneering low dose 3D imaging.....	22
2D and 3D imaging with one sensor	24
Extraoral bitewings.....	26
Quality cephalometry for orthodontics.....	28
Planmeca Romexis® – one software for all your needs.....	30
Reinventing 3D imaging	32
The complete implant workflow.....	34
3D tools for orthodontists and dental labs.....	36
Access to unique X-ray device data	38
Your mobile world of imaging.....	39
Share images and expertise online.....	40
Professionals proudly present the Planmeca ProMax® 3D family	42
Planmeca ProMax® 3D s.....	44
Planmeca ProMax® 3D Classic.....	46
Planmeca ProMax® 3D Plus.....	48
Planmeca ProMax® 3D Mid	50
Planmeca ProMax® 3D Max.....	52
Patient specific implants	54
Stand out with colour	56
Technical specifications	56

Fantastic five

Meet the Planmeca ProMax® 3D family

Mac and Windows
compatible

New!

Extended volume
Ø110 x 80 mm



Planmeca ProMax® 3D s



Planmeca ProMax® 3D Classic



Planmeca ProMax® 3D Plus



Planmeca ProMax® 3D Mid



Planmeca ProMax® 3D Max

*Planmeca ProMax® 3D is a product family consisting of exceptional all-in-one units. With three different types of three-dimensional imaging – as well as panoramic, extraoral bitewing and cephalometric imaging – **these intelligent products can meet all your maxillofacial imaging needs.***

True all-in-one units for all your imaging needs.

Unique 3D combination – an industry first

*We're the first company to combine three different types of 3D data with one X-ray unit. The **Planmeca ProMax® 3D** family brings together a Cone Beam Computed Tomography (CBCT) image, 3D face photo and 3D model scan into one 3D image – using the same advanced software. This 3D combination creates a virtual patient in 3D, helping you with all your clinical needs.*




CBCT

Why Cone Beam Computed Tomography?

Cone Beam Computed Tomography (CBCT) is an X-ray imaging technology where a large number of 2D images are taken of a patient from different angles. A 3D volumetric image is then calculated from these 2D projections. The resulting images can be viewed with our advanced imaging software from any angle, including the axial, coronal, sagittal and cross-sectional planes.

See more
than ever
before



Renowned dental implant surgeon Franck Renouard couldn't imagine working without his CBCT

"I acquired my **Planmeca ProMax® 3D Classic** in 2007 and was one of the very first users in France. The choice was quickly made, as Planmeca's unit was far ahead of its competitors."

All necessary diagnostic information from a single unit

"In implant cases, I usually start the analysis by taking a panoramic radiograph or a simple intraoral radiograph. As soon as I discover an ambiguity or low volume, I go for CBCT. For some indications, such as sinus lifts or onlay bone graftings, I always use CBCT."

I usually take large-volume studies right away. This enables diagnosis of endodontic or bone lesions that could otherwise go undetected in sites other than the implant site. The unit's resolution is more than sufficient for everyday examinations. In cases where a patient has a lot of metal restorations in their mouth or a problem with staying still, the artefact filter is very useful."

Full sinus visualisation

"I take a CBCT study systematically before a sinus lift. It provides visualisation of the sinus anatomy and allows me to see if there is an intrasinus pathology or anatomical features such as septa. It also allows detecting possible thick antral arteries, which are common in the bone wall. I need to be aware of these parameters before surgery."

Nowadays, I do not like to receive paper-based exams, as radiologists often do not provide the axial views which are essential in sinus study. When I take the radiographs myself, I can choose the slice that interests me."

Improved patient satisfaction

I then use CBCT to check the integration of my graft when I fill the tissue before implant placement. I also need CBCT when there are complications, which occur in 3–5% of cases. I can assure you that the patients are delighted to receive their diagnosis and treatment immediately, with no need to visit a specialist radiology centre."

Planmeca ProMax 3D Classic is a well-designed radiology unit. All in all, using 3D has become natural to me. I only wonder how we managed before 3D!"

Dr Franck Renouard, DDS, Paris, France



Dr Renouard specifies that he has not received any financial compensation or other benefit for this interview.

3D face photo

Planmeca ProFace® is an exclusive 3D face photo system available for all of our 3D X-ray units. This pioneering integrated system produces a realistic 3D face photo and CBCT image in a single imaging session. You can also take a separate 3D face photo without exposing your patient to any radiation.

The world's first
X-ray integrated
face camera



Planmeca ProFace® – the face in 3D

Designed to fulfil the most diverse diagnostic needs of today's maxillofacial and dental professionals, **Planmeca ProFace®** is a highly effective tool for pre-operative planning and treatment follow-up. It's also ideal for patient motivation and for sharing information with colleagues.

Safer and faster facial surgery

The 3D photo visualises soft tissue in relation to dentine and facial bones. As both a CBCT image and a 3D photo are generated in one imaging session, the patient position, facial expression, and muscle position remain unchanged – resulting in images that are perfectly compatible.


Careful pre-operative planning – where you can study the facial anatomy thoroughly using our **Planmeca Romexis®** software – facilitates accurate and detailed operations and enhances the aesthetic result.



Pre and
post-operative
comparisons



Measure distances
and relationships
between bone
and soft tissue



Superimpose images
for comparison



Shows deviation
to instantly see
any changes



Create a 2D photo
series automatically

3D model scanning

*You can use all X-ray units in the **Planmeca ProMax® 3D** family to scan both impressions and plaster casts – an exciting feature that was an industry first for our CBCT units. And with our advanced **Planmeca Romexis®** software, the digitised models are available immediately and stored for later use.*



Scanning a plaster cast to a digital model



Scanning an impression to a digital model

Advantages of 3D model scanning

Digital models save space


3D digital models are stored in the **Planmeca Romexis®** database in standard STL format, which reduces the need to make or maintain physical plaster casts.

Create your virtual patient


The scanned 3D model can be superimposed on to CBCT data, creating a virtual patient and helping you with all your clinical and treatment planning needs. The combined data set provides an artefact-free model of your patient's dentition including bone, crowns and soft tissue. This offers valuable new options for implant planning, surgical guide manufacturing, orthodontic purposes and orthognathic surgery.



Scanned impressions of upper and lower arches and bite index in 3D



Upper and lower arch models in occlusion. A useful tool for orthodontic treatment planning and patient progress follow-up



Superimposed CBCT and 3D model of upper jaw. Measure, compare, and track changes in teeth movements



Crown, impression scan, and CBCT for more accurate implant planning



Analyse the STL data further in the Planmeca Romexis® 3D Ortho Studio module, and then carry out a comprehensive dental cast analysis and create an orthodontic treatment plan

Real-time jaw movement – in 3D

Planmeca 4D™ Jaw Motion is the only true CBCT integrated solution for tracking, visualising, and analysing jaw movement in 3D. It offers incomparable visualisation of mandibular 3D movements in real-time – creating a fourth dimension in diagnostics.



Key components of Planmeca 4D™ Jaw Motion

Planmeca 4D® Jaw Motion adds a new dimension to 3D data by visualising a patient's jaw movement. First, a CBCT image (e.g. a **Planmeca Ultra Low Dose™** image) is acquired with a **Planmeca ProMax® 3D** unit with the patient wearing dedicated tracking devices. Integrated **Planmeca ProFace®** cameras are then used to track lower jaw movements in relation to the upper jaw. All movements are visualised, analysed, and stored to the **Planmeca Romexis®** imaging software in real time.

Applications:

Due to its capability to visualise mandibular jaw and condyle movement, Planmeca 4D Jaw Motion can be a supporting tool for:

- Temporomandibular (TMD) examinations
- Preoperative planning and postoperative treatment verifications
- Articulator programming

Key features:

- The only CBCT integrated jaw tracking solution
- Track, visualise, and record jaw movement in 3D
- Visualise movements in the Planmeca Romexis software in real time
- Record movements for later use and analysis
- Measure and visualise the movement paths of points of interest in frontal, sagittal, and axial movement graphs and in 3D
- Align digital dental models with a CBCT image for improved visualisation
- Export movement and measurement information to 3rd party software for analyses and treatment planning



Planmeca

ProMax® 3D family

Key features

Advanced technology:

- Ideal resolutions and optimal balance between image quality and patient dose – always complying with the ALARA (*As Low As Reasonably Achievable*) principle
- The pioneering **Planmeca Ultra Low Dose™** protocol enables CBCT imaging with an even lower dose than traditional 2D panoramic imaging
- Optimal volume size and location for every clinical need
- Special imaging protocols for dental and ENT applications
- Certified for use with the **suressmile** system for orthodontics

Effortless use:

- Effortless patient positioning and unmatched comfort
- True all-in-one X-ray units not only for 3D imaging, but 2D panoramic and cephalometric imaging as well
- Easy to use for a smooth workflow
- **ProTouch™ Desktop** for remote control panel operation on the imaging workstation
- **Planmeca Romexis®** software
- Mac and Windows support



Some modalities mentioned may not be compatible with the full range of X-ray units in the Planmeca ProMax® 3D family.
For full details on availability, see technical specifications at the end of the brochure.

Ease of operation

Our **Planmeca ProMax® 3D** units are known across the world for incredible ease of use and exceptional patient comfort. A relaxed patient means a smooth imaging workflow and the best quality images.

Open patient positioning


- Effortless positioning with open-face architecture
- Unrestricted view of your patient
- No claustrophobic feeling for your patient
- Fine adjustment using positioning lasers and joystick
- Verify correct positioning with a scout image
- Easy wheelchair accommodation with side-entry access

Unmatched
patient
support



Scout images for easy positioning

Scout images and 2D views help positioning and can even be used for preliminary diagnosis.



User-friendly Planmeca ProTouch™ control panel

- Clear and straightforward graphical user interface guides you smoothly through the work process
- Pre-programmed sites and exposure values for different image types and targets save you time and allow you to focus on your patients
- The control panel can also be operated remotely from the imaging workstation

Easy imaging with ready-designed protocols

- Imaging protocols designed for specific diagnostic tasks, areas, or target sizes
- Appropriate volume size, resolution, and exposure values
- Automatic selection and adjustment of the target position
- Reduced volume sizes for child patients to prevent unnecessary radiation



Advanced technology

Our intelligent high-tech solutions and algorithms guarantee an ideal imaging geometry, perfect usability, and crystal-clear images free from noise and artefacts.

Intelligent solutions for the best image quality

SCARA technology

The precise, free-flowing, computer-controlled SCARA (*Selectively Compliant Articulated Robot Arm*) arm construction can produce any movement pattern required. This enables accurate and reliable volume positioning and volume diameter adjustment, reducing the amount of radiation your patients are exposed to.

New 120 kV tube voltage


120 kV tube voltage enables optimised image quality for challenging targets – reducing artefacts and ensuring higher contrast images.

Optimised imaging modes for different needs

- **Low dose** mode captures an image with a minimal dose of radiation. Ideally suited for orthodontic, pediatric and sinus studies. Voxel size 400 or 600 μm
- **Normal** mode is the best choice for most common imaging needs. Voxel size 200 μm
- **High definition** mode is designed for imaging of small objects, such as ear bones. Voxel size 150 μm
- **Braces** protocol offers optimised exposure settings for imaging patients with brackets. Voxel size 150 μm
- **High resolution** provides more detail, when necessary. Voxel size 100 μm
- **Endodontic** mode offers the best resolution with the smallest size. Voxel size 75 μm

ROI for higher resolution images

The ROI (Region of Interest) reconstruction function can generate a new small voxel volume from the image data of a previously taken large voxel volume. This enables a more precise diagnosis without the need for an additional radiation dose for the patient.



Original volume,
400 μm



Reconstructed volume
with higher resolution, 200 μm

Never miss a shot with Planmeca CBCT units

Movement, metal artefacts, and small voxel sizes are generally recognised as challenges to CBCT image quality. With Planmeca CBCT units and their advanced image enhancement options, you can rise above these concerns and succeed every time. The options can either be selected preventively before imaging or utilised afterwards to achieve reliable results. The choice is yours!

Movement artefact correction with Planmeca CALM™

- Iterative movement correction algorithm
- Eliminates the need for retakes
- Cancels the effects of patient movement
- Excellent when imaging more lively patients

Metal artefact reduction with Planmeca ARA™

- Reliable algorithm for artefact-free images
- Removes shadows and streaks caused by metal restorations and root fillings
- Tried and tested – the results of extensive scientific research

Noise removal with Planmeca AINO™

- Noise-free images without losing valuable details
- Allows lower exposure values by reducing noise
- Improves image quality when using small voxel sizes (e.g. in the endodontic imaging mode)
- Enabled by default when using the Planmeca **Ultra Low Dose™** imaging protocol

Endodontic mode

The endodontic imaging mode provides perfect visualisation of even the finest anatomical details. This advanced imaging mode is an ideal choice for endodontics and other cases with small details.

- Extremely high resolution with 75 μm voxel size
- Enables precise diagnostics and treatment plans


Without movement artefact correction




With the Planmeca CALM™ movement removal algorithm




Without artefact removal





With the Planmeca ARA™ artefact removal algorithm



Without noise removal



With the Planmeca AINO™ noise filter



Pioneering low dose 3D imaging

Planmeca ProMax® 3D units offer a unique Planmeca Ultra Low Dose™ imaging protocol that enables CBCT imaging with an even lower patient radiation dose than standard 2D panoramic imaging.

More information, less radiation

Planmeca Ultra Low Dose™ can be used with all voxel sizes and in all imaging modes from Normal to Endodontic mode. Using the Planmeca Ultra Low Dose protocol reduces the effective patient dose by an average of 77% without a statistical reduction in image quality*.

The unique and pioneering imaging protocol is based on intelligent 3D algorithms developed by Planmeca. Our 3D imaging system always allows the clinician to choose the optimal balance between image quality and dose, based on the ALARA principle.



* Study of Orthodontic Diagnostic FOVs Using Low Dose CBCT protocol (Ludlow, John Barrett and Koivisto, Juha).

planmeca.com/ULD-poster

Ideal for many clinical cases


The Planmeca Ultra Low Dose protocol has proven to be ideal for many clinical cases.

- Orthodontics:
 - Defining the amount of bone around the root
 - Localising unerupted and impacted teeth before orthodontic treatment
 - Defining orthodontic landmarks for cephalometric analysis
- Post-operative and follow-up images in maxillofacial surgery
- Airway studies
- Sinus studies
- Implant planning




Planmeca ProMax® 3D Mid

- FOV Ø 200 x 170 mm / Voxel size 600 µm
- Effective patient dose 14.7 µSv




Planmeca ProMax® 3D Classic

- FOV Ø 40 x 50 mm / Voxel size 150 µm
- Effective patient dose 14.4 µSv



Planmeca ProMax® 3D Max

- FOV Ø 85 x 50 mm / Voxel size 400 µm
- Effective patient dose 4.0 µSv



Planmeca ProMax® 3D Mid

- FOV Ø 200 x 170 mm / Voxel size 600 µm
- Effective patient dose 29.2 µSv

The Planmeca Ultra Low Dose™ protocol has changed 3D imaging completely

We at MESANTIS® 3D DENTAL-RADIOLOGICUM produce about 7,500 CBCT images per year at eight locations in Germany.

Our main concern in X-ray imaging is to reduce the possible radiation dose as much as is reasonably achievable (ALARA principle). Traditional digital 2D X-rays at an orthodontist's clinic usually have an effective dose ranging between 26–35 µSv (ICRP 2007). Conventional CBCT images of the head with modern CBCT equipment show an effective dose ranging between 49–90 µSv.

The latest image protocol with a specific associated algorithm is called the Planmeca Ultra Low Dose™ protocol. In medical terms, it allows radiologists to adjust imaging parameters individually according to the clinical needs of each case. The mA-values, in particular, can be individually adjusted and reduced for each patient, as it is required according to all international scientific guidelines. Therefore, it is possible to

further reduce the effective dose significantly by using the Planmeca Ultra Low Dose protocol. Depending on the field of view, nowadays CBCT equipment with a Planmeca Ultra Low Dose algorithm has an effective dose between 4 to 22 or 10 to 36 µSv.

Our patients and referring colleagues are always happy to hear that the effective dose for certain indications is now even lower than in traditional 2D X-ray imaging. Since last year we have been able to replace the common CBCT protocols with the Planmeca Ultra Low Dose protocol.

At MESANTIS® 3D DENTAL-RADIOLOGICUM in Germany, the Planmeca Ultra Low Dose imaging protocol is used either with a small or large field of view. Using the new protocol, a lot of patients can benefit from improved 3D diagnostics without being exposed to a higher radiation dose.


Prof. Dr. Axel Bumann

Prof. Dr. Bumann states that he has not received any financial reward or other benefit for this interview.

Prof. Dr. Axel Bumann
DDS, PhD, Orthodontist,
Oral surgeon, Oral and
Maxillofacial Radiology,
MESANTIS® 3D
DENTAL-RADIOLOGICUM

2D and 3D imaging with one sensor

Our advanced **SmartPan™** imaging system uses the same 3D sensor also for 2D panoramic imaging.




2D SmartPan™ – unique panoramic imaging

- A unique system for 2D imaging
- Uses the same 3D sensor for 2D panoramic imaging, eliminating the need to change sensors
- Users can browse between panoramic images and select the most suitable one for diagnosis
- Same patient positioning and image processing parameters as in 2D imaging programs

2D programs

Standard: Basic panoramic programs	Standard panoramic Lateral TMJ (closed & open) PA TMJ (closed & open) PA sinus
Standard	Child (Paediatric) mode for each standard and optional program to reduce the dose
Optional	Horizontal and vertical segmenting for panoramic program
Optional	True Bitewing
Optional: Advanced panoramic programs	Interproximal panoramic Orthogonal (perio) panoramic Lateral-PA TMJ Lateral multiangle TMJ PA multiangle TMJ PA linear sinus Lateral sinus



*Normal **SmartPan™** produces 9 different parallel panoramic layers with about 2 mm shift and one autofocus layer*

***MultiView SmartPan™** calculates 9 different rotated panoramic layers. This allows adjusting the view angle for improved diagnosis.*



Extraoral bitewings


What if you could do all your routine diagnostic imaging extraorally?

Planmeca ProMax® extraoral bitewings are ideal for periodontics, elderly and child patients, claustrophobic patients, patients with a strong gag reflex, and patients in pain. Extraoral bitewings enhance clinical efficiency and take less time and effort than conventional intraoral bitewing imaging.

What are the advantages of extraoral bitewings?

- Ideal for all patients – no sensor positioning required
- Consistently opens interproximal contacts, giving better diagnostic value
- Larger diagnostic area than in intraoral modalities
- More clinical data: canine to third molar
- Enhanced clinical efficiency – takes less time and effort than conventional intraoral bitewings
- Enhanced patient experience and comfort – eliminates gagging

Better diagnostic value with extraoral bitewings




True Bitewing program, adult




Standard panoramic image of the same patient as the bitewing above

True bitewings only possible with our SCARA3 technology



True Bitewing program, 5-year-old child



True Bitewing program, 8-year-old child




Quality cephalometry for orthodontics

Our exceptional equipment and advanced software have been designed to meet all your orthodontic needs.

Cephalometric imaging with Planmeca ProMax® units

- The functional and easy-to-use head positioner ensures accurate positioning for all cephalometric projections
- The carbon fibre ear posts and nasal positioner are extremely stable, hygienic, and transparent to radiation
- The unit automatically aligns itself to take cephalometric exposures and then selects a corresponding collimator
- The rotating tube head in the 3D unit eliminates the need to remove the 3D sensor
- Dedicated collimation options for paediatric imaging



Two equipment options:

One-shot Planmeca ProCeph™ cephalostat

- Effective one-shot cephalostat
- Short exposure time – no motion artefacts, low patient dose
- Image sizes from 18 x 20 cm to 30 x 25 cm

Scanning Planmeca ProMax® cephalostat

- Digital cephalostat that scans your patient's head horizontally using a narrow X-ray beam with an extremely low effective dose of radiation
- Exceptional flexibility in image formats, with field sizes of up to 30 x 27 cm

Easier and
more accurate
than
ever before

Two options for cephalometric analyses:

Planmeca Romexis® Cephalometric Analysis module

Take advantage of the **Planmeca Romexis® Cephalometric Analysis** module's wide range orthodontic and orthognathic tools.

- Tools for creating cephalometric analyses, superimpositions, and surgical treatment plans (VTO) in minutes
- Fully customisable analyses, norms, and reports
- Microsoft Excel export and import function
- Compatible with the Windows operating system

Online automatic analysis service

Acquire cephalometric analyses regardless of time and place with the **Planmeca Romexis®** automatic cephalometric analysis service.

- Online automatic cephalometric tracing in a few seconds
- Over 50 analyses available for download immediately after tracing
- Direct link from the Planmeca Romexis 2D module for ordering analyses

New!

Planmeca Romexis® one software for all your needs


*We offer a revolutionary all-in-one software solution for clinics of all sizes. Our world-leading **Planmeca Romexis®** software is the brains behind all of our products, bringing together all the devices at your dental clinic from CAD/CAM to imaging devices and dental units. It supports the most versatile range of 2D and 3D imaging modalities.*



Imaging and CAD/CAM in one software
– an industry first

Reinventing 3D imaging

Our pioneering **Planmeca Romexis®** software offers specially designed tools for implantologists, endodontists, periodontists, prosthodontists, orthodontists, maxillofacial surgeons, and radiologists. You can also view your images wherever you are using our mobile apps, and enjoy unmatched compatibility with other systems.



Excellent tools for quality images

With a complete set of tools for image viewing, enhancement, measurement, drawing and annotations, **Planmeca Romexis®** improves the diagnostic value of radiographs. Versatile printing and image import and export functionalities are also included. The software consists of different modules – so you can choose those most suited to your needs.

Convenient 3D diagnosis


The Planmeca Romexis 3D rendering view gives an immediate overview of the anatomy and serves as an excellent patient education tool. The images can be instantly viewed from different projections or converted into panoramic images and cross-sectional slices. Measuring and annotation tools – such as nerve canal tracing – assist in safe and accurate treatment planning.

Free Planmeca
Romexis® Viewer
application
planmeca.com/Viewer


Full-featured viewer application
No installation required
Mac and Windows support
Distribute to specialists
or patients

Superimpose CBCT

New to Planmeca Romexis 3D, the module allows the superimposition of two CBCT images. It is a valuable tool for before-and-after comparisons and can be used for orthognathic surgery follow-ups, as well as orthodontic treatments, for example. The module also allows users to compare CBCT and MRI images side by side – providing a comprehensive view of a patient's anatomy.



Tampere University Hospital, Medical Imaging Center, Finland



Tooth segmentation


Planmeca Romexis provides a new, intuitive and efficient tool for segmenting a tooth and its root from a CBCT image. The guided process enables quick segmentation of a patient's full dentition. Surface models of segmented teeth can be visualised, measured and utilised e.g. in **Planmeca Romexis® 3D Ortho Studio** as part of orthodontic treatments.

Easy sharing of results

Studies can be quickly converted into multi-page printouts or handed out with the free **Planmeca Romexis® Viewer** media. Cases can be seamlessly transferred to mobile devices or partner clinics that also use Planmeca Romexis.

Best compatibility with other systems


Planmeca Romexis offers excellent compatibility with other systems, allowing you to freely use third-party products at your clinic. TWAIN support and DICOM standard compliance ensure that our flexible software can be used effortlessly with most systems.



Visualise and measure airways and sinus volumes before and after treatment for simplified diagnosis and treatment planning. Our advanced software tools allow accurate measurements in 3D space. Measurements can easily be reviewed using the saved views.

The complete implant workflow

Our **Planmeca Romexis® 3D Implant Planning** module offers all the necessary tools for fully digital implantology – from planning to guided surgery. The software's implant library includes realistic implant models and abutments, as well as collections of sleeves for guided surgery. After completing the implant plan, a surgical guide can be immediately designed in the same **Planmeca Romexis®** software with just a few clicks.




Realistic
implant models
from over 70
manufacturers

The **Planmeca Romexis®** software platform provides the perfect environment for top-down implant planning. By superimposing a crown and dental model onto CBCT data, users can create a complete virtual setup for optimally positioning the implant – taking prosthodontic and surgical perspectives into account.




Mark the nerve on the CBCT image


Superimpose the 3D model scan onto the CBCT image with the Planmeca Romexis® software



Use the Planmeca Romexis® crown library, or import a patient-specific crown from the CAD system to the software



Select the preferred implant (abutment) and sleeve from the extensive Planmeca Romexis® library and find the optimal position for it from a prosthetic and surgical perspective



Design the surgical implant guide with just a few clicks in Planmeca Romexis® – the software will create an open STL file of the design




Print the surgical guide with Planmeca Creo™ or any other suitable 3D printer

3D tools for orthodontists and dental labs



The **Planmeca Romexis® 3D Ortho Studio** module provides orthodontists and dental laboratories with several innovative tools. The advanced software is designed for the examination and analysis of digital dental models scanned e.g. with **Planmeca ProMax® 3D X-ray** units or the **Planmeca Emerald™** or **Planmeca PlanScan®** intraoral scanners. It includes an extensive set of premium tools for treatment planning in 3D.



Dental model analysis

Dental impressions and plaster casts scanned with the **Planmeca ProMax® 3D** model scanning mode can be aligned with the bite index using the **Planmeca Romexis®** software. Examination, analysis and treatment planning are then conveniently done in the **Planmeca Romexis® 3D Ortho Studio** module.

The module makes dental model analysis easier than ever by offering all the necessary tools for virtual base creation, occlusion inspection, and versatile teeth and arch measurements.





Planmeca Emerald™



Plaster cast in Planmeca ProMax® 3D



Impression scan in Planmeca ProMax® 3D



Treatment planning and verification in 3D

A staged treatment plan can be established in Planmeca Romexis 3D Ortho Studio by displacing the teeth in a virtual tooth setup while visualising intersections and contacts.

For improved visualisation, segmented roots and bone surfaces from CBCT images can be combined.

All the applied changes such as tooth movements, interproximal reductions, and tooth extractions are summarised in a detailed treatment plan report. The plan can be easily shared with others.

It is also possible to make 3D comparisons of treatment plan models and patient scans to verify treatment progress.



Export of digital dental models in STL format

Planmeca Romexis 3D Ortho Studio generates a series of digital dental models for each treatment stage. The models can be exported in STL format for 3D printing and custom appliance design and manufacturing.

The module is compatible with the Windows operating system.



The Braces imaging protocol for Planmeca ProMax® 3D units

Our special Braces* imaging protocol offers optimised exposure settings for imaging patients with brackets. The acquired images can be utilised when designing braces. The Braces imaging protocol is optimised for use with **suressmile**.


suressmile®
to be sure.

Suresmile® certification

Planmeca ProMax 3D units have been certified for use with **suressmile** by OraMetrix.

The comprehensive treatment management system **suressmile** is designed to enable orthodontists to visualise and simulate multiple diagnostic set-ups and design customised archwires for every patient. The accuracy of the patient scans plays a critical role in maximising the effectiveness of the system.

*Available for Planmeca ProMax® 3D Classic, Planmeca ProMax® 3D Mid and Planmeca ProMax® 3D Max.



Access to unique X-ray device data

Take the efficiency of your clinic to the next level with real-time information on networked equipment usage and events. Our **Romexis® Clinic Management** software offers several quality assurance and service benefits for local users, whereas **Romexis® Insights** allows you to remotely monitor your clinic from anywhere.

Romexis® Clinic Management – fluent and safe use of equipment

- See clear graphical overviews of a clinic – showing equipment status, occupancy, and users
- Enable local network access

Romexis® Insights – consolidated online monitoring of all equipment

- Monitor equipment from anywhere over the internet – also using mobile devices
- Utilise interactive dashboard views to see statistics from all clinics or individual locations and equipment
- Monitor trends and changes to clinic operations using informative graphics
- Allow stakeholders (such as service technicians) to securely access equipment information

Key benefits of networked equipment:

Planmeca equipment can be networked to gather valuable data on their use.

- Enhance operational planning – usage hours
- Use detailed event logs to improve quality assurance – including radiation hygiene
- Maximise equipment uptime with fast and accurate trouble-shooting

Real-time monitoring of day-to-day operations for clinic staff with Romexis® Clinic Management



All data in the clinic database

Local monitoring

Cloud reporting

Advanced operational data analytics for business stakeholders with Romexis® Insights



All data in a cloud database

Your mobile world of imaging

Our advanced **Planmeca mRomexis™** imaging application for iOS and Android allows you to flexibly view and capture images on mobile tablet devices. Remove the constraint of place – easily consult with colleagues and communicate with patients both in and outside your clinic.

Increased flexibility with Planmeca mRomexis™

Use our fast, easy, and light **Planmeca mRomexis™** mobile imaging application to view all your images in the **Planmeca Romexis®** database on a local network, or to carry images with you on your tablet device. You can also use the application to capture 2D X-ray images with Planmeca equipment, or to take photos with the tablet camera.

Expand the possibilities of Planmeca Romexis and experience the new level of freedom our mobile world can offer!

For
iOS and
Android



Key benefits:

- Available for both iOS and Android tablets
- Supports an extensive range of images – 2D and 3D X-ray images, 3D dental models, STL files, **Planmeca ProFace®** facial photos, and standard photos
- Direct connectivity with the **Planmeca Romexis®** server for retrieving or saving images
- Convenient acquisition of 2D X-ray images with Planmeca equipment
- Capturing photos with the camera of the mobile device
- Voice annotations to images can be recorded using the mobile device's microphone
- Flexible and secure retrieving of images via the **Planmeca Romexis® Cloud** image transfer service
- Excellent tool for patient education and communication



Download the **Planmeca mRomexis™** application for iOS and Android from the App Store or Google Play.



Share images and expertise online



Planmeca Romexis® Cloud is a secure image transfer service for Planmeca Romexis® users and their partners. Now you can easily share images and CAD/CAM cases with any specialist or patient.

Planmeca Romexis® Cloud

IMAGES
REFERRALS
INTERPRETATIONS
TREATMENT PLANS

Anybody, anywhere

- General practitioner
- Colleague
- Radiologist
- Specialist
- Dental lab
- Patient

Planmeca Romexis® user

- Radiology center
- General practice

Advantages

- Seamlessly integrated into **Planmeca Romexis®** ensuring an efficient workflow – no need for external applications or CDs and DVDs
- Automatic delivery of images and attachments
- Automatic notification to recipient of new cases

- Cases can be sent to any recipient who has an e-mail account
- Secure transfer and storage of information
- Streamline your communication with **Planmeca Romexis® Cloud**

Features

Sending images to recipient

- 2D images: panoramic, cephalometric, photos, intraoral X-ray images
- 3D images: CBCT, 3D photos, surface scans
- All annotations and other elements are included

Sending documents to recipient

- Attach one or more referrals, reports, or other documents

Versatile possibilities for communication

Recipients can download and view images at no cost using:

- Planmeca Romexis
- **Planmeca mRomexis™** imaging application for iOS and Android devices
- Free **Planmeca Romexis® Viewer**

.....
Planmeca Romexis® software and Planmeca Romexis® Cloud subscription are required for sending new cases.
 Visit <http://online.planmeca.com/> to subscribe and start sending images now.

Professionals proudly present the Planmeca ProMax® 3D family

Which one is right for you?

Planmeca ProMax® 3D s

Planmeca ProMax® 3D s is an ideal 3D unit for capturing small details. It is perfect for single implant, endodontic, and wisdom tooth cases.

Planmeca ProMax® 3D Classic

The **Planmeca ProMax® 3D Classic** imaging sensor covers the whole dentition area, so the unit gives a clear view of the mandible and maxilla.

Planmeca ProMax® 3D Plus

The newest member in our 3D family, **Planmeca ProMax® 3D Plus**, offers a wide variety of different volume sizes and is a great choice for any imaging need.

Planmeca ProMax® 3D Mid




Thanks to its wide volume size selection, **Planmeca ProMax® 3D Mid** handles a wide range of diagnostic tasks without compromising best practices.

Planmeca ProMax® 3D Max

Planmeca ProMax® 3D Max is a dedicated 3D imaging device that produces all required volume sizes when diagnosing the maxillofacial region – from the smallest special cases to images of the entire head.

Professionals proudly present

Planmeca ProMax® 3D s



Volume sizes

Ø50 x 80 mm
Ø50 x 50 mm
Stitched volume 90 x 60 x 80 mm



Long-term cooperation with Planmeca

"We purchased a **Planmeca ProMax® 3D s** for our dental clinic about four years ago or so. Before that, we had equipped our clinic with five Planmeca dental units, so it was only natural to continue the cooperation with Planmeca also on the X-ray side. Also, several radiologists recommended Planmeca's 3D units to us for their high quality.

We use the unit for implant cases, for lower third molar surgery, and for endodontic cases – particularly in difficult infection cases of teeth with multiple roots. Personally, I use the **Planmeca Romexis® 3D Implant Planning** module the most. It's very practical as I can virtually place the implants myself in the software.

The unit itself is very easy to use – our whole staff uses it, although mainly dentists take 3D images. Positioning is effortless and images are of high quality. And the unit's design is stylish and refined.

I would definitely recommend the unit to others. We have just taken the new sensor into use and I am very satisfied with the image quality. And the feedback from consulting radiologists has been good as well."

*Ari Mäkelä, Licentiate in Dentistry
Dental Care Center Janne, Järvenpää, Finland*



Chinese hospital chose Planmeca ProMax® 3D s

"I bought the **Planmeca ProMax® 3D s** system in September 2010, so I have been using it for over 2 years now. Factors influencing my decision were Planmeca's good reputation and quality-price ratio. For me, it is also important that everyday performance is excellent and if necessary, the after sales service works quickly.

I use my Planmeca 3D s system for various cases – for diagnosis in oral and maxillofacial surgery, for implantology, for diagnosis of periodontal and dental pulp diseases, and for orthodontics. The image quality is very clear, which makes diagnosis very easy with the excellent **Planmeca Romexis®** software.

In implant cases, Planmeca ProMax 3D s is very important for my preparation phase. The data I get from the image of the bone structure and thickness makes the operation easy and safe for the customer.

Planmeca ProMax 3D s really adds value to my work as I can perform many different kinds of tasks quickly and efficiently."


*Sun Zhizong, Dean
Donggang City Stomatology Hospital, Liaoning, China*

Professionals proudly present

Planmeca ProMax® 3D Classic


New extended volume size for the Teeth program

- Extended volume size that increases the diameter from Ø80 x 80 mm to Ø110 x 80 mm
- Captures a larger diagnostic area without increasing the patient dose
- Full dentition scans with the **Planmeca ProMax® 3D Classic** Teeth program
- Single rotation scans without stitching



Volume sizes

Ø80 x 80 mm
Ø80 x 50 mm
Ø50 x 80 mm
Ø50 x 50 mm
Extended volume Ø110 x 80 mm
Stitched volume 140 x 105 x 80 mm



*Dr Pekka Nissinen, GPD &
Dr Kim Lemberg, DDS,
PhD, Specialist in Oral and
Maxillofacial Radiology*

*West Vantaa Dental Clinic,
Finland*



Finnish dental clinic chooses Planmeca ProMax® 3D Classic

"We decided to purchase a **Planmeca ProMax® 3D Classic** 8x8 for our clinic as we wanted to start taking our own CBCT images and not have to send our patients elsewhere to have their 3D X-rays taken. In such cases, there is always the risk that the treatment process will suffer due to the patient's own lack of activity. Now we have our own radiologist and things have gone very smoothly. We also have two surgeons working with us, as we do a lot of implant treatments and treat also difficult endodontic cases."

Implant case acceptance has skyrocketed

"After acquiring the Planmeca ProMax 3D Classic, the amount of implant cases treated at our clinic has increased considerably. Patients are always amazed

when we offer to take their 3D images straight away. The unit is also especially suited to complicated endodontic cases, as you can notice everything in a 3D volume. It is also excellent for cases of wisdom teeth that have grown at a cumbersome angle.

The image quality produced by Planmeca ProMax 3D Classic is excellent. I think it is safe to say that we have the best 3D unit in Finland. This opinion is shared by our surgeons and many radiologists.

The **Planmeca Romexis®** software is a great working tool. It is logical, easy to use, and functions well – just a really good piece of software."

*Pekka Nissinen, GPD,
West Vantaa Dental Clinic, Finland*

Optimal image quality for every single field of dentistry




"I've been using Planmeca ProMax 3D Classic ever since its introduction to the market in 2007, and have used it for all imaging purposes. The image quality has proven to be reliable in every single field of dentistry, even in the most demanding imaging cases. The unit is very user-friendly, and all in all the imaging process can be carried out in an uncomplicated manner.

The Planmeca Romexis software is, in my opinion, the best software on the market when it comes to 3D imaging."

*Kim Lemberg, radiologist, West Vantaa
Dental Clinic, Finland*

Professionals proudly present

Planmeca ProMax® 3D Plus



Volume sizes

Ø160 x 90 mm
Ø160 x 50 mm
Ø90 x 90 mm
Ø90 x 50 mm
Ø70 x 70 mm
Ø70 x 50 mm
Ø40 x 70 mm
Ø40 x 50 mm



Dr. Dirk Ladig

*Oral surgery practice,
Hoyerswerda, Germany*



German oral surgery practice is impressed with the image quality of Planmeca ProMax® 3D Plus

"I have been using the **Planmeca ProMax® 3D Plus** unit in my oral surgery practice since 2013. Before that, I had good experience with Planmeca X-ray units. My panoramic X-ray unit ran smoothly for 19 years, the service was good and I was satisfied. Moreover, in 2000, I integrated cone beam computed tomography into my practice by adding a second unit. The decisive factor in purchasing the Planmeca ProMax 3D Plus unit was the radiographs of the new flat-panel devices shown to me by colleagues. The higher resolution of the images was very impressive! There was also a change in the physical layout of my practice. Instead of having two X-ray rooms, I wanted to have one. Planmeca ProMax 3D Plus combines two devices in one: OPG and CBCT. As a result, we need considerably less space.

More information in a single image

I use the device for different kinds of treatment planning; mainly implant cases, but also high-risk wisdom tooth surgery. In my view, a key benefit of the Planmeca ProMax 3D Plus is the possibility of displaying the entire mandible – including the ascending mandibular ramus and mandibular joint – in a single image. I also use the images for diagnosis of foreign body location, apical variances and inflammatory processes in the jaw area. CBCT provides much better diagnostic options for screening for infectious foci in patients with unclear symptoms or certain systemic diseases. Questions related to orthodontic treatments of impacted and displaced teeth, for example, can be easily solved on behalf of colleagues.

Low radiation exposure with adjustable volume sizes




What I really like about the unit is that I can select the volume according to the required image. The radiation exposure for patients is thus kept as low as possible. I use low-dose scans particularly with orthodontic diagnosis. The layer lights are especially useful when centring the image volume.

Operating and adjusting the unit is easy. What's more, the transition from analogue to digital control went well. Since the patients stand upright within the unit, positioning them is much easier than with the predecessor of the CBCT model (with patient bench), without having any problems with motion blur. The new device is also much more pleasant for the patients because there is no feeling of constriction."

*Dr. Dirk Ladig
Oral surgery practice,
Hoyerswerda, Germany*

Professionals proudly present

Planmeca ProMax® 3D Mid



Volume sizes

Ø200 x 170 mm
Ø200 x 100 mm
Ø200 x 60 mm
Ø160 x 170 mm
Ø160 x 100 mm
Ø160 x 60 mm
Ø100 x 100 mm
Ø100 x 60 mm
Ø80 x 80 mm
Ø80 x 50 mm
Ø40 x 80 mm
Ø40 x 50 mm



*Dr Carlo Pizzo, DDS &
Dr Gioia Amico, DDS*

A&P Clinic, Cittadella, Italy



Italian A&P Clinic opts for Planmeca ProMax® 3D Mid after a thorough market analysis

“In our new dental clinic, we have been using **Planmeca ProMax® 3D Mid** for six months now – and we are really satisfied with it.

We chose the unit after a thorough analysis of what the market was offering. We needed an imaging unit that could provide a wide range of FOV choices, the possibility to take panoramic images and cephalometric shots, and last but not least, software that could run natively on Mac OS, because our IT infrastructure was entirely built on Apple computers. The only unit that fulfilled all of these requirements was Planmeca ProMax 3D Mid.”

For every clinical application

“We love using it for taking panoramic images, preliminary treatment planning,

3D scans, wisdom teeth extractions and implant surgery. With **Planmeca Romexis®** – its dedicated software – we can virtually place the exact dental implants we are going to use by choosing them from the integrated 3D implant library. This feature works amazingly well.”

3D magic with the latest technology





“The machine and the software work seamlessly together: they are fast, reliable and easy to use. The 3D rendering is an incredibly powerful tool for us – for visualising the real bone morphology of the patients, and for the patients themselves to understand their clinical situation and the treatment we are offering them. So Planmeca Romexis can

become a really effective communication tool. For this reason, we adopted also the **Planmeca ProFace®** option. By superimposing a 3D scan of the patient's face and a CBCT X-ray image, we can show our clients an easy-to-understand image, in which they can really recognize themselves. Even today, this looks like magic for many of our patients!”

*Dr Carlo Pizzo & Dr Gioia Amico,
A&P Clinic, Cittadella, Italy*

Professionals proudly present

Planmeca ProMax® 3D Max



Volume sizes

Ø230 x 260 mm
Ø230 x 160 mm
Ø230 x 100 mm
Ø230 x 60 mm
Ø130 x 160 mm
Ø130 x 130 mm
Ø130 x 100 mm
Ø130 x 90 mm
Ø130 x 55 mm
Ø100 x 130 mm
Ø100 x 90 mm
Ø100 x 55 mm
Ø50 x 55 mm



Dr Corrado Gazzero

*MD, Specialist in Radiodiagnostics,
Qualified Expert in Radioprotection*

Studio Gazzero, Genoa, Italy



Radiologist praises the versatility of Planmeca ProMax® 3D Max

"I was the first **Planmeca ProMax® 3D Max** user in Italy and have been using it for about three years now. Before that, I used **Planmeca ProMax® 3D Classic 8x8** for 2 years. And I've been using Planmeca equipment since 1995 because of their image quality, their reliability, and the fast maintenance service.

I really enjoy working with Planmeca ProMax 3D Max. I have used it for every possible dental case, including all aspects of implantology, as well as endodontics, examining alterations of the bone structure, wisdom tooth extractions,

supernumerary teeth and more. In ENT cases, I have used the unit for the study of the paranasal sinuses and facial bone structures.

One of the most remarkable advantages is the possibility to choose the image quality and therefore to optimise the patient dose. The volume selection is complete, the imaging programs are easy to use and patient positioning is effortless."


Dr Gazzero, Studio Gazzero, Genoa, Italy

Patient specific implants

You imagine it. We make it.

***Planmeca ProModel™** offers patient-specific implants, surgical guides, and anatomical models – all individually designed for best possible surgical results.*

The implants are designed and manufactured to match any form, ensuring an exact fit to the patient's anatomy. The service also includes physical anatomical models and surgical guides for assisting in both pre-planning and the surgery itself.



The 3D design is created in an online meeting between surgeon and designer.



Anatomical model for surgical pre-planning and ready-to-use patient specific implant.

Planmeca ProModel™ service concept

- A unique service for creating patient specific implants, surgical guides and anatomical models from CBCT/CT images
- 3D implants are designed in an online session between the surgeon and Planmeca designer
- Ordering is quick and easy – from order to delivery in just a few business days
- Significantly lowers costs and reduces operation times by up to 4 hours
- Faster and more precise operations leading to better aesthetic results

Faster operations, precise fit and better aesthetic results

Stand out with colour

Complement the splendid design of your **Planmeca ProMax® 3D** X-ray unit by giving it a personal touch with your favourite colours. Select the perfectly matching shades from our exquisite and inspiring collection and create the looks of your dreams!



Technical specifications

Technical data

	3D s	3D Classic	3D Plus	3D Mid	3D Max
Anode voltage	60–90 kV	60–90 kV	60–90 kV	60–90 kV 60–120 kV	60–96 kV* 60–120 kV**
Anode current	1–14 mA	1–14 mA	1–14 mA	1–14 mA	1–12 mA
Focal spot	0.5 mm, fixed anode	0.5 mm, fixed anode	0.5 mm, fixed anode	0.5 mm, fixed anode	*0.6 mm, fixed anode **0.5 mm, fixed anode
Image detector	Flat panel	Flat panel	Flat panel	Flat panel	Flat panel
Image acquisition	Single 200 degree rotation	Single 200 degree rotation	200 / 360 degree rotation	200 / 360 degree rotation	210 / 360 degree rotation
Scan time	7.5–27 s	9–37 s	9–33 s	9–33 s	9–40
Typical reconstruction time	2–25 s	2–25 s	2–30 s	2–55 s	2–55 s

Comparison

	3D s	3D Classic	3D Plus	3D Mid	3D Max
3D dental programs	Yes	Yes	Yes	Yes	Yes
3D ENT programs	-	-	Yes	Yes	Yes
3D face photo	Yes	Yes	Yes	Yes	Yes
3D models scan	Yes	Yes	Yes	Yes	Yes
Suresmile certification	-	Yes	-	Yes	Yes
4D jaw motion	-	-	-	Yes	Yes
2D panoramic imaging	Yes	Yes	Yes	Yes	Yes
2D cephalometric imaging	Yes	Yes	Yes	Yes	-

Maximum volume sizes

	3D s	3D Classic	3D Plus	3D Mid	3D Max
Maximum volume without stitching	Ø50 x 80 mm	Ø80 x 80 mm	Ø160 x 90 mm	Ø200 x 100 mm	Ø230 x 160 mm
Extended volume without stitching		Ø110 x 80 mm			
Maximum volume with horizontal stitching	90 x 60 x 80 mm	140 x 105 x 80 mm			
Maximum volume with vertical stitching				Ø200 x 170 mm	Ø230 x 260 mm

Dental programs

Volume size (child mode) [mm]


	3D s	3D Classic	3D Plus	3D Mid	3D Max	Voxel size, isotropic
Tooth	Ø50 x 50 (Ø42 x 42)	Ø50 x 50 (Ø42 x 42)	Ø40 x 50 (Ø34 x 42)	Ø40 x 50 (Ø34 x 42)	Ø50 x 55 (Ø42 x 50)	75 µm*, 100 µm, 150 µm, 200 µm, 400 µm
	Ø50 x 80 (Ø42 x 68)	Ø50 x 80 (Ø42 x 68)	Ø40 x 70 (Ø34 x 60)	Ø40 x 80 (Ø34 x 68)		150 µm, 200 µm, 400 µm
Teeth		Ø80 x 80 (Ø68 x 68) Ø80 x 50 (Ø68 x 42) extended volume: Ø110 x 80	Ø70 x 50 (Ø60 x 42) Ø70 x 70 (Ø60 x 60) Ø90 x 50 (Ø75 x 42) Ø90 x 90 (Ø75 x 75)	Ø80 x 50 (Ø68 x 42) Ø80 x 80 (Ø68 x 68) Ø100 x 60 (Ø85 x 50) Ø100 x 100 (Ø85 x 85)	Ø100 x 55 (Ø85 x 50) Ø100 x 90 (Ø85 x 75)	150 µm, 200 µm, 400 µm
	triple scan: 90 x 60 x 80	triple scan: 140 x 105 x 80				200 µm, 400 µm
Jaw			Ø160 x 50 (Ø160 x 50) Ø160 x 90 (Ø160 x 90)	Ø160 x 60 Ø160 x 100 Ø200 x 60 (Ø200 x 60) Ø200 x 100 (Ø200 x 100)	Ø130 x 55 (Ø110 x 50) Ø130 x 90 (Ø110 x 75) Ø230 x 60 Ø230 x 100	200 µm, 400 µm, 600 µm
Face				Ø200 x 170 (Ø200 x 170)	Ø100 x 130 (Ø85 x 110) Ø130 x 130 (Ø110 x 110) Ø130 x 160 (Ø110 x 136)	200 µm, 400 µm
Skull					Ø230 x 160 Ø230 x 260	400 µm, 600 µm

ENT (Ear, Nose, Throat) programs


Volume size (child mode) [mm]

	3D Plus	3D Mid	3D Max	Voxel size, isotropic
Nose	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)	Ø100 x 90 (Ø85 x 75)	200 µm, 400 µm
Sinus	Ø90 x 90 (Ø90 x 90) Ø160 x 90 (Ø160 x 90)	Ø100 x 100 Ø100 x 170 Ø160 x 100 Ø160 x 170 Ø200 x 100 Ø200 x 170	Ø100 x 90 Ø100 x 130 Ø130 x 100 Ø130 x 130 Ø130 x 160	200 µm, 400 µm, 600 µm
Middle ear	Ø40 x 50 (Ø34 x 42)	Ø40 x 50 (Ø34 x 42)	Ø50 x 55 (Ø42 x 50)	75 µm*, 100 µm, 150 µm, 200 µm, 400 µm
	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)		150 µm, 200 µm, 400 µm
Temporal bone	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)	Ø100 x 90 (Ø85 x 75)	150 µm, 200 µm
Vertebrae	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)	Ø100 x 90 (Ø85 x 75) Ø100 x 130 (Ø85 x 110)	200 µm, 400 µm
Airways	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)	Ø100 x 90 (Ø85 x 75) Ø100 x 130 (Ø85 x 110) Ø130 x 130 (Ø110 x 110) Ø130 x 160 (Ø110 x 136)	200 µm, 400 µm

*Requires Endodontic imaging licence



Technical specifications



Dimensions

	3D s or 3D Classic	3D Plus or 3D Mid	3D Max
A	1298–2123 mm (51.1–83.5 in.)	1315–2095 mm (51.8–82.5 in.)	-
B	1560–2385 mm (61.4–93.8 in.)	1610–2390 mm (63.4–94.1 in.)	1582–2482 mm (62.3–97.7 in.)
C	1145 mm (45.1 in.)	1130 mm (44.6 in.)	-
D	850 mm (33.5 in.)	930 mm (36.6 in.)	930 mm (36.6 in.)
E	270 mm (10.6 in.)	247 mm (9.7 in.)	222 mm (8.7 in.)
F	698 mm (27.5 in.)	810 mm (32 in.)	788 mm (31 in.)
G	1250 mm (49.2 in.)	1366 mm (53.8 in.)	1351 mm (53.2 in.)
H	777 mm (30.6 in.)	756 mm (29.8 in.)	-
J	Ø820 mm (32.3 in.)	Ø1010 mm (39.8 in.)	Ø1010 mm (39.8 in.)

Physical space requirements

	3D s or 3D Classic	3D s or 3D Classic with cephalostat	3D Plus or 3D Mid	3D Plus or 3D Mid with cephalostat	3D Max
Width	115 cm (44 in.)	200 cm (79 in.)	118 cm (47 in.)	206 cm (82 in.)	116 cm (45.3 in.)
Depth	125 cm (49 in.)	125 cm (49 in.)	137 cm (54 in.)	137 cm (54 in.)	137 cm (54 in.)
Height*	153–243 cm (60–96 in.)	153–243 cm (60–96 in.)	161–239 cm (64–94 in.)	161–239 cm (64–94 in.)	161–239 cm (64–94 in.)
Weight	113 kg (lbs 248)	128 kg (lbs 282)	131 kg (lbs 289)	146 kg (lbs 322)	131 kg (lbs 289)

Minimum operational space requirements

	3D s or 3D Classic	3D s or 3D Classic with cephalostat	3D Plus or 3D Mid	3D Plus or 3D Mid with cephalostat	3D Max
Width	150 cm (59 in.)	215 cm (85 in.)	158 cm (63 in.)	225 cm (89 in.)	158 cm (63 in.)
Depth	163 cm (64 in.)	163 cm (64 in.)	175 cm (69 in.)	175 cm (69 in.)	175 cm (69 in.)
Height*	243 cm (96 in.)	243 cm (96 in.)	239 cm (94 in.)	239 cm (94 in.)	239 cm (94 in.)

*The maximum height of the unit can be adjusted for offices with limited ceiling space.

Example installation

Included in delivery	Planmeca ProMax 3D unit with 3D reconstruction server	
Minimum set up	Client workstation and database server <ul style="list-style-type: none">Planmeca Romexis 3D ExplorerDatabase serverPlanmeca Romexis Image Database The client workstation and database server can also be in separate computers.	
Additional equipment	Additional diagnostic workstations with different software configurations Planmeca Romexis tools: <ul style="list-style-type: none">3D Explorer3D Cross Sections module3D TMJ module3D Implant Planning moduleDICOM module	

Planmeca Romexis® imaging software

Supported 2D modalities	Intraoral Panoramic Cephalometric 2D linear tomography Photos Stack images (CBCT slices and panoramic slices)
Supported 3D modalities	3D CBCT 3D photo 3D surface scan
Supported photo sources	Intraoral camera Digital camera or scanner (import or TWAIN capture)
Operating systems	Win 7 Pro (64 bit) / Win 8.1 Pro (64 bit) / Win 10 Pro (64 bit) Win 2008 Server / Win 2012 Server Mac* (OS X or newer) For detailed information please see system requirements of Planmeca Romexis www.planmeca.com *Cephalometric Analysis module, 3D Ortho Studio module and Planmeca PlanCAD Easy are supported on Windows operating systems.
Image formats	JPEG or TIFF (2D image) DICOM (2D and 3D image) STL (3D image) TIFF, JPEG, PNG, BMP (import/export)
Image size	2D X-ray image: 1–9 MB 3D X-ray image: typically 50 MB–1 GB
Installation options	Client–Server Java Web Start deployment
DICOM 3.0 support	DICOM Import/Export DICOM DIR Media Storage DICOM Print SCU DICOM Storage SCU DICOM Worklist SCU DICOM Query/Retrieve DICOM Storage Commitment DICOM MPPS
Interfaces	TWAIN Client PMBridge (patient information and images) VDDS (patient information and images) InfoCarrier (patient information) Datagate (patient and user information)
3 rd party software integrations	Dolphin Imaging Nobel Clinician Materialise Dental Simplant Straumann coDiagnostiX Cybermed N-Liten

Find all the latest Planmeca news




www.facebook.com/PlanmecaOy



www.planmeca.com/newsroom

One software for all.



Planmeca Oy designs and manufactures a full line of industry-leading dental equipment, including 3D and 2D imaging devices, CAD/CAM solutions, dental care units and software. Planmeca Oy, the parent company of the Finnish Planmeca Group, is strongly committed to better care through innovation, and it is the largest privately held company in the field.

PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | fax +358 20 7795 555 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations.
Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, DentoVac, Digital perfection, Economat Plus, Elegant, Flexy, Mini-dent, Perio Fresh, PlanEasyMill, Planmeca 4D, Planmeca AINO, Planmeca ARA, Planmeca CAD/CAM, Planmeca CALM, Planmeca Chair, Planmeca Clarify, Planmeca Compact, Planmeca Creo, Planmeca Emerald, Planmeca FIT, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Lumo, Planmeca Maximity, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca mRomexis, Planmeca Noma, Planmeca Olo, Planmeca Online, Planmeca PlanCAD, Planmeca PlanCAM, Planmeca PlanClear, Planmeca PlanID, Planmeca PlanMill, Planmeca Planosil, Planmeca PlanPure, Planmeca PlanScan, Planmeca PlanView, Planmeca ProCeph, Planmeca ProFace, Planmeca ProlD, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca Serenus, Planmeca SingLED, Planmeca SmartGUI, Planmeca Solanna, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Viso, Planmeca Verity, Planmeca Waterline Cleaning System, Planmeca Xtremity, Proline Dental Stool, ProTouch, Saddle Stool, SmartPan, SmartTouch, Trendy, and Ultra Relax are registered or non-registered trademarks of Planmeca in various countries.

10032799/0518/en

