

# 2020 FALL PROMOTIONS

## EQUIPMENTSAVINGS

PAGES 31-43

## PROMOTIONS, REBATES AND SAVINGS

from A-dec, Midmark and Planmeca



Visit our website -  
[goetzedental.com/promotions](http://goetzedental.com/promotions)  
to see all of our current equipment promotions



## ADVANCEDENTAL

PAGES 03-22

### NOW AVAILABLE!

Advance Dental Products by Goetze Dental



KANSAS CITY

ST. LOUIS

DENVER

DES MOINES

OMAHA

SPRINGFIELD

WICHITA

since

**GOETZE**  
d e n t a l

# FALL PROMOTIONS

OCTOBER 01 - DECEMBER 31

# FOUR EASY WAYS TO ORDER

## CONTENTS

- 03-22 | ADVANCE DENTAL PRODUCTS
- 31-43 | EQUIPMENT SAVINGS

## VENDOR INDEX

- 23 | A-dec
- 38-42 | A-dec
- 24 | Crest & Oral-B
- 25 | Dentsply Sirona
- 34-35 | Dental EZ
- 28 | DMG
- 28 | GC America
- 26 | HuFriedyGroup
- 31-33 | Midmark
- 29 | Parkell
- 43 | Planmeca
- 30 | RevenueWell/EDDSON
- 27 | SDI
- 36-37 | UV Angel

### PLEASE NOTE:

Goetze Dental is not responsible for any typographical errors in the price or descriptions of items featured in this flyer. Prices subject to change due to manufacturers' price increases. Due to the COVID-19 pandemic, orders for all PPE and infection control products can not be cancelled, changed or returned. Visit [goetzedental.com/your-resources](http://goetzedental.com/your-resources) for the most up-to-date pricing.

**GAINSpeed™**  
desktop ordering application



**GoetzeOnline.com**  
e-commerce ordering option



**800.692.0804**  
8am - 5pm CST



Fax:  
**800.552.7791**





# Advance Advancing the Quality of Dental Care Products



## Topical Anesthetic Gel

Contains 20% Benzocaine for maximum strength pain relief while fast working to quickly relieve patient discomfort. Available in six great tasting flavors.

### Bottle 1oz

Cherry .....	<b>575018</b>
Mint .....	<b>575019</b>
Piña Colada .....	<b>575020</b>
Bubble Gum .....	<b>575021</b>
Strawberry .....	<b>575022</b>
Raspberry .....	<b>575023</b>
	<b>Ea \$6.75</b>



## Lidocaine HCl 2%

with Epinephrine 1:100,000

### Features:

- No-methylparaben formulation to avoid allergic reactions on sensitive patients.
- 100 percent latex free components to help reduce allergic responses.
- Terminal sterilization of product/cartridges for safe use.

Box/50 .....	<b>577673</b>
	<b>Ea \$41.87</b>



## Premium Non-Woven Sponges

- Unique rayon/polyester blend
- Provides maximum absorbency
- Less adhesion to wounds for faster healing
- Soft, Smooth, Strong, Virtually lint free
- Ideal for sustained patient use, disinfecting, general cleaning
- 4 ply non-sterile

2" x 2"	
Case/5000 .....	<b>097805</b>
	<b>Ea \$75.00</b>

4" x 4"	
Case/2000 .....	<b>097837</b>
	<b>Ea \$90.66</b>



- BLUE
- DUSTY ROSE
- SILVER GREY
- WHITE

## Plastic Cups

Horizontal embossed grip for easier dispensing and holding. Double rolled lip for easy separation. Can be used with any 5oz dispenser.

<b>5oz Cups Case/1000</b>	
Blue .....	<b>191885</b>
Dusty Rose .....	<b>191880</b>
Silver Gray .....	<b>191887</b>
White .....	<b>191880</b>
	<b>Ea \$28.45</b>

**GOETZE DENTAL** is not responsible for any typographical errors in the price or descriptions of items in this flyer. Prices subject to change due to manufacturer price increases and market fluctuations. Due to the COVID-19 pandemic, orders for all PPE and infection control products can not be guaranteed for delivery and are subject to backorder status, allocation of amounts and/or cancellation.

## Saliva Ejectors


Provides optimal suction without aspirating tissue using a unique wire-reinforced tube, which maintains any desired configuration. Latex Free

**Pkg/100**

White / White tip.....**197510**

Clear / Blue Tip.....**197508**

**Ea \$4.68**



## Air/Water Syringe Tips

- Rounded edge for maximum patient comfort
- Convenient dispenser bottle keeps work area clear and assists in preventing cross-contamination
- Rounded non-scratch end to minimize wear and tear of O-ring in hand piece

**Bag/250** .....**576916**

**Ea \$27.64**

## Surgical Aspirator Tips

- Designed to fit all high-volume evacuation systems
- 30 degree angle for maximum patient comfort
- Made of rigid opaque plastic



**Bag/25**

1/8" White ..... **576914**

1/4" Green ..... **576915**

**Ea \$5.99**

## Combo Tip Evacuators (HVE)

Double-ended design with both a vented and non-vented end on the same tip. Latex Free



**Pkg/100** .....**198174**

**Ea \$7.26**



Both Econoback and Polyback Towels come in a wide variety of colors!

## Econoback Towels 3-Ply

Highly absorbent tissue specially treated to enhance wet strength and tear resistance, and a bonded layer of poly film. Reinforced edge to provide strength for Towel clips.

### 2-ply tissue with 1-ply poly (3 ply bibs)

19" x 13" - Case/500

Aqua .....	<b>198252</b>
Beige .....	<b>198257</b>
Blue .....	<b>198260</b>
Dusty Rose .....	<b>198266</b>
Green .....	<b>198254</b>
Lavender .....	<b>198263</b>
Peach .....	<b>198264</b>
Silver Gray.....	<b>198268</b>
White .....	<b>198255</b>
Yellow .....	<b>198261</b>

**Ea \$30.71**

## Polyback Towels 4-Ply

Highly absorbent tissue specially treated to enhance wet strength and tear resistance, and a bonded later of poly film. Reinforced edge to provide strength for Towel clips.

### 3-Ply tissue with 1-ply poly (4 ply bibs)


19" x 13" - Case/500

Aqua .....	<b>198269</b>
Beige .....	<b>198270</b>
Blue .....	<b>198275</b>
Dusty Rose .....	<b>198292</b>
Green .....	<b>198280</b>
Lavender .....	<b>198282</b>
Peach .....	<b>198288</b>
Silver Gray.....	<b>198294</b>
White .....	<b>198285</b>
Yellow .....	<b>198290</b>

**Ea \$34.39**



**HEAVY weight paper!**



## Bracket Tray Covers

Heavy weight paper impedes the flow of moisture, minimizing liquid seep-through, protecting the covers and unit beneath the tray.

### Ritter- "Size B" (8.5" x 12.25")

**Box/1000**

Blue .....	<b>198350</b>
Dusty Rose .....	<b>198347</b>
Green .....	<b>198422</b>
Silver .....	<b>198358</b>
White .....	<b>198360</b>

**Ea \$29.32**



### General Purpose Cleaner

Non-ammoniated concentrated liquid formula employs surfactants and anti-corrosive additives to safely and thoroughly clean instruments prior to sterilization. Phosphate-free.

Gallon ..... **575655**  
Ea **\$19.75**

### Tartar & Stain Remover


Powerful, ready-to-use ultrasonic solution to remove stains from dentures, partials, appliances. Also removes temporary cement.

Gallon ..... **575656**  
Ea **\$19.75**

### Temporary Cement Remover

Strong ready-to-use ultrasonic solution to remove temporary cement.

Gallon ..... **575657**  
Ea **\$19.75**



### Self-Sealing Autoclave Pouches

Autoclavable self-sealing, blue tinted, see-through pouches. Easy-to-read color changing arrows are indicators for chemical vapor, steam or gas sterilization. Utilizes a triple seal plus a fourth "Secure Seal" at the chevron.

<b>Box/200</b>	
2.25" x 4" .....	<b>577685</b> Ea <b>\$9.24</b>
2.75" x 9" .....	<b>575073</b> Ea <b>\$11.20</b>
3.5" x 5.25" .....	<b>575074</b> Ea <b>\$11.72</b>
3.5" x 9" .....	<b>575075</b> Ea <b>\$11.14</b>
5.25" x 10" .....	<b>575077</b> Ea <b>\$16.14</b>

### Economy - Box/500

3.5" x 9" .....	<b>575076</b> Ea <b>\$27.17</b>
-----------------	------------------------------------







Easy removal  
with no residue  
left behind!



### Light Handle Cover

Specially designed barriers effectively cover equipment with near custom fit.

#### Box/500 T-style

4" x 5 3/4" ..... **575327**  
Ea **\$13.75**

### Syringe Sleeve

Syringe sleeves have pre-cut openings for easy application. Made of clear poly material that protects instruments from cross-contamination.

#### Box/500

2 1/2" x 10" ..... **575326**  
Ea **\$9.95**

### X-Ray Cover

Specially designed barriers effectively cover equipment with near custom fit.

#### Box/250

15" x 26" ..... **575328**  
Ea **\$21.95**

### Chair Cover

Tapered design conforms neatly to patient chairs, covering chair back.

#### Box/225

27 1/2" x 24" ..... **575324**  
Ea **\$23.48**

### Tray Sleeve

Back-flip closure ensures tray is completely covered. Made of a clear poly material that protects tray from cross-contamination.

#### Box/500

10 1/2" x 14" ..... **575325**  
Ea **\$16.95**

### Universal Cover Film

4" x 6" perforated sheets with non-adhesive side edges for easy removal. Specially formulated adhesive leaves no residue behind. Helps to prevent cross-contamination. Convenient dispenser box contains 1 roll of 1200 perforated 4" x 6" sheets

#### Roll/1200

Clear ..... **575330**  
Blue ..... **575329**  
Ea **\$14.95**

### Temporary C & B Material

10:1 Ratio

Syringeable bis-acryl composite for chairside provisional restorations. Available in 3 shades: A1, A2, A3. All are delivered in 76g (50mL) cartridges designed to fit on a 10:1 style automix gun. Ten waste saver tips are included per kit. 50mL Cartridge with 10 tips

A1 ..... **575295**  
A2 ..... **575296**  
A3 ..... **575297**  
Ea **\$83.07**



### Denture Boxes

- One-piece "hinged" lid provides easy access to contents inside
- Sturdy, colorful plastic helps protect dentures from damage
- Deep soaking well fits both uppers and lowers at the same time
- Flat bottom provides area for writing "Return in found" information
- Assorted colors 3 3/4" x 3" x 2 3/8"

Pkg/12 ..... **577907**  
Ea **\$13.19**

### Retainer Boxes

- Clamshell design with easy to open/close "snap-lock"
- Sturdy, colorful plastic helps protect retainer from dust and damage
- Vented top and bottom holes allow box to "breathe"
- Flat bottom provides area for writing "Return if found" information
- Assorted colors 3" x 3 1/8" x 1 1/8"

Pkg/12 ..... **577908**  
Ea **\$12.12**







Ideal For:

- Occlusal registration
  - Copying transfers
- Primary impressions with bite trays



## High Performance Bite Registration Material

Formulated to have a light, mousse-like consistency ideal for occlusal registration, copying transfers and primary impressions with bite trays. Our material has excellent durometer or hardness preventing distortion in disposable trays, specifically the "triple trays". Unique our Bite Registration approaches zero resistance allowing for unmodified natural movement of the jaw. The non-flowing, thixotropic formula stays in place and will not drip, run, or sag.

### Bite Registration Kit w/Tips - 50ml

Each pack contains (2) 50ml cartridges & (6) 6.5mm mixing tips

#### Regular Set

Work Time 2:00 Intraoral Set 1:00

Unflavored..... 575100  
Ea **\$26.20**

#### Fast Set

Work Time 0:30 Intraoral Set 1:00

Unflavored..... 575101  
Ea **\$26.20**

Vanilla..... 575102  
Ea **\$26.20**

#### Super Fast Set

Work Time 0:30 Intraoral Set 0:30

Unflavored..... 575099  
Ea **\$26.20**

Vanilla..... 575103  
Ea **\$26.20**

#### Rigid Bite - Fast Set

Unflavored..... 575104  
Ea **\$26.20**

#### Clear Bite - Fast Set

Unflavored..... 575105  
Ea **\$27.95**

## Multi-Color Disposable Impression Trays

Perforated

- Disposable premium impression trays made of high quality, rigid, transparent plastic
- Strategically perforated trays allow easy modification with dental bur to accommodate implants
- Can be heated and molded for a custom fit
- Bright color-coded body for quick and easy identification during procedure
- Easy to apply I.D. labels included in every bag
- Non-sterile

### Large Perforated Yellow - Bag/12

#1 - Upper..... 575279  
#2 - Lower..... 575280

### Medium Perforated Blue - Bag/12

#3 - Upper..... 575281  
#4 - Lower..... 575282

### Small Perforated Orange - Bag/12

#5 - Upper..... 575283  
#6 - Lower..... 575284

### Quadrant Perforated Green - Bag/12

#7 - UL/LR ..... 575285  
#8 - UR/LL ..... 575286

### Anterior Perforated Pink - Bag/12

#9 - Upper..... 575287  
#10 - Lower..... 575288

### Assorted Perforated

Bag/12..... 575289  
Ea **\$4.25**

## Premium 3 Way Impression Bite Trays

Color Coded

- Disposable Advance bite trays precisely capture the master impression, counter impression, and bite registration in one procedure
- Rigid frame with nylon webbing to ensure accurate bite
- Ergonomically designed thumb grip aids in placement and removal
- Patent Pending 1/4" D grooves provide superior material retention

### Box/12

Posterior (Blue)..... 575277  
Sideless (White) ..... 575275  
Anterior (Purple)..... 575276  
Quadrant (Yellow)..... 575274  
Assorted Sizes ..... 575278  
Full Arch (Pink) ..... 575273  
Ea **\$7.62**



Peppermint Flavor



## Pulp Vitality

Refrigerant Spray

For use of thermo testing tooth vitality.

Peppermint Flavor

6oz Spray..... 575267  
Ea **\$17.90**





## Disposable Prophyl Angles

Latex Free

Straight-body style in pink with a webbed cup that is soft and latex-free. This angle features a screw-type cup, allowing for stability. Fits comfortably and securely in the hand for hard-to-reach areas, while maintaining a favorable paste-holding capacity.

### Soft Gray Cup

Box/100 ..... 575078  
Ea **\$45.89**

Box/500 ..... 575079  
Ea **\$188.69**

### Firm White Cup

Box/100 ..... 575080  
Ea **\$45.89**

Box/500 ..... 575081  
Ea **\$188.69**

Contact your  
Goetze Dental  
Representative  
to order!  
Or Order at:  
**800.692.0804**



## Pre-bent Multipurpose Dispensing Tips

Dispensing Tips are pre-bent to the proper angle and have a blunt end for patient comfort. Great for placing sealants, etchants, cements and other flowable materials. Secure Luer lock, universal color coding for easy identification



### Tips - Pkg/100

18ga - Pink ..... 575290  
(Cavity liners, Bases & Composites)

20ga - Yellow ..... 575291  
(Composites & Sealants)

20ga - Black ..... 575292  
(Composites & Cements)

22ga - Gray ..... 575293  
(Sealants, Etchants & Cements)

25ga - Blue ..... 575294  
(Etchants)

Ea **\$8.12**

## Etching Gel 37% Phosphoric

Formulated for maximum working and handling characteristics.



### Syringes (1.2ml) with 12 Tips

Pkg/4 ..... 577676  
Ea **\$11.84**

### Syringes (1.2ml) with 25 Tips

Economy Pkg/12 ..... 577677  
Ea **\$30.00**

### Syringes (25ml) with 50 Tips

Jumbo Pkg/2 ..... 577679  
Ea **\$36.50**

### Pre-bent Tips (25ga) - Blue

Pkg/20 ..... 577681  
Ea **\$8.50**

### 3cc Empty Refillable Syringes

Pkg/5 ..... 577678  
Ea **\$12.50**



## Fluoride Trays

Fluoride trays are color-coded for easy identification and are disposable for your convenience. The Single Arch trays are made of soft, light foam for patient comfort and have a convenient handle for easy placement and removal. The Dual Arch trays have a fill line that helps reduce gel overflow and have locking handles for easy placement and removal.

### Single Arch - Pkg/100

Small ..... 575028  
Medium ..... 575029  
Large ..... 575030  
Ea **\$39.99**

### Dual Arch - Pkg/100

Pedo ..... 575024  
Small ..... 575025  
Medium ..... 575026  
Large ..... 575027  
Ea **\$39.99**



## Propyl Paste

### Mint w/Fluoride - Box/200

X-Coarse ..... 575114  
Coarse ..... 575113  
Medium ..... 575107

### BubbleGum w/Fluoride - Box/200

Coarse ..... 575115  
Medium ..... 575108  
Fine ..... 575112

### Raspberry w/Fluoride - Box/200

Coarse ..... 575117  
Medium ..... 575116

### Assorted w/Fluoride - Box/200

Coarse ..... 575110  
Medium ..... 575109

Ea **\$24.99**



**Advance Technical  
Adult 30 Tuft**

Assorted: 4 assorted colors  
Box/72 ..... 578066  
Ea **\$48.00**



**Compact Super Row Bristles  
Adult 35 Interdental Cut Tuft with  
Super Grip**

Assorted: Dark Purple, Dark Blue, Teal, Orange  
Box/72 ..... 578067  
Ea **\$34.35**



**FlexBrush Clean  
Adult 40 Tuft Interdental Cut**

Assorted: Dark Blue, Red, Dark Green, Yellow  
Box/72 ..... 578062  
Ea **\$32.85**



**Super Tip Bristle  
Adult 38 Tapered Tuft**

Assorted: Yellow, Red, Dark Blue, Dark Green  
Box/72 ..... 578063  
Ea **\$34.58**



**Advance Compact Micro Tip  
Adult 33 Tuft**

Assorted: Blue, Green, Pink, Purple  
Box/72 ..... 578068  
Ea **\$49.71**



**Dual Action Ergo Grip  
Gum Stimulator and Tongue Cleaner  
Adult 36 Tuft**

Assorted: Dark Pink, Bright Blue, Vibrant Green, Orange  
Box/72 ..... 578069  
Ea **\$34.35**



**Clear Grip with Indicator  
Adult 38 Tuft**

Assorted: Blue, Dark Green, Orange, Bright Purple  
Box/72 ..... 578057  
Ea **\$50.76**



**Compact Opaque Straight Handle  
Adult 35 Soft Rounded Tuft**

Assorted: Purple, Dark Blue, Dark Green, Orange  
Box/72 ..... 578065  
Ea **\$34.35**



**Triple Action Sport Grip  
Adult 32 Tuft Interdental with  
Tongue Cleaner and Power Tip**

Assorted: Green, Purple, Red, Blue  
Box/72 ..... 578058  
Ea **\$39.75**

**YOUTH BRUSHES**



**Ladybug Grip**  
Infant 24 Soft Tuft

Assorted: 4 assorted colors  
Box/72 ..... 578061  
Ea **\$31.15**



**Advance Color Youth**  
Age 8+ 32 Tuft

Assorted: Two Tone colors  
Box/72 ..... 578060  
Ea **\$34.45**



**Child's Sparkle Ergo Grip**  
Ages 5-7 27 Rounded Tuft

Assorted: 5 assorted two-tone colors with sparkles  
Box/72 ..... 578059  
Ea **\$33.45**







## Sutures

Super sharp siliconized needles for repeated passage through tissue. Needles are drilled-end to ensure smooth transition from the needle to the material and maximizes the needle attachment strength. Packaged in a double sterile pouch for sterile transfer, in color-coded suture boxes with a front dispenser. Packaged

12 per box

Silk = Light Blue  
 Chromic Gut = Beige  
 Plain Gut = Yellow

### 18" Plain Gut Sutures with Needle

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00071 591B	Absorbable	Plain Gut	3/8 circle -19mm reverse cut	4/0	18"	575058
	C-31	AD00065 594B	Absorbable	Plain Gut	1/2 circle -24mm reverse cut	3/0	18"	575056

Ea \$32.65

### 18" Chromic Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00558 558B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	4/0	18"	575061

Ea \$40.54

### 18" Silk Braided Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00066 781B	Non - Absorbable	Silk - Black Braided	3/8 circle -19mm reverse cut	4/0	18"	575057
	C-6	AD00072 784B	Non - Absorbable	Silk - Black Braided	3/8 circle -19mm reverse cut	3/0	18"	575059
	C26	AD00062 778B	Non - Absorbable	Silk - Black Braided	1/2 circle -15mm reverse cut	3/0	18"	575053
	C-31	AD00064 785B	Non - Absorbable	Silk - Black Braided	1/2 circle -24mm reverse cut	3/0	18"	575055

Ea \$27.64

### 27" Plain Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00822 554B	Absorbable	Plain Gut	3/8 circle -19mm reverse cut	3/0	27"	575065

Ea \$33.19

### 27" Chromic Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00636 559B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	3/0	27"	575064
	C-6	AD00635 560B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	4/0	27"	575063

Ea \$33.19

Contact your  
 Goetze Dental  
 Representative  
 to order!  
 Or Order at:  
**800.692.0804**



Standard Operatory: **FG** - Friction Grip **RA** - Right Angle **SS** - Short Shank **HP** - Hand Piece

	Round					Pear				Inverted Cone			Straight/Flat End Plain		Straight/Flat End Cross Cut		
SIZE	1/2	2	4	6	8	2P	4P	330	332	33 1/2	35	37	56	57	556	557	558
FG 10PK		575222	575223	575224	575225			575221			575219						575217
FG 100PK		575165						575159			575162		575163	575164	575160	575161	
RA 10PK		575178	575177	575176	575175												
RA 100PK		575154															
FG SS 10PK								575227									575226
FG SS 100PK								575215									

	Straight Round End Plain	Straight/Round End Cross Cut	Taper Flat End Plain	Taper Flat End Cross Cut	Amalgam Prep	Metal Cutter					Surgical Length						
SIZE	1157	1557	1558	169	701	245	2	330	557	1557	2	4	6	557	558	701	702
FG 5PK											575211	575212	575213	575214		575208	575210
FG 10PK		575155			575218	575220	575173		575174								
FG 100PK						575158	575172	575171		575170							
FG SS10PK																	
FG SS 100PK																	

	Trimming and Finishing																
SIZE	7006	7008	7204	7379	7404	7406	7408	7642	7664	7802	7803	7901	7902	7903	8404 20 Blade	9006 30 Blade	9008 30 Blade
FG 5PK	575202				575204	575205	575200			575199	575198	575197	575196				

	Trimming and Finishing			
SIZE	9406	9903	SE 9-10 10 Blade	SE 9-20 20 Blade
FG 5PK	575194	575195		

### Trimming & Finishing FG Burs

Twelve Bladed Burs - Pkg/5  
Ea **\$28.95**

### Trimming & Finishing FG Burs

Twelve Bladed Burs - Pkg/5  
Ea **\$21.50**

### FG Carbide Burs

Economy - Pkg/100  
Ea **\$104.99**

FG Carbide - Pkg/10  
Ea **\$13.75**

### Metal Cutting Carbide Burs

Economy - Pkg/100  
Ea **\$299.48**

Pkg/10  
Ea **\$30.77**

### RA Carbide Burs

Economy - Pkg/100  
Ea **\$124.98**

Pkg/10  
Ea **\$19.50**

### Surgical Length FG Carbide Burs

Pkg/5  
Ea **\$9.25**



## SonicTouch

Rechargeable Sonic Power Toothbrush



**Advanced Sonic Cleaning Technology**  
Produces 40,000 brush strokes per minute for a powerful dynamic cleaning.



**2 Minute Auto-Timer**  
Built-in timer to ensure dentist recommended brushing time.




**Base Charging Indicator**  
Built-in interval timer helps you brush each area of your mouth evenly.



**Waterproof Design**  
Protection from direct contact from water, spills or splashes.



**30 Second Interval**  
Built-in interval timer helps you brush each area of your mouth evenly.



**4 Brushing Modes:**  
Clean, Soft, Whiten, Massage  
Choose from Clean, Soft, Whiten, or Massage modes to cater all of your oral care needs.

### SonicTouch

Rechargeable Sonic Power Toothbrush

Advance SonicTouch technology provides a powerful dynamic cleaning action that helps to keep teeth and gums healthy by removing plaque.

This technology provides effective stain removal which helps to brighten teeth and restore them to their natural whiteness.

- 2 Minute Auto-Timer, 40,000 brush strokes per minute
- 30 second interval timer helps you brush each area evenly
- Base charging indicator
- Waterproof design
- 4 brushing modes: clean, soft, whiten, massage
- 3 brush heads
- 1 year warranty

SonicTouch Toothbrush ..... **578048**  
**Ea \$45.83**

### SonicTouch Refill Brush Heads

The rubberized brush head is designed to protect tooth enamel. The brush conforms to the contours of your teeth for effective, deep cleaning. High quality DuPont bristles promote naturally whiter teeth by gently polishing off surface stains, end rounded.



Pkg/2 ..... **578049**  
**Ea \$7.36**

### Advance Thick Flexible Interdental Picks

Cylinder interproximal picks. Perfect for cleaning between teeth and under orthodontic work.



Pack/72 (Thick) ..... **578185**  
**Ea \$53.08**

### Advance Thin Flexible Interdental Picks

Cylinder interproximal picks. Perfect for cleaning between teeth and under orthodontic work.



Pack/72 (Thin) ..... **578184**  
**Ea \$53.08**



### Advance Varnish

5% Sodium Fluoride

Easy non-messy single dose delivery system. Tooth-shaded color. Contains Xylitol. Tasteful and transparent 5% sodium fluoride.

#### Box/50

0.4ml unidose (Toffee) ..... **578143**  
0.4ml unidose (Mint) ..... **578144**  
**Ea \$50.42**



**Endodontic Irrigating Tips**

Endodontic irrigation needles used to irrigate canals, fistulas and pockets. Luer lock hub

**Safe End/Side Port - 100/bag**

23 Gauge .....	<b>579271</b>
27 Gauge .....	<b>579272</b>
30 Gauge .....	<b>579273</b>
	<b>Ea \$28.75</b>

**Slotted End/Side Vent - 100/bag**

23 Gauge .....	<b>579274</b>
27 Gauge .....	<b>579275</b>
30 Gauge .....	<b>579276</b>
	<b>Ea \$44.85</b>

**Endodontic Syringe with Needle**

- Used for the irrigation of pockets, fistulas and canals
- 3cc Luer syringe
- Slotted End/ Side Vent

**Endodontic Syringe with Needle - 100/bag**

23 Gauge .....	<b>579277</b>
27 Gauge .....	<b>579278</b>
30 Gauge .....	<b>579279</b>
	<b>Ea \$35.80</b>

**Viscous Brush Tips**

- Easy extrusion of viscous materials
- Photo sensitive opaque tip maintains the flow of light-cured material
- Luer lock hub



Pack/100 .....	<b>579288</b>
	<b>Ea \$52.19</b>

**Pre-bent Brush Tips**

- Curved cannula for optimal visibility
- Precision placement of hemostatic agents
- Padded-end brush tip controls flow
- Universal luer lock hub



19 Gauge - Pack/20.....	<b>579285</b>
	<b>Ea \$12.82</b>
19 Gauge - Pack/100.....	<b>579286</b>
	<b>Ea \$50.83</b>
19 Gauge - Pack/500.....	<b>579287</b>
	<b>Ea \$178.48</b>

**Low Viscosity Brush Tips**

- Curved cannula for optimal entry
- Opaque luer lock hub
- Nylon bristle end enables controlled delivery of solution



19 Gauge - Pack/100	<b>579289</b>
	<b>Ea \$42.16</b>

**Black Mini Brush Tips**

- Easy extrusion of viscous materials
- Photo sensitive opaque tip with luer lock hub



18 Gauge - Pack/100.....	<b>579291</b>
	<b>Ea \$64.86</b>


**Calcium Hydroxide**

- 12.5pH aqueous solution for temporary filling and cleansing root canals
- Contains no more than 35% Calcium Hydroxide paste

**Standard Kit**  
Contains: 4-1.2mL syringes with tips..... **579269**  
**Ea \$32.62**

**Value Pack**  
Contains: 20-1.2mL syringes with tips ..... **579270**  
**Ea \$75.51**

Contact your Goetze Dental Representative to order!  
Or Order at:  
**800.692.0804**



**Advance Interdental Brushes**

Great for removal of food and plaque on bridges, crowns, braces and between teeth.

- Wide Cone - Blue Grip**  
Box 72.....**579225**
- Moderate Cylinder - Red Grip**  
Box 72.....**579224**
- Ultra Tight - Yellow Grip**  
Box 72.....**579222**
- Tight Cylinder - Green Grip**  
Box 72.....**579223**  
**Ea \$54.37**



**Prepasted Disposable Adult 39 Tuft (Mint Flavor)**

Assorted: Yellow, Red, Green, Blue  
Box/100 .....**579708**  
**Ea \$16.30**



**Denture Brush**

Adult  
Assorted: Red, Green, Purple, Blue  
Pkg/12 .....**579496**  
**Ea \$7.95**



**Premium Disposable Impression Trays**

Disposable impression tray made from rigid, high-impact FDA approved styrene. Advance Premium disposable impression trays assist in generating a consistent and unsurpassed impression free of distortion or bite through. Smooth, rounded edges, interior retention slots, and vertical perforation enhance material retention.

- Perforated, 12/bag**
- #1 Large Upper.....**579517**
  - #2 Large Lower.....**579520**
  - #3 Medium Upper.....**579521**
  - #4 Medium Lower.....**579523**
  - #5 Small Upper.....**579524**
  - #6 Small Lower.....**579525**
  - #7S Upper Right/Upper Left.....**579532**
  - #8s Upper Left/Upper Right.....**579533**
  - #9S Lower Anterior .....**579534**
  - #10S Upper Anterior.....**579535**  
**Bag \$9.50**

- Solid, 12/bag**
- #7 Upper Right/Lower Left.....**579526**
  - #8 Upper Left/Upper Right.....**579528**
  - #9 Lower .....**579530**
  - #10 Upper .....**579531**  
**Bag \$9.50**

**Advance Ferric Sulfate**

- 15.5% Ferric Sulfate Hemostatic Gel
- Provides effective tissue management when completing impressions, restorations, crown & bridge procedures
- Will not flow onto surrounding tissue.

**Standard Kit**

Contains: 4-1.2ml prefilled syringes, 8 brush tips .....**579281**  
**Ea \$16.33**

**Bulk Syringe**

Contains: 30ml syringe .....**579460**  
**Ea \$35.17**



**Hemostatic Gel**

- Viscous gel will not flow onto surrounding tissue
- Surfactant lowers tension for better wound penetration
- Quickly controls bleeding

**Aluminum Chloride-25%**

4/pack with tips.....**579283**  
**Ea \$16.25**

30mL syringe .....**579284**  
**Ea \$28.41**



**Sodium Hypochlorite**

For the debridement and irrigation of root canals during and after instrumentation

16oz Bottle .....**579267**  
**Ea \$10.35**



# Advance BASIC Great Value on Everyday Products!



## BASIC Non-Woven Sponges

- Soft, viscose/poly blend, non-woven material for increased patient comfort.
- Economical alternative to cotton.
- No chemical additives to avoid allergic reactions.
- Non-sterile

**Non-Sterile - 2" x 2"**  
Case/5000 ..... **575070**  
Ea **\$39.99**

**Non-Sterile - 4" x 4"**  
Case/2000 ..... **575071**  
Ea **\$43.99**

## BASIC Bracket Tray Covers

High quality paper tray cover.



**Ritter - Size "B"**  
**(8.5" x 12.25") Case/1000**  
White ..... **561033**  
Ea **\$24.86**

## BASIC Cotton Rolls

Pliable, non-linting super absorbent rolls made from 100% cotton. They do not adhere to mucous membranes and will retain their shape and strength when wet.



**Non-Sterile #2 Medium**  
Case/2000 ..... **577297**  
Ea **\$18.75**



## BASIC Self-Sealing Sterilization Pouches

- Dual External Indicators, Built-in Internal Indicator
- Blue-Tinted, Transparent Film
- Triple Seal & Heat Processed
- Wide Self-Seal Adhesive Strip

**Box/200**  
2.25" x 4" ..... **576806**  
Ea **\$5.99**

2.75" x 9" ..... **576807**  
Ea **\$9.99**

3.5" x 5.25" ..... **576811**  
Ea **\$10.48**

3.5" x 9" ..... **576812**  
Ea **\$9.75**

5.25" x 10" ..... **576814**  
Ea **\$14.25**

7.5" x 13" ..... **576815**  
Ea **\$32.67**

## Economy Pkg/1000

2.75" x 9" ..... **576809**  
Ea **\$45.75**

3.5" x 9" ..... **576813**  
Ea **\$46.50**



## Basic Curved Handle

Adult 38 Soft Tuft Rounded

Assorted: Dark Blue, Dark Purple, Dark Green, Orange  
Box/72 ..... **578056**  
Ea **\$29.93**



## Basic Compact Head Straight Handle

Young Adult 26 Soft Tuft

Assorted: Orange, Red, Dark Blue, Aqua  
Box/72 ..... **578053**  
Ea **\$28.81**



## Basic Child's Curved Handle

Ages 3-5 23 Soft Tuft

Assorted: Red, Green, Blue, Purple, White  
Box/72 ..... **578055**  
Ea **\$19.15**



## Basic Straight Handle

Adult 39 Rounded Tuft

Assorted: Red, Dark Blue, Dark Purple, Dark Green  
Box/72 ..... **578054**  
Ea **\$27.36**



## Basic Twist

Compact Adult 27 Tuft

Assorted: Red, Green, Blue, Purple  
Box/72 ..... **578064**  
Ea **\$31.15**



**Advance BASIC GREAT VALUE on Everyday Products!**





# ENDODONTICS

Each Advance Endodontic NiTi File system includes matching:  
 - Gutta Percha  
 - Paper Points




**\$26.92** Pk/6  
 Compare to:  
 ProTaper® Gold and EdgeTaper Platinum™




**\$26.92** Pk/6  
 Compare to:  
 ProTaper® Universal and EdgeTaper™



**\$18.33** Pk/3   **\$24.43** Assorted Pk/4  
 Compare to:  
 Wave One Gold® and EdgeOne Fire™



**\$18.33** Pk/3   **\$24.43** Assorted Pk/4  
 Compare to:  
 Wave One® and EdgeOne Platinum™





- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with **ProTaper® Gold and Edge Taper Platinum™** file systems for shaping and finishing
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Multiple tapers and the Heat Matrix technology of the Advance Taper HT ensure added flexibility while easily cutting through dentin in specific canal zones
- Shaping files have multiple tapers that ensure flexibility and cut dentin in specific canal zones with greater efficiency



- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with **ProTaper® Universal and Edge Taper™** file systems for shaping and finishing
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Shaping files have multiple tapers that ensure flexibility and cut dentin in specific canal zones with greater efficiency



Heat Treated

Better Flexibility

Order#	Size	Taper	Length	Pk Size
578940	SX	04	19MM	6
578915	ASST	X	21MM	6
578921	F2 25	07	21MM	6
578294	F3 30	08	21MM	6
578927	F4 40	09	21MM	6
578930	F5 50	06	21MM	6
578933	S1 18	05	21MM	6
578936	S2 20	02	21MM	6
578916	ASST	X	25MM	6
578919	F1 20	07	25MM	6
578922	F2 25	08	25MM	6
578925	F3 30	09	25MM	6
578928	F4 40	06	25MM	6
578931	F5 50	05	25MM	6
578934	S1 18	02	25MM	6
578938	S2 20	04	25MM	6
578917	ASST	X	31MM	6
578920	F1 20	07	31MM	6
578923	F2 25	08	31MM	6
578926	F3 30	09	31MM	6
578929	F4 40	06	31MM	6
578932	F5 50	05	31MM	6
578935	S1 18	02	31MM	6
578939	S2 20	04	31MM	6


Order#	Size	Taper	Length	Pk Size
578981	SX	04	19MM	6
578957	ASST	X	21MM	6
578961	F2 25	07	21MM	6
578964	F3 30	08	21MM	6
578986	F4 40	09	21MM	6
578969	F5 50	06	21MM	6
578972	S1 18	05	21MM	6
578975	S2 20	02	21MM	6
578958	ASST	X	25MM	6
578962	F1 20	07	25MM	6
578965	F2 25	08	25MM	6
578967	F3 30	09	25MM	6
578970	F4 40	06	25MM	6
578973	F5 50	05	25MM	6
578976	S1 18	02	25MM	6
578979	S2 20	04	25MM	6
578960	ASST	X	31MM	6
578963	F1 20	07	31MM	6
578966	F2 25	08	31MM	6
578968	F3 30	09	31MM	6
578971	F4 40	06	31MM	6
578974	F5 50	05	31MM	6
578977	S1 18	02	31MM	6
578980	S2 20	04	31MM	6



- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with **Wave One® Gold and Edge One Fire™** file systems
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Multiple tapers and the Heat Matrix technology of the Advance One HT ensure added flexibility while the optimized tip promotes safe canal navigation
- Finishing files ensure the proper apical fit for the matching Advance One or your existing Gutta Percha points



- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with **Wave One®** file system
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Finishing files ensure the proper apical fit for the matching Advance One or your existing Gutta Percha points



Heat Treated Better Flexibility	Order#	Size	Length	Pk Size
	578823	ASST	21MM	4
	578833	SM	21MM	3
	578830	PR	21MM	3
	578827	MD	21MM	3
	578824	LG	21MM	3
	578822	ASST	25MM	4
	578834	SM	25MM	3
	578831	PR	25MM	3
	578828	MD	25MM	3
578825	LG	25MM	3	
578835	SM	31MM	3	
578832	PR	31MM	3	
578829	MD	31MM	3	
578826	LG	31MM	3	

Order#	Size	Length	Pk Size
578780	ASST	21MM	4
578847	SM	21MM	3
578844	PR	21MM	3
578840	MD	21MM	3
578837	LG	21MM	3
578781	ASST	25MM	4
578848	SM	25MM	3
578845	PR	25MM	3
578841	MD	25MM	3
578838	LG	25MM	3
578849	SM	31MM	3
578846	PR	31MM	3
578842	MD	31MM	3
578839	LG	31MM	3

## Advance ENDODONTICS

100% guarantee on all of our products!



## Advance TAPER ABSORBENT PAPER POINTS



Each Advance Taper NiTi File system includes matching:

- Paper Points
- Gutta Percha

## Advance TAPER GUTTA PERCHA POINTS




- Size-matched to the Advance Taper & Advance Taper HT File System
- Superior absorbency

- Size-matched to the Advance Taper & Advance Taper HT File System
- Excellent radiopacity

Paper Points	Size	Order #	PK Size
SX-F3 PP ASST	SX-F3/28	578897	60
F1 PP 28 MM	F1/28	578902	60
F2 PP 28 MM	F2/28	578905	60
F3 PP 28 MM	F3/28	578907	60
F4 PP 28 MM	F4/28	578908	60
F5 PP 28 MM	F5/28	578910	60

Gutta Percha	Size	Order #	PK Size
SX-F3 GP ASST	SX-F3/28	578878	60
F1 GP 28 MM	F1/28	578881	60
F2 GP 28 MM	F2/28	578884	60
F3 GP 28 MM	F3/28	578886	60
F4 GP 28 MM	F4/28	578888	60
F5 GP 28 MM	F5/28	578890	60

## Advance ONE ABSORBENT PAPER POINTS



Each Advance One NiTi File system includes matching:

- Paper Points
- Gutta Percha

## Advance ONE GUTTA PERCHA POINTS



- Size-matched to the Advance One & Advance One HT File System
- Millimeter-marked length indicator lines
- Superior absorbency

- Size-matched to the Advance One & Advance One HT File System
- Millimeter-marked length indicator lines
- Excellent radiopacity
- Color coded size identification

Paper Points	Size	Order #	PK Size
ADV ONE PP	ASST	578855	60
ADV ONE PP	SMALL	578866	60
ADV ONE PP	PRIMARY	578863	60
ADV ONE PP	MEDIUM	578861	60
ADV ONE PP	LARGE	578859	60

Gutta Percha	Size	Order #	PK Size
ADV ONE GP	ASST	578796	60
ADV ONE GP	SMALL	578808	60
ADV ONE GP	PRIMARY	578805	60
ADV ONE GP	MEDIUM	578803	60
ADV ONE GP	LARGE	578801	60

# 100% guarantee on all of our products!

# Optimize Tissue Recovery Time with an Anti-Bacterial Endodontic Sealer

Advance Dental Products **Bioceramic Endodontic Sealer** combines the proven advantages of **MTA (Mineral Trioxide Aggregate)** with a superior canal obturation product. Our bioceramic sealer is as biological as the actual pulp tissue, and is engineered to penetrate, fill, and seal main and lateral root canals. The 13% MTA based sealer maintains an optimal pH, which elicits anti-bacterial effects and encourages new tissue formation and recovery.

- Provides rapid recovery of tissue without causing inflammatory reaction
- Avoids tooth stain and discoloration; eugenol and bismuth oxide free
- Easy to use dual syringe automix
- Easy removal for retreatment
- No need for low temperature gutta percha cone use standard heat settings for warm vertical
- Working Time: 23 minutes, Set Time: 130 minutes

Up to  
**20% - 40%**  
savings over comparable  
leading brands  
per gram!



**\$119.94** Each  
#578778 - 1 Ea. 4g Dual Barrel Syringe  
**\$15.32** Package  
#578779 - 10 Pk. Auto Mix Tips

**Advance RCS**  
BIOCERAMIC ROOT CANAL SEALER

## Work with Advance Taper and Advance Taper HT NiTi Rotary Files

### Tri Auto ZX2

by Morita

- » Cordless micromotor with built-in apex locator
- » Apical patency by rotary instrumentation
- » Now possible with rotary files as small as 10#
- » Observe the file tip position with apex locator function

577194 - Tri Auto ZX2

EACH  
**\$2399.00**



### Endo-Mate TC2 Cordless Endodontic Micromotor


by NSK

- » Cordless endodontic handpiece with torque control
- » Large LCD panel
- » Access friendly Key, Auto Reverse and Alarm Function
- » Apex locator capability
- » Contains: MP-F16R Head, Motor, and Charger

571894 - Endo-Mate TC2 Cordless Endodontic Micromotor

EACH  
**\$1318.00**





## More Flexibility for Half The Cost

- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with Vortex™ Blue and Edge Sequel™ Sapphire file systems
- Creates reproducible access and predictable canal patency
- Increased torque strength for safe and predictable canal navigation while the unique design enhances tactile experience
- Guide files polish dentin surfaces for easier identification of hidden canals

Order#	Size	Length	PK Size
578810	15/.04	21MM	4
578812	15/.04	25MM	4
578814	15/.04	30MM	4
<hr/>			
578816	20/.04	21MM	4
578817	20/.04	25MM	4
578818	20/.04	30MM	4
<hr/>			
578819	25/.04	21MM	4
578820	25/.04	25MM	4
578821	25/.04	30MM	4
<hr/>			
578854	35/.04	21MM	4
578856	35/.04	25MM	4
578858	35/.04	30MM	4
<hr/>			
578860	40/.04	21MM	4
578862	40/.04	25MM	4
578864	40/.04	30MM	4
<hr/>			
578865	45/.04	21MM	4
578867	45/.04	25MM	4
578868	45/.04	30MM	4
<hr/>			
578869	50/.04	21MM	4
578871	50/.04	25MM	4
578872	50/.04	30MM	4
<hr/>			
578873	.04 (15-40)	25MM	4

Order#	Size	Length	PK Size
578874	15/.06	21MM	4
578875	15/.06	25MM	4
578876	15/.06	30MM	4
<hr/>			
578877	20/.06	21MM	4
578879	20/.06	25MM	4
578880	20/.06	30MM	4
<hr/>			
578883	25/.06	21MM	4
578885	25/.06	25MM	4
578887	25/.06	30MM	4
<hr/>			
578893	35/.06	21MM	4
578894	35/.06	25MM	4
578895	35/.06	30MM	4
<hr/>			
578896	40/.06	21MM	4
578898	40/.06	25MM	4
578900	40/.06	30MM	4
<hr/>			
578901	45/.06	21MM	4
578903	45/.06	25MM	4
578904	45/.06	30MM	4
<hr/>			
578906	50/.06	21MM	4
578909	50/.06	25MM	4
578911	50/.06	30MM	4
<hr/>			
578912	.06 (15-40)	25MM	4

\*Market research and test studies available upon request. 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, we will issue a full refund of the purchase price through your distributor. ProTaper® Gold, ProTaper® Universal, Vortex™ Blue, Wave One®, and Wave One® Gold are registered trademarks of Dentsply Sirona. Edge One Fire™ Edge Sequel™ Sapphire, Edge Taper™, and Edge Taper Platinum™ are registered trademarks of Edge Endo.



### Cobalt-TWIRL

ABSORBENT PAPER POINTS

- Size-matched to the Advance Cobalt File System
- Superior absorbency



### Cobalt-TWIRL

GUTTA PERCHA POINTS

- Size-matched to the Advance Cobalt File System
- Excellent radiopacity



Each Advance Cobalt Twirl NiTi File system includes matching:  
- Paper Points  
- Gutta Percha

Paper Points	Gutta Percha	Size	PK Size
578937	578787	15/.04	60
578941	578788	15/.06	60
578942	578789	20/.04	60
578943	578790	20/.06	60
578944	578791	25/.04	60
578945	578792	25/.06	60
578946	578793	30/.04	60
578947	578794	30/.06	60
578948	578795	35/.04	60
578949	578797	35/.06	60
578950	578798	40/.04	60
578951	578799	40/.06	60
578952	578800	45/.04	60
578953	578802	45/.06	60
578954	578804	50/.04	60
578955	578806	50/.06	60
578956	578807	ASST .04	60
578959	578809	ASST .06	60

### Dust Free Alginate

- Mint
- Color changing and dust free.
  - Easy to use, avoiding remakes.
  - Shorter duration in the mouth maximum patient comfort.
  - Pleasing mint flavor.



**1 lb. Bag - Color Changing**  
Fast Set Color Changing (Mint)..... **579373**  
Ea **\$11.38**

**1 lb. Bag**  
Regular Set (Mint)..... **579372**  
Fast Set (Mint)..... **579371**  
Ea **\$11.38**

**Alginate Substitute - 50ml 8/Pack, Fast Set**  
**579985** ..... **Pack \$53.92**

**Alginate Substitute - 380/mL, 2/pack, Fast Set**  
**579986** ..... **Pack \$84.44**

	Work Time	IO Set Time
Fast Set (50mL)	1:00	1:00
Fast Set (380mL)	1:15	1:00

### Advance VPS

Using state-of-the-art hydrophilic technology, Advance Vinyl Polysiloxane Impression Material has great "wetability". It won't "run away" from the prepared impression area. It will stay exactly where you want it.

- Uniform consistency assures crisply detailed impressions every time.
- Advance Vinyl Polysiloxane Impression Material retains its dimensional stability, even after disinfection by immersion.
- Pleasant, neutral taste and odor dramatically increase patient acceptance, helping to reduce gag reflexes during impressions.
- Advance Vinyl Polysiloxane Impression Material achieves its maximum physical properties upon setting. No "rest distortion" due to change over time.
- Available in three basic viscosities and a range of setting times and hardness values, you will find a combination of Advance Vinyl Polysiloxane Impression material to meet your every need.

FAST SET	Work Time	IO Set Time
Light Body (50mL)	1:00	1:00
Medium (50mL)	1:00	1:00
Heavy Body (50mL)	1:00	1:00
Heavy Body (380mL)	1:30	1:00

REGULAR SET	Work Time	IO Set Time
Light Body (50mL)	2:30	1:00
Medium (50mL)	2:30	1:00
Heavy Body (50mL)	2:30	1:00

PUTTY	Work Time	IO Set Time
Fast Set	1:00	2:00

### Advance VPS

(Contains: 4-50ml cartridges with tips)

#### Unflavored:

Fast Set, Light Body.....**579975**  
Regular Set, Light Body.....**579976**  
Fast Set, Medium Body.....**579979**  
Regular Set, Medium Body.....**579980**  
Fast Set, Heavy Body.....**579983**  
Regular Set, Heavy Body .....**579984**  
**4/Pack \$34.56**

#### Mint:

Fast Set, Light Body.....**579973**  
Fast Set, Medium Body.....**579977**  
Fast Set, Heavy Body.....**579981**  
**4/Pack \$34.56**

#### Vanilla:

Fast Set, Light Body.....**579974**  
Fast Set, Medium Body.....**579978**  
Fast Set, Heavy Body.....**579982**  
**4/Pack \$34.56**

Fast Set, Heavy Body,  
380ml Cartridges.....**579987**  
**4/Pack \$99.62**

Fast Set, Putty.....**579988**  
**4/Pack \$44.34**

Prices subject to change due to manufacturer price increases and market fluctuations. Orders for these supplies can not be guaranteed for delivery and are subject to backorder status, allocation of amounts and/or cancellation.

# A smart solution for easier compliance.



## lexa

You expect thorough sterilization. But what if a system offered complete traceability, too?

The Lexa Sterilizer delivers bone-dry loads using patented vacuum technology. It documents every cycle from initiation to release, and makes printing package labels almost effortless.

Lexa's a smart solution for easier compliance—and your busy dental practice.



To learn more about Lexa, contact your local authorized dealer.



# STAY WARM WITH HOT SAVINGS

from **Crest** + **Oral-B**

**BUY 12 OF SAME, GET 3 SAME FREE!**

NO MIX AND MATCH  
OFFER VALID 10/1/2020- 12/31/2020



**Implant System**

577119  
\$104.99/bundle



**Gingivitis System**

570811  
\$88.00/bundle



**OrthoEssentials System**

570810  
\$88.00/bundle



**Kids 6+ Electric Bundle**

566767  
\$29.86/bundle



**Kids 3+ Spiderman Electric Bundle**

577115  
\$29.86/bundle



**Kids 3+ Princess Electric Bundle**

577114  
\$29.86/bundle

**BUY 5 OF SAME, GET 1 SAME FREE!†**

NO MIX AND MATCH  
OFFER VALID 10/1/2020- 12/31/2020



**Oral-B® Glide™ PRO-HEALTH™  
Deep Clean Mint**  
4m 72/box  
213253 **\$43.26**



**Oral-B® Glide™ PRO-HEALTH™  
Floss Dispenser**  
200 meter vial-1/box  
213252 **\$8.29**



**Oral-B® Glide™ PRO-HEALTH™  
Floss Font Refills**  
200 meter font refills-2/box  
213251 **\$14.99**

Crest® + Oral-B® products are sold at dental professional pricing for limited use in your practice. Products purchased through professional dental distribution channels are not meant for resale outside your dental care practice. © Disney © Disney•Pixar © & TM Lucasfilm © 2020 P&G

continuing the care that starts in your chair





# CHALLENGE ACCEPTED.

Nobody thought a bulk fill flowable was possible. Until we invented Surefil® SDR flow+ material. Today, we **redefine restorative materials again with Surefil one™ restorative**, transforming how dentists manage challenging clinical situations where isolation is difficult or when patient compliance is an issue.

- **Self-adhesive**—no separate adhesive step needed
- **Unlimited depth of cure**—no layering of composite materials, making it faster than a traditional composite restoration, **especially in compromised cases.**

**Surefil one™ restorative. Progress where you need it most.**

**For more information, visit [www.dentsplysirona.com](http://www.dentsplysirona.com)**

Surefil one™  
Self-Adhesive Composite Hybrid



# INVEST IN YOUR PRACTICE

Now with Our Expanded Product Portfolio

October 1 – December 31, 2020

## PRACTICE ESSENTIALS


Buy 2 Enzymax™ Gallons,  
**RECEIVE 1 AT NO CHARGE<sup>1</sup>**

Buy 3 VistaTab™ Dental Waterline Cleaner Tablets,  
**RECEIVE 1 AT NO CHARGE**

Buy 1 bag of Safe-Flo™ SE Valves, **RECEIVE 1 BAG OF COMFORT PLUS™ SALIVA EJECTORS AT NO CHARGE**

Buy 4 of the same Hand Essentials™ Skin Repair Cream Products,  
**RECEIVE 1 AT NO CHARGE<sup>1</sup>**

Buy any 3 Hu-Friedy Sterilization Chemical Indicator or Integrator Monitoring Products,  
**RECEIVE 1 OF THE SAME AT NO CHARGE<sup>2</sup>**



## PERIODONTAL

Buy 3 boxes of twist™ plus Prophy Angles (soft cup DPA), **RECEIVE 1 AT NO CHARGE**

Buy any 5 Probes, Explorers, or Expros,  
**RECEIVE 1 AT NO CHARGE<sup>1</sup>**

Buy any 5 Scalers, **RECEIVE 1 AT NO CHARGE<sup>1</sup>**

Buy any 12 Scalers, **RECEIVE 3 AT NO CHARGE<sup>1</sup>**

Buy any 24 Scalers, **RECEIVE 6 AT NO CHARGE<sup>1</sup>**

Buy any 10 Diagnostic Instruments, **RECEIVE 3 AT NO CHARGE<sup>1</sup>**

Buy 2 Diamond Sharpening Cards, **RECEIVE 1 AT NO CHARGE**

Buy 1 Sidekick™ Sharpener, **RECEIVE 4 SCALERS AT NO CHARGE**



## ENDODONTIC

Buy 5 Endodontic Hand Instruments,  
**RECEIVE 2 RUBBER DAM CLAMPS AT NO CHARGE**

Buy 5 Endodontic Micro Instruments,  
**RECEIVE 1 AT NO CHARGE<sup>1</sup>**

Buy 6 Rubber Dam Clamps, **RECEIVE 2 AT NO CHARGE<sup>1</sup>**



## POWER

Buy 1 SWERV™ Unit, **RECEIVE 5 ULTRASONIC INSERTS AT NO CHARGE**

Buy 6 Ultrasonic Inserts, **RECEIVE 3 AT NO CHARGE<sup>1</sup>**

Buy 3 Ultrasonic Inserts, **RECEIVE 1 AT NO CHARGE<sup>1</sup>**

Buy 2 Boxes of AIRFLOW® Powder,  
**RECEIVE 1 BOX AT NO CHARGE<sup>1</sup>**

Buy 2 Piezon Tips, **RECEIVE 1 TIP AT NO CHARGE<sup>1</sup>**



## RESTORATIVE


Buy any 7 Stainless Steel Molar Crowns Refills,  
**RECEIVE 3 AT NO CHARGE**

Buy any 7 Operative Instruments, **RECEIVE 2 AT NO CHARGE<sup>1</sup>**

Buy any 4 XTS™ Instruments, **RECEIVE 1 AT NO CHARGE<sup>1</sup>**

Buy 3 Akro-Flex™ Instruments,  
**RECEIVE 1 XTS™ INSTRUMENT AT NO CHARGE**

Buy 5 XTS™ Instruments, **RECEIVE 2 BLACK LINE RUBBER DAM CLAMPS AT NO CHARGE**



## SURGICAL

Buy any 3 Boxes of Perma Sharp™ Sutures,  
**RECEIVE 1 BOX AT NO CHARGE<sup>1</sup>**

Buy any 3 Surgical Hand Instruments (mix and match),  
**RECEIVE 1 AT NO CHARGE<sup>1</sup>**

Buy any 3 Hinged Instruments (mix and match),  
**RECEIVE 1 AT NO CHARGE<sup>1</sup>**



## ORTHODONTIC

Buy 6 Cutters or Pliers, **RECEIVE 2 AT NO CHARGE<sup>1</sup>**

Buy 4 Cutters or Pliers,  
**RECEIVE 1 MATHIEU AT NO CHARGE<sup>3</sup>**

Buy 5 Orthodontic Hand Instruments,  
**RECEIVE 1 AT NO CHARGE<sup>1</sup>**



Visit us online at [Hu-Friedy.com/Spotlight](http://Hu-Friedy.com/Spotlight) | Promo code: 2863

Offer(s) valid in the 50 United States and District of Columbia. Participation is not permitted for customers who participate in any other pricing program, including government, group practices, community health centers, schools, institutions, and discount customers. Receipt of at no charge goods, after rebate, cannot be combined with other offers or promotions. Offer(s) not valid on Made-To-Order (MTO) products. To receive your at no charge product(s), please fax a copy of your distributor invoice dated October 1 – December 31, 2020 indicating the required purchases and your at no charge goods selection to 773-868-3560 or by email to [hufriedysales@usa-800.com](mailto:hufriedysales@usa-800.com). Proofs of purchases will not be returned and may be used for marketing purposes. Redemption must be postmarked by January 31, 2021. 1) Equal or lesser value. 2) Offer valid only on IMS™ Steam Sterilization Chemical Indicators and Integrators. 3) Only valid on 678-330 and 678-331.

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, or Crosstex International, Inc., its affiliates or related companies, unless otherwise noted. AIRFLOW is a registered trademark of EMS Electronic Medical Systems S.A. Marks not registered in all jurisdictions.

©2020 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-911/0820

Hu-Friedy and Crosstex are now proud members of  **HuFriedyGroup**  
The Best in Practice

# SDI SPECIALS

OCTOBER | NOVEMBER | DECEMBER 2020

COMPOSITES | GLASS IONOMERS  
TOOTH WHITENING | EQUIPMENT  
CEMENTS | ADHESIVES  
ETCHANTS | SEALANTS  
AMALGAMS | ACCESSORIES

## COMPOSITES

**3+1\*** BUY 3,  
GET 1 FREE

**10+5\*** BUY 10,  
GET 5 FREE

### AURA EASY

STOCK LESS, DELIVER MORE

- 4 shades cover 8 VITA classical A1-D4® shades
- Increased wear resistance and strength
- Smooth high gloss finish every time
- Non-sticky consistency makes it simple to sculpt
- 4 universal shades to minimize inventory



### AURA EASY FLOW

ULTRA UNIVERSAL RESTORATIVE MATERIAL

- Optimal flow for cavity adaptation
- High visibility on radiographs
- Strong mechanical properties
- Easy shade matching
- Optimized optical properties
- Stable aesthetics over time



Four unique Aura Easyflow shades cover 8 VITA® A1-B4® shades, enabling an easy switch from traditional shading systems.

**ae1** | **ae2** | **ae3** | **ae4**  
A1/B1 | A2/B2 | A3/B3 | A3.5/B4

Approximate equivalence to VITA® classical shade guide



**2+1\*** RADII CURING LIGHTS | BUY 2 GET 1 FREE



### RADII-CAL

- Pulse technology = High intensity 1200mW/cm<sup>2</sup> plus low heat emission
- Stays charged for 1200 ten second cures (equivalent to 3h 20 min continuous cure)



### RADII PLUS

- 5 year warranty~
- Pulse technology = High intensity 1500mW/cm<sup>2</sup> plus low heat emission
- Interchangeable LED attachments for curing, transillumination and whitening

\*5 year warranty on all parts except battery and full arch bleaching LED attachment which is 2 years.



### RADII XPERT

- Unique target assist for precision curing
- Optimally collimated beam
- Consistent light intensity over clinically relevant distances

## TRY THE RIVA RANGE TODAY!

**2+1\*** BUY 2 RIVA  
GET 1 FREE



## BUY 3 RIVA BOXES GET ULTRAMAT 2 FOR FREE

### ULTRAMAT 2

HIGH SPEED MULTI-USE TRITURATOR

- 5 year limited warranty
- Low vibration and extremely quiet
- Accurate and consistent trituration



COMPOSITES | GLASS IONOMERS | TOOTH WHITENING | EQUIPMENT | CEMENTS  
ADHESIVES | ETCHANTS | SEALANTS | AMALGAMS | ACCESSORIES

Redeem your free goods online using this link:  
<http://www.sdi.com.au/freegoods-us>

[facebook.com/sdiusa](https://www.facebook.com/sdiusa) [instagram.com/sdiusa](https://www.instagram.com/sdiusa)

\*Please allow 2-4 weeks for delivery. Please note offers cannot be combined. Offers cannot be used in conjunction with contract pricing. Valid on invoices dated October 1st - December 31st, 2020. All claims for free goods must be made within 60 days from the end of promotion date. Free products must be of equal or lesser value. Offers are valid in the US.



YOUR  
SMILE.  
OUR  
VISION.

SDI (North America) Inc.  
1279 Hamilton Parkway  
Itasca, IL 60143  
Toll Free 1 800 228 5166  
Fax: 630-361-9222  
[usa.canada@sdi.com.au](mailto:usa.canada@sdi.com.au)  
[www.sdi.com.au](http://www.sdi.com.au)





# EQUIA Forte® HT

Bulk Fill Glass Hybrid Restorative System

Building on the success of EQUIA Forte®

- Improved translucency and physical properties for esthetic, strong restorations
- Intelligent particle size distribution and interaction allowing superb handling characteristics
- Exceptional marginal seal for long-term resistance to microleakage
- High fluoride release at tooth-restorative interface
- Ideal restorative treatment for high-risk caries, pediatric and geriatric patients

### EQUIA Forte® HT Intro Kit: \$365.61

Contains: 48 capsules, one 4mL EQUIA Forte® Coat bottle and one 5.7mL CAVITY CONDITIONER.

579404 A1	579405 A2	579406 A3	579407 A3.5	579408 B1
579409 B2	579410 B3	579411 C4	579412 Assorted	

Buy 3,  
Get 1  
FREE!\*



### EQUIA Forte® HT Fil 48-Capsule Refills: \$335.16

Contains: 48 capsules.

579413 A1	579414 A2	579415 A3	579416 A3.5	579417 B1
579418 B2	579419 B3	579420 C4	579421 Assorted	

800.323.7063 • [www.gcamerica.com](http://www.gcamerica.com) • [www.gcamerica.com/training](http://www.gcamerica.com/training) • © 2020 GC America Inc.

\*Promo 2454. To receive your no charge goods, send a copy of your Goetze invoice dated 10/1/2020–12/31/2020 to GC America, 3737 West 127th Street, Alsip, IL 60803 or fax to 800-423-2963 Attn: Customer Care Group. Must be redeemed within 6 weeks of dealer invoice date. Offers cannot be combined with any other offers or promotions. No substitutions allowed. Offers are valid in the US and Canada only. No charge goods will be of equal or lesser value to the lowest selling unit price. Please allow 6-8 weeks for no charge goods delivery. No charge goods will be fulfilled by GC America unless noted otherwise. All promos expire 12/31/2020.

# LuxaCore® Z Dual

## The core material that cuts like dentin

Choose the material for core build-ups and single-step post cementation that cuts identical to dentin and provides unbeatable compressive and flexural strength.

Discover the Possibilities



BUY 3, GET 1 OF  
THE SAME FREE!



512340 (A3) 48gm Cartridge, 35 Automix & Intra-Oral Tips .....	\$237.00
551516 (Light Opaque) 48gm Cartridge, 35 Automix & Intra-Oral Tips .....	\$237.00
551515 (Blue) 48gm Cartridge, 35 Automix & Intra-Oral Tips .....	\$237.00
540593 (A3) 2-9gm Syringes, 10 Smartmix & Intra-Oral Tips .....	\$170.54
512341 (Light Opaque) 2-9gm Syringes, 10 Smartmix & Intra-Oral Tips .....	\$170.54
540593 (Blue) 2-9gm Syringes, 10 Smartmix & Intra-Oral Tips .....	\$170.54



Offer valid 10/01/20-12/31/20. To receive free goods, please fax your paid Goetze Dental invoice to DMG America at 201-894-0213, or email to [freegoods@dmg-america.com](mailto:freegoods@dmg-america.com) indicating shade of free good. The most popular shade will be shipped unless otherwise noted. Each qualifying purchase must be submitted on one invoice. Promotion cannot be combined with any other offer and may change or be discontinued at any time without notice. All offers must be redeemed within 30 days of purchase. Limit 5 offers per dental office. Promotions not available to dental supply dealers. Visit [www.dmg-america.com](http://www.dmg-america.com) or call 800.662.6383 for more information.

# Brush&Bond®

## Don't Leave Your Bonds Up to Chance

- Bonds direct and indirect composite restorations to all tooth surfaces.
- Permanently seals the tubules to prevent sensitivity.
- Protects the pulp against microleakage and bacterial invasion.
- Creates an acid and wear-resistant shield for the tooth.
- Doubles the surface hardness of the preparation.
- Cures with all currently available curing lights.
- Low film thickness won't affect crown fit.



Apply in Just Three Easy Steps! 1) Apply 2) Dry 3) Cure

**BUY 3, GET 1 FREE\***  
OF SAME PRODUCT



**Original Brush&Bond System 081300**  
**\$159.67, NET 3+1: \$119.75**

*Kit includes Brush&Bond liquid (3 ml) plus 100 standard activator brushes.*

**Original Brush&Bond Liquid 081302**  
**\$93.10, NET 3+1: \$69.83**

*Includes 1 bottle of Brush&Bond liquid (3 ml).*

**Brush&Bond Universal System 568399**  
**\$246.45, NET 3+1: \$184.84**

*Kit includes Brush&Bond Universal liquid (3 ml) plus 100 standard activator brushes.*

**Brush&Bond Universal Liquid 578416**  
**\$111.38, NET 3+1: \$83.54**

*Includes 1 bottle of Brush&Bond Universal liquid (3 ml).*

# Blu-Mousse®

VPS BITE REGISTRATION MATERIAL

**Bite registrations, matrices, model bases, bite plates and CAD scanning capabilities—the uses for Blu-Mousse just keep growing!**

- Exceptional accuracy.
- Choice of setting speeds: 30 seconds, 60 seconds, and 2 minutes.
- Plaster-like hardness (90-durometer).
- Fluffy, no-slump consistency.

Blu-Mousse remains one of the most accurate, most stable bite registration materials in dentistry. And even though the versatility of Blu-Mousse seems limitless, it is still the light, fluffy VPS material it was created to be. With its whipped cream texture, it offers no resistance when a patient bites down, so it won't trigger any reflexes that would affect the integrity of your bite registration.



**BUY 3, GET 1 FREE\***  
EQUAL OR LESSER VALUE



**Blu-Mousse Split Cartridge Kits**

**\$46.99, NET 3+1: \$35.24 ea.**

30-Second Set Time, Unscented 131220

60-Second Set Time, Unscented 574161

2-Minute Set Time, Unscented 131200

*Kits include (2) 50 ml cartridges (100 ml ttl.) and mixing tips.*

**Blu-Mousse Standard Cartridge Kits**

**\$41.11, NET 3+1: \$30.83 ea.**

30-Second Set Time, Unscented 571532

2-Minute Set Time, Unscented 571531

*Kits include (2) 50 ml cartridges (100 ml ttl.) and mixing tips.*

\*To receive your no additional charge goods, email a copy of your Goetze invoice to [dealerpromotions@parkell.com](mailto:dealerpromotions@parkell.com). Pricing and offer valid 10/01/20 through 12/31/20. Invoice must be received by 01/31/21. Please allow 3 – 6 weeks for delivery. Qualifying purchase must be on one invoice. Date of purchase must be between the dates on the offer.



# NOW INTEGRATED & AVAILABLE

THROUGH GOETZE DENTAL

**EDDSON**<sup>TM</sup>  
PRACTICE MANAGEMENT SOFTWARE



**RevenueWell**



**GOETZE**  
d e n t a l

Ready to get started?

VISIT  
[revenuwell.com/goetze](http://revenuwell.com/goetze)  
SCHEDULE A DEMONSTRATION  
800.692.0804





# SAFE. SIMPLE. EFFECTIVE.

Promote safe, effective sterilization with Midmark Steam Sterilizers. Plus, maintain more efficient and simpler compliance by going digital with your data logger. **Take advantage of this offer and put the power of simplicity to work for you.**



## Option 1: \$350 Rebate

Purchase any eligible\* Midmark M11, M9 or M3 Steam Sterilizer and receive a \$350 rebate.



## Option 2: \$550 Rebate

Purchase any eligible\* Midmark M11 or M9 Steam Sterilizer with a Sterilizer Data Logger and receive a \$550 rebate.

### TIMING

**Promotion order period:**  
September 1–December 31, 2020

**Last date to invoice and take shipment of products:**  
January 14, 2021

**Last date to claim rebates:**  
February 28, 2021

### REDEMPTION INSTRUCTIONS

1. Visit [midmark.com/simple2020](http://midmark.com/simple2020)
2. Upload invoice
3. Complete the information required on the form
4. Select **"Submit"**

### QUESTIONS

For product inquiries, please call Midmark Dental Customer Experience at 1.800.MIDMARK or visit [midmark.com](http://midmark.com).

For questions about the promotion, please call 1.800.MIDMARK, ext. 107575 or visit [midmark.com/simple2020](http://midmark.com/simple2020)

Eligible Model Numbers:  
Option 1: M9-040, M9-041, M11-040, M11-041, M3-040 and M3-041  
Option 2: M9-040-DL, M9-041-DL, M11-040-DL, M11-041-DL  
\*See website for full terms and conditions regarding this offer and any applicable federal regulatory conditions ([midmark.com/simple2020](http://midmark.com/simple2020)). Purchases must be made through Midmark authorized dealers to qualify ([midmark.com/distributor](http://midmark.com/distributor)).

009-10053-00 | © 2020 Midmark Corporation, Miamisburg, Ohio USA. All rights reserved.



USA

## MECHANICAL ROOM PROMOTION 2020



# SAVE NOW!

UP TO \$4,496 USD IN RETAIL SAVINGS

**FREE PRODUCT:** Rely on Midmark mechanical room equipment to keep your practice up and running so you can focus on what matters most—your patients.



## FREE PRODUCT

### Option 1

Purchase a single Midmark air compressor or vacuum, claim one free accessory:

- Buck Boost
- or
- Control Panel

### Option 2

Purchase a vacuum and air compressor combination, claim two free accessories:

- Buck Boost or Control Panel plus
- NXT Hg5 Amalgam Separator or Water Filter Bypass

May qualify for additional rebates with Midmark PLUS—visit [midmark.com/plus](http://midmark.com/plus) for details.

For a limited time, select free mechanical room accessories with your Midmark air compressor and vacuum purchases.



Buck boost



Control panel



NXT Hg5 amalgam separator



Water filter bypass separator

## TIMING

**Promotion order period:**  
May 11, 2020–December 31, 2020

**Last date to invoice and take shipment of products:**

January 14, 2021

## ELIGIBLE PRODUCTS

Buy any eligible\*\* compressor AND/OR vacuum to receive free product(s).

\*\*PowerAir®, ClassicSeries®, PowerVac®, PowerVac® G, PowerMax®

## QUESTIONS

For product and promotion inquiries, please call the Midmark dental customer experience department at 1.800.MIDMARK or visit [midmark.com/airvacSave](http://midmark.com/airvacSave).

Offer for dental practices only as sold by Midmark® authorized dental distributors in the U.S. All qualifying products must be entered on one purchase order. This offer constitutes a discount or reduction in price on the purchased product(s). Buyer/Lessee/Lessor is obligated to: (i) fully and accurately disclose the amount of this discount in Buyer/Lessee/Lessor's cost reports or claims for reimbursement to Medicare, Medicaid or other U.S. federal healthcare program that requires such disclosure; and (ii) provide documentation of the discount to representatives of the US Secretary of the Department of Health and Human Services and state agencies upon request. Buyer/Lessee/Lessor may request additional information from Midmark if necessary in order to meet its reporting or disclosure obligations. This offer may not be combined with any other offers, promotions or contract pricing except the Midmark PLUS Program.

[midmark.com/airvacSave](http://midmark.com/airvacSave)



## OPERATORY PACKAGE PROMOTION

While you focus on preventing the spread of the COVID-19 coronavirus with extra safety measures to keep your patients and team safe, you may be ready to equip your empty operatories to help manage your patient load. That's why we're offering savings on Midmark operatory equipment.



# Free product or discount

Purchase any eligible operatory package (chair and unit) and choose **ONE** of the following:



**Option 1: One Free Doctor's and Assistant's Stool Set**



**Option 2: Discount of 15% on Procenter Delivery Units**

### TIMING

**Promotion order period:**  
August 1 - December 31, 2020

**Last date to invoice and take shipment of products:**  
February 15, 2021

### ELIGIBLE PRODUCTS

- Elevance® Chair and Elevance Unit Packages
- Elevance Chair and Procenter Unit Packages
- Ultra Chair and Procenter Unit Packages
- Ultra Chair and Asepsis 21® Unit Packages

**NOTE: Package must include a chair and delivery unit.**

### QUESTIONS

For product inquiries, please call Midmark dental customer experience at 1.844.856.1235. For questions about the promotion, please contact 1.800.MIDMARK, ext. 108923.

See website for full terms and conditions regarding this offer and any applicable federal regulatory conditions ([midmark.com/OpSavings2020](http://midmark.com/OpSavings2020)).

[midmark.com/OpSavings2020](http://midmark.com/OpSavings2020)



# UTILITY Aeras Ramvac



## THE SMART UTILITY ROOM THAT MONITORS & COMMUNICATES



**TECHNICIAN:** The best part is, your service provider will get all the alerts and notifications you get, so you can stop worrying and leave the follow-up to them.



**MOBILE CONTROL:** Never worry about having forgotten to turn your system off at the end of the day. You can now flip the switch from your phone—anywhere, anytime.



**EMERGENCY ALERTS:** When something starts going wrong, the sensors will know, and you get alerted before it becomes a serious problem.



**QUICK FIXES:** The system not only tells you when trouble is brewing, it also troubleshoots for you, suggesting for example simple fixes like ventilating your back room to help avoid downtime.



**AUTOMATED SCHEDULING:** Automate daily activation and shutoff. Simply set it and forget it. When treatments run late, the equipment will sense activity and maintain power until you're done.



**MAINTENANCE NOTIFICATIONS:** Proactive service is critical to your system's health, but it's also easy to forget to do. Now your system will remind you well in advance.




**COMPREHENSIVE DASHBOARD:** A bird's eye view of your entire back office operation, including detailed technical outputs, configuration options, and service records.



**WARRANTY:** We believe in the durability of our smart system so much that we're increasing the industry's longest warranty period to 8 years. Your investment delivers peace of mind now and well into the future.\*

\*Terms apply

NOTE: All features above apply only if subscribing to monitoring service.



## THE INDUSTRY'S ONLY CONNECTED UTILITY ROOM

The Power and Reliability of  
Ramvac Vacuums and  
Compressors, Now in the Cloud.

# AERAS UTILITY WITH RAMVAC RELIABILITY

## COMPRESSORS

### DRYER, CLEANER AIR

The Aeras dual-desiccant setup produces an uninterrupted supply of clean air to power your practice. Air is stored at a pressure dew point of -40 C, preventing harmful bacteria growth, and particulates in the air are filtered to 0.01 microns, trapping most airborne bacteria and microbes.

### EFFICIENCY & DURABILITY

Our unique rocking-piston design doubles the life of our compression rings compared to standard piston design.

### FLEXIBLE DUAL-PRESSURE SETTINGS

A pressure toggle switch lets you easily switch between standard (85-115 PSI) and high-pressure (105-130 PSI), so you can add a milling machine without having to upgrade equipment.

## VACUUMS

### QUIET EFFICIENCY

Energy efficient Ramvac motors run at a minimum of 50% less BTUs than competitive technology, making them efficient —and as quiet as a conversation (66-79 dBA depending on model size).\*

### BUILT-IN SAVINGS

100% water-free operation means you spend less than with wet vacuum designs.

### UNRIVALED PERFORMANCE

The advanced design of our rotary vane pump ensures the vacuum is maintained even as additional HVEs are opened.

### MINIMAL MAINTENANCE

Ramvac vacuum systems are built with rotary vane pumps. Because they are lubricated, they require minimal maintenance and are designed to outlast non-lubricated systems.

### GENERAL COMPRESSOR INFORMATION

Model No.	Voltage	Dimensions W x D x H (in)	No. of Users	Peak Airflow Output @ 80 PSI (CFM)
C3-2	230V	29" x 21" x 29"	3	6.0
C4-1	115V	29" x 21" x 29"	4	6.6
C4-2	230V	29" x 21" x 29"	4	6.6
C6-2	230V	35" x 25" x 30"	6	9.3
C7-2	230V	35" x 25" x 30"	7	12.0
C11-2	230V	41" x 29" x 32"	11	18.1

Notes: Pressure is 85-115 PSI. High-Pressure Mode is 105-130 PSI. Pump Up Cycle Time, 85 to 115 PSI (Sec) is 40-45

\*Data on File.

### GENERAL VACUUM INFORMATION

Model No.	Input Power	Stacked Dimensions (Pump and Tank)	No. of Users	Breaker Size
RV4-12T	115v or 230v, 1 ph 60 Hz	25" x 19" x 44"	4	115v or 230v, 1 ph 60 Hz
RV5-2T	230v, 1 ph 60 Hz	32" x 21" x 58"	5	230v, 1 ph 60 Hz
RV7-2T	230v, 1 ph 60 Hz	32" x 21" x 58"	7	230v, 1 ph 60 Hz
RV10-2T	230v, 1 ph 60 Hz	32" x 22" x 58"	10	230v, 1 ph 60 Hz

Notes: All Aeras vacuum systems (pump and tank) except RV4-12T include a platform (575056) for stacking. Height dimensions do not include S-type filters

\*Data on File.

# DENTALEZ



UVANGEL®

## UV ANGEL AIR™

### Automated Continuous UV-C Air Treatment System

Our engineered air system uses patented UV-C light air purification technology to reduce levels of viruses, bacteria and fungi by automatically and continuously treating the air to create healthier environments.

**Simple** – UV Angel Air operates without interruption 24/7/365 and does not interfere with current staff workflow

**Seamless** – Designed with today's critical architectural and clinical considerations at the forefront, the system integrates a sealed UV-C air treatment chamber into existing or new construction in-ceiling lighting

**Effective** – Using the latest in advanced UV-C light purification technology, laboratory studies have shown effective removal of bacteria, fungus and viruses from the air

### How UV Angel Air Works


UV Angel uses years of advanced research and development in ultraviolet light and IOT enabled technologies to create a truly modern and effective air treatment system:

- Using patented UV-C treatment technology, air is quietly drawn into a sealed UV-C air chamber with a series of fans and filters
- Air is circulated through the UV-C air chamber where it is treated with an enclosed high intensity UV-C light to inactivate bacteria, fungus and viruses in the air
- Treated air is then returned to the room creating a healthier environment
- UV Angel Air is unobtrusive, works continuously, and with the in-ceiling design, maintains the valuable floor space in patient and staff areas



Research has shown that reducing contamination in the air also **REDUCES BACTERIA AND FUNGUS ON SURFACES.**

Hospital air samples, on average, are **UP TO 8 TIMES MORE CONTAMINATED THAN SURFACES.**<sup>1</sup>



<sup>1</sup>Lee, Linda D, DrPH, MBA, LV-17-C042, Can using active air UV-C technology reduce the amount of bacteria and/or fungus in the air and improve indoor air quality? ASHRAE Conference (2017)



# PREVENTATIVE CARE FOR THE AIR.

---

According to the WHO, some medical procedures performed on infectious patients can produce aerosols that are able to stay suspended in the air for longer periods of time. These aerosols can be inhaled by others if appropriate precautions aren't taken by wearing PPE and using isolation measures.

**INTRODUCING UV ANGEL AIR™ SERIES:** an unobtrusive, automated UV-C air treatment system that has been proven in laboratory testing to be up to 99.99% effective at eliminating harmful viruses, fungi, and bacteria.

For more information visit [uvangel.com/dental](http://uvangel.com/dental)



# FRESH START FOR 2021

## Save BIG plus get cash back with rebates.

Core Equipment Packages

Furniture Promotion

Sterilization Center Promotion

**ALL ORDERS MUST SHIP IN 2021**

(Must Ship January 1 – March 31, 2021)

Core Equipment Packages

**RECEIVE UP TO  
\$750 CASH REBATE\***

(Must be shipped January 1 – March 31, 2021)



**FREE STOOLS &  
SPECIAL PRICING ON SELECT  
CONFIGURED PACKAGES\***

(Must be shipped January 1 – March 31, 2021)

Furniture Promotion

**RECEIVE UP TO  
\$500 CASH REBATE\***

(Must be shipped January 1 – March 31, 2021)



**SAVE 10% ON ALL  
A-DEC FURNITURE WITH  
A CORE EQUIPMENT PACKAGE\***

(on orders of \$12,000 or more)

(Must be shipped January 1 – March 31, 2021)

Sterilization Center Promotion

**RECEIVE UP TO  
\$500 CASH REBATE\***

(Must be shipped January 1 – March 31, 2021)



**SAVE 10% ON  
STERILIZATION CENTERS**

(on orders of \$12,000 or more)

(Must be shipped January 1 – March 31, 2021)

To learn more about A-dec solutions visit [A-dec.com](http://A-dec.com)  
Promo runs October 1, 2020 - December 31, 2020. Must ship by March 31, 2021.  
All amounts shown in USD. See Terms & Conditions for complete details.  
©2020 A-dec Inc. All rights reserved.





**FREE STOOLS & SPECIAL PRICING ON CONFIGURED PACKAGES**



**RECEIVE UP TO \$750\* CASH REBATE**

**ALL ORDERS MUST SHIP IN 2021**

(Must Ship January 1 – March 31, 2021)

**500 PACKAGE**

FREE 521 & 522 Stools  
(\$2,965 Value)

500B Chair + 500B Delivery + Any 500 LED Light\*



**+\$750 REBATE**

Take delivery in 2021

**400 CONFIGURED PACKAGE**

\$25,525

**\$19,999**

+FREE 421 Stool  
(\$6,581 Combined Savings)

411 Chair + 332 Delivery + 572 LED Light package details on back



**+\$250 REBATE**

Take delivery in 2021

**300 FLEXIBLE PACKAGE**

+FREE 421 Stool  
(\$1,055 Value)

300 Delivery System  
+ 511A or 411 Chair + Any 500 LED Light\*



**+\$250 REBATE**

Take delivery in 2021

**SAVE 10% ON FURNITURE**

when you pair dental furniture with a core equipment package



**+\$500 REBATE**

Take delivery in 2021

**Take advantage of the 2021 Fresh Start  
Order October 1 through December 31 - Take delivery to 2021\*.**

Orders must ship by March 31, 2021.

To learn more about A-dec solutions visit [A-dec.com](http://A-dec.com)





## STACK UP THE SAVINGS.

**Save 10%**  
on Dry Vacuum Systems  
& Oil-Free Compressors  
by Dean Dental Systems

To learn more about A-dec solutions  
visit [A-dec.com](http://A-dec.com).

**Promo runs**

**October 1, 2020 – December 31, 2020**  
**AVAILABLE IN U.S. ONLY**

All amounts shown in USD.  
See Terms & Conditions for complete details.  
©2020 A-dec Inc. All rights reserved.



## The Next Generation-Mechanical Room

### Dry Vacuum System

- Variable Speed Technology
- Small Footprint
- Remote Monitoring Enabled

### Oil-Free Compressor

- Compact and Quiet
- Low Maintenance Dry Air
- Intelligent Controls



**SAVE 10% ON A-DEC  
STERILIZATION CENTERS**



**RECEIVE \$500 CASH REBATE**

**ALL ORDERS MUST SHIP IN 2021**

(Must Ship January 1 - March 31, 2021)



**Take advantage of the 2021 Fresh Start  
Order October 1 through December 31 - Take delivery to 2021\*.**

Orders must ship by March 31, 2021.

To learn more about A-dec solutions visit [A-dec.com](http://A-dec.com)

All amounts shown in USD.  
See Terms & Conditions for complete details. ©2020 A-dec Inc. All rights reserved.



THE  
CHOICE IS  
YOURS.

A-dec 400



## THERE'S AN A-DEC FOR EVERYONE.

Every dental practice is unique. That's why A-dec offers a full line of dental chairs, with a wide variety of options and price points, all designed for comfort, performance, and reliability.

From the entry-level A-dec 300, to the classic A-dec 400, to the premium A-dec 500, there's an A-dec chair to fit your practice, your personal style, and your budget.


VISIT

[a-dec.com/choice](https://a-dec.com/choice)

to see our award-winning family of A-dec dental chairs—and find the one that's right for you.



A-dec 300



A-dec 500

**a-dec**



# What is your formula for **RECOVERY?**



3D  
IMAGING

+

INTRAORAL  
SCANNING

+

CHAIRSIDE  
MILLING

=

MORE DENTISTRY,  
MORE EFFICIENTLY



It takes unique problem solving to succeed in this new era. Making wise investments in technology unlocks the potential for exponential growth.

**PLANMECA**  
Better care through innovation

Let us show you the solution for success...  
[info.planmecausa.com/formula](http://info.planmecausa.com/formula)



## ABOUT US

Goetze Dental is a family-owned business serving private dental practices just like yours since 1884.

We work to understand your goals and dreams and strive to provide the best dental supplies, equipment, technology and service.

### We believe the best patient care is delivered in the private practice.

100% of our energy and resources will continuously center on the products and services we offer to the privately held dental practice through integrity, education, innovation and customer service.

We stand united with the artists, the scientists and the explorers of the dental community, whose interest in providing exceptional patient care never fades.

### Notice about PPE and Infection Control Supply Availability and Purchases

During these unprecedented times, we want you to know that we are committed to continuing to do what's best for you and your team, and the entire community of private dental practices that we serve, as well as for our employees.

We are monitoring this quickly changing situation and working closely with our manufacturer partners to continue to serve you. Unfortunately, due to global demand for PPE and other infection control supplies, orders for these supplies will not be guaranteed for delivery and are subject to backorder status, allocation of amounts and/or cancellation.


Orders placed cannot be cancelled or returned.

We have compiled additional resources and information relevant to the dental community at [www.goetzedental.com/covid19](http://www.goetzedental.com/covid19).

**Please let us know how we can help. Reach out to your private practice specialist or call Customer Service at 800-692-0804.**

# FOUR EASY WAYS TO ORDER

**GAINSpeed™**  
desktop ordering application



**GoetzeOnline.com**  
e-commerce ordering option



**800.692.0804**  
8am - 5pm CST



Fax:  
**800.552.7791**

