



Periodontics



Endodontics



Implantology



Implant  
Maintenance



Prophylaxis



Caries  
Detection

Practical clinical solutions,  
**effective workflows.**



Follow us on social media  
[@acteon\\_north\\_america](https://twitter.com/acteon_north_america)  
[www.acteongroup.com/us](https://www.acteongroup.com/us)



**Q1 2024  
PROMOTIONS**

# About Implantology workflow

ACTEON®'s implantology workflow leverages its expertise in different areas of the dental care pathway to help practitioners address some of the challenges associated with implantology to facilitate treatment, reduce treatment costs, increase surgical success, reduce chair time and improve the patient experience.



## Image Acquisition

I need to have access the high quality images to ensure my diagnosis is accurate




X-MIND® prime makes your image acquisition and diagnosis as easy as possible. High definition imaging combined with exceptional ease of use and ergonomics make it a powerful asset to any dental practice.



## Diagnosis

To succeed in my implant procedures, I need tools to properly assess my patient's morphology and perform effective surgical treatment planning



AIS 5 is a versatile software that includes many options to facilitate your procedures, such as bone density assessment, virtual planning tools and surgical guided design.



## Treatment

I need to make my treatment as efficient and minimally invasive and as possible for my patients



With Piezotome® Cube, you perform extractions at the same speed as rotary instruments with more safety and less trauma. You can ensure better socket preservation and minimize the need of bone grafting\*.

## IMPLANT SUCCESS

Bundle Package Includes  
**XMIND PRIME 3D/3D CEPH**  
**AIS**  
**PIEZOTOME CUBE**  
**PLUS EXTRACTION KIT**


**SAVINGS UP TO 33%**

**SCAN**

**1**



**x-mind  
prime 2**



**acteon  
imaging suite**

**2**

**PLAN**

**3**

**PREP**



**piezotome  
cube**

PROMOTION REDEMPTION INSTRUCTIONS IMPLANT SUCCESS SOLUTION includes all the technology described together: X-MIND Prime 3D, AIS Design Module Software, Piezotome CUBE, and Extraction tip kit. Offers cannot be combined with any other offers & substitutions are not permitted / one promotional offer per package. Promotional offers are good for purchases made between the dates of January 1, 2024 and March, 31 2024. Dealer invoice must be received by ACTEON prior to May 1, 2024. Promotion is void if purchased product is returned. Promotions only apply to new products and do not include demonstration or refurbished products.



# Image Acquisition

## PRIME + 1

Buy An X-Mind Prime 3D

Select

1 INTRAORAL IMAGING DEVICE

PRIME + 1  
**SPECIAL SAVINGS**



x-mind  
prime 2



1 Intraoral  
Imaging Device

## PRIME UPGRADE PROGRAM



x-mind  
prime 2

### Advanced Image Quality

High Resolution Sensor, up to 75µm resolution

### ACTEON Imaging Suite 5.0

Access All Of The Latest Features

### Enhanced Noise Reduction Filter & Reconstruction Algorithm

Easier Diagnosis

**Existing Customers Save Up To 50% off retail**

**Competitive Customers Making A Switch Save Up To 35% off retail**

## IMAGING ESSENTIALS

Bundle Package Includes

XMIND PRIME 2D or 3D

XM-Lite

SOPIX2

SOPRO 717 FIRST

**SPECIAL SAVINGS**



xm-lite



sopix 2



sopro717first



x-mind  
prime 2



# Image Acquisition

## SOPIX2 TRADE IN PROGRAM

# SAVE 40%

on a new **sopix2** sensor with trade in of any old or damaged Intraoral Sensor



**sopix 2**

**ASK YOUR  
ACTEON REP  
FOR MORE  
DETAILS**

### XM -Lite Sopix 2, PSPIX SAVINGS BUNDLE

Purchase 2 or more products for one low price.



**xm-lite**




**sopix 2**



**pspix2**

### Intraoral 3 Camera SAVINGS BUNDLE

Purchase 3 cameras & docks for the price of two



**soprocare**



**sopro717first**



**C50**

PROMOTION REDEMPTION INSTRUCTIONS Offers cannot be combined with any other offers & substitutions are not permitted / one promotional offer per package. Promotional offers are good for purchases made between the dates of January 1, 2024 and March 31, 2024. To redeem free goods, visit [acteongroup.com/us/promotions](http://acteongroup.com/us/promotions) and submit the free goods request form with dealer invoice prior to May 1, 2024. Extended warranty does not cover workstation PC. Promotion is void if purchased product is returned. Allow 4-6 weeks for delivery of product and free goods.





Treatment

# newtron


faster, more comfortable treatments

newtron  
booster



\$100 Rebate

newtron  
p5



\$150 Rebate

# newtron p5XS

\$200 Rebate

PERFECT FOR



Prophylaxis



Periodontics



Endodontics



Implant  
Maintenance

# newtron

total protect package

SAVINGS UP TO 30%

## BEST VALUE & CONVENIENCE

back to back procedures, no time to autoclave

NEWTRON B.LED

2 INTRA SCALING KITS

2 NEWTRON B.LED HANDPIECES

Contact Your ACTEON Rep for more information



x2



# greater variety of procedures

because there are over 60 unique tips

6 tip pack  
special savings

Purchase 6 tips  
for the price of 5



ACTEON tips



Tips	Newtron® Devices	POWER	IRRIGATION
------	---------------------	-------	------------

### PROPHYLAXIS

1 / 2 / 3 / 1S	14	
10P	14	
10X / 10Z	12	
EX1' / EX2' / EX3' / EXL' / EXR'	12	

### PERIODONTICS

H1 / H2L / H2R / H3 / H4L / H4R	2	
P2L / P2R	3	
TK1-1S	2	
TK1-1L / TK2-1L / TK2-1R	2	

### IMPLANT CARE

PH1 / PH2L / PH2R	2	
IP1	3	
IP2L / IP2R / IP3L / IP3R	5	

### ENDODONTICS

CAP1	10	
CAP2 / CAP3	10	
ET18D	10	
ET20 / ET25 / ET25S / ETBD	7	

Tips	Newtron® Devices	POWER	IRRIGATION
------	---------------------	-------	------------

### ENDODONTICS

ET20D / ET25L / ET40 / ET40D	7	
IRR20-21 / IRR20-25 / IRR25-21 / IRR25-25	6	
K10 / K15 / K25 / K30	6	
AS3D / AS6D	7	
AS9D	6	
ASLD / ASRD	7	
P14D / S12-70D	7	
P15LD / P15RD	7	
SO4	7	

### PROSTHESIS & ESTHETICS

PM1 / PMS1'	15	
PM2 / PMS2'	10	
PM3 / PMS3'	8	
PM4	15	
PMV1' / PMV2' / PMV3'	15	
PMV4' / PMV5' / PMV6'	10	
5AE / ETPR	20	
C20	11	



drop-by-drop



medium spray



high spray



no irrigation



# piezotome

## PRO PACKAGE cube

### Pro Tip Kits

Extraction, Bone Surgery, Sinus Lift, Intralift



#### PRO PACKAGE

##### PIEZOTOME CUBE & PRO TIP KITS

(Included Tip Kits Pictured, No Substitutions)

Extraction Kit, Bone Surgery Kit,  
Sinus Lift Kit, Intralift Kit

Contact Your ACTEON  
Rep for more information





# LESS TRAUMA, LESS DRAMA

cutting selectivity, minimal bone loss, power(60 watts vs 30 or less),  
range of tips, **superior outcomes**

**6 tip pack**  
special savings

purchase 6 tips  
for the price of 5

**extraction**  
**essential**  
**intralift**  
**sinuslift**

**bonesurgery**  
**piezocision**  
**crestsplitting**  
**crownextension**

**1 tip kit + 1 tip**  
special savings

purchase 1 tip kit & 1 tip  
for the price of 1 tip kit



II SL5 6



expasyl  
exact



**EXPASYL  
EXACT Pack**

4 box pack  
(20/box or 50/box)  
Buy 4 boxes for the price of 3

Promotion also includes traditional boxes of Expasyl  
(Regular 261030, Strawberry 261001)

**33% SAVINGS**



**pain free, effective gingival  
retraction that is easy & exact**

# riskontrol

**4 box bundle  
special savings**

purchase 4 boxes of tips for the price of 3

**adapter / box bundle  
special savings**

purchase 1 adapter & box of tips  
for the price of 1 adapter



**with single-use riskontrol, 20 times less  
chance of direct or cross-contamination risks.**

# air-n-go easy

**SPECIAL SAVINGS**

**air-n-go easy  
package**

air-n-go easy kit  
1 air-n-go perio powder

