

There's an A-dec for everyone.

A-dec offers creative, reliable solutions for practices of every size and specialty. Including yours. No matter which dental chair and delivery system you choose, you'll have peace of mind, knowing that they're designed to work the way you do—and to work without fail. Engineered with legendary quality and reliability, A-dec equipment is tested for 20 years of life, and backed with superior service and support. So you can practice with confidence.



A-dec 500

Premier comfort, performance, and connectivity

Bring intelligent simplicity and infallible peace of mind to your operatory, and superior comfort to your patients.

A-dec 500 Chair

- Exquisitely smooth and quiet, cradles patients with whole body support
- Impressive range of motion, including our lowest base-down position
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for optimal access and healthy, natural posture

A-dec 500 Pro Delivery

- Six adjustable control head positions for customisable placement
- DS7 touchscreen offers a simple, dynamic interface for intuitive, connected control
- Multi-user settings can be individually programmed and personalised
- Enabled with A-dec+ software, to easily add more services and features as technology evolves



A-dec 400

Contemporary luxury with proven reliability

Sleek design, patient comfort, and uninhibited dental team access—all combined with reliable performance.

A-dec 400 Chair

- Distinctive contemporary style with contoured patient comfort
- Smooth, efficient positioning, plus 30° rotation on each side of center
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for excellent access and healthy, natural posture

Delivery Options

- Pairs with A-dec 300 Pro delivery
- Choose from Continental® or Traditional deliveries
- Radius® option rotates easily around the chair to support both left- and right-handed dentistry



A-dec 300

A small footprint with full-size performance
Accessible, flexible, and evolved to bring you rich
feature sets and the option of digital connectivity.

A-dec 300 Chair

- Compact design fits a variety of operatory spaces, with minimal maintenance
- Contoured, pressure-mapped surfaces support patients in comfort
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for excellent access and healthy, natural posture

A-dec 300 Pro Delivery

- Open platform design is easy to configure, customise, and upgrade in the future
- Pivoting control head reduces hand strain and fatigue
- Intuitive, responsive user interface shows just what you need to know in the moment
- Add the optional A-dec+ Gateway to access the A-dec+ software platform

Chair	A-dec 500	A-dec 400	A-dec 300
Ultra-soft start/stop hydraulic drive	•	_	
Hybrid drive (hydraulic base/electric tilt)		•	•
Swivel (30° each side of center)	•	•	_
Virtual pivot for added patient comfort	•	•	•
Thin backrest 1.25" (32mm)	•	•	•
Dual-articulating gliding headrest	Lever release	Locking knob or lever release	Locking knob
Multi-position armrest settings	2	2	2
Chair at lowest point	13.5" (343mm)	13.75" (349mm)	13.75" (349mm)
Chair at highest point	31.5" (800mm)*	31.5" (800mm)	31.5" (800mm)
Integrated 300 watt power supply	•	•	•
Load rating	500 lb. (227 kg)	400 lb. (181 kg)	400 lb. (181 kg)
Choice of sewn or formed upholstery	•	•	•

Delivery System	A-dec 500 Pro**	A-dec 500	A-dec 300 Pro**
Chairs	A-dec 500	A-dec 500	A-dec 500, A-dec 400, A-dec 300
Mounting options	Radius	Radius and 12 o'clock	Radius, side, and support center
Control pad	DS7	Deluxe Plus	CP5, CP5i
A-dec+ connectivity	•	_	Optional (requires A-dec+ Gateway)
Four-position control block	•	•	Optional
Three-position control block	_	_	•
Holder positions	6	6	5
Clinical device integration	5	5	2
USB pass-through port	•	•	(on deliveries with CP5i)
Support center	A-dec 400	A-dec 400	A-dec 400, A-dec 300
Self-contained water system with quick-disconnect 2-litre bottle	•	•	•
Accommodates amalgam separator	•	•	•
Cuspidor swivel	•	•	A-dec 400 only
Programmable timing for cup and bowl	•	•	•
Dental light	A-dec 500 LED	A-dec 500 LED	A-dec 500 LED, A-dec 300 LED

^{*}Using maximum upper limit. The factory set upper limit is 29" (737mm). ** Availability varies by market

Reliable Creative Solutions

FOR DENTAL PROFESSIONALS WORLDWIDE













More than a warranty. It's a commitment.

At A-dec, our focus is on making dentistry better—which includes superior customer service along with superior dental products. Your A-dec 500 equipment is backed by a five-year warranty, and a lifetime of customer support. We'll be here today, tomorrow, and well into the future. That's our commitment to you.

Any colours and products shown may vary and are subject to change without notice. Contact your authorised A-dec dealer for the most current product information. All trademarks in this document are subject to the trademark rights of A-dec, its affiliates, and their collective suppliers in the United States and foreign countries.

©2023 A-dec Inc. All rights reserved. 85.0394.25/E-only/11-21(9-23)/Rev C3

A-dec Headquarters

2601 Crestview Drive Newberg, Oregon 97132 USA P: 1.800.547.1883 within USA/CAN P: +1.503.538.7478 outside USA/CAN a-dec.com

A-dec Australia

P: 1.800.225.010 within AUS **P:** +61.(0)2.8332.4000 outside AUS

A-dec China

P: 400.600.5434 within CN **P:** +86.571.89026088 outside CN

A-dec United Kingdom

P: 0800.233285 within UK **P:** +44.(0).24.7635.0901 outside UK

