

# PRO▶VECTA® HD SPECIAL



PN A6350-15S  
PN A6350-24S  
PN A6350-33S  
PN A6350-15D  
PN A6350-24D  
PN A6350-33D

**Purchase a  
ProVecta HD  
Intraoral  
X-ray  
for just  
\$4,100**

**Save  
Up To \$1,100**  
off retail

Promo Code: Q119PVHD

## PERFECT IN FUNCTION & DESIGN

- Factory programmed dosage specifically set for each view when using ScanX® flexible phosphor plates (PSPs)
- DC technology allows constant radiation levels to ensure reliable image quality
- Choose from short (66"), medium (72") and long (84") extended arm lengths
- Right-angle/parallel technology for precise alignment
- Adjustable mA and kV levels
- Optional rectangular radiation field collimator
- Generator head handle allows precise positioning that stays in place