

Invest in the Air Techniques Advantage

600%

MORE CAPACITY THAN THE LEADING BRAND

mL by volume

1500
1200
900
600
300
0

Leading Competitor

Acadia Plus



- Save time and money by replacing your filter less often
- Ideal for busy dental practices

BETTER FLOW=BETTER PERFORMANCE

- Increase your vacuum's performance and get better flow from your operatories
- Using the larger 2" connections with existing 2" plumbing lines can increase your performance (CFM) by 57% over the maximum 1.5" connection on competitor models
- Don't let your amalgam separator limit your utility equipment

57%
MORE
CFM

ACADIA PLUS

LEADING
COMPETITOR

30%
LESS
EXPENSIVE

SUPERIOR VALUE

- When comparing the overall cost of the filter to its capacity, the Acadia Plus is 30% less expensive than the leading competitor.
- If you have a large dental practice, this can result in a substantial yearly cost savings!



For more information, to schedule a demo, or to place an order, contact your local dealer representative. To learn more, including current promotions, scan the QR code, or visit www.airtechniques.com

**AIR
TECHNIQUES** equipped for life®