

Operatory Seating

Selection and Specification Guide

9500 Series

The evolution of ergonomic comfort continues with the state-of-the-art in the fight against back pain.

- Exclusive dynamic motion technology in the seat and backrest
- New, low-profile, wrap-around lumbar support
- Easily adjust desired level of synchronized seat with HybriGel™ and backrest tilt to help maintain neutral posture
- Independently adjust backrest to allow greater range of movement
- Standard height adjustment of 19.5" - 25.5" (lower & higher heights available)
- Assistant model features 25" - 33" pneumatic height adjustment, ratcheted body support (left or right) with HybriGel™ and adjustable foot ring



Model 9500B

Model 9520BL

Notes

9200 Series

Complete ergonomic support and adjustability with standard size seat and smaller backrest.

- Sculpted backrest and large (19" W x 16.5" D) seat with HybriGel™
- Continuous, synchronized seat and back tilt
- Independently adjust backrest to allow greater range of movement
- Standard height adjustment of 19" - 25" (lower & higher heights available)
- Assistant model features 24.5" - 32.5" pneumatic height adjustment, ratcheted body support (left or right) with HybriGel™ and adjustable foot ring



Model 9200BST

Model 9220BLST

9100 Series

Our number one seller! Designed for superior ergonomic support. Great for smaller operatories.

- Sculpted backrest and smaller (19.5" W x 15" D) seat
- Continuous, synchronized seat with HybriGel™ and back tilt
- Independently adjust backrest to allow greater range of movement
- Standard height adjustment of 19" - 25" (lower & higher heights available)
- Assistant model features 24.5" - 32.5" pneumatic height adjustment, ratcheted body support (left or right) with HybriGel™ and adjustable foot ring



Model 9100B

Model 9120BL

9000 Series

Plush design with superior ergonomic features.

- Large (19.5" W x 15" D) backrest and cut out seat with HybriGel™
- Continuous, synchronized seat and back tilt
- Independently adjust backrest to allow greater range of movement
- Standard height adjustment of 19" - 25" (lower & higher heights available)
- Assistant model features 24.5" - 32.5" pneumatic height adjustment, ratcheted body support (left or right) with HybriGel™ and adjustable foot ring



Model 9000BST

Model 9020BLST

Operatory Seating

Selection and Specification Guide

New **135AS & 135JS Saddle Series**

Introducing two extensions of our 360° Dynamic Motion technology – the Amazone Series & the Jumper Series.

- Exclusive 360° Dynamic Motion technology helps strengthen back muscles, reduce back pain & minimize sitting fatigue
- Easily adjust seat height and level of dynamic motion desired with convenient control ring for optimal personal comfort
- Amazone and Jumper Ultraleather upholstery available in Dove Gray, Raven Wing & Papyrus
- Amazone is available in lower and standard heights (21"-27" & 24"-32") with a narrower upholstered seat design
- Jumper is available in standard and higher heights (24"-32" & 27"-37") with a wider upholstered seat design



*Model 135AS
Amazone*

*Model 135JS
Jumper*

Notes

135DSS Split Saddle Seat

Innovative, dynamic split seat technology promotes back and circulatory health.

- Patented split seat design offers superior support and cooling comfort
- Unique 360° circular movement can be customized to personal preference
- Standard height adjustment of 22.5"-29.75" (lower & higher heights available)



Model 135DSS

3300 Series

Traditional style seating offers ergonomic support.

- Large (19" W x 17" D) contoured, waterfall seat and fully adjustable, contoured backrest
- Synchronized seat and back tilt
- Standard height adjustment of 18.5"-24.5" (higher height available)
- Assistant model features 23"-31" pneumatic height adjustment, ratcheted body support (left or right) and adjustable foot ring



Model 3335B

Model 3345BL

3100 Series

Traditional style seating offers basic ergonomic support.

- Large (19" W x 17" D) contoured, waterfall seat and contoured, dual adjustable, tilting backrest
- Standard height adjustment of 18"-24" (higher height with foot ring available)
- Assistant model features 23"-31" pneumatic height adjustment, ratcheted body support (left or right) and adjustable foot ring



Model 3125B

Model 3145L

The Power to Advance. Safely. Financially. Competitively.

www.brewercompany.com | P: 1.888.273.9371

The Brewer Company is an FDA registered & ISO 13485:2003 certified company.

LIT-D-006 Rev B