

*DCI*EDGE



DCI EDGE SERIES 4 OPERATORY PACKAGE

BLENDING STYLE AND VALUE

The Series 4 operator package combines design simplicity with style at a price that fits any budget. At DCI, we are committed to delivering easy-to-use, high-quality, high-value products you can count on.

Swing Mount Delivery System

Delivery system and light easily convert from left to right handed operation.





Post Mount Delivery System

The post mount delivery system features an integrated vacuum trap, water system, as well as air and water quick connects.




Optional chair control touchpad

RELIABLE INSIDE & OUT

For over 30 years, DCI has been the industry's trusted source for the most reliable dental parts on the market. Creating long lasting, easy to maintain dental equipment defines our culture. As a result, DCI is the preferred brand among service technicians, setting the standard for reliability.

Low Maintenance

Simple design with high quality components improve reliability and ease of service.



Cuspidor

Durable metal construction with porcelain bowl.

Rear Assistant's Package

Mounts to rear of chair, includes chair touchpad.



Telescoping Assistant's Package

Featuring the most reliable, ergonomic syringes and vacuum valves on the market.



Hygiene Delivery System

Everything needed for hygiene in one place. Integrated vacuum lines provide a sleek, modern look. Available in multiple configurations.

CONFIGURED YOUR WAY

The Series 4 delivery system is offered in multiple mounting options to match your clinical needs today and as your practice grows.



Side Delivery System

The DCI Edge Side Delivery tucks cleanly out of the way. Available in cabinet or wall mount. Water bottle system included.



Series 4 Cart System

Sturdy metal frame on adjustable cart base comes standard with DCI self-contained water system and junction box with utility center.



Operatory Support Cart

Featuring a durable wood top and high quality assistant's instrumentation. Available with optional storage organizer.

GET COMFORTABLE, STAY HEALTHY

The DCI Edge operatory chair provides comfort, style and efficiency. Our engineers have designed an ergonomic chair that combines a thin, narrow back with a low chair height. This allows you to be in the proper working position regardless of patient size.



ErgoBack
Plush



Narrow Back
Asepsis



Narrow Back
w/ Slings

Doctor's and Assistant's Stools



Headrest Scoop

Headrest and back design support patients of all statures



Thin/Narrow Back

Thin, narrow back design provides comfortable, ergonomic access to the oral cavity




Foldout Armrests

Armrest rotation facilitates easy entry and exit of patients



Heavy-duty Construction

High quality steel and aluminum castings stand the test of time



FRESH, NEW DESIGN

DCI's Series 4 operator light was engineered utilizing reflective LED technology to deliver next level lighting performance. The brand new design blends innovation with beauty on every level, redefining what you should expect from an operator light.



Advanced Features


3rd Axis Rotation for improved positioning.

Dual Membrane Integrated Touch Controls allow you to easily access and control three intensity settings of 17,500, 26,000 and 35,000 lux.



Durable Design

DCI's heavy-duty metal structure and handle design along with high quality components demonstrate a strong commitment to building the highest quality dental equipment.



Mounting Configurations



Track Mount Light



Cabinet Mount Light



Post Mount Light



Ceiling Mount Light



Started over 35 years ago, DCI was founded on three guiding principles: making our customers' lives easier, providing dental care for those in need, and enhancing the lives of our employees. These principles give us our edge. It's an edge we want to pass along to you.

DCI is a family owned company for two generations which is based in Newberg, Oregon, the heart of dental equipment manufacturing in America. Over the years, John & Jason Spencer redefined the dental landscape by designing dental equipment that challenged the status quo. Today, the family brings their design expertise and commitment to value and service to the forefront with the DCI Edge Series 5 Operator Package.



Jason, Amy and John Spencer

Upholstery Chart

Nauga Soft



Ultra Leather

