

DCI **EDGE**



DCI EDGE SERIES 5 OPERATORY EQUIPMENT

CELEBRATING
40 YEARS
OF PASSION
AND PROMISE

DCI is a multigenerational family-run company based in Newberg, Oregon. We believe in managing every aspect of the business around one simple goal: **making our customer's life easier.**

This focus has made DCI Edge the fastest growing dental equipment brand in North America. Our passion is designing and building the best dental equipment in the industry. Our promise is to manufacture with an unwavering focus on reliability and optimal value.



CRAFTED WITH
DECADES OF
EXPERIENCE

DCI edge



CHAIR CITRUS ULTRALEATHER™ UPHOLSTERY

844-854-8555 | dciedge.com




FINANCIAL FREEDOM

It's your practice. It's your budget. Keep your options open. DCI Edge delivers quality equipment at a price point that allows you to optimize the design of your operatory without compromise. This cost savings empowers your financial freedom.



REINVEST

Upgrade your practice with the latest and greatest technology, and secure your edge in the ever-changing tech world.



SAVE

There is peace of mind that comes with money in the bank. Save for a rainy day or for your retirement plan.



CHOICE


You choose how to balance your funds to make your life easier.

GET COMFORTABLE STAY HEALTHY



ERGONOMICS

Thin, narrow back and dual articulating headrest allow for ideal positioning for the provider and patient. This leads to enhanced comfort and efficiency for the entire dental team.



SIT OR STAND?

The Series 5 chair height ranges from 16" to 31" at the seat so you can comfortably practice sitting or standing.

ACCESS TO THE ORAL CAVITY

Unique contoured back design ensures full access between the 3 and 9 o'clock working positions while comfortably supporting the patient.



SERIES 5 CHAIR TIN MAN PEARLIZED ULTRALEATHER
UPHOLSTERY



ARE YOU IN THE RIGHT POSITION?

Optimal positioning near the patient ensures a healthy posture, reducing pain and fatigue.



COMFORT THAT KEEPS YOUR PATIENTS SMILING



CHAIR & STOOL APRICOT ULTRALEATHER UPHOLSTERY



ERGOBACK PLUSH

PATIENT SUPPORT

The innovative headrest scoop and back design provide patients with back support regardless of height.

SOFT START AND STOP

Adjust your chair without disrupting your patient with smooth start and stop technology.

DESIGN APPEAL

Available in ergoback or narrow back, plush or asepsis, and in a large variety of colors. Upholstered armrests add to the comfort and style.

MULTI-LAYER SYSTEM

Four layers of specialized foam with varying levels of density deliver premium comfort and support. Adapts to all body types.




NARROW BACK ASEPSIS



NARROW BACK PLUSH



ERGOBACK ASEPSIS



BEAUTIFUL INSIDE AND OUT




 CHAIR CITRUS ULTRALEATHER UPHOLSTERY



RELIABLE AND BUILT TO LAST

Our design philosophy blends innovation with high quality construction. We are proud to show you what's under the hood.



PERFORMANCE MEETS DESIGN



EXCLUSIVE EXTENDED WARRANTY


Ask about the DCI Edge exclusive NSK pricing and extended warranty on a variety of air and electric attachment bundles.



INTEGRATED AND READY TO GO

Easily access essential clinical and diagnostic devices such as the Dentsply Sirona Cavitron system, the ACTEON® SOPRO family of intraoral cameras, piezoelectric scalers, and curing lights.

DUAL INTEGRATED TOUCHPADS AND BRAKE HANDLES



ERGONOMIC ACCESS TO INSTRUMENTS

Additional head pivot swivels to improve access and reduce unnecessary movement.



Lights On/Off



Return To Last Position



01 Programmable Buttons for Preset Positions



Backrest Incline and Recline
Chair Up and Down

INTERNAL COMPONENTS

Reliable pneumatic components and ample room for high tech integration. As the largest manufacturer of dental parts in the world, DCI is known for quality parts and components.



INTEGRATED TECHNOLOGY

Experience award-winning design featuring WAVEONE™ technology. The NSK NLZE Electric Motor is engineered with superior power, greatly expanding dental treatment possibilities.



INNOVATIVE REFLECTIVE TECHNOLOGY INCREASED VISIBILITY



INCREASED VISIBILITY

DCI delivers cutting edge illumination technology by utilizing over 500 reflective facets to bounce light back to the oral cavity. The result is a uniform light pattern.

LESS EYE FATIGUE

Optimal light performance reduces shadows, improving illumination of the oral cavity.

ACCURATE SHADE MATCHING

5000 Kelvin color temperature with a color rendering index (CRI) of 90 produces a neutral white light, which is ideal for color matching at any intensity level.

NO-CURE SETTING

Filters out UV light to prolong the working time of photo-initiated resins.



PROGRAMMABLE NO-TOUCH SENSORS

Control the light without contact.



INTEGRATED TOUCHPADS



No-Cure




Intensity



On/Off

ACCURATE SHADE MATCHING 5000 K




MOUNTING CONFIGURATIONS



SYSTEM SOLUTIONS

MAXIMUM FUNCTIONALITY



SWING MOUNT DELIVERY SYSTEM

Delivery system and light easily convert from left to right handed operation. The compact swing radius optimizes space making it perfect for operatories of any size.





SHOWN WITH CUSPIDOR



POST MOUNT DELIVERY SYSTEM

The post mount delivery system features an integrated vacuum trap and water system, along with air and water quick connects.

SLEEK CABINETS DESIGN TO COMPLETE YOUR OPERATORY

Designed for optimal organization and efficiency, DCI Edge Cabinetry combines the work and storage space you need with showroom aesthetic appeal. Purposeful and streamlined construction at a price point intended to offer the best value for your budget.



CHAIR & STOOL
SOLID SURFACE
CABINETS

SKY ULTRALEATHER UPHOLSTERY

MORNING ICE

UPPER : ITALIAN SILVER ASH, LOWER: STERLING ASH





DCI EDGE CABINETS

QUALITY, RELIABILITY, AFFORDABILITY

HIGH QUALITY CONSTRUCTION

Structural steel support with a 6-inch subbase provides durability while lifting wooden parts off the floor and away from moisture. Thermal-fused edge-banding make doors and drawers water resistant and protects edges from damage.

PRECISION MANUFACTURING

Designed with cutting edge software, and cut with state-of-the-art CNC Machinery, DCI Edge cabinetry represents quality at the highest level.

RELIABLE COMPONENTS

Engineered with top-of-the-line components that stand the test of time.

FUNCTIONAL AND AFFORDABLE

Built to support all operatory needs, including units and lights, at a competitive price point.

ATTENTION TO DETAIL

Enjoy detailed features that make routine tasks more efficient: soft-close doors and drawers, undermount sinks that are easily cleaned after each patient, powder-coated steel sub-bases that stand up to repeated cleaning.


COMPLETE SOLUTIONS

From side cabinets to 12 o'clock workstations, center island consoles to sterilization centers and lab casework, DCI Edge cabinetry offers complete solutions to meet your individual needs.



CHAIR & STOOL
SOLID SURFACE
CABINETS

CITRUS ULTRALEATHER UPHOLSTERY
SILVER SMOKE
UPPER & LOWER: MAGNOLIA, ACCENT: NIGHTFALL





BEAUTY MEETS FUNCTIONALITY

With DCI Edge cabinetry, aesthetics are not compromised for utility. With a designer's attention to visual appeal, our cabinetry provides ergonomic access to materials and instrumentation, ample workspace to accommodate any clinical procedure, and optimal storage space. Additionally, the units are easily adjustable for left or right capability.



36" 12 O'CLOCK CABINET
WITH UNDER COUNTER MOUNT
DELIVERY SYSTEM



42" 12 O'CLOCK CABINET
WITH CABINET MOUNT
DELIVERY SYSTEM



42" 12 O'CLOCK CABINET
WITH REAR FLOOR-MOUNTED
DELIVERY SYSTEM



DESIGNED FOR DENTAL ORGANIZATION AND EFFICIENCY

FLEXIBLE DESIGN

DCI Edge provides modular customization to accommodate even the most challenging spaces for optimal operatory organization.

TECHNOLOGY FRIENDLY

Tech-friendly modules allow for easy and efficient computer and technology integration. Designed with paperless offices in mind.

SOLUTIONS FOR PRODUCTIVITY

Adjustable pull-out ergo shelves, tip-out bins, flexible monitor mount, integrated electrical and USB ports, and colored drawer organizers ensure that everything you need is efficiently placed within reach.



INTEGRATED DESIGN YOUR CHOICE OF OPTIONS

Do it your way! Choose the configuration and carefully considered options that empower you and your Assistant to organize efficiently without leaving the chair.

REAR FLOOR-MOUNTED DELIVERY UNIT FEATURES

- Accommodates ancillary integration, now or in the future
- Mobile surface has an adjustable height range of 5"
- Unique design accommodates easy positioning for left or right-handed practitioners
- 6" side to side movement without encroaching on operator space
- Alignment of collection trap and water bottle enable convenient service position for Assistant
- Assistant touchpad pivots 180 degrees
- Work surface tray holder provides additional space for supplies



12 O'CLOCK FREE STANDING CABINET FEATURES

- Assistant's work surface can increase effortlessly to accommodate any clinical procedure
- Colored drawer organizers optimize storage space
- Ergonomic access to tubs, trays, cassettes, materials, and instrumentation
- Upper doors feature textured glass with LED lighting
- Pass-through upper with textured glass inserts
- Mid-section power station with tip out bins. Two access doors with hidden duplex outlets for clean cord management
- Swivel steel powder-coated ergo shelves to hold tubs and trays, or act as a work surface. Robust, die cast arms, that hold up to 10 pounds
- Standard 24" or extended 30" upper cabinets
- LED task lighting and clock timer
- Dual glove dispensers
- Removable back panel with free standing 12 o'clock
- Lower center module with CPU storage, bi-fold door and adjustable shelf
- Lower Assistant's module side opening bi-fold door and pull-out drawers
- Easy-glide monitor mount can be installed at top or bottom of upper cabinet for convenient viewing
- Optional nitrous flow meter cut-outs available
- Accent colors on mid-section



CABINETS STREAMLINED FOR STERILIZATION



FLIP DOWN WRAP STATION

Two conveniently spaced drawers designed to pull out and flip down, ensuring ample space to wrap cassettes of any size.



SYSTEMATIC DECONTAMINATION

Optimize workflow with specialized waste containers, stainless steel undermount sink, dual integrated drain boards and conveniently placed glove dispensers.



FACTORY INTEGRATED ACCESSORIES

Our designers and engineers have curated a list of the best options in dental sterilization. Factory installation provides optimal placement and smooth workflow.



STERILIZATION TOWERS BUILT TO LAST

Built with waterproof materials that guarantee protection from water damage. Multiple tower options accommodate your specific sterilizer selections.



SIMPLE STERILE ORGANIZATION

Pull-out storage with wire racking holds numerous cassettes, tubs, trays and wrapped instruments. Adjustable shelves allow for fast and easy management.



SPECIALIZED LIGHTING FOR VISUAL WORKFLOW

Red lights indicate contaminated storage while blue lights illuminate sterile storage. Additional task lighting brightens the entire work surface for maximum visibility.



CABINETS STREAMLINED FOR WORKFLOW



RECEIVING, CLEANING & DECONTAMINATION




PREPARATION & PACKAGING



STERILIZATION



STERILE STORAGE & MAKE READY



A WORKSPACE BUILT JUST FOR YOU

Sterilization is one of the most important steps in a practice, and one size never fits all. With DCI Edge, you can create a Sterilization Center that fits your unique needs. All options feature a process-based flow that keeps productivity high and risk low.

STERILIZATION WITH STANDARDS

Designed to follow CDC guidelines with exceptional focus on infection control and patient safety.


STERILIZE WITH STYLE

Fully functional clock timer enhances productivity, while stainless steel valance adds to the beautiful design aesthetic.


STRAIGHT RUN




GALLEY STYLE



L SHAPE



U SHAPE





EXPERIENCE DCI EDGE

3 SHOWROOMS

2 MANUFACTURING FACTORIES



DCI Edge Showroom & Manufacturing Facility Oregon, USA
705 S Springbrook Road
Newberg, OR 97132

DCI Edge Showroom and NSK America Corporate Headquarters
1800 Global Parkway
Hoffman Estates, IL 60192

DCI Edge Showroom & Manufacturing Facility North Carolina, USA
9201 Stockport Place
Charlotte, NC 28273

[dciedge.com](https://www.dciedge.com)
equipment@dcionline.com
1-844-854-8555



50169/February2023-001
©2023 DCI

