

## DentalEZ® Delivery Units: Reliable and Simple

Our delivery units help you create an operatory that is organized and functional, placing everything you need at your fingertips for maximum productivity. Their sleek, modern designs enhance the appearance of the dental office.

### CORE™ Delivery Unit

- Secure tray placement
- Air gauge for pressure monitoring
- Integrated holders
- Easy access control block and internal components for simple maintenance reducing downtime
- Controlled by flex arm for stability and ease of positioning
- Optional integrated touchpads
- Magellan or Fixed Mount configurations available



### Rear Assistant's Instrumentation


- Centrally located solids collector
- Optional integrated touchpad
- Simple footprint and aesthetically pleasing
- Redesigned for 2017

## Delivery Unit Mounting Options



### Fixed Mount

CORE Chair, CMU with Assistant's Instrumentation & Cuspidor





### Magellan Mount

CORE Chair with Magellan Unit

## Galaxy® Delivery Units / Simplicity® Delivery Units

- Units available in fixed chair-mounted or Magellan left/right mounted
- 3 styles of control heads designed for up to 4 handpieces (3 standard) and 2 static hangers (1 standard)
- Galaxy® heads standard with die cast aluminum pinch valve system providing years of reliable performance
- Clean water system is standard, city water kit available
- Galaxy with pinch valve system (A)
- Galaxy Euro style with pinch valve system (B)
- Simplicity® style with control block system (C)
- Simplicity heads offer a sleek modern style to enhance the appearance of any dental office



## Assistant Instrumentation/Cuspidors:

- Galaxy rear-mounted assistant's instrumentation & cuspidor (D)
- Simplicity cuspidor with assistant's instrumentation (E)
- Simplicity telescoping assistant's instrumentation (F)
- Simplicity rear-mounted assistant's instrumentation (G)



## Delivery Unit Add-Ons:

- The Ultrasonic Integrated Scaler has an electronic module that is fully integrated into the DentalEZ dental treatment unit. This saves valuable countertop space in treatment room (1)
- Optional Integrated Sopro Camera handpieces can be connected and disconnected very easily (2)
- Our LED Curing Light has a durable anodized aluminum housing that is self-cooling, silent, shock-resistant and easy to clean (3)
- Piezo Scaler is a new generation multi-purpose ultrasonic module that offers a wide range of power for flexibility (4)



For more information, please visit [www.DentalEZ.com](http://www.DentalEZ.com)