

DentalEZ Lights:

LED Lighting Made Easy

Providing ideal lighting for every dental procedure is a critical aspect of dentistry. Our dental lights enhance visual acuity in the operator by providing the precise type, color and intensity of light required for every procedure. In addition, our lighting is designed to provide dependable performance, simple operation, easy cleaning and maintenance, and exceptional value in the dental practice.



CORE™ LED Light

- Rectangular light focus
- Magnetic removable and autoclavable handles that can be adjusted 360°
- Hands-free adjustable light intensity to 30,000, 20,000 and 10,000 lux
- Programmable sensor for three brightness settings
- 5500° – 6300° Kelvin color temperature
- Removable lens for easy cleaning
- Available in a variety of mounting options



EverLight®

- Improved CRI versus the original model (+23% better CRI)*
- Nine color and intensity settings, including color matching
- Available in six mounting configurations
- Three Kelvin color temperatures, three lux intensity settings
- 12 Philips® LEDs maximize light output
- Focused light pattern provides proper color-corrective light for most procedures
- Reduces energy costs up to 70% versus halogen lights
- LED lights last 10 times longer than halogen lights

**Based on DentalEZ Testing*



Simplicity® Halogen Light

- Full focus from 20 to 40 inches from the shield allows the dental team to reposition the chair without adjusting the light
- Light pattern integrity provides consistent illumination even with obstruction by the dental team
- 3-axis adjustment allows for a wide range of positioning
- Aluminum casting light frame for better aesthetics and durability
- Ergonomically designed handles for ease of movement



For more information, please visit www.DentalEZ.com