

Simplicity™ LED

LED Lighting Made Easy

- Smooth rotation and 3 axis adjustment allow for a wide range of light positioning to better view the oral cavity
- Control manually or with hands free sensor
- Ergonomic handles for comfort of operator
- Cast aluminum light head for better aesthetics and durability
- Low, medium and high settings for flexibility from 18,000, 27,000 and 36,000 lux
- Rectangular light pattern (6½" x 3½") lights the oral cavity with minimal shadows
- Color temperature ranges from 5,000 to 5,700 Kelvin
- 5 Year Limited Warranty



Upgrade your existing
Simplicity halogen light to LED
with our retrofit kit!



Ordering Information:

Simplicity LED Retrofit Kit	Part #3802-581
Chair/Post Mounted	Part #3544-234
Ceiling Mounted 8'-8'6"	Part #3544-238
Ceiling Mounted 8'6"-10'	Part #3544-237
Wall Mounted	Part #3544-239
Cabinet Mounted	Part #3544-240
Track Mounted 8'-8'6"	Part #3544-241
Track Mounted 8'6"-10'	Part #3544-242
Track Mounted Dual Heads/Trolley 8'-8'6"	Part #3544-243
Track Mounted Dual Heads/Trolley 8'6"-10'	Part #3544-244



Please visit DentalEZ.com/SimplicityLED for more information.