

# HIGH DEFINITION IRIS HD

Introducing  
**IRIS HD USB 3.0**  
Intraoral Camera.

**Experience the Clarity.**

1280 x 1024  
High Definition (HD)

640 x 480  
Standard

The new IRIS **HD** USB 3.0 brings the smallest details into brilliant full-screen high-resolution focus. With dual capture buttons and one-touch focus, control is always within reach. IRIS HD even knows to turn itself on and off when you pick it up or set it down.

Extreme clarity in HD is obtained with the IRIS HD precision optical lens array and Sony's advanced HD sensor for the finest picture quality.



A touch of the button and the newly designed electronic focus will advance automatically to the next focus mode, all five positions are shown onscreen.

Equipped with a motion sensor, it will Auto-On when lifted and Auto-Off after ten minutes of inactivity. The streamlined soft-tip has brilliant 8-point LED lighting and a new tapered USB 3.0 connector that fits in a standard delivery unit.

See why Digital Doc is the #1 choice in intraoral cameras. Call 1-800-518-1102 or visit us at [www.digi-doc.com](http://www.digi-doc.com).



© Digital Doc LLC, 4789 Golden Foothill Parkway,  
El Dorado Hills, CA 95762. The Digital Doc logo  
and IRIS HD logo are trademarks of Digital Doc LLC.  
All rights reserved.

