


**AT FOREST,  
OUR MISSION IS TO CREATE  
STREAMLINED, INNOVATIVE  
AND IMAGINATIVE  
TREATMENT ROOMS.  
OUR SOLUTIONS PRIORITIZE  
THE HAPPINESS  
AND COMFORT  
OF PATIENT, STAFF  
AND DENTIST.  
WE ALWAYS LOOK FOR  
A BETTER WAY TO DO THINGS,  
EVEN IF IT MEANS BREAKING INDUSTRY NORMS.  
WE DELIVER CHOICES,  
COLORS AND EVERY DETAIL  
DONE JUST RIGHT.**







Rear Mount Side Delivery




## Forest Dental 7000 Series

We are very proud of our 7000 Rear Mount, Side Delivery Series because we believe when compared feature for feature with any other manufacturer, it is the highest quality delivery system made in North America. Once you realize the benefits of the Rear Mount, Side Delivery concept relative to chair-mounted deliveries, you will quickly see why Forest's 7000 Delivery System is rapidly increasing in popularity.



## Key Reasons Behind the Growing Popularity of Rear Mount, Side Deliveries

**Historical Perspective:** The handpiece tubing limitations that originally inspired the design of chair-mounted delivery systems have been virtually eliminated. Handpiece lubricant used to cause tubing to stiffen over time, creating an irritating torque on the back of the handpiece. This problem could only be mitigated by keeping the tubing between the back of the handpiece and the control head as straight as possible. As a result, the control head was positioned over the patient, close to the mouth, when using a chair mounted whip arm system or beside the patient with either a chair or wall mounted delivery. Because tubing no longer hardens in reaction to lubricants and swivels eliminate the remaining torque, handpieces can now effortlessly come from beside the patient's head in what is called a Rear Mount, Side Delivery. We strongly feel that our Forest 7000 Series is, feature for feature, the best available on the market today.



Lifetime Warranty



Easy Access



Level Foundation



Handpiece Lockouts



Fine Hardware



Metal Plugs

## The 7000 is Our Fastest Growing Delivery Concept Because of Three Key Reasons:


- 1. Efficient.** Handpieces, vacuum devices and consumables are more accessible regardless of where you're positioned around the chair. Dentists new to the Rear Mount, Side Delivery concept appreciate the efficiency of not having to reposition a chair-mounted control for patient entry and exit or as they move around the mouth. To meet individual styles of instrument and material management, the Forest 7000 can be ordered with up to four cassette/tub holders, three storage organizer bins and a monitor mount. Unique to the 7000 are handpiece lock-outs, making it the only control on the market that allows a burr change while another handpiece is in use—the key feature enabling true four-handed dentistry.
- 2. Cost Effective.** Because the 7000 has one mount and one set of plumbing, it costs less to install than a chair-mounted handpiece delivery with assistant's work surface, which requires two mounts and plumbing locations.
- 3. Patient Experience.** By eliminating the cluttered "Christmas tree" effect of over-the-patient controls, the 7000 looks less overwhelming, which in turn creates a more patient-friendly environment.

# INTRODUCING THE 7020 PRO SERIES.


NO MATTER HOW YOU SPIN IT, IT'S ALL  
ABOUT REAR-MOUNT SIDE DELIVERY.

“The ergonomics of the rear-mount side-delivery system are especially important because improper sitting and reaching is what shortens our careers. With this system, you’re sitting in the proper position at 11 o’clock, bent from the hip, and you don’t have to take your eyes away from the patient or reach. Everything is where it needs to be at all times.”

- Brad Boeke, DDS




XX






**CONSIDER THE SPINE**

**DENTISTRY IS DESIGN CENTRIC.  
THE 7020 PRO SERIES EMBRACES  
ERGONOMICS AND BEAUTY...**









**WHY FOREST?**





**EVERY DETAIL  
DONE RIGHT.**



WE KEEP IT TIGHT



# 7020 PRO SERIES TOP DOWN FEATURES

Many little things done well. It's all in the details more metal and less plastic for a long beautiful life.

## 1. Upper Shelf

- Procedure Shelf in 3/16" rust-free aluminum 10 x 32"
- Close to work surface for less reach
- Aluminum radius-ed 3/16" rail support tubs
- Secondary work surface at bottom of shelf

## 2. Shelf Tower

- 1" round aluminum stock
- Optional tray holders for additional work surface
- Optional monitor mount

## 3. Bin Rack & Bins

- 3/16" aluminum plate bin rack
- ABS molded plastic bins
- Customizable dividers and shelves
- Brush-anodized aluminum bin lids with quick-release stainless steel spring latch

## 4. Work Surface

- Choice of Corian® solid surfaces
- Large work surface, finished size 480sq." (16.5" X 30")
- Large aluminum reinforcing plate
- Balanced with 60lb air-spring
- Rounded corners for easy pass-by

## 5. Chassis

- Life-time warranty
- Aluminum construction
- Large capacity for ancillary options (full 492 sq inches of usable space)
- Powder coat paint or electroless nickel plating on all surfaces
- 4 TIG-welded gusseted corners

## 6. Key-board tray

- Movable left to right
- Stores under control head
- Structural billet aluminum arm
- Aluminum construction

## 7. Control

- Life-time warranty on all metal parts
- No water in master control block for long trouble-free life
- Standard hp lock-outs accommodating 4-handed dentistry (allows burr removal with another hp in use)
- Very fine water coolant adjustment
- Tool kit for functionality adjustments
- Room for ancillaries including electric hp, Piezo & magneto strictive scalar, hp light source

## 8. Burr-Buddy

- Extra work surface for mirror, explorer

## 9. Chair Control Touchpad (Optional)

- Mounted in middle of hp holder bar



**DRESSED UP**

# DRESSED DOWN




**MAKE IT SIMPLE**






**MOUNT IT**



## Beneath Counter Mount


The very popular beneath counter mount is anchored by a half-inch steel plate that can be integrated into your custom-made rear cabinet, or as part of a Modular & Custom Cabinetry rear cabinet (PFR8144), which includes the steel plate pre-installed.



## Wall Mount

The wall mount is often used with an upper cabinet module for the storage of tubs and cassettes. The 7020 can replace the lower cabinet because it contains the industry's largest storage space "under the hood" for ancillaries that would otherwise be installed in cabinetry. The money saved by eliminating the lower cabinet goes a long way toward purchasing other practice-enhancing items and is an excellent solution for shorter treatment rooms.







## Under Cabinet Mount

The under cabinet mount will free up space for a drawer or nitrous mount by moving the mount underneath the cabinet. Designed to be compatible with all big-brand rear cabinets, the 7020 under cabinet mount is an excellent alternative to their delivery systems.

OR PUT IT ON A CART





## Cart Mount

The cart mount can be ordered with your desired length of umbilical, allowing for a greater range of movement left to right.





**TIME TO PLAY**



## Forest 3900 Chair

---





## Forest 3900 Plush

The 3900 Plush upholstery is designed with features ideal for any practice. Wide at the elbow and tapering in directly at the shoulder, the chair's backrest allows for close access to the patient while accommodating people of all sizes. Patients will enjoy the comfort of ample lumbar support, high-quality foam and standard glove-soft Ultraleather™ upholstery.



Aluminum Back



Thin Back



Sleek Design





## Forest 3900 Memory Comfort

The Memory Comfort upholstery option is an exceptional selection for long procedures or for those who are especially concerned with the patient experience. Simply stated, this upgrade makes it the most comfortable patient chair available today. With an unprecedented three inches of Aerus memory foam layered atop a five inch base of high density foam, Forest's Memory Comfort offers support where it's needed most, particularly the lumbar area. Patients will enjoy pillow-like comfort comparable to the finest luxury memory foam mattress. Opt for the Soft-touch Massage and Knee Pillow—your patients won't stop talking about his or her dental experience, a very powerful marketing tool indeed.



Knee Pillow



Back Massage



Memory Foam Headrest



## Forest 3900 Narrow Back With or Without Slings

If you need more space when working around the patient, the 3900 Narrow Back—with or without slings—is a smart choice. The optional slings lend a feeling of extra support to patients of any size, especially children. They also provide a contact barrier between you and your patient to help them feel relaxed and at ease. The 3900 slings are enhanced with a unique stiffening web at the elbow to prevent collapsing. The Narrow Back chair, with or without slings, could be the upholstery alternative that best suits your needs.



Pedo Wedge



Narrow Back Option



Supportive Slings



# Every Detail Done Well

Affordable Excellence signifies Forest's commitment to providing the same quality of material and workmanship as the big name brands for less. We use the highest quality bearings, hydraulics, surface treatments and electronics available, while placing an even greater focus on gaps, tolerances, fit and finish than our competition.



Easy Release



Solid Stainless



Drop Down Design



Pinch Release Headrest




Contrast Stitching



Cast Aluminum Back and Base



Touch Pad



Robust Engineering



Durable Pressure Formed Plastic



**Forest Dental**  \* sales videos available online at [forestdental.com](http://forestdental.com)

Phone

Fax

Digital

• 800 . 423 . 3555  
• 503 . 640 . 3012

• 503 . 693 . 9715

• [sales@forestdental.com](mailto:sales@forestdental.com)  
• [forestdental.com](http://forestdental.com)

01.15.18 • 98-906 REV. -



Comfort Stools



## Forest Comfort Stools

Finally, a stool you can get excited about. Forest offers North America's only doctor, assistant and office stools designed for forward-leaning practitioners. This exclusive feature allows for hours of pain-free work.

Every aspect of the stool has been designed to:

- *Remove stress on lumbar disks.* Because the DuPont® mesh seat prevents slipping, it supports legs in a decline position, keeping your back in the proper 90 degree angle relative to the seat.
- *Allow for natural movement while supporting your spine and releasing tension.* Because the backrest is rounded on the horizontal plane rather than the industry standard concave, it rolls up and down as you move forward and back, massaging the tissue and increasing circulation in the lumbar area.
- *Maintain proper alignment.* The Active-Tilt™ seat plate mechanism allows you to move and lock into a decline, neutral and recline position, eliminating pressure points in any posture.
- *Eliminate noise.* The all-surface 60mm polyurethane casters are virtually silent on any surface.

In addition, the Assistant's Comfort Stool offers ratcheting torso support and a foot ring with height adjustment. The Office Comfort Stool features adjustable, removable armrests, making it a smart investment for everyone at your practice. Thanks to intelligent engineering you and your staff will enjoy the comfort and back health you deserve.

## Affordable Excellence

Forest Dental led the industry in growth through the Great Recession because doctors fundamentally appreciate the value of many small details done well. With metal air and water Ts, chrome delivery arm end caps, oil impregnated bearings, premium hydraulics, tight gaps and tolerances, the quality fit and finish of Forest equipment is truly unparalleled. Committed to providing big brand features for less, we use the highest quality materials to ensure your equipment will perform for years to come.

## Designer Friendly®

Forest Dental is the only North American dental equipment manufacturer to offer equipment finished in a paint color other than white; we call this customization Designer Friendly®. Your dental equipment will enhance the interior design of your practice when you choose from one of 18 standard Ultraleather™ upholstery colors, nine powder coat colors, three carbon fiber look trims and a chrome option. Patients won't stop talking about their experience in your beautiful treatment rooms because Designer Friendly® is patient-friendly.



Shock Absorbers



Adjustable Armrests



No Pressure Points



Designed to Massage Lumbar



6187 Assistant's Comfort Stool




6176 Doctor's Comfort Stool



6198 Office Comfort Stool





Control Heads





## Forest Standard IC

Equipped with an ultra reliable pneumatic control block, the Standard IC is acclaimed for its simple design and straight-forward technology. We use one Clippard anti-retraction valve per handpiece enabling each one to operate independently. If a water drip were to occur in one position due to a stuck valve, two back-ups are available, unlike handpieces from other manufacturers with only one valve for all three positions. Because no water enters the main control block, you won't experience the damaging effects of water deposits and contamination. Affordable Excellence is achieved when every detail is done well.



Aluminum Lid



Metal Air/Water T's



Appropriate Technology (Built-in Tool Kit)



## Forest Euro

The Euro control head pivots at the center of the chassis, allowing it to be positioned both next to and over the patient in the transthoracic position. Designed for maximum reach to the patient, the whips are mounted at the front of the control. The true beauty of Forest's European style delivery is how smoothly the whips operate. They pull forward effortlessly and move side to side with just enough resistance to retract back. This control is top-rated in the eyes of the discerning European dentists who have tested it against the vast array of whip system controls available to them.



Cast Iron Chair Mounts



Feather-Touch Movement



Center Pivot





## Forest Value

The Value Series has earned a reputation for remarkable reliability at an affordable price. The durable pinch valve block system makes it ideal for high-production clinics, teaching facilities and the hygiene room. After seven to ten years of use, typical maintenance involves no more than moving the existing tubing a quarter of an inch away from the old pinch-point. This series is another great example of Affordable Excellence, simple yet appropriate technology.



Adjustable Holders



Simple Effective Technology




Designer Friendlytm

## Control Head Mounting Options


Whether mounting your Forest control on the chair, wall, floor or a cart, you can be certain that with our simple and appropriate technology you will enjoy peace of mind for years to come.

### Pivot Chair Mount




Pivot Chair Mounts easily move left to right on large roller bearings to accommodate left and right-handed practitioners. Unlike a Fixed Chair Mount delivery system, the mounting arm comes from underneath the patient rather than over the top, making it less obtrusive and more patient-friendly. Pivot Assistant's Instrumentation, Operatory Lights and Monitor Mounts are popular options to support your style of practice.

### Fixed Chair Mount With Sidebox




Fixed Chair Mounts With Sidebox allow for a Cuspidor, Assistant's Instrumentation a Monitor and Light to mount on a sidebox, which is large enough to accommodate a host of ancillary support products, finish the aluminum side covers in one of our nine powder coat paint colors and choose a complementary carbon fiber trim to finish the top cap. If you use a cuspidor, select our beautiful high-gloss ceramic black for a customized and beautiful patient-friendly delivery system.

# Fixed Chair Mount Without Sidebox




Fixed Chair Mounts Without Sidebox are the most cost effective solution of our three chair-mounted configurations. With the ability to mount a cuspidor, Assistant's Instrumentation, a Light and Monitor, this mount not only eliminates the cost of the sidebox, but it also has a smaller footprint, perfect for the narrow treatment room.


## Designer Friendly™ Surface Options




Forest Grey




Shadow




White




Arctic Silver




Pewter




Midnight



Copper




Charcoal




Bronze

Powder Coat

## Carbon Fiber



Laser Leather



Lemans



Greywood



## Forest 9080 Halogen Light

---

Affordable Excellence + Designer Friendly®







## Forest Halogen Light (Models 9080, 9090)

Sales of the 9080 halogen light grew 120 percent during the last four years of the recession thanks to Forest's exclusive concepts of Affordable Excellence and Designer Friendly®. Choose from one of nine lustrous, two-coat powder coat paint colors, accented with chromed billet aluminum flex arm end caps that won't crack, fade or fall off. A third axis of rotation allows you to place the light over your shoulder while keeping the long axis of the light pattern parallel with the mouth. The reflector is embedded into the glass preventing it from peeling or cracking over time. An adjustable flex arm on/off actuator, a long-life halogen bulb, a scratch-resistant Lexan light shield and more make this light an outstanding value.



Adjustable Light Pattern



Scratch-Resistant Light Shield



Third Axis of Rotation





9078 LED Light





## Forest LED Light (Model 9072, 9077)


The 9072 brings you all the benefits associated with the latest LED bulbs in an all-aluminum light head housing. Aluminum is superior to plastic as it can be powder coated in one of our nine colors to complement your interior design; the tough finish easily withstands the rigors of the treatment room and will not fade, discolor or crack.

The optics of the 3X5 inch pattern are designed to keep light out of the patient's eyes. The color of the light can be adjusted from 4200 to 6000 Kelvin for tissue examination and shade/color matching respectively. A much appreciated feature of the LED bulb is that the same color is maintained regardless of the intensity level, enabling you to turn down the light to eliminate glare, without affecting color. Brightness or LUX can be adjusted from 8000 to 35000. A third axis of rotation allows the light head to be positioned over the shoulder while keeping the pattern parallel with the long axis of the mouth so you don't have to compete with the light for space above the patient's head. The 9072 is the only light of its kind with autoclavable handles small enough to fit in a standard sterilizer, especially important for patients that require hospital grade infection-control.


# Forest Light Mounting Options

Forest offers five different mounting options for our lights: wall, cabinet, ceiling, track and chair. All mounts are heavy-duty in construction and designed with leveling adjustments to eliminate drift. The bearings and bushings in key pivot points are of the highest quality, designed for a long, trouble-free life.


## Wall & Cabinet Mount




## Ceiling Mount



## Track Mount



## Chair Mount



## Forest Dental Cuspidors

---

*Affordable Excellence + Designer Friendly\**



# Forest Deluxe Cuspidor

The Forest Deluxe Cuspidor is constructed of earthenware ceramic clay and finished with a premium high-gloss glaze. Available in black and white, all of our cuspidors are the most stain and fade resistant on the market today. A cup fill lever moves to either side of the bowl giving you the ability to control its activation with the option to add a unit-mounted activator as well. The Deluxe Cuspidor swings in toward the patient or completely out of the way for convenience. It can be mounted on our pivot left-to-right system or positioned on our chair mounted sidebox to meet your exact needs. The sidebox is then finished with a complementary carbon fiber look cover and painted in one of our nine powder coat paint colors. The cuspidor is equipped with a triple chrome-plated cup fill and bowl rinse spout. You, your staff and customers will enjoy the best looking and the only Designer Friendly® cuspidor on the market today.

## Forest Basic Cuspidor

Our Basic Cuspidor employs simple, easy-to-service and maintain mechanics with a timed bowl rinse and manual cup fill. The ceramic bowl is available in both white and black, which can be easily removed for cleaning. The cover comes in three carbon fiber look finishes or the industry standard Forest Grey and white. The Basic Cuspidor can be unit-mounted with or without a sidebox. Simply stated, it is an expensive-looking Designer Friendly® cuspidor, at an affordable price.



**Forest Dental**

#### Phone

- 800 . 423 . 3555
- 503 . 640 . 3012

#### Fax

- 503 . 693 . 9715

#### Digital

- sales@forestdental.com
- forestdental.com

#### Address

- 6200 ne cherry drive hillsboro oregon 97124



\* sales videos available online at forestdental.com



White Deluxe Cuspidor




Black Basic Cuspidor




White Basic Cuspidor







Forest Grey




Shadow




White




Arctic Silver




Pewter




Midnight



Copper




Charcoal



Bronze

Powder Coat


## Formica®




Surf




Folkstone



Rocky Road Granite




Oyster Grey




White Blanc



Moonstone




Silver Mist



Platinum Granite



Sanibel Granite




Goldrush Gemstone


## Gibraltar®



Arctic Melange




Baja Melange




Antarctica

## Corian®



Cinder Granite




Starry Starry Night





Raven Wing



Cave



Stone



Doe




Papyrus




Charcoal




Diplomat Blue



Dove Grey



Eucalyptus



Riviera



Pearl Midnight



Pearl Pewter



Pearl Abalone


Ultraleather®




Baltic




Parrot



Lemon



Red




Apricot



Laser Leather



Lemans



Greywood

Carbon Fiber



# UPDATE

## YOUR WAY

Updating your practice over a period of time, as you do with your home, following a plan that works with your budget, is one of your most important long-term financial considerations. You'll minimize the financial impact while improving curb appeal and the enjoyment of patients and staff.

# Why Update?

---

Hello,

Ten years ago, after 25 years selling dental equipment, I bought 38-year-old Forest Dental. I am proud to say that during that time we have been the North American growth leader in our industry and a top choice for treatment room updates.

The history of our success is the result of many factors, including a commitment to innovation, teamwork, and attentively listening to customer needs. But in great part it is also due to a single guiding principle: dare to be different. Along with producing state-of-the-art dental equipment, we are the only manufacturer that offers paint colors other than white. Have fun re-creating the colors of your favorite sports team or school. Or match the paint and upholstery colors of your living room. How about black upholstery with red contrast stitching on a charcoal gray metallic patient chair? Unleash your imagination, and, yes, dare to be different.

Forest customers are quick to point out that their Forest Treatment Rooms are a key internal marketing component because they create a fun atmosphere that speaks volumes about their aesthetics and attention to detail. Insist on being shown the beauty of Forest.

Thank you,  
Hank Barton

## Practice

---

Incremental upgrades provide immediate rewards as well as long-term returns.

Resale Value

---

Efficiency

---

Productivity

---

Patient Loyalty



Treatment Speed

---


New Services

---

Less Waste

---

Staff Retention



## Patients

---

Value. Comfort. Trust.

---

Your patients will appreciate your services more when they are performed in a beautiful, up-to-date environment.

## Doctor

---

Modern treatment room layout, equipment placement and interior design work together to enhance your joy in dentistry.

## Productivity

---

Plan a custom configuration that best suits your style of dentistry.

## Satisfaction


---

Enjoy your work environment and reduce stress.

## Ease

---

Adopt better ergonomics for more comfort and choice.





Only Forest offers complete color integration of dental equipment with your interior design. Choose from nine powder coat paint colors, three carbon fiber trim choices as well as a chrome option. Create a personalized look that will make you, your staff and your patients smile.


## Forest Fundamentals

---



### WOW FACTOR

Color is important, especially in a dental treatment room. Just as a tooth needs to match the one next to it, your treatment room surfaces can blend seamlessly together: a big step toward total color integration.



### INNOVATION IMPERATIVE

Habits are hard to change, but the upside of change is innovation. Forest embodies the latest thinking in ergonomics, materials, instrument management and infection control.



### FUN FACTOR

More and more people are realizing FUN is the most important WHY in life. Forest embraces fun. Whether it's red stitching on black upholstery or choosing the colors of your favorite sports team, all you have to do is set your imagination free.



### BUILT TO LAST

New doesn't have to mean plastic. While our competitors continue to move toward more and more plastic parts, Forest uses metal whenever possible so you get real value from your investment.

## Affordable Excellence


---

Big brand features for less is the brand strategy that helped us become an industry growth leader. We pay extra attention to detail and quality, from tight tolerances between parts to durable finishes and metal components instead of plastic, so you get maximum value from what you pay for.



### WORD OF MOUTH

We feel good when being in a warm and friendly environment. We tell our friends about a restaurant with great food at a reasonable price in the perfect setting. The same is true with dentistry. We've seen that within a year of upgrading, practices experience up to a 40% increase in new patients.



### FORM FACTOR

Modern equipment design and layout removes clutter, frees up space and adds new levels of flexibility so everyone in your office can work at their optimum.




### GO FURTHER

New materials and instrument management systems greatly eliminate inventory levels and shorten patient turnaround times, so you can make the most of your resources and your time.



### MORE CHOICE

Forest offers North America's greatest array of delivery options so your upgrade can include no more or less than exactly what your style of dentistry demands.



Update Your Practice One Step (or One Room) at a Time.

## PARK PLACE

### Treatment Room

- 3900C Hydraulic Operatory Chair
- 3907-999-DC Dual Upholstery
- 6176 Doctor Comfort Stool
- 6187 Assistants Comfort Stool
- 9072P LED Unit Light
- 4485PI Standard IC Delivery System  
Black Tubings & Holders
- 1112-644 Light Post
- 1700-081-CC Unit Mounted Touch Pad

## FENWAY

### Treatment Room

- 3900C Hydraulic Operatory Chair
- 3907-999-DC Dual Upholstery
- 3907-999-DC Deviated/Contrast Stitching
- 6176 Doctor Comfort Stool
- 6187 Assistants Comfort Stool
- 9072-C LED Ceiling Light
- 7000BC Duo Swing Beneath Counter Mount DU
- Solid Work Surface : Black Tubings & Holders
- 1700-081-CS Unit Mounted Touch Pad




Pass Go and Collect \$200

We spend 60% of our time at work. Shouldn't we make sure that this environment is as fun and welcoming as possible? Forest's Designer Friendly™ lets you create offices and treatment rooms that will delight your patients and staff—and reflect your personality and attitude. Choose from:

- 9 powder coat paint colors
- 3 carbon fiber-like trim options
- Chrome finishes
- 24 standard ultra-leather upholstery colors
- 95 optional ultra-leather upholstery colors

No longer are you restricted to white dental equipment. If you are a sports fan, each room could represent a favorite sports team. How about a nautical theme, or the same color palette you use in your living room? Can you picture your clinic and its treatment rooms featured in Architectural Digest? We can, thanks to total color integration.



# WHY

## CHOOSE FOREST

Our Brand Strategy Affordable Excellence. Simply put, with Forest you get big brand features (and more!), for less.

## A Few Others to Consider

---

- Life-time warranty for all the parts “under the hood”
- Life-time chassis warranty on our best-in-class 7000 series controls
- Five-year warranty for product in general
- Durable chairs, lights and delivery systems with greater resistance to UV and cleaners thanks to more metal, less plastic
- Aluminum LED light head
- Aluminum 7000 series controls
- Cast iron chair mount adapter
- Triple chrome-plated aluminum light and delivery system flex arm end caps
- Aluminum control head covers
- Aluminum utility center base (to protect against damage from cleaning)
- Aluminum patient chair back and base and utility center base
- Brass air and water Ts
- Pressure- (not vacuum-) formed chair and utility center covers that will maintain their shape over time
- Two-coat powder coat paint with high-quality clear coat that produces show-car shine with maximum durability
- 3900 patient chair with independent chair speed adjustments for base and back, up and down, so you can tailor chair speed to your practice style
- Optional soft-touch back massage with integrated compressor (eliminating impact on your central compressor)
- Exclusive hand piece lockouts on 7000 series controls, opening the door to 4-handed dentistry
- Complete toolkit with everything you need to adjust controls for pressure, hp angles and more



## What People Say About Forest

---

- An investment in Forest equipment enhances the interior of any practice, offering customized paint colors, upholstery, and surface treatments—surpassing big brand innovation, for less.
- I love that Forest's equipment blends seamlessly with my treatment room design. My staff and patients love our new look! Making my patients feel at ease is a priority. With Forest, it's as if they're walking into their favorite coffee shop. I could see a difference in my patients' mood right away.
- The color choices drew me in. We turned choosing colors for our new upholstery into a team-building activity. I couldn't be happier with the results.
- With dentistry, the devil's in the details. Forest customized our look down the contrast stitching. It was a very cool detail that I never would have picked up on without Forest's helpful guidance.
- Frankly, plastic just doesn't hold up like metal to the disinfectants we use. More metal, less plastic. That's what I saw, and it won me over.
- The handpiece lockouts on the Forest 7000 are a great time saver. My assistant can change a burr while I'm working. It suits us perfectly.
- The lifetime parts warranty put me at ease and made the decision to replace my master control block much easier.
- Practicing dentistry for the past 15 years has taken a toll on my body. The forest comfort stool has been a great relief for my back! Thank you, Forest, for considering the health of the the doctor.
- My new chair isn't just the most comfortable one I've ever had; the Ultraleather™ upholstery makes it luxurious.
- The money I saved on my 10 rooms of Forest helped me buy a new cone beam.
- Forest customer service is just as good, even better, than stores like Nordsrtom. Complete toolkit with everything you need to adjust controls for pressure, hp angles and more.



Forest Dental

---

### Phone

- 800 . 423 . 3555
- 503 . 640 . 3012

---

### Fax

- 503 . 693 . 9715

---

### Digital

- sales@forestdental.com
- forestdental.com

---

### Address

- 6200 ne cherry drive hillsboro oregon 97124