

FX-FCD Dental Aerosol Extractor



Features

- ✓ 3" x 5' Arm Reach
- ✓ 8" Clear Fume Hood Included
- ✓ 4-Stage Filtration + UV Light
- ✓ 10 Speed Control
- ✓ Control Panel & Remote
- ✓ Audible Filter Change Alarm
- ✓ Sturdy Metal Construction
- ✓ Long-Lasting, Low-Maintenance DC Brushless Motor w/ 40,000 hour life
- ✓ Casters
- ✓ Germicidal UV Light
- ✓ 3-Year Warranty

Tech Data

Airflow:	200 CFM @ Hood 270 CFM @ Inlet
Arm Diameter:	3"
Dimensions:	17" L x 10" W x 18" H
Input Power:	110V / 3 AMPs / 350W
Net Weight:	37 lbs. (55 lbs. w/ Filters)
Noise Level:	30-62 dB
Filter Stages:	1. Washable Pre-Filter 2. Water Impingement Filter 3. HEPA H13 Filter (99.95% Efficient at 0.3μ) 4. 7 lb. Carbon Filter 5. Germicidal UV Light