

# 2021 SUMMER PROMOTIONS

## FEATURED PRODUCTS

PAGES 30-34



PAGE 32

ORCA™ FOAM  
Absorbable Gelatin  
Sponge Hemostat, USP  
By ORCA™ Professional



PAGE 30

Crest® Whitening  
Emulsions with LED  
accelerator light  
By Crest® + Oral-B®



PAGE 31

Crest® + Oral-B® iO  
Transformational  
Gum Health Electric  
Rechargeable System  
By Crest® + Oral-B®

3 Questions to  
Consider When  
Evaluating Cone  
Beam 3D CBCT

PAGES 35-37

## ADVANCEDENTAL

PAGES 06-28



PAGE 28

LC Nano Hybrid  
Flowable Composite



PAGE 28

Universal Bond Plus  
& Universal Bond  
Plus DCA



PAGE 28

Temporary Luting  
Cement



NOW AVAILABLE!

Advance Dental Products by Goetze Dental

KANSAS CITY

ST. LOUIS

DENVER

DES MOINES

OMAHA

SPRINGFIELD

WICHITA

GOETZE  
d e n t a l

18  
since

# SUMMER PROMOTIONS

JULY 01 - SEPTEMBER 30

# FOUR EASY WAYS TO ORDER

## CONTENTS

- 06-28 ADVANCE DENTAL PRODUCTS
- 30-34 FEATURED PRODUCTS
- 35-37 3 QUESTIONS TO CONSIDER WHEN EVALUATING CONE BEAM 3D CBCT

## SUMMER PROMOTIONS

- |       |           |       |                    |
|-------|-----------|-------|--------------------|
| 38-40 | 3M        | 56-59 | SDI                |
| 42-43 | DMG       | 60-61 | Septodont          |
| 66    | Hu-Friedy | 47-50 | Solometex          |
| 44-46 | Kerr      | 51    | SS White           |
| 52-53 | Morita    | 64    | Surgical Esthetics |
| 54-55 | NSK       | 62-63 | Young              |
| 41    | Parkell   |       |                    |
| 65    | Pulpdent  |       |                    |

## ADVERTISER INDEX

- |       |                                     |       |                   |
|-------|-------------------------------------|-------|-------------------|
| 29    | A-dec   KaVo                        | 73    | Ivoclar Vivadent  |
| 77    | Awaken2Sleep                        | 71    | Kulzer            |
| 67    | 3M                                  | 74    | NSK               |
| 80-81 | Clear Comfort Aligners              | 72    | Orabloc           |
| 76    | Crest & Oral-B                      | 78-79 | ORCA Professional |
| 72    | DMG                                 | 82    | PacDent           |
| 83    | EDDSON Practice Management Software | 70    | Planmeca          |
| 71    | Integra Miltex                      | 68    | SDI               |
| 69    | Isodam                              | 69    | SS White          |
|       |                                     | 75    | Triology          |

PLEASE NOTE: Goetze Dental is not responsible for any typographical errors in the price or descriptions of items featured in this flyer. Prices subject to change due to manufacturers' price increases. Due to the COVID-19 pandemic, orders for all PPE and infection control products can not be cancelled, changed or returned. Visit [goetzedental.com/your-resources](http://goetzedental.com/your-resources) for the most up-to-date pricing.

**GAINSpeed™**  
desktop ordering application



**GoetzeOnline.com**  
e-commerce ordering option



**800.692.0804**  
8am - 5pm CST



Fax:  
**800.552.7791**



MEET EDDSON

# EDDSON™

PRACTICE MANAGEMENT SOFTWARE

Simplify the way you work.

As a dental team, you juggle a variety of tasks in an effort to provide care to your patients. EDDSON™ practice management software simplifies the way you work by streamlining and anticipating dental practice workflows.

→ Learn more or schedule a free demo at [meeteddson.com](http://meeteddson.com) or call 800.692.0804



## TECHNOLOGY



### INFORMATION TECHNOLOGY FOR THE PRIVATE DENTAL PRACTICE

Dental practices face unique challenges with information and technology. Making sure dental practice equipment and technology are working in harmony is a crucial part of keeping your practice thriving.

Goetze Dental's Information Technology Team specializes in integrating dental equipment and technology and supporting the dental IT environment. Goetze HelpDesk Agents are experts in Dental Technology HIPAA compliance. As a full-solution dental practice IT provider, we provide both stand-alone solutions and managed services for private dental practices in the Midwest.

→ Learn more at [goetzedental.com](http://goetzedental.com) or call 800.692.0804, option 7 for a FREE IT Practice Evaluation.



**4** EASY WAYS TO ORDER

**01** **GAINSpeed™**  
Our proprietary desktop ordering app

**02** **GoetzeOnline.com**  
E-commerce ordering platform

**03** **800.692.0804**  
Phone

**04** **800.552.7731**  
Fax



## ABOUT

Goetze Dental is a family-owned business serving private dental practices just like yours since 1884.

We work to understand your goals and dreams and strive to provide the best dental supplies, equipment, technology and service.

We believe the best patient care is delivered in the private practice.

### OUR COMMITMENT TO PRIVATE PRACTICES

100% of our energy and resources will continuously center on the products and services we offer to the privately held dental practice through integrity, education, innovation and customer service.

We stand united with the artists, the scientists and the explorers of the dental community, whose interest in providing exceptional patient care never fades.

## SUPPLIES



### DENTAL SUPPLIES

As a clinician, you are driven to provide the best care possible for your patients. Choosing dental supplies is a vital component in that care. We partner with the top manufacturers in dental so you have the products you need, when you need them.



**AUTHORIZED  
DEALER OF  
THESE BRANDS  
AND MORE**

3M  
Advance  
Ansell Microflex  
Aurelia Supermax  
Butler  
// GUM  
// Sunstar  
Coltene

Crest Oral-B (P&G)  
Crosstex  
Dentsply Sirona  
DMG  
Glaxosmithkline  
GC America  
Halyard  
Kulzer

Hu-Friedy  
IDS Megagen  
Ivoclar Vivadent  
Integra-Miltex  
Johnson & Johnson  
KaVo Kerr  
Parkell  
Premier

Pulpdent  
Quala  
SDI  
Septodont  
Sempermed  
SS White  
Sultan  
Surgical Esthetics

Young  
// Biotrol  
// Denticator  
Vista  
Zirc


# EQUIPMENT

Goetze Dental Equipment Specialists and Technology Technicians are masters at integrating dental equipment, software, hardware and dental technology.

## AUTHORIZED DEALER OF THESE DENTAL EQUIPMENT BRANDS AND MORE

A-dec	DCI Edge	KaVo Kerr
Accutron	DentalEZ	// Dexis
Acteon	Dentsply Sirona	// Gendex
Air Techniques	Digital Doc	// iCat
Belmont	Forest Dental	// KaVo Imaging
Brewer	Hu-Friedy	// Nomad
Coltene	IDS Dental	// Pelton & Crane
Crown Seating	Morita	Ivoclar Vivadent
		Midmark

NSK  
Planmeca  
Porter Instrument  
Royal  
SciCan  
Vatech



# SERVICE

## EQUIPMENT REPAIR

Dental equipment and technology malfunctions and it's usually at the most inconvenient time. Just a few hours of lost productivity can cost you thousands in lost production.

Goetze Dental's Service Team is committed to providing the solutions you need to keep your equipment and technology functioning, so you can continue to provide exceptional patient care.

## SCHEDULE A SERVICE CALL

at [goetzedental.com](http://goetzedental.com) or call **800.692.0804**, option 5



# DENTAL OFFICE DESIGN


Look at your practice through the eyes of a patient... Is it time for a change?

Whether a remodel is in store or you're ready to start from scratch, we are quite passionate about creating DREAM practices and we have a bit of experience doing so.

Our Design-Build Team has years of experience with all types of dental office design-build projects. They are experts at helping keep your project on time AND within your budget.

**DREAM IT. DESIGN IT. BUILD IT.**

Let Goetze Design Services help you create the perfect office.



Schedule a complimentary design consultation at [goetzedental.com](http://goetzedental.com) or call 800.692.0804

**GOETZE**  
dental



## Topical Anesthetic Gel

Contains 20% Benzocaine for maximum strength pain relief while fast working to quickly relieve patient discomfort. Available in six great tasting flavors.

**Bottle 1oz**

Cherry .....	<b>575018</b>
Mint .....	<b>575019</b>
Piña Colada .....	<b>575020</b>
Bubble Gum .....	<b>575021</b>
Strawberry .....	<b>575022</b>
Raspberry.....	<b>575023</b>
	<b>Ea \$6.97</b>



## Lidocaine HCl 2% with Epinephrine 1:100,000

Features:

- No-methylparaben formulation to avoid allergic reactions on sensitive patients.
- 100 percent latex free components to help reduce allergic responses.
- Terminal sterilization of product/cartridges for safe use.

Box/50 .....	<b>577673</b>
	<b>Ea \$37.79</b>



## Premium Non-Woven Sponges

- Unique rayon/polyester blend
- Provides maximum absorbency
- Less adhesion to wounds for faster healing
- Soft, Smooth, Strong, Virtually lint free
- Ideal for sustained patient use, disinfecting, general cleaning
- 4 ply non-sterile

<b>2" x 2"</b>	
Case/5000 .....	<b>097805</b>
	<b>Ea \$75.00</b>

<b>4" x 4"</b>	
Case/2000 .....	<b>097837</b>
	<b>Ea \$90.65</b>



## Plastic Cups

Horizontal embossed grip for easier dispensing and holding. Double rolled lip for easy separation. Can be used with any 5oz dispenser.

**5oz Cups Case/1000**

Blue .....	<b>191885</b>
Dusty Rose .....	<b>191889</b>
Silver Gray.....	<b>191887</b>
White .....	<b>191880</b>
	<b>Ea \$28.45</b>

**GOETZE DENTAL** is not responsible for any typographical errors in the price or descriptions of items in this flyer. Prices subject to change due to manufacturer price increases and market fluctuations. Due to the COVID-19 pandemic, orders for all PPE and infection control products can not be guaranteed for delivery and are subject to backorder status, allocation of amounts and/or cancellation.



Both Econoback and Polyback Towels come in a wide variety of colors!

## Econoback Towels 3-Ply

Highly absorbent tissue specially treated to enhance wet strength and tear resistance, and a bonded layer of poly film. Reinforced edge to provide strength for Towel clips.

2-ply tissue with 1-ply poly  
(3 ply bibs)  
19" x 13" - Case/500

Aqua .....	198252
Beige .....	198257
Blue .....	198260
Dusty Rose .....	198266
Green .....	198254
Lavender .....	198263
Peach .....	198264
Silver Gray .....	198268
White .....	198255
Yellow .....	198261
<b>Ea</b>	<b>\$30.71</b>

## Polyback Towels 4-Ply

Highly absorbent tissue specially treated to enhance wet strength and tear resistance, and a bonded later of poly film. Reinforced edge to provide strength for Towel clips.

3-Ply tissue with 1-ply poly  
(4 ply bibs)  
19" x 13" - Case/500

Aqua .....	198269
Beige .....	198270
Blue .....	198275
Dusty Rose .....	198292
Green .....	198280
Lavender .....	198282
Peach .....	198288
Silver Gray .....	198294
White .....	198285
Yellow .....	198290
<b>Ea</b>	<b>\$34.39</b>



HEAVY weight paper!



## Bracket Tray Covers

Heavy weight paper impedes the flow of moisture, minimizing liquid seep-through, protecting the covers and unit beneath the tray.

Ritter- "Size B" (8.5" x 12.25")  
Box/1000

Blue .....	198350
Dusty Rose .....	198347
Green .....	198422
Silver .....	198358
<b>Ea</b>	<b>\$29.32</b>
White .....	198360
<b>Ea</b>	<b>\$28.25</b>

## Saliva Ejectors


Provides optimal suction without aspirating tissue using a unique wire-reinforced tube, which maintains any desired configuration. Latex Free

Pkg/100

White / White tip.....	197510
Clear / Blue Tip .....	197508
<b>Ea</b>	<b>\$4.68</b>

## Air/Water Syringe Tips


- Rounded edge for maximum patient comfort
- Convenient dispenser bottle keeps work area clear and assists in preventing cross-contamination
- Rounded non-scratch end to minimize wear and tear of O-ring in hand piece



Bag/250 .....	576916
<b>Ea</b>	<b>\$27.64</b>

## Surgical Aspirator Tips

- Designed to fit all high-volume evacuation systems
- 30 degree angle for maximum patient comfort
- Made of rigid opaque plastic



Bag/25

1/8" White .....	576914
1/4" Green .....	576915
<b>Ea</b>	<b>\$5.99</b>

## Combo Tip Evacuators (HVE)

Double-ended design with both a vented and non-vented end on the same tip. Latex Free

Pkg/100 .....	198174
<b>Ea</b>	<b>\$7.26</b>





### General Purpose Cleaner

Non-ammoniated concentrated liquid formula employs surfactants and anti-corrosive additives to safely and thoroughly clean instruments prior to sterilization. Phosphate-free.

Gallon ..... **575655**  
Ea **\$19.75**

### Tartar & Stain Remover

Powerful, ready-to-use ultrasonic solution to remove stains from dentures, partials, appliances. Also removes temporary cement.

Gallon ..... **575656**  
Ea **\$19.75**

### Temporary Cement Remover

Strong ready-to-use ultrasonic solution to remove temporary cement.

Gallon ..... **575657**  
Ea **\$19.75**



### Self-Sealing Autoclave Pouches

Autoclavable self-sealing, blue tinted, see-through pouches. Easy-to-read color changing arrows are indicators for chemical vapor, steam or gas sterilization. Utilizes a triple seal plus a fourth "Secure Seal" at the chevron.

#### Box/200

2.25" x 4" ..... **577685**  
Ea **\$9.24**

2.75" x 9" ..... **575073**  
Ea **\$11.20**

3.5" x 5.25" ..... **575074**  
Ea **\$11.72**

3.5" x 9" ..... **575075**  
Ea **\$11.14**

5.25" x 10" ..... **575077**  
Ea **\$16.14**

#### Economy - Box/500

3.5" x 9" ..... **575076**  
Ea **\$27.17**







### Light Handle Cover

Specially designed barriers effectively cover equipment with near custom fit.

**Box/500 T-style**  
4" x 5 3/4" ..... **575327**  
**Ea \$16.67**

### Syringe Sleeve

Syringe sleeves have pre-cut openings for easy application. Made of clear poly material that protects instruments from cross-contamination.

**Box/500**  
2 1/2" x 10" ..... **575326**  
**Ea \$11.08**

### X-Ray Cover

Specially designed barriers effectively cover equipment with near custom fit.

**Box/250**  
15" x 26" ..... **575328**  
**Ea \$23.34**

### Chair Cover

Tapered design conforms neatly to patient chairs, covering chair back.

**Box/225**  
27 1/2" x 24" ..... **575324**  
**Ea \$24.86**

### Tray Sleeve

Back-flip closure ensures tray is completely covered. Made of a clear poly material that protects tray from cross-contamination.

**Box/500**  
10 1/2" x 14" ..... **575325**  
**Ea \$17.97**



Easy removal with no residue left behind!

### Universal Cover Film

4" x 6" perforated sheets with non-adhesive side edges for easy removal. Specially formulated adhesive leaves no residue behind. Helps to prevent cross-contamination. Convenient dispenser box contains 1 roll of 1200 perforated 4" x 6" sheets

**Roll/1200**  
Clear ..... **575330**  
Blue ..... **575329**  
**Ea \$16.57**



### Temporary C & B Material

10:1 Ratio

Syringeable bis-acryl composite for chairside provisional restorations. Available in 3 shades: A1, A2, A3. All are delivered in 76g (50mL) cartridges designed to fit on a 10:1 style automix gun. Ten waste saver tips are included per kit. 50mL Cartridge with 10 tips

A1 ..... **575295**  
A2 ..... **575296**  
A3 ..... **575297**  
**Ea \$84.73**



### Retainer Boxes

- Clamshell design with easy to open/close "snap-lock"
- Sturdy, colorful plastic helps protect retainer from dust and damage
- Vented top and bottom holes allow box to "breathe"
- Flat bottom provides area for writing "Return if found" information
- Assorted colors 3" x 3 1/8" x 1 1/8"

**Pkg/12** ..... **577908**  
**Ea \$12.12**

### Denture Boxes

- One-piece "hinged" lid provides easy access to contents inside
- Sturdy, colorful plastic helps protect dentures from damage
- Deep soaking well fits both uppers and lowers at the same time
- Flat bottom provides area for writing "Return if found" information
- Assorted colors 3 3/4" x 3" x 2 3/8"

**Pkg/12** ..... **577907**  
**Ea \$13.19**





## Disposable Prophyl Angles

Latex Free

Straight-body style in pink with a webbed cup that is soft and latex-free. This angle features a screw-type cup, allowing for stability. Fits comfortably and securely in the hand for hard-to-reach areas, while maintaining a favorable paste-holding capacity.

### Soft Gray Cup

Box/100 ..... **575078**  
Ea **\$47.97**

Box/500 ..... **575079**  
Ea **\$197.26**

### Firm White Cup

Box/100 ..... **575080**  
Ea **\$47.97**

Box/500 ..... **575081**  
Ea **\$197.26**

Contact your  
Goetze Dental  
Representative  
to order!  
Or Order at:  
**800.692.0804**



## Pre-bent Multipurpose Dispensing Tips

Dispensing Tips are pre-bent to the proper angle and have a blunt end for patient comfort. Great for placing sealants, etchants, cements and other flowable materials. Secure Luer lock, universal color coding for easy identification

### Tips - Pkg/100

18ga - Pink ..... **575290**  
(Cavity liners, Bases & Composites)

20ga - Yellow ..... **575291**  
(Composites & Sealants)

20ga - Black ..... **575292**  
(Composites & Cements)

22ga - Gray ..... **575293**  
(Sealants, Etchants & Cements)

25ga - Blue ..... **575294**  
(Etchants)

Ea **\$8.12**

## Etching Gel 37% Phosphoric

Formulated for maximum working and handling characteristics.



### Syringes (1.2ml) with 12 Tips

Pkg/4 ..... **577676**  
Ea **\$11.84**

### Syringes (1.2ml) with 25 Tips

Economy Pkg/12 ..... **577677**  
Ea **\$30.00**

### Syringes (25ml) with 50 Tips

Jumbo Pkg/2 ..... **577679**  
Ea **\$36.50**

### Pre-bent Tips (25ga) - Blue

Pkg/20 ..... **577681**  
Ea **\$8.50**

### 3cc Empty Refillable Syringes

Pkg/5 ..... **577678**  
Ea **\$12.50**



## Fluoride Trays


Fluoride trays are color-coded for easy identification and are disposable for your convenience. The Single Arch trays are made of soft, light foam for patient comfort and have a convenient handle for easy placement and removal. The Dual Arch trays have a fill line that helps reduce gel overflow and have locking handles for easy placement and removal.

### Single Arch - Pkg/100

Small ..... **575028**  
Medium ..... **575029**  
Large ..... **575030**  
Ea **\$18.16**

### Dual Arch - Pkg/100

Pedo ..... **575024**  
Small ..... **575025**  
Medium ..... **575026**  
Large ..... **575027**  
Ea **\$41.32**



## Propyl Paste

**Mint w/Fluoride - Box/200**

X-Coarse ..... **575114**  
Coarse ..... **575113**  
Medium ..... **575107**

**BubbleGum w/Fluoride - Box/200**

Coarse ..... **575115**  
Medium ..... **575108**  
Fine ..... **575112**

**Raspberry w/Fluoride - Box/200**

Coarse ..... **575117**  
Medium ..... **575116**

**Assorted w/Fluoride - Box/200**

Coarse ..... **575110**  
Medium ..... **575109**  
Ea **\$24.99**



Ideal For:

- Occlusal registration
  - Copying transfers
- Primary impressions with bite trays



## High Performance Bite Registration Material

Formulated to have a light, mousse-like consistency ideal for occlusal registration, copying transfers and primary impressions with bite trays. Our material has excellent durometer or hardness preventing distortion in disposable trays, specifically the triple trays. Unique our Bite Registration approaches zero resistance allowing for unmodified natural movement of the jaw. The non-flowing, thixotropic formula stays in place and will not drip, run, or sag.

### Bite Registration Kit w/Tips - 50ml

Each pack contains (2) 50ml cartridges & (6) 6.5mm mixing tips

#### Regular Set

Work Time 2:00 Intraoral Set 1:00

Unflavored..... 575100  
Ea **\$26.20**

#### Fast Set

Work Time 0:30 Intraoral Set 1:00

Unflavored..... 575101  
Ea **\$26.20**

Vanilla..... 575102  
Ea **\$26.20**

#### Super Fast Set

Work Time 0:30 Intraoral Set 0:30

Unflavored..... 575099  
Ea **\$26.20**

Vanilla..... 575103  
Ea **\$26.20**

#### Rigid Bite - Fast Set

Unflavored..... 575104  
Ea **\$26.20**

#### Clear Bite - Fast Set

Unflavored..... 575105  
Ea **\$27.95**

## Multi-Color Disposable Impression Trays

Perforated

- Disposable premium impression trays made of high quality, rigid, transparent plastic
- Strategically perforated trays allow easy modification with dental bur to accommodate implants
- Can be heated and molded for a custom fit
- Bright color-coded body for quick and easy identification during procedure
- Easy to apply I.D. labels included in every bag
- Non-sterile

### Large Perforated Yellow - Bag/12

#1 - Upper..... 575279  
#2 - Lower..... 575280

### Medium Perforated Blue - Bag/12

#3 - Upper..... 575281  
#4 - Lower..... 575282

### Small Perforated Orange - Bag/12

#5 - Upper..... 575283  
#6 - Lower..... 575284

### Quadrant Perforated Green - Bag/12

#7 - UL/LR..... 575285  
#8 - UR/LL..... 575286

### Anterior Perforated Pink - Bag/12

#9 - Upper..... 575287  
#10 - Lower..... 575288

### Assorted Perforated

Bag/12..... 575289  
Ea **\$4.44**

## Premium 3 Way Impression Bite Trays

Color Coded

- Disposable Advance bite trays precisely capture the master impression, counter impression, and bite registration in one procedure
- Rigid frame with nylon webbing to ensure accurate bite
- Ergonomically designed thumb grip aids in placement and removal
- Patent Pending grooves provide superior material retention

### Box/12

Posterior (Blue)..... 575277  
Sideless (White)..... 575275  
Anterior (Purple)..... 575276  
Quadrant (Yellow)..... 575274  
Assorted Sizes..... 575278  
Full Arch (Pink)..... 575273  
Ea **\$8.10**



Peppermint Flavor



## Pulp Vitality

Refrigerant Spray

For use of thermo testing tooth vitality.

Peppermint Flavor

6oz Spray..... 575267  
Ea **\$18.26**





**Advance Technical  
Adult 30 Tuft**

Assorted: 4 assorted colors

Box/72 ..... **578066**  
Ea **\$48.00**



**Compact Super Row Bristles  
Adult 35 Interdental Cut Tuft with  
Super Grip**

Assorted: Dark Purple, Dark Blue, Teal, Orange

Box/72 ..... **578067**  
Ea **\$34.35**



**FlexBrush Clean  
Adult 40 Tuft Interdental Cut**

Assorted: Dark Blue, Red, Dark Green, Yellow

Box/72 ..... **578062**  
Ea **\$32.85**



**Super Tip Bristle  
Adult 38 Tapered Tuft**

Assorted: Yellow, Red, Dark Blue, Dark Green

Box/72 ..... **578063**  
Ea **\$34.58**



**Advance Compact Micro Tip  
Adult 33 Tuft**

Assorted: Blue, Green, Pink, Purple

Box/72 ..... **578068**  
Ea **\$49.71**



**Dual Action Ergo Grip  
Gum Stimulator and Tongue Cleaner  
Adult 36 Tuft**

Assorted: Dark Pink, Bright Blue, Vibrant Green, Orange

Box/72 ..... **578069**  
Ea **\$34.35**



**Clear Grip with Indicator  
Adult 38 Tuft**

Assorted: Blue, Dark Green, Orange, Bright Purple

Box/72 ..... **578057**  
Ea **\$50.76**



**Compact Opaque Straight Handle  
Adult 35 Soft Rounded Tuft**

Assorted: Purple, Dark Blue, Dark Green, Orange

Box/72 ..... **578065**  
Ea **\$34.35**



**Triple Action Sport Grip  
Adult 32 Tuft Interdental with  
Tongue Cleaner and Power Tip**

Assorted: Green, Purple, Red, Blue

Box/72 ..... **578058**  
Ea **\$39.75**

**YOUTH BRUSHES**



**Ladybug Grip**

Infant 24 Soft Tuft



Assorted: 4 assorted colors

Box/72 ..... **578061**  
Ea **\$31.15**



**Advance Color Youth**

Age 8+ 32 Tuft



Assorted: Two Tone colors

Box/72 ..... **578060**  
Ea **\$34.45**



**Child's Sparkle Ergo Grip**

Ages 5-7 27 Rounded Tuft



Assorted: 5 assorted two-tone colors with sparkles

Box/72 ..... **578059**  
Ea **\$33.45**





## SonicTouch

Rechargeable Sonic Power Toothbrush

- Advanced Sonic Cleaning Technology**  
Produces 40,000 brush strokes per minute for a powerful dynamic cleaning.
- 2 Minute Auto-Timer**  
Built-in timer to ensure dentist recommended brushing time.
- Base Charging Indicator**  
Built-in interval timer helps you brush each area of your mouth evenly.
- Waterproof Design**  
Protection from direct contact from water, spills or splashes.
- 30 Second Interval**  
Built-in interval timer helps you brush each area of your mouth evenly.
- 4 Brushing Modes: Clean, Soft, Whiten, Massage**  
Choose from Clean, Soft, Whiten, or Massage modes to cater all of your oral care needs.

### SonicTouch

Rechargeable Sonic Power Toothbrush

Advance SonicTouch technology provides a powerful dynamic cleaning action that helps to keep teeth and gums healthy by removing plaque.

This technology provides effective stain removal which helps to brighten teeth and restore them to their natural whiteness.

- 2 Minute Auto-Timer, 40,000 brush strokes per minute
- 30 second interval timer helps you brush each area evenly
- Base charging indicator
- Waterproof design
- 4 brushing modes: clean, soft, whiten, massage
- 3 brush heads
- 1 year warranty

SonicTouch Toothbrush ..... **578048**  
Ea **\$45.83**

### SonicTouch Refill Brush Heads

The rubberized brush head is designed to protect tooth enamel. The brush conforms to the contours of your teeth for effective, deep cleaning. High quality DuPont bristles promote naturally whiter teeth by gently polishing off surface stains, end rounded.



Pkg/2 ..... **578049**  
Ea **\$7.36**

### Advance Thick Flexible Interdental Picks

Cylinder interproximal picks. Perfect for cleaning between teeth and under orthodontic work.



Pack/72 (Thick) ..... **578185**  
Ea **\$53.08**

### Advance Thin Flexible Interdental Picks

Cylinder interproximal picks. Perfect for cleaning between teeth and under orthodontic work.



Pack/72 (Thin) ..... **578184**  
Ea **\$53.08**



### Advance Varnish

5% Sodium Fluoride

Easy non-messy single dose delivery system. Tooth-shaded color. Contains Xylitol. Tasteful and transparent 5% sodium fluoride.

**Box/50**

0.4ml unidose (Toffee) ..... **578143**

0.4ml unidose (Mint) ..... **578144**

Ea **\$50.42**



## Sutures

Super sharp siliconized needles for repeated passage through tissue. Needles are drilled-end to ensure smooth transition from the needle to the material and maximizes the needle attachment strength. Packaged in a double sterile pouch for sterile transfer, in color-coded suture boxes with a front dispenser. Packaged

12 per box

Silk = Light Blue  
Chromic Gut = Beige  
Plain Gut = Yellow

### 18" Plain Gut Sutures with Needle

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00071 591B	Absorbable	Plain Gut	3/8 circle -19mm reverse cut	4/0	18"	575058
	C-31	AD00065 594B	Absorbable	Plain Gut	1/2 circle -24mm reverse cut	3/0	18"	575056

Ea \$33.29

### 18" Chromic Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00558 558B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	4/0	18"	575061

Ea \$41.34

### 18" Silk Braided Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00066 781B	Non - Absorbable	Silk - Black Braided	3/8 circle -19mm reverse cut	4/0	18"	575057
	C-6	AD00072 784B	Non - Absorbable	Silk - Black Braided	3/8 circle -19mm reverse cut	3/0	18"	575059
	C-31	AD00064 785B	Non - Absorbable	Silk - Black Braided	1/2 circle -24mm reverse cut	3/0	18"	575055

Ea \$28.19

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C26	AD00062 778B	Non - Absorbable	Silk - Black Braided	1/2 circle -15mm reverse cut	3/0	18"	575053

Ea \$28.96

### 27" Plain Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00822 554B	Absorbable	Plain Gut	3/8 circle -19mm reverse cut	3/0	27"	575065

Ea \$33.86

### 27" Chromic Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	AD00636 559B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	3/0	27"	575064
	C-6	AD00635 560B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	4/0	27"	575063

Ea \$33.86

Contact your  
Goetze Dental  
Representative  
to order!  
Or Order at:  
**800.692.0804**

Standard Operatory: **FG** - Friction Grip **RA** - Right Angle **SS** - Short Shank **HP** - Hand Piece

	Round					Pear				Inverted Cone			Straight/Flat End Plain		Straight/Flat End Cross Cut		
SIZE	1/2	2	4	6	8	2P	4P	330	332	33 1/2	35	37	56	57	556	557	558
FG 10PK		575222	575223	575224	575225			575221			575219						575217
FG 100PK		575165						575159			575162		575163	575164	575160	575161	
RA 10PK		575178	575177	575176	575175												
RA 100PK		575154															
FG SS 10PK								575227									575226
FG SS 100PK								575215									

	Straight Round End Plain		Straight/Round End Cross Cut		Taper Flat End Plain		Taper Flat End Cross Cut		Amalgam Prep	Metal Cutter				Surgical Length					
SIZE	1157	1557	1558	169	701	245	2	330	557	1557	2	4	6	557	558	701	702		
FG 5PK											575211	575212	575213	575214		575208	575210		
FG 10PK					575218	575220	575173		575174										
FG 100PK		575155				575158	575172	575171		575170									
FG SS10PK																			
FG SS 100PK																			

	Trimming and Finishing																
SIZE	7006	7008	7204	7379	7404	7406	7408	7642	7664	7802	7803	7901	7902	7903	8404 20 Blade	9006 30 Blade	9008 30 Blade
FG 5PK	575202				575204	575205	575200			575199	575198	575197	575196				

	Trimming and Finishing			
			Safe End	
SIZE	9406	9903	SE 9-10 10 Blade	SE 9-20 20 Blade
FG 5PK	575194	575195		

### Trimming & Finishing FG Burs

Twelve Bladed Burs - Pkg/5  
Ea **\$28.95**

### Trimming & Finishing FG Burs

Twelve Bladed Burs - Pkg/5  
Ea **\$21.50**

### FG Carbide Burs

Economy - Pkg/100  
Ea **\$104.99**

FG Carbide - Pkg/10  
Ea **\$13.75**

### Metal Cutting Carbide Burs

Economy - Pkg/100  
Ea **\$299.48**

Pkg/10  
Ea **\$30.77**

### RA Carbide Burs

Economy - Pkg/100  
Ea **\$124.98**

Pkg/10  
Ea **\$19.50**

### Surgical Length FG Carbide Burs

Pkg/5  
Ea **\$9.25**



**Endodontic Irrigating Tips**

Endodontic irrigation needles used to irrigate canals, fistulas and pockets. Luer lock hub

**Slotted End/Side Vent - 100/bag**

- 23 Gauge ..... 579271
  - 27 Gauge ..... 579272
  - 30 Gauge ..... 579273
- Ea \$31.62**

**Safe End/Side Port - 100/bag**

- 23 Gauge ..... 579274
  - 27 Gauge ..... 579275
  - 30 Gauge ..... 579276
- Ea \$44.85**

**Endodontic Syringe with Needle**

- Used for the irrigation of pockets, fistulas and canals
- 3cc Luer syringe
- Slotted End/ Side Vent

**Endodontic Syringe with Needle - 100/bag**

- 23 Gauge ..... 579277
  - 27 Gauge ..... 579278
  - 30 Gauge ..... 579279
- Ea \$39.37**

**Viscous Brush Tips**

- Easy extrusion of viscous materials
- Photo sensitive opaque tip maintains the flow of light-cured material
- Luer lock hub



Pack/100 ..... 579288  
**Ea \$57.41**

**Pre-bent Brush Tips**

- Curved cannula for optimal visibility
- Precision placement of hemostatic agents
- Padded-end brush tip controls flow
- Universal luer lock hub



- 19 Gauge - Pack/20..... 579285  
**Ea \$12.82**
- 19 Gauge - Pack/100..... 579286  
**Ea \$50.83**
- 19 Gauge - Pack/500..... 579287  
**Ea \$178.48**

**Low Viscosity Brush Tips**

- Curved cannula for optimal entry
- Opaque luer lock hub
- Nylon bristle end enables controlled delivery of solution



19 Gauge - Pack/100 579289  
**Ea \$42.16**

**Black Mini Brush Tips**

- Easy extrusion of viscous materials
- Photo sensitive opaque tip with luer lock hub



18 Gauge - Pack/100..... 579291  
**Ea \$71.35**

**Calcium Hydroxide**

- 12.5pH aqueous solution for temporary filling and cleansing root canals
- Contains no more than 35% Calcium Hydroxide paste

**Standard Kit**  
Contains: 4-1.2mL syringes with tips..... 579269  
**Ea \$32.62**

**Value Pack**  
Contains: 20-1.2mL syringes with tips ..... 579270  
**Ea \$83.05**

Contact your  
Goetze Dental  
Representative  
to order!  
Or Order at:  
**800.692.0804**






**Advance Interdental Brushes**

Great for removal of food and plaque on bridges, crowns, braces and between teeth.

- Wide Cone - Blue Grip**  
Box 72..... **579225**
- Moderate Cylinder - Red Grip**  
Box 72..... **579224**
- Ultra Tight - Yellow Grip**  
Box 72..... **579222**
- Tight Cylinder - Green Grip**  
Box 72..... **579223**  
**Ea \$54.37**



**Prepasted Disposable Adult 39 Tuft (Mint Flavor)**

Assorted: Yellow, Red, Green, Blue  
Box/100 ..... **579708**  
**Ea \$16.30**



**Denture Brush**

Adult  
Assorted: Red, Green, Purple, Blue  
Pkg/12 ..... **579496**  
**Ea \$7.95**



**Premium Disposable Impression Trays**

Disposable impression tray made from rigid, high-impact FDA approved styrene, Advance Premium disposable impression trays assist in generating a consistent and unsurpassed impression free of distortion or bite through. Smooth, rounded edges, interior retention slots, and vertical perforation enhance material retention.

- Perforated, 12/bag**
- #1 Large Upper..... **579517**
- #2 Large Lower..... **579520**
- #3 Medium Upper..... **579521**
- #4 Medium Lower..... **579523**
- #5 Small Upper..... **579524**
- #6 Small Lower..... **579525**
- #7S Upper Right/Upper Left..... **579532**
- #8s Upper Left/Upper Right..... **579533**
- #9S Lower Anterior ..... **579534**
- #10S Upper Anterior..... **579535**  
**Bag \$9.50**

- Solid, 12/bag**
- #7 Upper Right/Lower Left..... **579526**
- #8 Upper Left/Upper Right..... **579528**
- #9 Lower ..... **579530**
- #10 Upper ..... **579531**  
**Bag \$9.50**

**Advance Ferric Sulfate**

- 15.5% Ferric Sulfate Hemostatic Gel
- Provides effective tissue management when completing impressions, restorations, crown & bridge procedures
- Will not flow onto surrounding tissue.

**Standard Kit**  
Contains: 4-1.2ml prefilled syringes, 8 brush tips ..... **579281**  
**Ea \$16.33**

**Bulk Syringe**  
Contains: 30ml syringe ..... **579460**  
**Ea \$35.17**



**Hemostatic Gel**

- Viscous gel will not flow onto surrounding tissue
- Surfactant lowers tension for better wound penetration
- Quickly controls bleeding

**Aluminum Chloride-25%**  
4/pack with tips..... **579283**  
**Ea \$16.25**

30mL syringe ..... **579284**  
**Ea \$28.41**



**Sodium Hypochlorite**

For the debridement and irrigation of root canals during and after instrumentation

16oz Bottle ..... **579267**  
**Ea \$11.38**

# Advance BASIC Great Value on Everyday Products!



## BASIC Level 1 Earloop Face Mask

Level 1 protection for low-risk procedures. Ideal for patient appointments and exams.

- Inner and outer layer constructed for fluid resistance, with ultrasonically sealed, soft, stretchy, latex-free earloops for superior comfort.
- Adjustable nosepiece for custom fit and easy breathing.

### Level 1 Mask - Box/50

Blue .....	575007
Pink .....	575006
Lavender .....	575008
<b>Ea \$13.16</b>	

## BASIC Level 3 Earloop Face Mask

### Level 3 Mask - Box/50

White .....	575009
<b>Ea \$15.60</b>	



## BASIC Non-Woven Sponges

- Soft, viscose/poly blend, non-woven material for increased patient comfort.
- Economical alternative to cotton.
- No chemical additives to avoid allergic reactions.
- Non-sterile

### Non-Sterile - 2" x 2"

Case/5000 .....	575070
<b>Ea \$39.99</b>	

### Non-Sterile - 4" x 4"

Case/2000 .....	575071
<b>Ea \$43.99</b>	



## BASIC Cotton Rolls

Pliable, non-linting super absorbent rolls made from 100% cotton. They do not adhere to mucous membranes and will retain their shape and strength when wet.

### Non-Sterile #2 Medium

Case/2000 .....	577297
<b>Ea \$18.75</b>	



## BASIC Self-Sealing Sterilization Pouches

- Dual External Indicators, Built-in Internal Indicator
- Blue-Tinted, Transparent Film
- Triple Seal & Heat Processed
- Wide Self-Seal Adhesive Strip

### Box/200

2.25" x 4" .....	576806
<b>Ea \$5.99</b>	
2.75" x 9" .....	576807
<b>Ea \$9.99</b>	
3.5" x 5.25" .....	576811
<b>Ea \$10.48</b>	
3.5" x 9" .....	576812
<b>Ea \$9.75</b>	
5.25" x 10" .....	576814
<b>Ea \$14.25</b>	
7.5" x 13" .....	576815
<b>Ea \$32.67</b>	

### Economy Pkg/1000

2.75" x 9" .....	576809
<b>Ea \$45.75</b>	
3.5" x 9" .....	576813
<b>Ea \$46.50</b>	

## BASIC Bracket Tray Covers

High quality paper tray cover.



### Ritter - Size "B"

(8.5" x 12.25") Case/1000

White .....	561033
<b>Ea \$24.86</b>	



## Basic Curved Handle

Adult 38 Soft Tuft Rounded

Assorted: Dark Blue, Dark Purple, Dark Green, Orange

Box/72 .....	578056
<b>Ea \$29.93</b>	



## Basic Compact Head Straight Handle

Young Adult 26 Soft Tuft

Assorted: Orange, Red, Dark Blue, Aqua

Box/72 .....	578053
<b>Ea \$28.81</b>	



## Basic Child's Curved Handle

Ages 3-5 23 Soft Tuft

Assorted: Red, Green, Blue, Purple, White

Box/72 .....	578055
<b>Ea \$19.15</b>	




## Basic Straight Handle

Adult 39 Rounded Tuft

Assorted: Red, Dark Blue, Dark Purple, Dark Green

Box/72 .....	578054
<b>Ea \$27.36</b>	



## Basic Twist

Compact Adult 27 Tuft

Assorted: Red, Green, Blue, Purple

Box/72 .....	578064
<b>Ea \$31.15</b>	

Advance BASIC  
**GREAT VALUE**  
on  
Everyday Products!



# Advance

# ENDODONTICS

Each Advance Endodontic NiTi File system includes matching:

- Gutta Percha
- Paper Points

**Advance TAPER HT**  
HEAT TREATED NITI ROTARY FILES

**\$26.92** Pk/6

Compare to:  
ProTaper® Gold and EdgeTaper Platinum™



**Advance TAPER**  
NITI ROTARY FILES

**\$26.92** Pk/6


Compare to:  
ProTaper® Universal and EdgeTaper™



**Advance ONE HT**  
HEAT TREATED NITI ROTARY FILES

**\$18.33** Pk/3    **\$24.43** Assorted Pk/4


Compare to:  
Wave One Gold® and EdgeOne Fire™



**Advance ONE**  
NITI ROTARY FILES

**\$18.33** Pk/3    **\$24.43** Assorted Pk/4

Compare to:  
Wave One® and EdgeOne Platinum™





- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with **ProTaper® Gold and Edge Taper Platinum™** file systems for shaping and finishing
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Multiple tapers and the Heat Matrix technology of the Advance Taper HT ensure added flexibility while easily cutting through dentin in specific canal zones
- Shaping files have multiple tapers that ensure flexibility and cut dentin in specific canal zones with greater efficiency



- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with **ProTaper® Universal and Edge Taper™** file systems for shaping and finishing
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Shaping files have multiple tapers that ensure flexibility and cut dentin in specific canal zones with greater efficiency



**Heat Treated**  
**Better Flexibility**

Order#	Size	Taper	Length	Pk Size
578940	SX	04	19MM	6
578915	ASST	X	21MM	6
578921	F2 25	07	21MM	6
578294	F3 30	08	21MM	6
578927	F4 40	09	21MM	6
578930	F5 50	06	21MM	6
578933	S1 18	05	21MM	6
578936	S2 20	02	21MM	6
578916	ASST	X	25MM	6
578919	F1 20	07	25MM	6
578922	F2 25	08	25MM	6
578925	F3 30	09	25MM	6
578928	F4 40	06	25MM	6
578931	F5 50	05	25MM	6
578934	S1 18	02	25MM	6
578938	S2 20	04	25MM	6
578917	ASST	X	31MM	6
578920	F1 20	07	31MM	6
578923	F2 25	08	31MM	6
578926	F3 30	09	31MM	6
578929	F4 40	06	31MM	6
578932	F5 50	05	31MM	6
578935	S1 18	02	31MM	6
578939	S2 20	04	31MM	6

Order#	Size	Taper	Length	Pk Size
578981	SX	04	19MM	6
578957	ASST	X	21MM	6
578961	F2 25	07	21MM	6
578964	F3 30	08	21MM	6
578986	F4 40	09	21MM	6
578969	F5 50	06	21MM	6
578972	S1 18	05	21MM	6
578975	S2 20	02	21MM	6
578958	ASST	X	25MM	6
578962	F1 20	07	25MM	6
578965	F2 25	08	25MM	6
578967	F3 30	09	25MM	6
578970	F4 40	06	25MM	6
578973	F5 50	05	25MM	6
578976	S1 18	02	25MM	6
578979	S2 20	04	25MM	6
578960	ASST	X	31MM	6
578963	F1 20	07	31MM	6
578966	F2 25	08	31MM	6
578968	F3 30	09	31MM	6
578971	F4 40	06	31MM	6
578974	F5 50	05	31MM	6
578977	S1 18	02	31MM	6
578980	S2 20	04	31MM	6



- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with **Wave One® Gold and Edge One Fire™** file systems
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Multiple tapers and the Heat Matrix technology of the Advance One HT ensure added flexibility while the optimized tip promotes safe canal navigation
- Finishing files ensure the proper apical fit for the matching Advance One or your existing Gutta Percha points



- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with **Wave One®** file system
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Finishing files ensure the proper apical fit for the matching Advance One or your existing Gutta Percha points



Heat Treated	Order#	Size	Length	Pk Size	
	578823	ASST	21MM	4	
	578833	SM	21MM	3	
	578830	PR	21MM	3	
	578827	MD	21MM	3	
	578824	LG	21MM	3	
	<hr/>				
	578822	ASST	25MM	4	
	578834	SM	25MM	3	
	578831	PR	25MM	3	
578828	MD	25MM	3		
578825	LG	25MM	3		
<hr/>					
578835	SM	31MM	3		
578832	PR	31MM	3		
578829	MD	31MM	3		
578826	LG	31MM	3		

Better Flexibility	Order#	Size	Length	Pk Size	
	578780	ASST	21MM	4	
	578847	SM	21MM	3	
	578844	PR	21MM	3	
	578840	MD	21MM	3	
	578837	LG	21MM	3	
	<hr/>				
	578781	ASST	25MM	4	
	578848	SM	25MM	3	
	578845	PR	25MM	3	
578841	MD	25MM	3		
578838	LG	25MM	3		
<hr/>					
578849	SM	31MM	3		
578846	PR	31MM	3		
578842	MD	31MM	3		
578839	LG	31MM	3		

## Advance ENDODONTICS

100% guarantee on all of our products!

## Advance TAPER ABSORBENT PAPER POINTS



## Advance TAPER GUTTA PERCHA POINTS



Each Advance Taper NiTi File system includes matching:

- Paper Points
- Gutta Percha

- Size-matched to the Advance Taper & Advance Taper HT File System
- Superior absorbency

- Size-matched to the Advance Taper & Advance Taper HT File System
- Excellent radiopacity


Paper Points	Size	Order #	PK Size
SX-F3 PP ASST	SX-F3/28	578897	60
F1 PP 28 MM	F1/28	578902	60
F2 PP 28 MM	F2/28	578905	60
F3 PP 28 MM	F3/28	578907	60
F4 PP 28 MM	F4/28	578908	60
F5 PP 28 MM	F5/28	578910	60

Gutta Percha	Size	Order #	PK Size
SX-F3 GP ASST	SX-F3/28	578878	60
F1 GP 28 MM	F1/28	578881	60
F2 GP 28 MM	F2/28	578884	60
F3 GP 28 MM	F3/28	578886	60
F4 GP 28 MM	F4/28	578888	60
F5 GP 28 MM	F5/28	578890	60

## Advance ONE ABSORBENT PAPER POINTS



## Advance ONE GUTTA PERCHA POINTS



Each Advance One NiTi File system includes matching:

- Paper Points
- Gutta Percha

- Size-matched to the Advance One & Advance One HT File System
- Millimeter-marked length indicator lines
- Superior absorbency

- Size-matched to the Advance One & Advance One HT File System
- Millimeter-marked length indicator lines
- Excellent radiopacity
- Color coded size identification

Paper Points	Size	Order #	PK Size
ADV ONE PP	ASST	578855	60
ADV ONE PP	SMALL	578866	60
ADV ONE PP	PRIMARY	578863	60
ADV ONE PP	MEDIUM	578861	60
ADV ONE PP	LARGE	578859	60

Gutta Percha	Size	Order #	PK Size
ADV ONE GP	ASST	578796	60
ADV ONE GP	SMALL	578808	60
ADV ONE GP	PRIMARY	578805	60
ADV ONE GP	MEDIUM	578803	60
ADV ONE GP	LARGE	578801	60

# 100% guarantee on all of our products!

Advance dental products carry our 100% no-questions-asked, money-back guarantee. We want you to be completely satisfied with your purchase of Advance dental products. If for any reason you are not completely satisfied with an Advance dental product, we will completely refund your purchase price.

# Optimize Tissue Recovery Time with an Anti-Bacterial Endodontic Sealer

Advance Dental Products **Bioceramic Endodontic Sealer** combines the proven advantages of **MTA (Mineral Trioxide Aggregate)** with a superior canal obturation product. Our bioceramic sealer is as biological as the actual pulp tissue, and is engineered to penetrate, fill, and seal main and lateral root canals. The 13% MTA based sealer maintains an optimal pH, which elicits anti-bacterial effects and encourages new tissue formation and recovery.

- Provides rapid recovery of tissue without causing inflammatory reaction
- Avoids tooth stain and discoloration; eugenol and bismuth oxide free
- Easy to use dual syringe automix
- Easy removal for retreatment
- No need for low temperature gutta percha cone use standard heat settings for warm vertical
- Working Time: 23 minutes, Set Time: 130 minutes

Up to  
**20% - 40%**  
savings over comparable  
leading brands  
per gram!



**\$119.94** Each  
578778 - 1 Ea. 4g Dual Barrel Syringe  
**\$15.32** Package  
578779 - 10 Pk. Auto Mix Tips

**Advance RCS**  
BIOCERAMIC ROOT CANAL SEALER

## Work with Advance Taper and Advance Taper HT NiTi Rotary Files

### Tri Auto ZX2

by Morita

- » Cordless micromotor with built-in apex locator
- » Apical patency by rotary instrumentation
- » Now possible with rotary files as small as 10#
- » Observe the file tip position with apex locator function

577194 - Tri Auto ZX2

EACH  
**\$2399.00**



### Endo-Mate TC2 Cordless Endodontic Micromotor


by NSK

- » Cordless endodontic handpiece with torque control
- » Large LCD panel
- » Access friendly Key, Auto Reverse and Alarm Function
- » Apex locator capability
- » Contains: MP-F16R Head, Motor, and Charger

571894 - Endo-Mate TC2 Cordless Endodontic Micromotor

EACH  
**\$1338.00**





## More Flexibility for Half The Cost

- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with Vortex™ Blue and Edge Sequel™ Sapphire file systems
- Creates reproducible access and predictable canal patency
- Increased torque strength for safe and predictable canal navigation while the unique design enhances tactile experience
- Guide files polish dentin surfaces for easier identification of hidden canals

Order#	Size	Length	Pk Size
578810	15/.04	21MM	4
578812	15/.04	25MM	4
578814	15/.04	30MM	4
578816	20/.04	21MM	4
578817	20/.04	25MM	4
578818	20/.04	30MM	4
578819	25/.04	21MM	4
578820	25/.04	25MM	4
578821	25/.04	30MM	4
578851	30/.04	21MM	4
578852	30/.04	25MM	4
578853	30/.04	30MM	4
578854	35/.04	21MM	4
578856	35/.04	25MM	4
578858	35/.04	30MM	4
578860	40/.04	21MM	4
578862	40/.04	25MM	4
578864	40/.04	30MM	4
578865	45/.04	21MM	4
578867	45/.04	25MM	4
578868	45/.04	30MM	4
578869	50/.04	21MM	4
578871	50/.04	25MM	4
578872	50/.04	30MM	4
578873	.04 (15-40)	25MM	4

Order#	Size	Length	Pk Size
578874	15/.06	21MM	4
578875	15/.06	25MM	4
578876	15/.06	30MM	4
578877	20/.06	21MM	4
578879	20/.06	25MM	4
578877	20/.06	30MM	4
578883	25/.06	21MM	4
578885	25/.06	25MM	4
578887	25/.06	30MM	4
578889	30/.06	21MM	4
578891	30/.06	25MM	4
578892	30/.06	30MM	4
578893	35/.06	21MM	4
578894	35/.06	25MM	4
578895	35/.06	30MM	4
578896	40/.06	21MM	4
578896	40/.06	25MM	4
578900	40/.06	30MM	4
578901	45/.06	21MM	4
578903	45/.06	25MM	4
578904	45/.06	30MM	4
578906	50/.06	21MM	4
578909	50/.06	25MM	4
578911	50/.06	30MM	4
578912	.06 (15-40)	25MM	4

\*Market research and test studies available upon request. 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, we will issue a full refund of the purchase price through your distributor. ProTaper® Gold, ProTaper® Universal, Vortex™ Blue, Wave One®, and Wave One® Gold are registered trademarks of Dentsply Sirona. Edge One Fire™ Edge Sequel™ Sapphire, Edge Taper™, and Edge Taper Platinum™ are registered trademarks of Edge Endo.



### Cobalt TWIRL

ABSORBENT PAPER POINTS

- Size-matched to the Advance Cobalt File System
- Superior absorbency



### Cobalt TWIRL

GUTTA PERCHA POINTS

- Size-matched to the Advance Cobalt File System
- Excellent radiopacity



Each Advance Cobalt Twirl NiTi File system includes matching:  
 - Paper Points  
 - Gutta Percha

Paper Points	Gutta Percha	Size	PK Size
578937	578787	15/.04	60
578941	578788	15/.06	60
578942	578789	20/.04	60
578943	578790	20/.06	60
578944	578791	25/.04	60
578945	578792	25/.06	60
578946	578793	30/.04	60
578947	578794	30/.06	60
578948	578795	35/.04	60
578949	578797	35/.06	60
578950	578798	40/.04	60
578951	578799	40/.06	60
578952	578800	45/.04	60
578953	578802	45/.06	60
578954	578804	50/.04	60
578955	578806	50/.06	60
578956	578807	ASST .04	60
578959	578809	ASST .06	60



**Advance**

Advance dental products carry our 100% no-questions-asked, money-back guarantee.



### Dust Free Alginate

- Mint
- Color changing and dust free.
  - Easy to use, avoiding remakes.
  - Shorter duration in the mouth offers maximum patient comfort.
  - Pleasing mint flavor.

**1 lb. Bag - Color Changing**  
Fast Set Color Changing (Mint)..... **579373**  
Ea **\$11.38**

**1 lb. Bag**  
Regular Set (Mint)..... **579372**  
Fast Set (Mint)..... **579371**  
Ea **\$11.38**

### Alginate Substitute

An A-silicone replacement impression material that is ready to use and requires no messy mixing of powders and liquids. Ideal for preliminary impressions, orthodontic and study models and provisional and temporary impressions. Dimensionally stable when cured, eliminates the need to pour immediately and able to do multiple pours without loss of accuracy. Can be disinfected using standard methods..

**50mL - Pkg/8**  
Mint (Fast Set)..... **579985**  
Ea **\$53.92**

**380mL - Pkg/2**  
Mint (Fast Set)..... **579986**  
Ea **\$84.44**

	Work Time	IO Set Time
Fast Set (50mL)	1:00	1:00
Fast Set (380mL)	1:15	1:00



### VPS Impression Material - Putty

This hand kneadable addition cure silicone material is an ideal tray material for a 2-step putty/wash technique. This easy to mix, non-oily material when cured hardens to a high durometer which helps resist distortion and provides excellent hydraulic force when combined with a light body material.

- Putty Kit contains:
- 300 mL Base Jar
  - 300 mL Catalyst Jar.

**300mL**  
Putty Kit (Fast Set) ..... **579988**  
Pack **\$44.34**

	Work Time	IO Set Time
Fast Set	1:00	2:00

### VPS Impression Material

A hydrophilic impression material with super wetting agents that improve wettability resulting in a highly accurate impression. Excellent tear strength and dimensional stability also help to reduce the chances for distortion. Available in multiple viscosities, flavored and non-flavored in both fast and regular set.

**Light Body - 50mL - Pkg/4**  
Mint (Fast Set)..... **579973**  
Vanilla (Fast Set) ..... **579974**  
Unflavored (Fast Set) ..... **579975**  
Unflavored (Regular Set). ..... **579976**

**Regular Body - 50mL - Pkg/4**  
Mint (Fast Set)..... **579977**  
Vanilla (Fast Set) ..... **579978**  
Unflavored (Fast Set) ..... **579979**  
Unflavored (Regular Set). ..... **579980**

**Heavy Body - 50mL - Pkg/4**  
Mint (Fast Set)..... **579981**  
Vanilla (Fast Set) ..... **579982**  
Unflavored (Fast Set) ..... **579983**  
Unflavored (Regular Set). ..... **579984**  
Pack **\$34.56**

**Heavy Body - 380mL - Pkg/2**  
Unflavored (Fast Set) ..... **579987**  
Pack **\$99.62**

FAST SET	Work Time	IO Set Time
Light Body (50mL)	1:00	1:00
Medium (50mL)	1:00	1:00
Heavy Body (50mL)	1:00	1:00
Heavy Body (380mL)	1:30	1:00

REGULAR SET	Work Time	IO Set Time
Light Body (50mL)	2:30	1:00
Medium (50mL)	2:30	1:00
Heavy Body (50mL)	2:30	1:00

Contact your  
Goetze Dental  
Representative  
to order!  
Or Order at:  
**800.692.0804**

**NEW**  
**Advance**  
INSTRUMENTS

For more information about Advance Dental Products contact your Goetze Dental Representative or call 1.800.692.0804

- Single Ended - #23 Shepherd's Hook 581395** **Ea \$9.88**
- Double Ended - #2 581396** **Ea \$10.50**
- Double Ended - #5 581397** **Ea \$10.50**

**Explorers**

Explorers are used to detect caries and calculus. The knurled handle provides a solid, ergonomic grip for greater tactile sensitivity.

- #CP11 Color Coded - 3,6,8,11 581398**
- #CP12 Color Coded - 3,6,9,12 581399**
- Williams Color Coded - 1,2,3,5,7,8,9,10 581400**

**Expos**

Expos have an explorer on one end and a probe on the other. This ergonomic design allows the practitioner to examine and evaluate the oral cavity with more precision without putting the instrument down.

**Ea \$17.38**

- #CP12 Color Coded - 3,6,9,12 581401**
- Williams Round - 1,2,3,5,7,8,9,10 581402**

**Probes**

Probes are used to assess pocket depth and measure dental anomalies and lesions during diagnostic examinations.

**Ea \$15.44**



- Front Surface Cone Socket**
- #4 - 581407**
- #5 - 581409**

- Front Surface Simple Stem**
- #4 - 581408**
- #5 - 581410**

**Mirrors**

Front surface mirrors give a true reflection, eliminating the ghost or double image seen on plane front mirrors. The stainless steel housing on this mirror offers a longer life than typical brass housings.

**12/Pack \$24.75**

- Octagonal Mirror Handle - Cone Socket 581405**
- Octagonal Mirror Handle - Simple Stem 581406**
- Round Mirror Handle - Cone Socket 581403**
- Round Mirror Handle - Simple Stem 581404**

**Mirror Handles**

Our solid construction knurled mirror handles are the perfect complement to our front surface mirrors, offering a tight fit and tactile connection to the instrument.

**Ea \$2.75**

- Serrated 6" College Pliers - Serrated 6" 581392**
- Self-Locking 6" College Pliers - #317A X-Fine Serrated Tip 581393**
- Narrow Tip 6" College Pliers - #317B 581394**

**Ea \$8.14**   **Ea \$8.11**   **Ea \$8.25**

**College Pliers**

Our college pliers can be used to grasp gauze or other materials from the tray or oral cavity. Precision tips and our two-stage guide pin ensure accuracy and prevent glove punctures when closed.

**NEW**  
**Advance**  
**PRODUCTS!**



**LC Nano Hybrid Flowable Composite**

Light-cured, flowable, 80% nanohybrid composite with a 3.5mm depth of cure. Natural aesthetics with lower shrinkage, the wear resistance of a hybrid, and excellent color retention. Radiopaque.\*

\*Application requires layering technique. Light-cure layer by layer in 2mm increments.

- Shade A1 - 1.8gr Syringe with 5 Tips..... **581809**
- Shade A2 - 1.8gr Syringe with 5 Tips..... **581810**
- Shade A3 - 1.8gr Syringe with 5 Tips..... **581811**  
**Ea \$17.78**

**Universal Bond Plus & Universal Bond Plus DCA**

Light-curing universal adhesive. Can be applied with or without phosphoric acid etching. Suitable for direct or indirect restorations. Provides secure adhesion, without the need for additional primer, for various materials such as metals, zirconium, or aluminum oxide, as well as silicate ceramics. Dualcure activator for self-curing or dual curing. Applied in one layer: only 35 seconds total working time.

- Universal Bond Plus - 5ml Bottle ..... **581812**  
**Ea \$99.15**
- Universal Bond Plus DCA - 2ml..... **581813**  
**Ea \$48.68**

**Advance Glutaraldehyde 2.6% Sterilant**

- Gallon ..... **581866**  
**Ea \$33.57**
- Indicator Test Strips 40/bottle ..... **581867**  
**Ea \$52.35**

**Temporary Luting Cement**

High-retentive, eugenol-free, zinc-oxide based temporary luting cement containing calcium hydroxide. Indicated for luting temporary provisional restorations and temporary definitive restorations.

- 5ml Syringe with 10 Tips..... **581814**  
**Ea \$33.58**



## Elevate your practice with A-dec® and KaVo®.

Now, two of the most respected names in the dental industry are collaborating to bring you unparalleled performance with effortless integration.

KaVo's impressive range of air-driven, electrical and hygiene handpieces are the ideal match for the legendary quality and reliability of A-dec dental chairs, delivery systems, dental lights and dental furniture. Enjoy an exceptional practice experience, with premier equipment that reflects your commitment to superior patient care.

**That's Excellence Amplified.**

Visit [a-dec.com/ExcellenceAmplified](https://a-dec.com/ExcellenceAmplified) for details, then contact your authorized A-dec representative to order.

# FEATURED PRODUCTS

## Crest® Whitening Emulsions with LED accelerator light

by Crest® + Oral-B®

### You've never seen teeth whitening like this!

The all-new Crest Whitening Emulsions is a breakthrough formula featuring a thin layer of 5x more active peroxide droplets\* suspended in a hydrating base to whiten teeth with virtually NO sensitivity. This means it starts working instantly and stays on teeth 10X longer, with no need to rinse or brush off. Just squeeze onto the wand, apply on your teeth, and boost with the LED Accelerator Light for even better teeth whitening results. Crest Whitening Emulsions is the best way to whiten. It removes stains, prevents new ones and leaves your teeth 100% whiter.

#### Crest Whitening Emulsions Features:

- » **BEST IN CLASS RESULTS:** Using first-of-its-kind technology, Crest Whitening Emulsions delivers 5x more active hydrogen peroxide compared to other whitening gels and pens.
- » **VIRTUALLY NO SENSITIVITY:** Our breakthrough formula features powerful peroxide droplets and a water-resistant hydrating base to protect teeth from sensitivity pain.
- » **SWIPE AND GO CONVENIENCE:** 1 Minute Apply & Go - Crest Whitening Emulsions can be applied in seconds and starts working instantly for a teeth whitening experience that's easier than ever. Get a better, faster, 100% whiter smile by including Crest Whitening Emulsions in your beauty and grooming routine
- » **BOOST WITH LIGHT:** Our LED accelerator light safely penetrates the enamel surface to weaken tough, yellow stains, for even better teeth whitening results

Returns for all Crest Whitening Emulsions products are not accepted due to temperature control requirements.



582400 - Crest® Whitening Emulsions with LED accelerator light

EACH  
**\$69.22**

**BUY 6  
GET 1**

CREST® WHITENING EMULSIONS WITH LED ACCELERATOR LIGHT  
**BUY 6 CASES, GET 1 CASE OF CREST® ORAL-B® SENSITIVE MANUAL SOLUTION FREE!**  
3 Kits/Case. See page 76 for details.



P&G Professional Oral Health

# Crest® + Oral-B® iO Transformational Gum Health Electric Rechargeable System

by Crest® + Oral-B®

## Transform Patients' Gum Health, Starting Day 1

Clinically proven to transform your patients' gum health starting day 1. In a clinical study, 100% of patients moved to a state of healthy gums after just 12 weeks of use. It delivers Superior Gum Health\* with an Irresistible, WOW Brushing Experience. Let your patients see for themselves what sets Oral-B® iO™ apart. With its breakthrough features, the Oral-B® iO™ Electric Rechargeable Toothbrush delivers a WOW experience your patients can't resist, and an exceptional clean for healthier teeth and gums\*\*. The innovative magnetic drive delivers concentrated energy where it's needed most, allowing the brush head to glide tooth by tooth for a smooth, quiet cleaning experience. The Oral-B® iO™ brush head combines oscillation rotations with the gentle energy of micro-vibrations for an enjoyable experience and a sensational clean. A Smart Pressure Sensor lights up green when optimal pressure is applied and red when the pressure is too hard. It includes Artificial intelligence technology with an easy-to-navigate interactive display with 5 brushing modes for your best clean in all areas of the mouth complete with a pro-timer to encourage brushing for at least 2 minutes. Powers a long-lasting charge in just 3 hours.

\*vs. a regular manual toothbrush and regular anti-cavity toothpaste  
 \*\*vs. a regular manual toothbrush



### Bundle Includes

- » Oral-B® iO Electric Rechargeable Toothbrush  
 5 Modes: Daily Clean, Sensitive, Whitening, Gum Care & Tongue Clean, Only Oral-B® iO™ Refills will fit on the Oral-B® iO™ toothbrush.
- » Oral-B® Glide® PRO-HEALTH Advanced Floss (15m)  
 Our most advanced floss is designed for patients with mild to moderate gingivitis
- » Crest® PRO-HEALTH Gum Detoxify Deep Clean Toothpaste (4.1oz)  
 Elevate expectations even around the gumline with Crest Gum Detoxify toothpaste. Get clinically proven healthier gums\* with stannous fluoride. \*helps prevent gingivitis
- » Crest® PRO-HEALTH™ Multi-Protection Mouthwash (250ml)  
 Provides 12-hour protection against plaque and gingivitis. Kills germs without the burn of alcohol
- » Replacement Brush Head Refills: 1 Oral-B® iO Series Ultimate Clean Refill and 1 Oral-B® iO Series Gentle Care Refill
- » Reusable Tote Bag

582397 - Bundle

BUNDLE  
**\$190.00**

**BUY 6 GET 1**  
 CREST® + ORAL-B® iO™ TRANSFORMATIONAL GUM HEALTH SYSTEM  
**BUY 6 CASES, GET 1 CASE OF CREST® ORAL-B® SENSITIVE MANUAL SOLUTION FREE!**  
 3 Bundles/Case. See page 76 for details.





582399  
 Oral-B Gentle Care iO  
 refill brush head

EACH  
**\$11.06**



582398  
 Oral-B Ultimate Clean  
 iO refill brush head

EACH  
**\$11.06**

# FEATURED PRODUCTS



## ORCA™ FOAM

by ORCA™ Professional

Absorbable Gelatin Sponge Hemostat, USP (Standard)  
ORCA Foam is an absorbable gelatin sponge hemostat. It's natural, safe, and stops bleeding effectively. ORCA Foam works great as a carrier medium for ORCA Dry Socket.

- » Natural & Safe
- » Stops bleeding effectively
- » Fast, complete biodegradation

Foam Hemostatic Sponges  
**582056** - 10MM 09, 50 individually blister packaged, sterile cubes (1cm x 1cm x 1cm each)

50/BOX  
**\$300.00**

## ORCA™ Dry Socket

by ORCA™ Professional

ORCA Dry Socket helps reduce pain from dry sockets. As with all ORCA products, it is non-medicated and made with natural essential oils. Clinically proven. Superior results.



Dry Socket Liquid  
**582058** - 15ML 10993

BOTTLE  
**\$135.82**



# Get More From Less

## Goodbye complexity with 3M™ RelyX™ Universal Resin Cement & 3M™ Scotchbond™ Universal Plus Adhesive

3M™ RelyX™ Universal Resin Cement works both as a standalone, self-adhesive cement and as an adhesive cement when combined with 3M™ Scotchbond™ Universal Plus Adhesive. Together they eliminate the hassle and confusion of multiple resin cements, primers and adhesives, simplifying your direct and indirect restorative workflows. With fewer products to stock and replenish, you save space, time and money.

### RelyX™ Universal Resin Cement

by 3M™

Goodbye complexity

3M™ RelyX™ Universal Resin Cement comes in a unique self-sealing automix syringe with minimized paste waste. It covers virtually all adhesive and self-adhesive resin cement indications with excellent bond strength, high esthetics and easy excess clean-up.

3M™ RelyX™ Universal Resin Cement works both as a standalone, self-adhesive cement and as an adhesive cement when combined with 3M™ Scotchbond™ Universal Plus Adhesive. Together they eliminate the hassle and confusion of multiple resin cements, primers and adhesives, simplifying your direct and indirect restorative workflows. With fewer products to stock and replenish, you save space, time and money. Efficiently treat routine cases with RelyX Universal Resin Cement offering outstanding self-adhesive bond strength

3M™ RelyX™ Universal Resin Cement is indicated for the final cementation of:

- » Veneers and occlusal veneers (tabletops)
- » Inlays and onlays
- » Crowns
- » Bridges
- » Maryland- and inlay/onlay-bridges
- » Posts and screws



Features and Benefits:

- » Simplifies workflows - saves space, time and money
- » Unique self-sealing, hygienic automix syringe with minimized waste
- » Easy excess clean-up
- » Excellent bond strength in self-cure and light-cure
- » High esthetics by 4 color stable, fluorescent shades

#### RelyX™ Universal Resin Cement Trial Kit

- 1 - 3.4g syringe
- 15 - 3M™ RelyX™ Universal Micro Mixing Tips,
- 5 - Elongation Tips
- 1 - 1.5ml vial 3M™ Scotchbond™ Universal Plus Adhesive
- 1 - 3ml 3M™ Scotchbond™ Universal Etchant,
- 25 - blue etchant tips
- 50 - orange microbrushes,
- Technique guides, Instructions for Use

581348 - Translucent  
581349 - A1

TRIAL KIT  
**\$230.77**

#### RelyX™ Universal Resin Cement Intro Kit

- 4 - 3.4g syringes shades TR, A1, A30, WO
- 60 - 3M™ RelyX™ Universal Micro Mixing Tips
- 20 - Elongation Tips
- 1 - 5ml vial 3M™ Scotchbond™ Universal Plus Adhesive
- 1 - 3ml 3M™ Scotchbond™ Universal Etchant
- 25 - blue etchant tips
- 50 - orange microbrushes
- Technique guides and Instructions for Use

581347 - Kit

INTRO KIT  
**\$669.23**

Value Pack ..... \$537.92

Contains:  
3 - 3.4g syringe shade, 45 - 3M™ RelyX™ Universal Micro Mixing Tips, 15 - Elongation Tips

581354 - Translucent  
581355 - A1

Refill ..... \$199.23

Contains:  
1 - 3.4g syringe shade, 15 - 3M™ RelyX™ Universal Micro Mixing Tips, 5 - Elongation Tips

581350 - Translucent  
581351 - A1  
581352 - A30  
581353 - WO



### Scotchbond™ Universal Plus Adhesive

by 3M™

The original, now better than ever!

3M™ Scotchbond™ Universal Plus Adhesive is a radiopaque, single-component, light cured universal adhesive for all direct and indirect bonding indications and all etching techniques. It can also serve as a universal primer for all restoration materials.



3M™ Scotchbond™ Universal Plus is the

first radiopaque universal adhesive, minimizing the risk of X-ray misdiagnosis and invasive overtreatment. It provides gold standard adhesion to all dental substrates, including glass ceramic, without the need for a separate primer, silane or activator. It supports minimally invasive dentistry by bonding and sealing caries affected dentin. The formulation is free of BPA derivatives like BisGMA. And there is virtually no post-operative sensitivity.

- » Bonding all composite filling, core build-up, sealant materials and resin cements
- » System adhesive for 3M™ RelyX™ Universal Resin Cement
- » Repair of all direct and indirect restorations
- » Primer for all restoration materials
- » Desensitization and seal of dentin

Features and Benefits:

- » Dentin-like radiopacity
- » Gold standard adhesion to all tooth and restoration substrates
- » No need for a separate primer or activator
- » Bonds and seals caries affected dentin
- » Free of BPA derivatives like BisGMA
- » Virtually no post-operative sensitivity



581358 - Intro Kit - Vial ..... \$198.77  
Includes 1 Vial 5ml, 50 Disposable Brushes, 1 Syringe Etchant 3ml, 25 Dispensing Tips, Step-by-Step Cards

581359 - Refill - Vial ..... \$168.14  
Includes 1 Vial 5ml

581360 - Value 3 Pack - Vial ..... \$453.98  
Includes 3 Vials (5ml)

581361 - Introductory Kit - 50 Unit Dose ..... \$152.62  
50 L-Pop, 1 Syringe Etchant 3ml, 25 Dispensing Tips, Step-by-Step Cards

581363 - Unit Dose 100 L-Pop ..... \$337.85  
100 L-Pop, 1 Syringe Etchant 3ml, 25 Dispensing Tips, Step-by-Step Cards

581364 - Unit Dose 200 L-Pop ..... \$615.23  
200 L-Pop, 1 Syringe Etchant 3ml, 25 Dispensing Tips, Step-by-Step Card

581365 - Unit Dose 400 L-Pop ..... \$1,138.46  
400 L-Pop, 1 Syringe Etchant 3ml, 25 Dispensing Tips, Step-by-Step Cards

# FEATURED PRODUCTS

## ACERO CROWNS



### REDUCED COSTS

Our crowns are crafted using 316L stainless steel are identical in size to the competitors yet cost only HALF THEIR PRICES!

### COMPATIBILITY

Our entire line of crowns is completely interchangeable and compatible with other brands.

#### Acero Crowns Advantage

- » Simplified universal numbering system, crown label corresponds to tooth number
- » Surgical-grade 316 U.S. stainless steel
- » User-friendly tray layout
- » Secure tray case



Contact your Goetze Dental Private Practice Specialist to see a demo kit.

### ACERO CROWNS<sup>XT</sup>

#### AceroXT Crowns

- » Half the price of other brands
- » Preformed and pre-contoured
- » Surgical-grade 316 stainless steel
- » Snaps into place
- » Slightly shorter crown height

XT Refill - 5/Pack.....\$23.00 Each



### ACERO CROWNS<sup>ES</sup>

#### AceroES White Crowns

- » Cost 25% less than other leading brands
- » Preformed and pre-contoured to prevent microgaps
- » Lifelike translucency
- » 100% compatible and interchangeable with leading brands

ES Refill.....\$23.00 Each  
ESL Refill.....\$23.00 Each



We use the universal numbering system for primary teeth to identify our Crowns, the system all dentists and dental assistants are familiar with.

# 3 QUESTIONS TO CONSIDER WHEN EVALUATING CONE BEAM 3D CBCT

BY ROB HAYS  
REGIONAL SALES DEVELOPER,  
GOETZE DENTAL

The pandemic has brought change to many areas of the dental practice. With change comes the opportunity to explore new protocols, systems, and workflows. Some practices are looking to leverage technology to optimize the patient experience by expanding the scope of procedures offered, reducing the number of visits required and streamlining procedural times.

3D CBCT can play a significant role in achieving those goals, but with so many options on the market, it can be difficult to decide which manufacturer and model will be the best fit for your practice. In this article we're going to share a few questions and thoughts to get you started or help you continue down your path to choosing a 3D CBCT.

Many dentists initially consider 3D when looking to implement implantology in their practice. And certainly, 3D imaging, working in conjunction with implant planning software



ProMax® 3D Mid  
by Planmeca

provides a complete view of the bone and surgical site before you begin. From this, proper implant sizes can be determined, and surgical guides can be fabricated to consistently provide accurate placement. This has a profound impact on streamlining workflow and reducing procedure time.

But implants are just the tip of the iceberg when it comes to the impact that 3D can have on your patients and on your practice. Below are the top questions our equipment team asks every dentist who is evaluating CBCT.

# 1

## What procedures would you like to add to your practice?

How would 3D CBCT enable that and how do you calculate the impact that would have on production?

Here are a few examples of procedures that are enhanced with the use of 3D:

- Implantology
- Endodontics
- Orthodontics
- Sleep and Airway
- TMJ Treatment

### BONUS

- Enhanced diagnosis
- Wow factor with patients for improved case acceptancs



Orthopantomograph™ OP 3D Pro  
by Kavo



Veraviewepocs 3D R100  
by J. Morita

# 2

## What volumetric size would best support the scope of procedures you provide in the practice?

If you hope to use 3D primarily for single-site implants, you may not need a full volume. Have you entertained adding sleep treatment to your practice? If so, you might consider exploring a larger field of view to accommodate airway study.

Volumetric sizes (called Field of View, or FOV) range from 3x3 to 30x30 in some systems, providing tremendous flexibility in how you use your CBCT. If you're at the stage of evaluating 3D systems with this level of detail, download our 3D Buyer's Checklist to help compare your options side-by-side.

Having a range of FOV on your 3D can also help expand or limit the scope of your inquiry (and responsibility) and help increase or reduce radiation dosages as needed as well so that you can be confident your patients are getting the the top quality of care you deliver with the highest levels of safety.

# 3

## How often do you refer a patient out to a specialist for a 3D? Are you taking 2D pans in your practice? If so, how often?

This is a fairly clear pathway to understanding the current imaging rhythm in your practice, and the number of patients accepting your recommendation for these diagnostic imaging protocols.

As an example, if you currently refer four patients per month to your specialist to take a 3D, and that specialist charges \$300 per scan, that represents \$1,200 you could retain and/or capture in your own practice. Some CBCT units allow you to switch back and forth between traditional 2D and full 3D volume, and even those systems that are exclusively 3D, you are able to export the scan.

Many practices find that they take more scans once they have a 3D on site, as it so profoundly impacts their diagnostic ability. The marginal cost of each additional rotation of the machine is \$0, so if diagnostic or treatment outcomes



**Green CT 2**  
by Vatech

will improve as a result, some dentists will opt to take an image at no charge.

If you can see more pathology more accurately, thanks to 3D CBCT, that should directly correlate to treatment that you recommend and complete.

# 3.5

## How do I evaluate all of the bells and whistles, and which features are important in my practice?

I felt the need to add a final question here, as 3D CBCT is a significant investment in your practice, and you have a lot of options.

Options like 3D facial renderings that overlay the scan for a completely digital patient could be important if you plan to take Ortho to the next level. Jaw movement tracking features have a number of applications for restorative dentistry and TMJ treatment. CAD/CAM integration is increasingly utilized for implant planning and design and fabrication of surgical guides.

Other questions to consider are those surrounding training, warranty, upgrade pathways and serviceability. Once you have a broad picture of your goals using the questions above, an expert Goetze Dental Equipment Specialist can help you filter through the details and work collaboratively to pair you with the manufacturer and model that best fits your needs.

## Ready to take the next step?

Contact your Goetze Dental Private Practice Specialist or call 800-692-0804.

Your specialist is happy to meet in person or virtually.

VISIT [GOETZEDENTAL.COM/3D](https://www.goetzedental.com/3d) TO VIEW DETAILS ON 3D CBCT MANUFACTURERS, BRANDS AND MODELS



**COMPARE YOUR OPTIONS.** Scan this QR code to Download our 3D CBCT Comparison Chart worksheet. Use it to capture and compare, side-by-side, the features and benefits of the brands and models you are considering.

# 3M™ Paradigm™ Brand: An exceptional value – from the company that stands for quality.


## 3M™ Paradigm™ Nano Hybrid Universal Restorative

- Easy to polish
- 12 shades that cover most dentists' daily needs
- Fluoresces
- Great handling without slumping or sticking to instruments
- Easy to use in posterior and anterior restorations
- Good fracture toughness, strength and wear resistance
- Easy to adapt to cavity
- Radiopaque



## 3M™ Paradigm™ VPS Impression Material

- Excellent detail reproduction
- Hydrophilic material
- Excellent flow of the wash material
- Optimal working and setting times
- Good tear strength for both regular and low viscosity washes
- Excellent dimensional stability



## 3M™ Paradigm™ Temporization Material

- Three of the most popular shades: A1, A2, A3 match 3M™ Filtek™ Supreme Flowable Restorative
- Strong, fracture resistant, and high toughness
- Smooth glossy surface requires minimal polishing
- Nanofiller technology provides strength and esthetics
- Temporaries that look natural in most conditions



**3+1 Buy Any 3,  
Get 1 Free!\***

**Offer Valid July 1 – September 30, 2021**

\*Free goods must be of equal or lesser value of lowest price of purchased product. See back for details.


## Adhesive Dentistry

### 3M™ Composites

- 3M™ Filtek™ Bulk Fill Flowable Restorative
- 3M™ Filtek™ One Bulk Fill Restorative
- 3M™ Filtek™ Supreme Flowable Restorative
- 3M™ Filtek™ Supreme Ultra Universal Restorative
- 3M™ Filtek™ Universal Restorative
- 3M™ Paradigm™ Nano Hybrid Universal Restorative

### 3M™ Adhesives

- 3M™ Scotchbond™ Universal Plus Adhesive
- 3M™ Scotchbond™ Universal Adhesive



**2x2+1x1**

Buy 2 Composites and Buy 2 Adhesives,  
Get 1 of Each Free!\*

## CAD/CAM and Cements

- 3M™ Chairside Zirconia
- 3M™ RelyX™ Dental Cements



**1x1+1**

Buy Any 1 RelyX Dental Cement and Buy Any 1 Chairside Zirconia, Get 1 of the Same Chairside Zirconia Free!\*

## Curing Lights

- 3M™ Elipar™ DeepCure-S LED Curing Light
- 3M™ Paradigm™ DeepCure LED Curing Light



**2+1**

Buy 2 of the Same, Get 1 of the Same Free!

## Impressioning

### 3M™ Penta™ Impression Material

- 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material
- 3M™ Impregum™ Penta™ Soft Polyether Impression Material
- 3M™ Impregum™ Penta™ Polyether Impression Material
- 3M™ Imprint™ 3 Penta™ VPS Impression Material
- 3M™ Imprint™ 4 Penta™ VPS Impression Material
- 3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material



**3+1 Penta™**

Buy Any 3, Get 1 Free!\*

### 3M™ Garant™ Impression Material

- 3M™ Impregum™ Super Quick Polyether Impression Material
- 3M™ Impregum™ Soft Polyether Impression Material
- 3M™ Imprint™ 3 VPS Impression Material
- 3M™ Imprint™ 4 VPS Impression Material
- 3M™ Imprint™ 4 Preliminary VPS Impression Material
- 3M™ Imprint™ 4 Bite VPS Registration Material



**3+1 Garant™**

Buy Any 3, Get 1 Free!\*

## Preventive

- 3M™ Clinpro™ 5000 Anti-Cavity Toothpaste



**3+1**

Buy Any 3 Cases, Get 1 Case Free!

- 3M™ Vanish™ White Varnish
- 3M™ Fast Release Varnish



**2+1**

Buy Any 2 Same Count Packs, Get 1 of the Same Free!

\*Bonus good of equal or lesser value of lowest price of purchased products.





# Predicta® Bioactive Desensitizer

**Buy 1 Predicta® Bioactive Desensitizer Three-Syringe Kit, Get 1 Predicta® Bioactive Desensitizer Single-Syringe Kit Free\***

- S625
- S626

**NEW Predicta Bioactive Desensitizer** is changing the way clinicians treat dentin hypersensitivity by using the elements that are naturally found in the tooth, (calcium and phosphate) and restoring them to the tooth in the form of hydroxyapatite plugs that quickly block the open tubules. Then, as calcium and phosphate ions continue to release, mineralization occurs to further seal the tubules.



## SMALL EQUIPMENT



### Buy 1 Parkell Ultrasonic Scaler

(TurboSENSOR®, TurboVue®, Integra™)

#### Get 1 Free Insert\*

- D560-110
- D570-110
- D680-110

### Buy 1 Digitest® 3 Pulp Vitality Tester, Get 1 Predicta Bioactive Desensitizer Single-Syringe Kit Free\*

- D655



### Buy Any 5 Inserts, Get 2 Free\*

(of equal or lesser value - mix & match)

- DLH30Tri
- DUI30
- DUI25
- DTI30
- DTI25
- DPR30
- DPL30
- DPI30
- DPI25
- DBI30
- DSG30T
- DSG30R
- DCG30Uni-PF
- DCG30Tri
- DCG30Slim
- DSG30P
- DSG30L
- DLH30U
- DLH30T
- DLH30P



## MATERIALS



### Predicta® Bioactive Bulk & Core— Buy 3, Get 1 Free\*

(of equal or lesser value - mix & match)

#### PLUS 1 Free Predicta Desensitizer Single Syringe Kit\*

- S600
- S601
- S605
- S606
- S610
- S611
- S615
- S616

### Blu-Mousse® — Buy 3, Get 1 Free\*

(mix & match)

- S440S
- S438S
- S437S
- S448S
- S457
- S436
- S457S
- S438

### Brush&Bond® — Buy 3, Get 1 Free\*

(Original & Universal kits and refills)

- S284
- S288
- S285
- S240
- S286
- S242
- S287

### Cinch® — Buy Any 3, Get 1 Free\*

(of equal or lesser value - mix & match)

- S435
- S461S
- S453S
- S463S
- S454S
- S464
- S456S
- S465
- S460S
- S466



\*To redeem your no charge goods, email a copy of your qualifying invoice dated 7/1/21 – 9/30/21 to Parkell at [dealerpromotions@parkell.com](mailto:dealerpromotions@parkell.com). No Charge Goods must be of equal or lesser value than product purchased. Please include the desired product SKU and/or description in the email. Invoices must be received by 10/31/21. Please allow 4-6 weeks for delivery of no charge goods. Offer valid in U.S. only. This offer cannot be combined with any other offer and is subject to change without notice. Excludes Tip Refills and Gun.

# Q3 2021 END USER SPECIALS

## LuxaCrown® The chairside way to long-lasting crowns Buy 3 LuxaCrown®, Get 1 FREE

- › Emergency cases
- › Implant cases
- › Pediatric cases
- › Palliative cases
- › Healthy aging adults
- › Phased treatments
- › Bruxism
- › Cosmetic cases



**LuxaCore® Z Dual**  
Buy 3 LuxaCore® Z Dual Automix or Smartmix, Get 1 of the same FREE



**Luxatemp®**  
Buy 3 Automix or Smartmix, Get 1 of the same FREE



**StatusBlue**  
Buy 3 StatusBlue MixStar or Automix Canisters, Get 1 of the same FREE



**TempoCem®ID**  
Buy 3 TempoCem®ID, Get 1 FREE



**LuxaFlow® or LuxaFlow® Ultra**  
Buy 3 LuxaFlow® or LuxaFlow Ultra Refills, Get 1 FREE



**Honigum or Honigum Pro**  
Buy 4 Honigum or Honigum Pro Heavy, Mono, Light, Heavy QuadFast, Light QuadFast, MixStar or Automix Cartridges, Get 1 of the same FREE



**Icon® Caries Infiltrant**  
Buy 1 Icon® Cube (Smooth Surface or Proximal), Get 1 Icon® Etch Refill FREE



**Kolorz ClearShield™**  
Buy 3 Kolorz® ClearShield™ Varnish (35 or 200 pack), Get 1 of the same FREE




**Kolorz Prophy Paste: Automatic 1 + 1, all year, all orders!**

To redeem free goods, fax your paid dealer invoice to DMG America at 201-894-0213 or email to [freegoods@dmg-america.com](mailto:freegoods@dmg-america.com), indicating shade of free good. Each qualifying purchase must be submitted on one invoice within 30 days. Cannot be combined with other offers and may change or be discontinued at any time without notice. Limit 5 per dental office. Not available to dental supply dealers. Valid 7/1/21 – 9/30/21.

DMG America, LLC  
65 Challenger Road, Suite 340, Ridgefield Park, NJ 07660  
tel: 800.662.6383 email: [info@dmg-america.com](mailto:info@dmg-america.com)  
web: [www.dmg-america.com](http://www.dmg-america.com)





## The quality is in the details.



Discover the possibilities with Honigum Pro

[www.dmg-america.com/honigumpro](http://www.dmg-america.com/honigumpro)

### Honigum Pro

Honigum Pro is the only VPS impression material system that:

- ↗ Goes where you want it
- ↗ Stays where you place it
- ↗ Has the unique ability to thin under load

Thanks to its true thixotropic formula, the material won't slump or run but remains stable where you place it, making it the ideal material for taking precise impressions, especially implant impressions.

Plus, with unmatched flowability, it makes its way into tiny crevices and hard-to-reach areas that others can miss, and its remarkable strength saves valuable time by preventing the need for retakes.

The better the impression, the better the restoration.

1 (800) 662-6383 | [info@dmg-america.com](mailto:info@dmg-america.com) | [dmg-america.com](http://dmg-america.com) | [dmg-connect.com](http://dmg-connect.com)



# KERR DENTAL Q3 SPECIALS

Effective July 1 - September 30, 2021

## Nexus™ RMGI

*Resin-Modified Glass Ionomer Luting Cement*

Nexus RMGI combines One-Peel™ cleanup, tack cure capability, and outstanding bond strength for consistent predictable outcomes. As the first RMGI available with advanced Nexus Technology, Nexus RMGI provides an optimal gel state and 2-3 second tack cure capability to ensure easier cleanup each and every time. Nexus RMGI utilizes a unique dual adhesive system with the highest bond strengths for an RMGI.



**BUY 2 refills of Nexus RMGI,  
GET 1 of the same AT NO CHARGE!\***

## ZenFlex™

*NiTi Rotary Shaping File*

Designed with proprietary heat treatment, ZenFlex delivers the qualities you need to handle even the most complex root anatomies.

- High cutting efficiency
- Minimally invasive (1mm maximum flute diameter)
- Ideal balance of strength & flexibility
- Controlled memory



**BUY 15 ZenFlex file packs,  
GET 5 of the same AT NO CHARGE!\***

\*Promotional goods must be of equal or lesser value and may not be combined with any other offers. Purchase must be made on one invoice between 7/1/21 and 9/30/21, unless otherwise noted. Redeem by 10/31/21. To receive your promotional goods, email or fax a copy of your authorized Kerr dealer invoice noting promo code **NQ321** and product(s) desired to: [kerrpromo@kerrdental.com](mailto:kerrpromo@kerrdental.com) or FAX: 888.727.2614. Limit to 3 redemptions. Incomplete submissions will not be processed. Allow 8 weeks for delivery. Offers valid in the 50 United States only. Promotions are subject to change or cancellation without notice. Offer void if purchased product is returned. Note that you may have an obligation under federal, state or local law to reflect discounts on product given pursuant to this promotion on any cost report forms submitted to a federal or state government or private payer who provides reimbursement for that product.

1-800-KERR-123 | [kerrdental.com](http://kerrdental.com)

MKT-21-0473\_Rev1 ©2021 Kerr Dental

**Kerr™**

# KERR DENTAL Q3 SPECIALS

Effective July 1 - September 30, 2021

## KaVo Kerr Carbides

Kerr offers a full line of over 300 high quality Operative, Trimming and Finishing, & Specialty Carbides to meet all of your needs. Our carbides go through extensive quality validation to ensure they meet or exceed international standards.

**BUY 3 packs of KaVo Kerr Carbides,  
GET 1 of the same AT NO CHARGE!\***



## CaviCide™

*Surface Disinfectant*

CaviCide is a convenient, ready-to-use, intermediate-level disinfectant which is effective against TB, HBV, HCV, HIV, viruses (hydrophilic and lipophilic), bacteria (including MRSA and VRE) and fungi. It can be used in all areas of the facility including ICU, operating rooms, isolation rooms, patient care areas and laboratories. When used as directed, it will also effectively clean and decontaminate critical and semi-critical instrumentation. CaviCide does the job of three products, thus saving money by reducing repetitive inventory and time.

1. List N: Disinfectants for Use Against SARS-CoV-2. Feb 9, 2021, [www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2](http://www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2).



**BUY 11 CaviCide 24 oz sprayer bottles,  
GET 1 of the same AT NO CHARGE!\***

\*Promotional goods must be of equal or lesser value and may not be combined with any other offers. Purchase must be made on one invoice between 7/1/21 and 9/30/21, unless otherwise noted. Redeem by 10/31/21. To receive your promotional goods, email or fax a copy of your authorized Kerr dealer invoice noting promo code **NQ321** and product(s) desired to: [kerrpromo@kerrdental.com](mailto:kerrpromo@kerrdental.com) or FAX: 888.727.2614. Limit to 3 redemptions. Incomplete submissions will not be processed. Allow 8 weeks for delivery. Offers valid in the 50 United States only. Promotions are subject to change or cancellation without notice. Offer void if purchased product is returned. Note that you may have an obligation under federal, state or local law to reflect discounts on product given pursuant to this promotion on any cost report forms submitted to a federal or state government or private payer who provides reimbursement for that product.

1-800-KERR-123 | [kerrdental.com](http://kerrdental.com)

MKT-21-0473\_Rev1 ©2021 Kerr Dental

**Kerr™**

# KERR DENTAL Q3 SPECIALS

Effective July 1 - September 30, 2021

## Restoratives



BUY 3 OptiBond™ Solo Plus,

**GET 1 of the same AT NO CHARGE!\***

(Kits or Refills)



BUY 4 refills of Harmonize™ Universal Composite,

**GET 1 of the same AT NO CHARGE!\***

## Endodontics



BUY 15 packs of NiTi Endodontic Files,

**GET 5 of the same AT NO CHARGE!\***



BUY 1 elements™ e-motion Endodontic Motor,

**GET 5 Endo NiTi and 5 Stainless Steel Handfiles AT NO CHARGE!\***

## Rotary



BUY 10 packs of NTI™ Standard or Specialty Diamond Burs,

**GET 5 of the same AT NO CHARGE!\***



BUY 3 packs NTI™ Trimming & Finishing, Operative or Operative Big Packs (100-ct), or Specialty Carbides

**GET 1 of the same AT NO CHARGE!\***

## TotalCare



BUY 4 gallons Empower™ Enzymatic Cleaning Solution,

**GET 1 of the same AT NO CHARGE!\***



BUY 4 boxes of Seal-Tight™ Air/Water Syringe Tips,

**GET 1 of the same AT NO CHARGE!\***

\*Promotional goods must be of equal or lesser value and may not be combined with any other offers. Purchase must be made on one invoice between 7/1/21 and 9/30/21, unless otherwise noted. Redeem by 10/31/21. To receive your promotional goods, email or fax a copy of your authorized Kerr dealer invoice noting promo code **NG321** and product(s) desired to: [kerrpromo@kerrdental.com](mailto:kerrpromo@kerrdental.com) or FAX: 888.727.2614. Limit to 3 redemptions. Incomplete submissions will not be processed. Allow 8 weeks for delivery. Offers valid in the 50 United States only. Promotions are subject to change or cancellation without notice. Offer void if purchased product is returned. Note that you may have an obligation under federal, state or local law to reflect discounts on product given pursuant to this promotion on any cost report forms submitted to a federal or state government or private payer who provides reimbursement for that product.

KAVO Kerr

PREFERRED LOYALTY PROGRAM

Get **DOUBLE** loyalty reward points on these featured products:



SimpliShade™ Universal Composite



Maxcem Elite™ OR Maxcem Elite™ Chroma Universal Resin Cement



ZenFlex™ NiTi Rotary Shaping File



VioNexus™ No-Rinse Spray

To view a list of all products that qualify for double reward points visit [kavokerrpreferred.com/memberoffers](http://kavokerrpreferred.com/memberoffers)

To participate in points-based offer, customers must be registered for KaVo Kerr Preferred and enter promo code "More Points" prior to being invoiced. New participants can enter the code at registration while returning members should sign into their account and enter the code on our home page. Customers can enter the promo code and register at [www.kavokerrpreferred.com](http://www.kavokerrpreferred.com). This offer cannot be combined with other similar offers unless otherwise stated. Purchase must be made from a qualifying dealer. Participating customers will receive points automatically after purchase. Offer good until 9/30/21.

1-800-KERR-123 | [kerrdental.com](http://kerrdental.com)

MKT-21-0473\_Rev1 ©2021 Kerr Dental

Kerr™

# 2021 Q3 Specials

July 1 - September 30, 2021



## OFFER #1

New Customers Only

Purchase the NXT Hg5 or NXT Hg5 mini Amalgam Separator for only \$346 (the cost of a collection container) **AND** receive a complimentary case of PowerScrub™ Vacuum Line Cleaner\* (\$782 value!)



Item#	Description	Suggested Retail	Promo Price
NXT-HG5-BDL1	NXT Hg5 Amalgam Separator	\$910.00	<b>\$346.00</b>
NXT-HG5-MINIBDL1	NXT Hg5 mini Amalgam Separator	\$910.00	<b>\$346.00</b>
NXT-HG5-HV	NXT Hg5 High Volume Amalgam Separator	\$3124.00	\$3124.00
PCS-VLCK	PowerScrub Vacuum Line Cleaner Case Kit	\$207.00	<b>\$0.00</b>

\*Rebate form submission required. NXT Hg5 must be drop-shipped from Solmetex.

## OFFER #2

Existing Customers Only

Upgrade your Hg5 Amalgam Separator to the new NXT Hg5 Amalgam Separator for only \$346 **AND** receive a complimentary case of PowerScrub™ Vacuum Line Cleaner\*



Item#	Description	Suggested Retail	Promo Price
NXT-HG5-001RPL	NXT Hg5 Amalgam Separator	\$910.00	<b>\$346.00</b>
NXT-HG5-MINIRPL	NXT Hg5 mini Amalgam Separator	\$910.00	<b>\$346.00</b>
PCS-VLCK	PowerScrub Vacuum Line Cleaner Case Kit	\$218.00	<b>\$0.00</b>

\*Rebate form submission required. NXT Hg5 must be drop-shipped from Solmetex.

## OFFER #3

Buy a case of PowerScrub™ Vacuum Line Cleaner and get a Sidekick Dispenser **FREE!**\*



PowerScrub meets all EPA requirements for pH levels between 6 and 8 and is the **only** line cleaner recommended for use with the NXT Hg5 Amalgam Separator

Item#	Description	Suggested Retail	Promo Price
PCS-VLCK	PowerScrub Vacuum Line Cleaner Case Kit	\$207.00	\$207.00
PCS-VLSK	PowerScrub Vacuum Line Cleaner Sidekick	\$53.99	<b>\$0.00</b>

\*Rebate form submission required for free product.

## OFFER #4

Get a \$25 rebate on the purchase of a Compliance Kit\*

Kit contains one collection container and one 1.25-Gallon amalgam bucket



Item#	Description	Suggested Retail
NXT-HG5-CK	NXT Hg5 Compliance Kit	\$425.00
HG5-CK	HG5 Compliance Kit	\$425.00

\*Rebate form submission required for free product.

To redeem offers, complete rebate form and submit to Solmetex. See reverse side for details. Questions? Call 877-207-1551.

## Rebate Form

**Offer #1 Qualifying products for complimentary case of PowerScrub™ Vacuum Line Cleaner with purchase of an NXT Hg5 series Amalgam Separator**

Quantity	Item#	Description
_____	NXT-HG5-BDL1	NXT Hg5 Amalgam Separator
_____	NXT-HG5-MINIBDL1	NXT Hg5 mini Amalgam Separator
_____	NXT-HG5-HV	NXT Hg5 High Volume Amalgam Separator

Yes! Please send me a complimentary case of PowerScrub (Item #PCS-VLCK).  
\*Be sure to include invoice with the form below. Enrollment in Solmetex Maintenance Program required.

**Offer #2 Qualifying products for complimentary case of PowerScrub™ Vacuum Line Cleaner when upgrading to the NXT Hg5 series Amalgam Separator**

Quantity	Item#	Description
_____	NXT-HG5-001RPL	NXT Hg5 Amalgam Separator Upgrade
_____	NXT-HG5-MINIRPL	NXT Hg5 mini Amalgam Separator Upgrade

Yes! Please send me a complimentary case of PowerScrub (Item #PCS-VLCK).  
\*Be sure to include invoice with the form below. Enrollment in Solmetex Maintenance Program required. Upgrade will not be processed without the serial number of the old Hg5 system. Locate serial number on the white PVC pipe of the old Hg5 Amalgam Separator.

**Offer #3 Qualifying products for FREE PowerScrub Sidekick Dispenser**

Quantity	Item#	Description
_____	PCS-VLCK	PowerScrub Vacuum Line Cleaner Case Kit

Yes! Please send me a complimentary PowerScrub Sidekick (Item #PCS-VLSK).  
\*Be sure to include invoice with the form below.

**Offer #4 Qualifying products for \$25 End User Rebate on the purchase of a Solmetex Compliance Kit**

Quantity	Item#	Description
_____	NXT-HG5-CK	NXT Hg5 Compliance Kit
_____	HG5-CK	HG5 Compliance Kit

\*Be sure to include invoice with the form below.

Form must be completed to process your rebate. Rebates require enrollment in the Solmetex Maintenance Program. **Dealer invoice must be submitted along with this rebate form for processing.**

Doctor name: \_\_\_\_\_ Practice name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip code: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_ Dealer: \_\_\_\_\_

**OFFER #1 ONLY:** Serial No. of NXT Hg5 Amalgam Separator (Example "NXTHG5-K-12345") required:  
**NXTHG5-** \_\_\_\_\_ - \_\_\_\_\_ Invoice No.: \_\_\_\_\_

*(Rebate will not be processed without the serial number which is located on the upper chamber of the NXT Hg5 Amalgam Separator)*

Make rebate check payable to: \_\_\_\_\_

Authorizing signature: \_\_\_\_\_

**Email completed form to [sales@solmetex.com](mailto:sales@solmetex.com)**

**PLEASE ALLOW 4-6 WEEKS TO PROCESS**

Limited time offer, specials only valid on purchases dated between July 1, 2021 – September 30, 2021. All 'Q3 Specials' promotions require enrollment in the Solmetex Maintenance Program. **Rebate forms must be received by October 15, 2021.** Offers not valid with any other promotions. Offers may be discontinued or changed without notice and void where not valid. Offer only applies to USA dealers and customers.





# SOLMETEX MAINTENANCE PROGRAM™

Set your  
schedule  
and we do  
the rest

Alleviate the guesswork around EPA compliance and protect your office from unnecessary fines. Enroll in the Solmetex Maintenance Program.

Our new program includes auto-delivery of all the products required by the EPA's "Dental Rule" regulation on amalgam waste. You set the custom schedule that works best for your practice and Solmetex handles the rest. *No Contracts. No Fees.*



COMPLIANCE  
CENTER



AMALGAM  
BUCKET



COLLECTION  
CONTAINER



VACUUM LINE  
CLEANER

(877) 207-1551 | [www.solmetex.com](http://www.solmetex.com)

# Solmetex™ Proudly Introduces the



## SOLMETEX MAINTENANCE PROGRAM™

Enroll today, be compliant tomorrow

In response to the EPA Dental Rule, Solmetex is excited to introduce its proprietary **Solmetex Maintenance Program**. A program designed as a turnkey solution to alleviate guesswork associated with the regulation to ensure your practice is always in compliance. Enrollment enables real-time access to all your practice's documentation and recycling activity online 24/7/365, access to documents that a local inspector may request at any time.

The Solmetex Maintenance Program covers the **four main components of the EPA Dental Rule** that offices are required to have to maintain compliance:

### Hg5 Series of Collection Containers with Recycle Kit



**NXT Hg5  
Collection Container**

with Recycle Kit

Item # **NXT-HG5-002CR**

- The only collection container approved for use with the NXT Hg5 Amalgam Separator
- Clear by design for easy visual inspection
- Simple change-out process when solids reach the full line or at 12 months, whichever occurs first.
- Packaging and shipping approved by the US DOT.
- Certificates of Recycling housed and available 24/7 in the Solmetex Compliance Center



**Hg5 [Legacy]  
Collection Container**

with Recycle Kit

Item # **HG5-002CR**

### PowerScrub™ Vacuum Line Cleaner



#### PowerScrub™ Vacuum Line Cleaner Case Kit

Item #

**PCS-VLCK**

- Per EPA regulation requirements, formulation consists of a neutral pH level between 6 and 8.
- Case kit includes four 64oz bottles.
- Recommended for use with the NXT Hg5 Amalgam Separator
- Maximizes capacity of the NXT Hg5 collection container.
- Compatible with all competing amalgam separators
- Compatible with all dry/wet vacuums
- Works at all water temperatures

### Practice Waste Solutions Amalgam Bucket



#### PWS Amalgam Waste Buckets

1.25 Gallon Bucket	Item # <b>PWS-AB-1</b>
2.5 Gallon Bucket	Item # <b>PWS-AB-2</b>
3.5 Gallon Bucket	Item # <b>PWS-AB-3.5</b>
5 Gallon Bucket	Item # <b>PWS-AB-5</b>

- Designed to address the "dry waste" requirements of the EPA Dental Rule
- Per EPA regulation requirements, clearly labeled for amalgam waste.
- Available in four (4) sizes
- Cost of PWS amalgam bucket includes cost of shipping and recycling.
- Certificates of recycling available 24/7 in the Solmetex Compliance Center

### Solmetex Compliance Center



- Available 24/7/365
- Recycling History
- Certificates of Recycling
- Certificates of Installation
- Certificates of Warranty

Participation in the Solmetex Maintenance Program is accompanied by a **non-contractual**, automatic fulfillment program for all the products and technologies necessary to ensure compliance with the EPA's regulation on the proper handling of amalgam waste within the dental practice.

# Beautiful Restorations to last a lifetime.

Sterile burs available that can positively impact your practice by decreasing sterilization costs so you can focus on producing better patient outcomes.

Effective July 1, 2021 - September 30, 2021

ALL STERILE and NON-STERILE CARBIDES  
**BUY 5 PACKS  
 GET 2 AT NO CHARGE\*\***

CHOOSE FROM: STANDARD OPERATIVE CARBIDES, FG, FGSL, RA, RASL, HP, GREAT WHITE® GOLD, TRIMMING & FINISHING, FISSUROTOMY®

CARBIDE BULK PACK INCLUDING: FG, RA, HP, GREAT WHITE®, FISSUROTOMY®, TRIMMING AND FINISHING

\*\* Can mix and match shapes to equal or lesser value. Free product to equal or lesser value Valid for 5, 10 and 25 Packs and Bulk Packs. Offer excludes Great White® Ultra, EndoGuide® Burs, Endo Safe End and Oral Surgery Burs.

STERILE R AVAILABLE IN ROLLS OF 25



SmartBurs® II - No anesthesia needed\*  
 Happier patients means more referrals  
\*University clinical study available upon request.

Fissurotomy® - Conservative exploration of fissures with a patient-friendly technique, often without the use of anesthesia.\*

## SMARTBURS® II AND FISSUROTOMY® BURS BUY 5 GET 2 AT NO CHARGE\*\*

\*\*CAN MIX AND MATCH TO EQUAL OR LESSER VALUE • NO CHARGE PRODUCT TO EQUAL OR LESSER VALUE

STERILE R AVAILABLE IN ROLLS OF 25

100% Satisfaction Guaranteed†

**SS White Dental®**

Better Patient Outcomes  
 Improved Efficiency  
 Faster Practice Growth

1145 Towbin Avenue Lakewood, New Jersey 08701  
 Visit us on the web at [www.sswhitedental.com](http://www.sswhitedental.com)

©2021 SS White® All Rights Reserved. SS White®, SmartBurs®II, Fissurotomy®, Revelation® and Great White® are registered trademarks of SS White Dental.  
 \*Data on file and available upon request. \*\*Free goods must be of equal or lesser value. All offers are limited to stock on hand and are not valid on prior purchases. Offers may not be combined with any other sale, promotion, discount, code, coupon and/or offer. Sales promotion choice must be indicated on invoice sent to SS White® on applicable order. SS White® reserves the right to discontinue a promotion at any time. †SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, SS White® will issue a full refund of the purchase price through your distributor. To claim free goods: On or before October 31, 2021 mail a copy of your invoice to: SS White Burs, Inc., Attn: Promotions Department, 1145 Towbin Ave., Lakewood, NJ 08701 or email [freegoods@sswhitedental.com](mailto:freegoods@sswhitedental.com). Please note, these are end-user doctor promotions and must be redeemed by the qualifying practice. Please allow 4-6 weeks for free delivery.

Thinking ahead. Focused on life.

# Q3 Promotions

J. MORITA USA | July 1 – September 30, 2021

## TwinPower Turbine:

### Buy 1

MORITA Type Optic Handpiece  
(Excludes Non-Optic & Basic Series)  
4HX-O, 4HEX-O, 4HUEX-O, 4HUMX-O,  
4HEX-O-45: (MORITA)



### Get 1

CP4-LD, 4-Hole  
LED Coupler

### Buy 1

MORITA Type Non-Optic Handpiece  
(Excludes Optic & Basic Series)  
4HX, 4HEX, 4HUEX, 4HEX-45: (MORITA)



### Get 1

CP4, 4-Hole Coupler  
Without Light

### Buy 2

Optic Handpieces  
(Excludes Non-Optic & Basic Series)  
4HX-O, 4HEX-O, 4HUEX-O, 4HUMX-O,  
4HEX-O-45: (MORITA, KV, NK, SR & WH\*)



### Get 1

Additional Handpiece

### Buy 3

Non-Optic or Basic Series Handpieces  
(Excludes all Optic Handpieces)  
4HX, 4HEX, 4HUEX, 4HEX-45, 4HEX-B,  
4HEX-B-45: (MORITA, KV, NK, SR, WH\*)



### Get 2

Additional Handpieces  
(Free handpieces must be of  
equal or lesser value)

### Buy 6

Handpieces  
(Excludes Basic Series)  
4HX-O, 4HX, 4HEX-O, 4HEX, 4HUEX-O,  
4HUEX, 4HUMX-O, 4HEX-O-45, 4HEX-45:  
(MORITA, KV, NK, SR, WH\*)



### Get 1

Lubrina 2 Handpiece  
Maintenance System

### Buy 9

Handpieces  
(Excludes Basic Series)  
4HX-O, 4HX, 4HEX-O, 4HEX, 4HUEX-O, 4HUEX,  
4HUMX-O, 4HEX-O-45, 4HEX-45: (MORITA, KV,  
NK, SR, WH\*)



### Get 1 + 4

(1) Lubrina 2 Handpiece  
Maintenance System &  
(4) maintenance couplers

### Buy 3

TwinPower Couplers



### Get 1

(Free coupler must be of equal  
or lesser value)



Add a Lubrina 2 maintenance unit onto your TwinPower purchase & receive 6 additional months of handpiece warranty.<sup>1</sup>

## Lubrina 2 Handpiece Maintenance System:

### Buy 1

Lubrina 2



### Get 1 + 4

(1) TwinPower Turbine Handpiece &  
(4) Maintenance Couplers **FREE\*\***



## Tri Auto ZX2 Endodontic Handpiece:

### Buy 1

Tri Auto ZX2



### Get \$100 Rebate

### And Buy 1

Extra Tri Auto ZX2  
Contra Head



### Get Additional \$100 Rebate



## Root ZX II & Root ZX mini Apex Locators:

### Buy 1

Root ZX II or Root ZX mini  
Apex Locator



### Get \$50 Rebate

### Buy 1

OTR Low Speed  
Handpiece Module



### Get \$100 Rebate



## Foundation Bone Augmentation Material:

### Buy 3

Boxes of Foundation  
(any size)



### Get 3 Medium Bullets **FREE**

OR

### Get 6 Small Bullets **FREE**



## Microsurgical Blades:

### Buy 4

Packs of Microsurgical  
Blades



### Get 1 Pack **FREE**

(Offer applies to blades 391, 390C, 390,  
370 & 350 ONLY)



## PerfectIM Systems Impression Material:

### Buy 5

Boxes (any viscosity)



### Get 1 Box **FREE**



To order, contact your dealer or call 888-JMORITA (566-7482).

Rebates listed in USD. Offers cannot be combined. Offers expire September 30, 2021. To receive your free gift or promotional offer, fax a copy of your dealer invoice dated July 1, 2021 through September 30, 2021 by October 15, 2021, to J. MORITA USA at (949) 465-1009 Attn: Q3 Promotions, or email to: promotions@jmoritau.com. J. MORITA USA reserves the right to change or cancel these promotions at any time. Students are not eligible for quarterly specials and discounts.

\*Not all handpieces available in all types. TwinPower connection options: TwinPower Morita, KaVo® MULTiflex® LUX, SIRONA® R/F, W&H® Roto Quick, and NSK® Mach/Phatelus®/FlexiQuick Coupling. KaVo MULTiflex LUX is a registered trademark of Kaltenbach & Voigt GmbH. Sirona is a registered trademark of Sirona Dental Systems GmbH. W&H and Roto Quick are registered trademarks of W&H Dental Bürmoos GmbH. NSK and Phatelus are registered trademarks of NAKANISHI INC. Name of coupling varies by country.

\*\*Lubrina 2 Coupling Options: TwinPower Morita, KaVo® MULTiflex® LUX, SIRONA® R/F, W&H® Roto Quick, and E-Type (4-Hole, 5-Hole, and 6-Pin compatibility).

†Invoice for both Lubrina 2 and Handpiece(s) must be submitted together. Call for details.

www.morita.com/usa

# TIME TO SWITCH

Aerosol reduction with the flip of a switch.



## NEW Ti-Max Z95L

### Easy On/Off Switch

Switch between Water Jet (water only) and Water Spray (water+air) in an instant.

Instantly switch between Water Spray or Water Jet **with the flip of a switch** with the world's first two-way spray technology. When extra safety is needed, just the Water Jet provides maximum protection while continuing to offer optimum cooling.

### Handpiece Offers

**Buy 2** Electric handpieces  
Z95L, Z85L, Z25L, Z45L,  
nano95LS, nano25LS



**Get 1** Electric handpiece  
(of equal or lesser value)



**Buy 2** Air-Driven handpieces  
Z800L, Z800KL, Z900L, Z900KL



**Get 1** Air-Driven handpiece  
(of equal or lesser value)



**Buy 3** Z95L or nano95LS handpieces



**Get 1** NLZ Electric Micromotor System



**Buy 3** Z95L or nano95LS handpieces AND 1 Endo handpiece



**Get 1** NLZ E Electric Micromotor System with Endodontic Function



**Buy 2** FX Series handpieces and motors  
FX25, FX57, FX65, FX204



**Get 1** FX Series handpiece or motor  
(of equal or lesser value)



# Buy 1 Operator Bundle, Get 1 of the same FREE

Choose between an Electric or an Air Bundle!

Air Operator Bundle		
Product	Quantity	Description
Z800L/Z800KL Z900L/Z900KL	3	Ti-Max Z Titanium Miniature Head (Z800L/Z800KL) or Standard Head (Z900L/Z900KL) High Speed Handpieces
PTL-CL-LED or KCL-LED	1	PTL Coupler - 6 pin LED <b>or</b> Kavo-Type Multiflex Coupler - 6 pin LED
FX-204	1	4 Hole Air Motor
FX25	2	FX Series, Non-Optic 1:1 Direct Drive Attachment
FX65	1	FX Series, Non-Optic 1:1 Straight Attachment
iProphy Air	1	iProphy Air Prophylaxis Air Motor

Retail Value:  
\$5,827

Electric Operator Bundle		
Product	Quantity	Description
Z95L	3	Ti-Max Z 1:5 Increasing Attachment
Z25L	2	Ti-Max Z 1:1 Push Button Latch Attachment
M65	1	S-Max M 1:1 Straight Attachment

Retail Value:  
\$8,556

## Maintenance Offer

Buy 1 iCare  
Maintenance System



Get 1 Z Series  
handpiece

Z95L, Z85L, Z45L, Z25L, nano95LS,  
nano25LS, Z800L, Z800KL, Z900L,  
Z900KL



Ti-Max Z Series

## Oral Hygiene Offers

Buy 2 iProphy Air  
Prophylaxis Air Motors



Get 2 iProphy Air  
Prophylaxis Air Motors



iProphy Air

Buy 1 iProphy Mobile  
Complete Set



Get 500 Handpiece Sleeves  
144 Preventech – Pivot Angles



iProphy Mobile

## Surgical Offers

Buy 1 Surgic Pro Set  
Set includes one 20:1  
surgical handpiece



Get 1 Surgical handpiece  
X-SG20L, SG20

OR

Get 1 iSD900 Cordless  
Prosthodontic Screwdriver



iSD900

Buy 2 Implant  
handpieces  
X-SG20L, SG20, X-DSG20L



Get 1 Implant  
handpiece  
(of equal or lesser value)



Ti-Max X-DSG20L

\*All NSK no charge goods requests must be submitted to NSK. Please send a copy of your distributor invoice to NSK promotions, 1800 Global Parkway, Hoffman Estates, IL 60192; or fax to: 1-800-838-9328 or email to: orders@nskdenal.com. iProphy Mobile offer includes one box of Pivot Plus Soft Spiral Cup DPAs (1133343). Send a copy of the dealer invoice to: Preventech, 4430-C Matthews- Indian Trail Rd., Indian Trail, NC 28079 or fax to 704-849-2417, email to: customerservice@preventech.com. To be eligible, you must submit by October 31, 2021.



Scan to Watch  
our Video!



NSK America Corp. [www.nskdenal.com](http://www.nskdenal.com)

1800 Global Parkway, Hoffman Estates, IL. 60192 • Tel: (888) 675-1675 • Fax: (800) 838-9328



# SPECIALS

JULY | AUGUST | SEPTEMBER 2021

COMPOSITES | GLASS IONOMERS  
TOOTH WHITENING | EQUIPMENT  
CEMENTS | ADHESIVES  
ETCHANTS | SEALANTS  
AMALGAMS | ACCESSORIES

**POLA**  
ADVANCED TOOTH  
WHITENING SYSTEM

## ALL TAKE HOME WHITENING

# 3+1\* OR 10+5\*

PROMO CODE: Q3POLA10R5



## POLA DISPENSER KITS

# 2+1 BUY 2 GET 1 FREE\*

POLA LIGHT  
REFILL  
SYRINGES



### POLA DAY

HYDROGEN PEROXIDE BASED  
ADVANCED TRAY TOOTH  
WHITENING SYSTEM

- Available in 7.5% and 9.5% HP, 35% CP

### POLA NIGHT

CARBAMIDE PEROXIDE BASED  
ADVANCED TRAY TOOTH  
WHITENING SYSTEM

- Available in 16% and 22% CP

PROMO CODE: Q3PDB2G1





# POLA IN CHAIR WHITENING BUNDLE

BUY RADII XPERT CURING LIGHT, RADII XPERT BLEACH ARCH, 1 POLA RAPID (3 PATIENT KIT)

CREATED BY SDI

# POLA

ADVANCED TOOTH WHITENING SYSTEM

OFFICIAL SMILE OF  
MISS \* USA  
MISSTEEN \* USA



RADII XPERT

+



RADII XPERT BLEACH ARCH

+



POLA RAPID (3 PATIENT KIT)



## GET \$1000 WORTH OF PRODUCTS FREE\*

x1



1 X POLA STAND

+

x1



1 X POLA RAPID (3 PATIENT KIT)

+

x3



3 X POLA DAY OR POLA NIGHT 10 SYRINGE KITS

+

x3




3 X POLA LUMINATE

PROMO CODE: POLABUNDLE2021

DROPSHIP FROM SDI AND BILLED THROUGH YOUR PREFERRED DENTAL DISTRIBUTOR.

## BONUS

Free Pola marketing package and Virtual CE or in person education program.



# ALL COMPOSITES 3+1\* OR 10+5\*

## AURA EASY

**STOCK LESS, DELIVER MORE**

- 4 shades cover 8 VITA classical
- Increased wear resistance and strength



## LUNA

**ANTERIOR/POSTERIOR NANOHYBRID COMPOSITE**

- Nanohybrid technology for superior strength and aesthetics
- Chameleon effect



## AURA EASYFLOW

**COMBINATION OF STRENGTH, VISIBILITY AND AESTHETICS IN A FLOWABLE**

- Established nanohybrid filler system
- Optimal flow for cavity adaptation
- High visibility on radiographs



## AURA BULK FILL

**THE SIMPLICITY OF BULK FILLING**

- One bulk fill placement up to 5mm depth of cure
- 20 sec. LED light cure
- Natural-looking finish with just one universal shade
- Easy placement in a matter of seconds



## WAVE

**FLUORIDE RELEASING, FLOWABLE RADIOPAQUE COMPOSITES**

- High compressive strength
- Nano technology
- Three viscosities
- Natural fluorescence
- No after cure shade change

## ALLOY PROMO CODE: Q3ALLOY

**BUY 2 ALLOY GET 1 FREE\***



**OR BUY 3 JARS OF ALLOY GET 1 FREE TRITURATOR + BONUS RIVA GLASS SC STARTER KIT\***



## ADHESIVES PROMO CODE: Q3ZIPB2G1

**BUY 2 ZIPBOND GET 1 FREE + BONUS SUPER ETCH\***

## ZIPBOND

UNIVERSAL ADHESIVE

- Consistent bond strengths to dentin and enamel
- Wet and dry dentin bonding
- Proprietary combination of adhesive monomers including MDP
- Predictable bond strengths to indirect substrates
- Effective hybridization



SPEND **\$1500** ON  
SDI'S RESTORATIVE PRODUCTS  
**GET YOUR CHOICE OF  
RADII PLUS OR RADII XPERT\***

[Excludes Alloy]



PROMO CODE: Q3RSBUNDLE

BUY ANY 3 BOXES FROM THE RIVA RANGE **GET 1 FREE\***

[Excludes RIVA STAR]



**AND A BONUS RIVA STAR  
BOTTLE KIT FREE\***



PROMO CODE: Q3RIVAB3G1

COMPOSITES | GLASS IONOMERS | TOOTH WHITENING | EQUIPMENT | CEMENTS | ADHESIVES | ETCHANTS | SEALANTS | AMALGAMS | ACCESSORIES



Redeem your free goods online using this link:  
<http://www.sdi.com.au/freegoods-us>

\*Please allow 2-4 weeks for delivery. Please note offers cannot be combined. Offers cannot be used in conjunction with contract pricing. Valid on invoices dated July 1st - September 30th, 2021. All claims for free goods must be made within 60 days from the end of promotion date. Free products must be of equal or lesser value. Offers are valid in the US only.

**SDI (North America) Inc.**  
1279 Hamilton Parkway  
Itasca, IL 60143

Toll Free 1 800 228 5166  
Fax: 630-361-9222  
[usa.canada@sdi.com.au](mailto:usa.canada@sdi.com.au)  
[www.sdi.com.au](http://www.sdi.com.au)

facebook.com/sdiusa | instagram.com/sdiusa

**KENDRA QUINN**  
Regional Sales Manager - USA  
Cell: (614) 558-6850  
[Kendra.Quinn@sdi.com.au](mailto:Kendra.Quinn@sdi.com.au)

**PHILLIP ZHOU**  
PA, DE, Southern NJ, MD  
Cell: (717) 395-0399  
[Phillip.Zhou@sdi.com.au](mailto:Phillip.Zhou@sdi.com.au)

**STEVE WISNIEWSKI**  
Michigan, Northern Indiana  
Cell: (313) 520-4114  
[Steve.Wisniewski@sdi.com.au](mailto:Steve.Wisniewski@sdi.com.au)

**MARIANNA OSIPOVA**  
NYC Metro and Northern New Jersey  
Cell: (718) 308-5588  
[Marianna.Osipova@sdi.com.au](mailto:Marianna.Osipova@sdi.com.au)

**MARK GUZIOR**  
MA, NH, NJ, CT, VT, RI  
Cell: (603) 533-3121  
[Mark.Guzior@sdi.com.au](mailto:Mark.Guzior@sdi.com.au)

**CARRIE COOMER**  
OH, KY, Southern IN  
Cell: (937) 214-7410  
[Carrie.Coomer@sdi.com.au](mailto:Carrie.Coomer@sdi.com.au)

**KATIE BLEVINS**  
Florida  
Cell: (561) 358-9291  
[Katie.Blevins@sdi.com.au](mailto:Katie.Blevins@sdi.com.au)

**COURTNEY BUCZEK**  
Central/Northern IL, WI  
Cell: (630) 542-7969  
[Courtney.Buczec@sdi.com.au](mailto:Courtney.Buczec@sdi.com.au)

**JIM POTEAT**  
GA, SC, NC, TN  
Cell: (704) 995-7303  
[Jim@Poteatdentalsolutions.com](mailto:Jim@Poteatdentalsolutions.com)

**MIKE MATUCH**  
Upstate NY, Erie/Pittsburg, PA, WV  
Cell: (716) 913-8534  
[MichaelMatuch@yahoo.com](mailto:MichaelMatuch@yahoo.com)

**LOUIE KALOUDIS**  
MS, AL, DC, VA  
Cell: (630) 363-9224  
[Louie.Kaloudis@sdi.com.au](mailto:Louie.Kaloudis@sdi.com.au)

**CYNTHIA CRALLIE**  
Los Angeles, Ventura, Kern Counties  
Cell: (630) 306-4509  
[Cynthia.Crallie@sdi.com.au](mailto:Cynthia.Crallie@sdi.com.au)

**LAVETTE THOMAS**  
No. Cal, Reno, NV  
Cell: (408) 893-8928  
[Lavette.Thomas@sdi.com.au](mailto:Lavette.Thomas@sdi.com.au)

**CYNTHIA CRALLIE**  
LA, AR, MS  
Cell: (630) 306-4509  
[Cynthia.Crallie@sdi.com.au](mailto:Cynthia.Crallie@sdi.com.au)

**CYNTHIA CRALLIE**  
DC, San Diego, Las Vegas, HI  
Cell: (630) 306-4509  
[Cynthia.Crallie@sdi.com.au](mailto:Cynthia.Crallie@sdi.com.au)

**LORI KELLAR**  
WA, OR, AK, ID, MT  
Cell: (206) 696-1485  
[Lori.Kellar@sdi.com.au](mailto:Lori.Kellar@sdi.com.au)

**KARLYN RASMUSSEN-DODD**  
TX, OK  
Cell: (832) 928-6938  
[Karlyn.Rasmussen-dodd@sdi.com.au](mailto:Karlyn.Rasmussen-dodd@sdi.com.au)

**KAREN ZAFFINO**  
CO, UT, WY, NM, El Paso, AZ  
Cell: (720) 635-3400  
[Zaffinoenterprise@gmail.com](mailto:Zaffinoenterprise@gmail.com)

**TIM VANWYCHEN**  
MN, ND, IA, SD, NE  
Cell: (651) 231-2801  
[VanWychen@hotmail.com](mailto:VanWychen@hotmail.com)

**DOUG SCHROEDER**  
MO, KS  
Cell: (314) 422-0061  
[ds1td@charter.net](mailto:ds1td@charter.net)

# 2021 Spotlight Products

JULY 1 - SEPTEMBER 30

**Buy 1 + Get 1** (5 pack only)

## Biodentine®



Bioactive Dentin Substitute • Learn more at [biodentine.com](http://biodentine.com)

- Resin-free = no shrinkage and micro-leakage resistant
- Non-toxic, preserves pulp vitality, use in direct/indirect pulp caps
- High pH of 12, high alkalinity, has antimicrobial properties and helps arrest residual caries

#01C0605 5 capsules & 5 single-dose pipettes • #01C0600 15 capsules & 15 single-dose pipettes



**Buy 1 BioRoot + Get a Biodentine 5 pack**

## BioRoot™ RCS



Bioceramic Root Canal Sealer

- Excellent results utilizing cold obturation techniques
- Outstanding adhesion to dentin and gutta-percha points
- High biocompatibility, very pure mineral tri-calcium silicate formulation

#01E0300 35 applications per pack



**Buy 1 R.T.R.+ + Get 1 Membrane** (15x20mm only)

## R.T.R.+ (Resorbable Tissue Replacement)

Biphasic Synthetic Bone Substitute

Two ratios of tri-calcium phosphate to hydroxyapatite: 40/60 and 80/20

- Infrabony defects • Peri-implant defects • Sinus lift • Ridge augmentation • Cystic cavities

#01S0540 40%  $\beta$ -TCP | 60% Hydroxyapatite • #01S0530 80%  $\beta$ -TCP | 20% Hydroxyapatite



## R.T.R Membrane (Resorbable Tissue Replacement)

Absorbable Collagen Membrane

- Protect periodontal ligament development from fast growing soft tissues
- Maintain bone grafting material to the region to be reconstructed
- Support wound healing & guide bone regeneration

#01S0510 (15 x 20mm), 1 per box • #01S0520 (20 x 30mm), 1 per box




## Aspirating Syringes

For injection of local anesthetic solution

- Ergonomic - Available in 3 distinct sizes to provide better control and added comfort when aspirating
- Petite Aspirating Syringe is ideal for smaller hands - equipped with a smaller thumb ring, finger grip and a shorter plunger length to reduce hand stress during aspiration

#01N2100 - 1 Petite Syringe (blue)    #01N2105 - 1 Fusion Syringe (titanium)  
#01N2120 - 1 Petite Syringe (pink)    #01N2110 - 1 Standard Syringe (bronze)



\* To receive rebate, submit dealer invoice to [marketing@septodont.com](mailto:marketing@septodont.com). Please allow 4-6 weeks for processing.


# 2021 Spotlight Products

JULY 1 - SEPTEMBER 30

## Septocaine®

(Articaine Hydrochloride 4% and Epinephrine)

- Predictable duration of anesthesia: up to 60 minutes for infiltration injections and 120 minutes for nerve blocks
- Fast onset (1 to 9 minutes) to avoid delay before effective analgesia is achieved



#01A1400 1:100,000 Injection (Articaine Hydrochloride 4% and Epinephrine)

#01A1200 1:200,000 Injection (Articaine Hydrochloride 4% and Epinephrine)

**Buy 1 + Get 1**

## DycloPro™ Topical Anesthetic

(Dyclonine Hydrochloride Topical Solution 0.5%)

A topical anesthetic useful to suppress the gag reflex, can also be applied to the gums prior to scaling (prophylaxis)

#99510 30g: Multiple doses in each bottle



## Septoject® Evolution

Sterile single use dental needle

- Scalpel-designed bevel requires less force for the needle insertion
- Less tissue is displaced resulting in reduced discomfort

#01N1600 27g Short • #01N1630 30g Short • #01N1650 30g Extra Short



**Buy 1 + Get 1**

## OraVerse™ (Phentolamine Mesylate)

The First and Only Local Anesthetic Reversal Agent

Allows you to bring patients “back to normal” sensation twice as fast

#99210 1 blister of 10 x 1.7ml glass cartridges



Available in the  
Apple Store & Google Play

**Download**  
**The Septodont App**

gives instant access to  
product literature,  
videos, CE webinars,  
testimonials, and more!

Offers valid July 1, 2021 through September 30, 2021.

Not to be combined with any other offers and are subject to change. Promotional offers are stocked and shipped with your order by your dealer. Questions about Septodont products? Call our friendly customer service team at 800-872-8305 or visit septodontusa.com to learn more.



**DISPOSABLE PROPHY ANGLES**

**YOUNG®**

Buy 4 - 200ct boxes of Young Splatter Guard Disposable Prophylaxis Angles, **GET 1 FREE!**  
Use promo code Q3YSG200

Buy 2 - 1200ct boxes of Young Splatter Guard Disposable Prophylaxis Angles, **GET 2 YOUNG HYGIENE HANDPIECES FREE!**  
Use promo code Q3YSG1200

Buy 3 - 200ct boxes of Young Classic, Contra or Slim Disposable Prophylaxis Angles, **GET 1 FREE!\***  
Use promo code Q3YDPA200

Buy 2 - 1200ct boxes of Young Classic, Contra or Slim Disposable Prophylaxis Angles, **GET 2 YOUNG HYGIENE HANDPIECES FREE!\***  
Use promo code Q3Y1200

**VERA®**

Buy 4 - 125ct boxes of Vera Splatter Guard Disposable Prophylaxis Angles, **GET 1 FREE!**  
Use promo code Q3VSG125

Buy 3 - 750ct boxes of Vera Splatter Guard Disposable Prophylaxis Angles, **GET 2 YOUNG HYGIENE HANDPIECES FREE!**  
Use promo code Q3VSG750

Buy 3 - 125ct boxes of Vera Classic or Contra Disposable Prophylaxis Angles, **GET 1 - 50CT BOX OF D-LISH VARNISH FREE!\***  
Use promo code Q3VDPA125

Buy 3 - 750ct boxes of Vera Classic or Contra Disposable Prophylaxis Angles, **GET 2 YOUNG HYGIENE HANDPIECES FREE!\***  
Use promo code Q3VDPA750



**DENTICATOR®**

Buy 1 - 1000ct boxes of eco™ Slim Disposable Prophylaxis Angles, **GET 1 FREE!\***  
Use promo code Q3DECO

Buy 3 - 125ct boxes of eco™ Slim Disposable Prophylaxis Angles, **GET 1 FREE!\***  
Use promo code Q3ECO125

**MIX & MATCH - ORIGINAL GREEN™, DIAMOND OR DENTICOLORS™**

Buy 3 - 144ct boxes, **GET 1 FREE!\***  
Use promo code Q3D144

Buy 3 - 500ct boxes, **GET A PROPHY PAL® HYGIENE HANDPIECE FREE!\***  
Use promo code Q3D500

Buy 2 - 500ct boxes, **GET 1 OF THE SAME FREE!\***  
Use promo code Q3D500

**ZOOPY®**

Buy 2 - 500ct boxes of Zooby Disposable Prophylaxis Angles, **GET 1 YOUNG HYGIENE HANDPIECE FREE!\***  
Use promo code Q3ZDPA500

BUY 2 - 100ct boxes of Zooby Disposable Prophylaxis Angles, **GET 1 - FREE!\***  
Use promo code Q3ZDPA100

**HYGIENE HANDPIECES**

**YOUNG®**

Buy 3 Young™ Hygiene Handpieces, **GET 1 FREE!**  
Available in silver or purple.  
Use promo code Q3YHP

Buy 1 Young™ Infinity Cordless Handpiece, **GET 1 - 125CT COUNT BOX OF VERA SPLATTER GUARD PROPHY ANGLES FREE!\***  
Use promo code Q3INFINITY

Buy 2 Young™ Infinity Cordless Handpieces, **GET 1 FREE!\***  
With the Purchase of the a 750ct Box of Vera Prophylaxis Angles or a 1200ct Box of Young Prophylaxis Angles  
Use promo code INFINITY2+1

**DENTICATOR®**

Buy 3 ProphylaxisPal® Hygiene Handpieces, **GET 1 FREE!\***  
Available in classic silver, prophylaxis pink or wild thing green.  
Use promo code Q3PPAL



**VARNISH**

**D-LISH®**

Buy 2 - 50ct boxes of D-Lish® Fluoride Varnish, **GET 1 FREE!**  
Use promo code Q3DV50

Buy 2 - 200ct boxes of D-Lish® Fluoride Varnish, **GET 1 FREE!**  
Use promo code Q3DV200

**ZOOPY®**

Buy 2 - 50ct boxes of Zooby® Fluoride Varnish, **GET 1 FREE!**  
Use promo code Q3ZV50



Zooby Fluoride Varnish Available in Growlin' Grrrape®, Turtle Melon® & Happy Hippo Cake®.



**HAND INSTRUMENTS**

**AMERICAN EAGLE INSTRUMENTS**

Buy 7 Stainless Steel or XP Sharpen-Free American Eagle Instruments, **GET 3 FREE!**  
Offer includes Quik-Tips  
Use promo code Q3AEI7+3

Buy 25 Stainless Steel or XP Sharpen-Free American Eagle Instruments, **GET 15 FREE!**  
Offer includes Quik-Tips  
Use promo code Q3AEI25+15




American Eagle Instruments

**ASPIRATION INSTRUMENTS**

**CRYSTAL TIP®**

Buy 2 - 250ct boxes of Crystal Tips®, **GET 1 ASSORTED FREE!**  
Available in assorted or pacific blue.  
Use promo code Q3CT250



Crystal Tip® Aspiration Instruments

**INFECTION CONTROL**

**BIREX QUAT**

Buy 7 - 24oz bottles of Birex Quat Ready-to-Use Spray, **GET 5 FREE!\***  
Use promo code Q3BRTU

Buy 2 - 1 gal. bottles of Birex Quat Ready-to-Use Spray Refill **GET 2 FREE!\***  
Use promo code Q3BRTU

**OPTI-CIDE3**

Buy 12 - 24oz bottles of Opti-Cide3 Spray, **GET 2 FREE!\***  
Use promo code Q3OPTIC1

Buy 2 - 1 gal. bottles of Opti-Cide3, **GET 2 FREE!\***  
Use promo code Q3OPTIC2

Buy 1 - 2.5 gal. bottles of Opti-Cide3 with Spigot, **GET 1 FREE!\***  
Use promo code Q3OPTIC3

Buy 6 canisters of Opti-Cide3 Wipes, **GET 3 FREE!\***  
Use promo code Q3OPTIC5





# SPECIALS!

July-Sept 2021



Mix n Match okay  
Cannot combine invoices  
Allow 4-6 weeks for delivery


email invoices to [promos@surgicalesthetics.com](mailto:promos@surgicalesthetics.com)

**New Customers**



**Buy 2 Bone Grafts  
Get 1 Free!**  
Limit 3 free

**Returning Customers**



**Buy 5 Bone Grafts, Get 1 Free!**  
or  
**Buy 12 Bone Grafts, Get 3 Free!**  
or  
**Buy 30 Bone Grafts, Get 10 Free!**



**Buy 4 CollaShield membranes,  
Get 1 Free!**



**ACTIVA™ Presto™  
Universal Stackable Light Cure**

**BUY 2 ACTIVA Presto Kits  
GET 1 FREE (Limit 2 kits FREE, A1 or A2 Shade)**

- Mimics nature
- Mineral enriched
- Moisture friendly
- No BPA, BIS-GMA, Bis-Phenol A derivatives



**ACTIVA™ BioACTIVE-Restorative  
and ACTIVA™ Kids Starter Kit**

**BUY 1 ACTIVA BioACTIVE-Restorative and/or Kids Starter Kit  
GET 1 Single Refill FREE (A1 or A2 shade)**



**ACTIVA™ BioACTIVE  
Value Refills Mix n Match**

**BUY Any 3 ACTIVA BioACTIVE  
Value Refills  
GET 1 Single Refill FREE  
(Limit 2 Kits—A1 or A2 FREE)**



**Embrace™ Wetbond™  
Pit & Fissure Sealant**

**BUY 3 Embrace Pit & Fissure Kits  
GET 1 Embrace Pit & Fissure Sealant Kit FREE  
(Limit 3 Kits Free)**



**Tuff-Temp Plus™**

*Save more on Provisional Materials*

**BUY Any 3 TTP 50mL Refill Kits  
(A1, A2 or A3 Shade)  
GET 1 Refill Kit FREE  
(A2 Shade)**



**Embrace™ Varnish**

**BUY 3 Embrace Varnish FV50 GET 1 FREE  
BUY 3 Embrace Varnish FV200 GET 1 FREE  
BUY 3 Embrace Varnish FVT GET 1 FREE  
BUY 3 Embrace Varnish FVX100 GET 1 FREE**



**1. BUY 2 ACTIVA™ Presto™ Kit\***  
2 x 1.2mL/2gm syringes + 20 applicator tips  
GET 1 ACTIVA™ Presto™ Kit FREE (A1 or A2 Shade)  
VPP1\* - 2 x 1.2mL/2gm syringes + 20 applicator tips  
\*Available shade—A1, A2, A3, A35, A4, A6, B1, BW

**5. BUY Any 3 Tuff-Temp Plus 50mL Refill Kits  
(A1, A2, A3)**  
GET 1 Refill Kit FREE  
TTP50A1, TTP50A2, TTP50A3

**2. BUY 1 ACTIVA™ BioACTIVE or ACTIVA Kids Starter Kit**  
**VR\*** - 1 x 5mL/8gm syringe, ACTIVA-Spenser, 20 mix tips  
GET 1 ACTIVA™ BioACTIVE Restorative Single Refill FREE  
**VR1\*** - 1 x 5mL/8gm syringe + 20 mix tips  
\*Shades Available: A1, A2, A3, A35, Pedo

**6. BUY 3 Embrace Varnish FV50  
GET 1 FV50 FREE**  
**FV50**—Box of 50 x 0.4mL Packets

**BUY 3 Embrace Varnish FV200  
GET 1 FV200 FREE**  
**FV200**—Box of 200 x0.4mL Packets

**BUY 3 Embrace Varnish FVT  
GET 1 FVT FREE**  
**FVT**— 12mL Tube

**BUY 3 Embrace Varnish FVX100  
GET 1 Embrace Varnish FVX100 FREE**  
**FVX100**—Economy Pack:100 x 0.4mL Packets (no brush)

**3. BUY 3 ACTIVA™ BioACTIVE Value Refills Mix n Match**  
GET1 ACTIVA Single Refill Free (Restorative)

**VR1** - 1 x 5mL syringe + 20 mix tips (Restorative)  
**VC1** - 1 x 5mL syringe + 20 mix tips (Cement)  
**VB1**— 1 x 5mL syringe + 20 mix tips (Base-Liner)  
**VK1P**— 1 x 5mL syringe + 20 mix tips (Kids)

**4. BUY 3 Embrace™ Pit & Fissure Sealant Kits**  
GET 1 Embrace Pit & Fissure Sealant FREE

**EMS** (Natural Shade) – 4x1.2mL syringes + tips  
**EMSW** (Off-White Shade) - 4x1.2mL syringes + tips



To redeem free goods, fax, mail or email your dealer invoice to Pulpdent by 10/31/2021 with proof of purchase from 7/1/2021-9/30/2021. Special offer purchases must appear on one invoice to be eligible. Free goods are of equal or lesser value than purchased goods. Offers valid in U.S. ONLY. Please allow 4-6 weeks for receipt of free goods.

## CROSSTEX™ FACE MASKS WITH SECURE FIT™ TECHNOLOGY

**Better fit. Better filtration. Better protection.**

- **Patented design features dual nose and chin strips** for a snug, custom fit that minimizes dangerous gaps and provides superior protection over standard ear loop masks.<sup>2</sup>
- **Made-in-USA** surgical masks available in 3 levels of protection—so you'll always have the right mask for the task.



To learn more, visit  
[Crosstex.com/SecureFit](https://Crosstex.com/SecureFit)

**BUY 3** Boxes of Crosstex™ Isofluid™, Procedural or Ultra™ Surgical Face Masks with Secure Fit™ Technology **GET 1** of the same at no additional charge!



## HU-FRIEDY HAND ESSENTIALS™ SOAP, SANITIZER AND LOTION PRODUCTS



- Our moisturizers, soaps and sanitizers will be there for you with **every glove change, disinfecting hand wash** and day-ending skin routine.
- **Non-irritating, cost effective** and designed for repeated use, our Hand Essentials are the perfect choice to care for your two most precious instruments.



To learn more, visit  
[Hu-Friedy.com/HandEssentials](https://Hu-Friedy.com/HandEssentials)

**BUY 4** of the same Hu-Friedy Hand Essentials™ Soap, Sanitizer and Lotion Products **GET 1** at no additional charge!

## DENTAPURE™ DP365B INDEPENDENT WATER BOTTLE CARTRIDGE

- The DentaPure™ Cartridge is **EPA registered** to provide **compliant treatment water for one year or 240L of water**, ensuring that your practice meets current **EPA microbiological water quality** standards of <500 CFU/mL.
- **EPA registered** Liquid Ultra™ Solution is a powerful and versatile dental unit waterline **cleaner/shock treatment** for external bottle systems that **kills biofilm bacteria**; removes existing biofilm and prevents/suppresses formation of biofilm in dental unit waterlines.



To learn more, visit  
[Crosstex.com/DentaPure](https://Crosstex.com/DentaPure)

**BUY 3** DentaPure™ Bottle Cartridges **GET 1** box (10 sets) of Liquid Ultra™ Solution at no additional charge



Hu-Friedy and Crosstex are now proud members of **HuFriedyGroup**  
The Best In Practice

**Buy 2 – 3M™ Composites and  
Buy 2 – 3M™ Adhesives,  
Get 1 – of Each Free!\***



# A full line of restorative solutions to fit your unique needs.

3M™ Filtek™ Supreme Ultra Universal Restorative	3M™ Filtek™ Supreme Flowable Restorative	3M™ Filtek™ Universal Restorative	3M™ Filtek™ One Bulk Fill Restorative	3M™ Filtek™ Bulk Fill Flowable Restorative
Exceptional esthetics and polish retention	Virtually no bubbles or material “run-on”	Excellent, creamy handling	Increased opacity for improved esthetics	Virtually no bubbles or material “run-on”
Exceptional handling	Easy-to-bend cannula for improved access	Low shrinkage stress/ stress relief	One-step placement up to 5mm	Easy-to-bend cannula for improved access
Outstanding strength for anterior and posterior	Ergonomic design is easy to hold and inject	Excellent wear resistance	Excellent adaptation	Ergonomic design is easy to hold and inject
Wear resistance similar to enamel <sup>1</sup>	Excellent adaptation, polish retention and wear resistance	8 designer shades plus Extra White (XW) cover the 19 VITA classical and bleach shades <sup>2</sup>	Excellent wear resistance	Injectable, one-step placement up to 4mm
36 shades, 4 opacities	Matches shades for Filtek Supreme Ultra Universal Restorative and Filtek Universal Restorative	Pink Opaquer (PO) with excellent masking ability	Excellent polish retention	Excellent adaptation and low shrinkage

Ask your Goetze Dental representative for complete ordering information.

1. In vitro data on file. 2. 3M internal data.

\*Select from 3M™ Filtek™ Restoratives, 3M™ Paradigm™ Nano Hybrid Universal Restoratives and 3M™ Scotchbond™ Universal Adhesives. To obtain bonus goods, you can redeem by email. Proof of purchase on one Goetze Dental invoice dated July 1 – September 30, 2021 is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). Redeem by October 31, 2021. Redeem by email: [3MOraCarePromos@mmm.com](mailto:3MOraCarePromos@mmm.com). Specify bonus product selections and include your business phone. Allow 3–5 weeks for delivery. Offers valid in the United States. Offer void if purchased product is returned. 3M, Filtek, Paradigm and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. © 3M 2021. All rights reserved.

**See page 40 for 3M product redemption coupon**

# SDI SPECIALS

JULY | AUGUST | SEPTEMBER 2021

COMPOSITES | GLASS IONOMERS  
TOOTH WHITENING | SDF RIVA STAR  
EQUIPMENT | CEMENTS | ADHESIVES  
ETCHANTS | SEALANTS  
AMALGAMS | ACCESSORIES

## ALL TAKE HOME WHITENING 3+1<sup>x</sup> OR 10+5<sup>x</sup>

PROMO CODE: Q3POLA10RS



**POLA MINI KITS**  
Regular Price = \$22.28  
Net 3+1 = \$16.71  
Net 10+5 = \$14.85



**POLA 10 SYRINGE KITS**  
Regular Price = \$52.07  
Net 3+1 = \$39.05  
Net 10+5 = \$34.71



**POLA LIGHTS**  
Regular Price = \$98.28  
Net 3+1 = \$73.71  
Net 10+5 = \$65.52

## POLA DISPENSER KITS

### 2+1<sup>x</sup>

POLA LIGHT  
REFILL  
SYRINGES



### POLA DAY

• Available in 7.5% and 9.5% HP, 35% CP



### POLA NIGHT

• Available in 16% and 22% CP

**POLA DISPENSER KITS**  
Regular Price = \$336.53  
Net 2+1 = 224.35

PROMO CODE: Q3PDB2G1

## COMPOSITES

PROMO CODE: Q3COMP10RS

# ALL COMPOSITES 3+1<sup>x</sup> OR 10+5<sup>x</sup>

### AURA EASY

**STOCK LESS, DELIVER MORE**  
Regular Price = \$103.45  
Net 3+1 = \$77.59  
Net 10+5 = \$68.97



### AURA BULK FILL

**THE SIMPLICITY OF BULK FILLING**  
Regular Price = \$103.45  
Net 3+1 = \$77.59  
Net 10+5 = \$68.97



### LUNA

**ANTERIOR/POSTERIOR NANOHYBRID  
COMPOSITE**  
Regular Price = \$103.40  
Net 3+1 = \$77.55  
Net 10+5 = \$68.93



### WAVE

**FLUORIDE RELEASING, FLOWABLE RADIOPAQUE COMPOSITES**  
Regular Price = \$26.27  
Net 3+1 = \$19.70  
Net 10+5 = \$17.51



### AURA EASYFLOW

**COMBINATION OF STRENGTH, VISIBILITY AND  
AESTHETICS IN A FLOWABLE**  
Regular Price = \$51.71  
Net 3+1 = \$38.78  
Net 10+5 = \$34.47



COMPOSITES | GLASS IONOMERS | TOOTH WHITENING | SDF RIVA STAR | EQUIPMENT | CEMENTS  
ADHESIVES | ETCHANTS | SEALANTS | AMALGAMS | ACCESSORIES

Redeem your free goods online using this link:  
<http://www.sdi.com.au/freegoods-us>

facebook.com/sdiusa instagram.com/sdiusa

\*Please allow 2-4 weeks for delivery. Please note offers cannot be combined.  
Offers cannot be used in conjunction with contract pricing. Valid on invoices dated  
July 1st - September 30th, 2021. All claims for free goods must be made within 60  
days from the end of promotion date. Free products must be of equal or lesser value.  
Offers are valid in the US only.



YOUR  
SMILE.  
OUR  
VISION.

SDI (North America) Inc.  
1279 Hamilton Parkway  
Itasca, IL 60143  
Toll Free 1 800 228 5166  
Fax: 630-361-9222  
usa.canada@sdi.com.au  
www.sdi.com.au

# 4D Isodam®


## Blue Polyisoprene Dental Dam



WHEN FAILURE IS NOT AN OPTION,  
CHOOSE ISODAM®.

### DENTAL ADVISOR

Non-Latex Dental Dam Tear Testing Comparison



For clinical procedures where dental dam stretch over 30 mm is required, the 4D Rubber **ISODAM®** has best-in-class tear resistance among the non-latex dental dams tested.

- Dental Advisor

**3+1**  
Standard package

STANDARD PACKAGE				\$26.67	ECONOMY PACKAGE				\$95.49
5" x 5"	Medium	20/pack	ISO01800520		5" x 5"	Medium	100/pack	ISO01800510	
5" x 5"	Heavy	20/pack	ISO02300520		5" x 5"	Heavy	100/pack	ISO02300510	
6" x 6"	Medium	15/pack	ISO01800615		6" x 6"	Medium	75/pack	ISO01800675	
6" x 6"	Heavy	15/pack	ISO01800615		6" x 6"	Heavy	75/pack	ISO02300675	

To read the full Dental Advisor research report, visit the ISODAM® product page on our website and click "data sheets".

[www.fourrubber.com](http://www.fourrubber.com)

Promotion valid July 1 to September 30, 2021



Not made with natural rubber latex

# Beautiful Restorations to last a lifetime.

Sterile burs available that can positively impact your practice by decreasing sterilization costs so you can focus on producing better patient outcomes.

Effective July 1, 2021 - September 30, 2021

ALL STERILE and NON-STERILE CARBIDES

**BUY 5 PACKS  
GET 2 AT NO CHARGE\*\***

CHOOSE FROM: STANDARD OPERATIVE CARBIDES, FG, FGSL, RA, RASL, HP, GREAT WHITE® GOLD, TRIMMING & FINISHING, FISSURATOMY®

CARBIDE BULK PACK INCLUDING: FG, RA, HP, GREAT WHITE®, FISSURATOMY®, TRIMMING AND FINISHING

SPECIALTY BURS: SMARTBURS® II

\*\*Can mix and match shapes to equal or lesser value. Free product to equal or lesser value Valid for 5, 10 and 25 Packs and Bulk Packs. Offer excludes Great White® Ultra, EndoGuide® Burs, Endo Safe End and Oral Surgery Burs.



STERILE | AVAILABLE IN ROLLS OF 25

©2021 SS White® All Rights Reserved. SS White®, SmartBurs® II, Fissurotomy® and Great White® are registered trademarks of SS White Dental. \*Data on file and available upon request. \*\*Free goods must be of equal or lesser value. All offers are limited to stock on hand and are not valid on prior purchases. Offers may not be combined with any other sale, promotion, discount, code, coupon and/or offer. Sales promotion choice must be indicated on invoice sent to SS White® on applicable order. SS White® reserves the right to discontinue a promotion at any time. †SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, SS White® will issue a full refund of the purchase price through your distributor. To claim free goods: On or before October 31, 2021 mail a copy of your Goetze invoice to: SS White Burs, Inc., Attn: Promotions Department, 1145 Towbin Ave., Lakewood, NJ 08701 or email freegoods@sswhitedental.com. Please note, these are end-user doctor promotions and must be redeemed by the qualifying practice. Please allow 4-6 weeks for free delivery.

**SWITE Dental®**

Better Patient Outcomes  
Improved Efficiency  
Faster Practice Growth


1145 Towbin Avenue Lakewood, New Jersey 08701  
Visit us on the web at [sswhitedental.com](http://sswhitedental.com)



**Bold**

**Flexible**

**Powerful**



SCAN CODE TO LEARN MORE



Discover a new  
CAD/CAM Solution

**PLANMECA**



1+1

**iBOND® Universal**

Every surface, every indication, every technique—every time.

iBOND Universal is the **key to superior bonds** in every situation.

- ✓ Easy to dispense
- ✓ Ultra strong
- ✓ No need for a DC activator

**BUY 1 iBOND Universal Kits or Refills,**

**GET 1 of the same at no charge!** Campaign number: 1000806125



**KULZER**  
MITSUI CHEMICALS GROUP

Giving a hand to oral health.

TO REDEEM ADDITIONAL PRODUCTS: Purchase products from your Goetze representative between 07/01/2021 and 09/30/21. Write campaign number on Goetze invoice and fax to Kulzer at (877) 271-5211. Invoices may be mailed to: Kulzer, ATTN: Free Goods Fulfillment, 4315 S. Lafayette Blvd., South Bend, IN 46614, or emailed to: freegoods.NA@kulzer-dental.com. Invoice must be received by 10/31/20. All offers subject to change. Not to be combined with any other offers. Please allow 6–8 weeks from receipt of invoice for delivery. OFFER EXPIRES 09/30/21. © 2021 Kulzer, LLC. All Rights Reserved. iBOND® is a registered trademark of Kulzer, GmbH. Distributed by: Kulzer, LLC. 4315 S. Lafayette Blvd., South Bend, IN 46614. Phone: (800) 431-1785 Fax: (877) 271-5211 kulzerUS.com. Participation is not permitted for customers who participate in any other pricing programs, including government, group practices, special markets, community health centers, schools, institutions, and discount customers. Receipt of free goods, after rebate, cannot be combined with other offers or promotions.


## Integra® Miltex® Products Limit Uncertainty



### Collagen Special ■ July 1 - September 30, 2021

#### Buy 4 Boxes of Collagen, Get 1 (Mix & Match)



HeliTape®, HeliCote® and HeliPlug®, Collagen Wound Dressing, HeliMend® and HeliMend® Advanced Collagen Membrane



Goetze Dental	Description	Quantity	Price
221390	HeliTape® Absorbable Collagen Wound Dressing, 1" x 3" (2.5cm x 7.5cm)	10/box	\$195.70
552270	HeliCote® Absorbable Collagen Wound Dressing, 3/4" x 1-1/2" (2cm x 4cm)	10/box	\$161.26
540526	HeliPlug® Absorbable Collagen Wound Dressing, 3/8" x 3/4" (1cm x 2cm)	10/box	\$128.53

Goetze Dental	Description	Quantity	Price
553286	HeliMend® Absorbable Collagen Membrane, 15mm x 20mm	1/box	\$111.68
553287	HeliMend Absorbable Collagen Membrane, 20mm x 30mm	1/box	\$137.28
553288	HeliMend Absorbable Collagen Membrane, 30mm x 40mm	1/box	\$224.05

Goetze Dental	Description	Quantity	Price
552271	HeliMend Advanced Collagen Membrane, 15mm x 20mm	1/box	\$130.74
553289	HeliMend Advanced Collagen Membrane, 20mm x 30mm	1/box	\$158.00
553290	HeliMend Advanced Collagen Membrane, 30mm x 40mm	1/box	\$230.58



To receive your promotional goods, send an email to [IntegraPromos@integralife.com](mailto:IntegraPromos@integralife.com). Include your complete shipping/contact information, product/quantity desired, along with your copy of your qualifying invoices. Include promo code **GOETZEQ32021** in the subject line.

Offer valid **07/01/2021 - 09/30/2021** in the USA for end-users only and does not apply to government, institutional or dealer stocking orders. By participating in this offer, Customer agrees, as applicable, to properly disclose the discounts or reductions in price in the costs claimed or charges made under any federal or state health care program which provides cost or charge based reimbursement to Customer for the items covered under this offer. Promotion void if a copy of the invoice is not included with promotion redemption request. Offers cannot be combined unless invoice totals meet promotion criteria. All Redemption requests must be submitted within 30 days of expiration date of offer. Please allow 4-6 weeks for goods to be delivered.

- \* Redeemed goods must be of equal or lesser value of lowest priced purchased item.
- \* Redemption invoices subject to verification by Integra Miltex.
- \* Buy/Get promotion applies only to Integra Miltex products advertised in this ad unless otherwise noted.

Integra, the Integra logo, HeliCote, HeliPlug, HeliMend, HeliTape, and Miltex are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. ©2021 Integra LifeSciences Corporation. All rights reserved. 0488744-1-EN 2021-07



# Orabloc®



ARTICAINE HCl 4% and Epinephrine 1:100,000 or 1:200,000



Product Codes:  
1:100,000 566898  
1:200,000 566899

The Quality that you need,  
the Price that you want!

**BUY 10 GET 1 FREE** or  
**BUY 18 GET 2 FREE!**

www.orabloc.com



Give your patients a taste of  
something better with Kolorz.



- Gluten free
- Allergy friendly
- Natural sweeteners
- Single doses
- Great-tasting flavors

**ClearShield™ 5% Sodium Fluoride Varnish**  
**Buy 3 Boxes, Get 1 Box FREE!\***

**Kolorz® Prophylaxis Paste**  
**Buy 1 Box, Get 1 Box FREE!\***

Kolorz® products offer the benefits of advanced GLUTEN-FREE formulas, convenient application options, non-caloric natural sweeteners and XYLITOL. Our proprietary flavorings were developed in consultation with food industry professionals that are guaranteed to taste great.

\*Free Kolorz Prophy Paste is automatically shrink-wrapped with the purchased product. \*\*Offer valid 7/1/21-9/30/21. To receive free goods, fax your paid Goetze invoice to DMG America at 201-894-0213, or email to [freegoods@dmg-america.com](mailto:freegoods@dmg-america.com), indicating type of free good. Each qualifying purchase must be submitted within 30 days. Promotion cannot be combined with any other offer and may change or be discontinued at any time without notice. Limit 5 offers per dental office. Not available to dental supply dealers.



[dmg-connect.com/offers](http://dmg-connect.com/offers) | 1.800.662.6383



# Ivoclar Vivadent 2021 Q3 Specials

## CURING LIGHTS

Bluephase® PowerCure™ System  
Bluephase® G4  
Bluephase® Style

**BUY:** \$3,500 of composite refills\*  
OR \$4,500 of clinical product\*

**GET:** 1 Curing Light FREE!

\*The composite material or clinical products purchase must be on a single invoice.



## COMPOSITES

Tetric® Prime  
Tetric® PowerFill  
Tetric® PowerFlow

High performance, universal composite system featuring Prime handling



**SAVE 25%  
ON 4 REFILLS**

(purchase 3 refills and receive 1 refill free of same brand and delivery form)

Tetric EvoCeram®  
Tetric EvoFlow®  
IPS Empress® Direct



**SAVE 20%  
ON 5 REFILLS**

(purchase 4 refills and receive 1 refill free of same brand and delivery form)

**SAVE 30%  
ON 13 REFILLS**

(purchase 10 refills and receive 3 refills free of same brand and delivery form)

## CEMENTS

Variolink® Esthetic  
Adhesive Cementation System

Recommended for Lithium Disilicate and Composite

Multilink® Automix  
Adhesive Cementation System

Recommended for Lithium Disilicate, Zirconia, Metal Ceramic or Composite



**SAVE 25%  
ON 4 REFILLS**

(purchase 3 refills and receive 1 refill free of same brand and delivery form)

## ADHESIVE

Adhese® Universal  
Dental Adhesive

- Efficient VivaPen® delivery
- Universal application
- Predictable results



**SAVE 25%  
ON 4 REFILLS**

(purchase 3 refills and receive 1 refill free of same brand and delivery form)

Qualifying purchases must appear on a single invoice and must be made between 07/01/21 – 09/30/21. Request for Free Goods must be received no later than 10/31/21. To redeem your promotional goods, please visit: [www.ivoclarvivadent.us/shop/promo](http://www.ivoclarvivadent.us/shop/promo) Ivoclar Vivadent reserves the right to modify or cancel these promotions at any time. These specials cannot be combined with any other offer or rebate program.

[ivoclarvivadent.com](http://ivoclarvivadent.com)  
**Making People Smile**

**ivoclar  
vivadent®**



# TIME TO SWITCH

Aerosol reduction with the flip of a switch.

**Water Jet**  
(Water Only)

**Water Spray**  
(Water + Air)



## NEW Ti-Max Z95L

### Easy On/Off Switch

Switch between Water Jet (water only) and Water Spray (water+air) in an instant.

Instantly switch between Water Spray or Water Jet **with the flip of a switch** with the world's first two-way spray technology. When extra safety is needed, just the Water Jet provides maximum protection while continuing to offer optimum cooling. Switching back to Water Spray is a breeze – just twist the switch!

## Valid July 1 - September 30

**Buy 2** Z Series Electric handpieces

Z95L, Z85L, Z25L, Z45L, nano95LS, nano25LS



**Get 1** Z Series Electric handpiece

(of equal or lesser value)

**Buy 3** Z95L or nano95LS handpieces



**Get 1** NLZ Electric Micromotor System


\*All NSK no charge goods requests must be submitted to NSK. Please send a copy of your Goetze Dental invoice to NSK promotions, 1800 Global Parkway, Hoffman Estates, IL 60192; or fax to: 1-800-838-9328 or email to: orders@nskdenal.com. To be eligible, you must submit by October 31, 2021.



NSK America Corp. [www.nskdenal.com](http://www.nskdenal.com)

1800 Global Parkway, Hoffman Estates, IL 60192, USA • TEL: +1 888 675 1675 • FAX: +1 800 838 9328

# Does your staff know why patients prefer **TRI-OLGY®** for *Mouth Pain & Sores?*



TV Commercial



Brochures



Contact your GOETZE Rep to order these **TRI-OLGY®** Brochures for your office - **FREE!**

# WE'RE INVESTED IN YOU

**100%** of our energy and resources will continuously center on the products and services we offer to the **privately held dental practice** through —

INTEGRITY  
EDUCATION & INNOVATION  
& customer service

**GOETZE**  
d e n t a l

EXCITING INNOVATIONS FROM **Crest** + **Oral-B**

**BUY 6 CASES, GET 1 CASE OF CREST® ORAL-B® SENSITIVE MANUAL SOLUTION FREE!†**

SHIPPED WITH ORDER

SAY HELLO TO **Oral-B** iO™



**CREST® + ORAL-B® iO™  
TRANSFORMATIONAL GUM HEALTH SYSTEM**

- Oral-B® iO™ Electric Toothbrush
- Crest® PRO-HEALTH™ Gum Detoxify™ Toothpaste (4.1 oz)
- Crest® PRO-HEALTH™ Multi-Protection Clean Mint Mouthwash (250 mL)
- Oral-B® Glide™ PRO-HEALTH™ Advanced Floss (15 m)

Brush head refills include: Ultimate Clean, Gentle Care

3 units/case (582397)

\$570.00/case - just \$190.00 per patient


**TRANSFORM PATIENTS' GUM HEALTH, STARTING DAY 1.**

In a clinical study, **100%** of patients moved to a state of healthy gums after just **12 weeks** of use with the Oral-B® iO™ Electric Toothbrush.



**CONCENTRATED ENERGY**

Innovative magnetic drive delivers concentrated energy to the tips of the bristles where it's needed most, allowing the brush head to glide tooth by tooth for a smooth, quiet cleaning experience.



**OSCILLATIONS + MICRO-VIBRATIONS**

Combines the oscillation rotations from Oral-B®'s iconic round brush head with the gentle energy of micro-vibrations for an enjoyable experience and a sensational clean.



**SMART PRESSURE SENSOR**

The Smart Pressure Sensor provides guidance on how to brush with the right amount of pressure to help protect gums. The ring lights up green when optimal pressure is applied, and red when the pressure is too hard.

**BUY 6 CASES, GET 1 CASE OF CREST® ORAL-B® SENSITIVE MANUAL SOLUTION FREE!†**

SHIPPED WITH ORDER

**NEW!**

**Crest**®  
**WHITENING EMULSIONS** | with LED ACCELERATOR LIGHT

AT HOME PROFESSIONAL LEVEL WHITENING TREATMENT

+

3 MINUTE BOOST LIGHT WITH VIRTUALLY NO SENSITIVITY

(582400) 3/CASE

\$207.66/CASE - JUST \$69.22 PER PATIENT



\*vs. a leading whitening pen brand

† Offer valid 7/1/21-9/30/21. Limit 1 of each of the above deals per office. No mix and match. One case of Crest Oral-B Sensitive Manual Solution contains 72 bundles.

Crest® + Oral-B® products are sold at dental professional pricing for limited use in your practice. Products purchased through professional dental distribution channels are not meant for resale outside your dental care practice. © 2021 P&G ORAL-28449

continuing the care that starts in your chair

**Crest** + **Oral-B**

# THE HOME SLEEP TEST SOLUTION YOU NEED

SIMPLE | FAST | ACCURATE



Simplify the Process



Get Access To Sleep Physicians



Treat More Patients

## WE'VE GOT YOU COVERED

25% Of Your Patients Need A Sleep Test

HAVE YOUR OWN	USE OURS
Hand out your test	Submit the order
Upload results	A2S ships the test
Use our sleep physicians	Receive HST report

**GET STARTED TODAY!**

(866) 631-2232 or [Info@awaken2sleep.com](mailto:Info@awaken2sleep.com)

# BRIDGING THE SUPPLY CHAIN GAP\*

# ORCA™ FOAM

50 individually blister packed, sterile cubes - 1cm x 1cm x 1cm each

## Absorbable Gelatin Sponge Hemostat, USP (Standard)



- Ideal size for easy application
- Absorbs up to 35x its weight
- Rapid hemostasis
- Aids in control of hemorrhage in surgical or extraction sites
- No cutting required
- Helps prevent cross contamination
- Fast, complete resorption
- Excellent carrier medium for fluids such as antibiotics or ORCA Professional Dry Socket
- Designed and packaged for oral surgery and dental application



Patient Safety and Quality Control are our top priorities. ORCA Foam™ is a natural and safe FDA approved product made from 100% pharmaceutical-grade gelatin and manufactured in a cGMP certified facility.



### About ORCA Professional

The founder of ORCA Professional is an Oral and Maxillofacial Surgeon with 25+ years of clinical practice experience who clearly understands the needs of doctors in his field of practice to successfully treat patients. He personally uses ORCA products daily and stands behind their quality and integrity. He is always generous in offering his time to provide product knowledge and information to assist sales representatives, distributors and practitioners.



*Innovation. Integrity. Trust.*

\*ORCA understands the deficiency regarding lack of product availability and will do everything possible to insure product supplies remain available

# A GAME CHANGING APPLICATION FOR DRY SOCKETS\*

## ORCA™ DRY SOCKET

for dry sockets and oral pain

Non-Medicated • Only Natural Ingredients • Clinically Proven



Designed and Packaged for Oral Surgery and Dental Applications

Highly Effective, Non-Medicated Alternative, Consists only of Essential Oils, Aloe Vera and other Proprietary Natural Ingredients

Outperforms Other Products Hands Down!

Clinically used in Oral and Maxillofacial Surgery practice for over 16 years with thousands of applications. Over 95% of patients experience complete resolution of dry socket symptoms with one single treatment

Simple and Effective Delivery

Recommended to be applied to socket or extraction site with ORCA Foam™ which will dissolve within a week without the need for removal. Delivery methods may also include an alternative carrier that will dissolve or gauze strip (removal required)



**Directions:** For adults and children 12 years of age and older. In a sterile fashion, place 8 to 10 drops onto ORCA Foam™, similar resorbable gelatin sponge or gauze strip. Apply directly to dry socket or areas of discomfort in oral region.

FOR PROFESSIONAL USE ONLY

### About ORCA Professional

The founder of ORCA Professional is an Oral and Maxillofacial Surgeon with 25+ years of clinical practice experience who clearly understands the needs of doctors in his field of practice to successfully treat patients. He personally uses ORCA products daily and stands behind their quality and integrity. He is always generous in offering his time to provide product knowledge and information to assist sales representatives, distributors and practitioners.

**ORCA™**  
**PROFESSIONAL**  
*Innovation. Integrity. Trust.*

\*The founder of ORCA has had years of referrals from colleagues that have had dry sockets which were refractory to other products and ORCA Professional has always risen to the occasion to resolve discomfort!

# Clear Comfort ALIGNERS



Invisible  
Comfortable  
Convenient



Advancing the Brilliance of Your Smile  
[ClearComfortAligners.com](http://ClearComfortAligners.com)



# Why Choose Clear Comfort Aligners?



## More Comfort, Less Time

Efficiency driven treatment planning combined with our patented bio-progressive aligners (soft, medium, hard) result in superior comfort for your patient and less time in treatment. Clear Comfort achieves up to 1.0mm of movement every three weeks while other aligners are only able to accomplish approximately 0.45mm every four weeks.



## Half the Cost of Other Brands

Only pay for what you use! Clear Comfort is priced according to the actual number of steps that are required instead of arbitrary fixed tiers. Cases consisting of one to twelve steps include one refinement at no additional charge and cases requiring thirteen steps or more include unlimited refinements. All cases include retainers.



## Attachments Not Required

The extended gingival coverage of Clear Comfort Aligners adds comfort and generates the leverage, force, and control needed for achieving challenging movements like rotations and angulation without attachments.



## Refinements the Exception, Not the Rule

True orthodontic principles, better fitting aligners, and advanced technology means that tooth movement is more precise and predictable. The Clear Comfort Aligner System delivers great results and the need for refinements is significantly reduced.



## Always Your Partner, Never Your Competition

Other companies seem determined to go direct to the consumer and eliminate you, the professional. Clear Comfort is committed to preserving your indispensable role as the provider as well as to helping you take advantage of the exploding consumer demand for clear aligner treatment.



**Clear Comfort**  
ALIGNERS

# Introducing **Pac-Cide XT™** Disinfectant Spray, Wipes and Dispenser



70.5% Ethanol and 0.2%  
Chlorhexidine Gluconate



Safe on Surface  
and Skin



Effective in  
One Minute



Potent  
Disinfectant



## Order Information

ITEM#	DESCRIPTION
PCXT-3L	Pac-Cide XT Surface Disinfectant Cleaner, 3 L solution
PCD-180	Pac-Cide XT Dry Wipes, 6.29" x 7.08", 180 wipes
PCXT-02	Pac-Cide XT Surface Disinfectant Cleaner Starting Kit

\* To qualify and receive your promotional goods, please submit a valid proof of purchase to [freegoods@pac-dent.com](mailto:freegoods@pac-dent.com). Cannot be combined with any other offer. No Cash Value. Equal or lesser value items(s) only. Offers are valid only at participating dealers in the continental US, excluding AK, and HI. Promotions are subject to change or cancellation without notice. Please allow 6-8 weeks for delivery.

📍 670 Endeavor Circle, Brea, CA 92821

☎ 909.839.0888

✉ [info@pac-dent.com](mailto:info@pac-dent.com)

🌐 [www.pac-dent.com](http://www.pac-dent.com)

MEET EDDSON



Introducing all of the tools you need to run a successful dental practice.



Intelligent Workflow for you and your team.



Enhanced Business Performance with dynamic and interactive Accounting and Reporting tools.



PEAKS Dashboard and Analytics built directly into EDDSON.



Native Web Forms to simplify capturing and digitizing patient information.



Plus many other features and modules engineered around the way you actually work to delivery efficiency and effectiveness in your busy day.

	BEST VALUE		
	Core	Plus	Pro
<b>Monthly Fees</b>	<b>\$250</b> /mo	<b>\$375</b> /mo	<b>\$495</b> /mo
<b>Updates &amp; Support</b>			
Software Updates	✓	✓	✓
Email & Realtime Support	✓	✓	✓
Training & Expert Services Discounts	—	10%	15%
<b>EDS Services</b>			
E-Claim Unit Price	35¢	30¢	25¢
Insurance Eligibility & Electronic Remittance Service	\$40	\$20	✓
<b>E-Prescribe by DoseSpot</b>			
E-Prescribe Per Provider	\$50	\$45	\$40
<b>Bill Pay</b>			
DocPay™ Integration Module*	\$20	✓	✓
<b>EDDSON Optional Modules</b>			
Time Clock	✓	✓	✓
Patient Check In/Update Module	\$10	✓	✓
Membership Module	\$30	✓	✓
Web Forms	\$40	\$20	✓
Performance Dashboard	✓	✓	✓
PEAKS Analytics	\$80	\$40	✓
QuickBooks® Integration	\$50	\$30	✓
VOIP Integration**	\$70	\$40	✓
Text Messaging & Appointment Confirmation***	\$40	\$20	✓
Patient Engagement & Review Engine - Powered by Swell™	\$219	\$199	\$169



Visit [eddsondental.com](http://eddsondental.com) or call 800.692.0804 to learn more.

Additional charges apply for training and software conversion.

\*Requires service with DocPay™

\*\*Requires phone system and service sourced by EDDSON from AVID Communications

\*\*\*2,500 messages per month included + \$0.03 for each additional message



## ABOUT US

Goetze Dental is a family-owned business serving private dental practices just like yours since 1884.

We work to understand your goals and dreams and strive to provide the best dental supplies, equipment, technology and service.

### We believe the best patient care is delivered in the private practice.

100% of our energy and resources will continuously center on the products and services we offer to the privately held dental practice through integrity, education, innovation and customer service.

We stand united with the artists, the scientists and the explorers of the dental community, whose interest in providing exceptional patient care never fades.

### Notice about PPE and Infection Control Supply Availability and Purchases

As we work through the changes brought on by the pandemic, we are committed to doing what's best for you and your team, and the entire community of privately-owned dental practices we serve, as well as for our employees.

We continue to work closely with our manufacturing partners to provide accurate and current information.

Due to the continued global demand for PPE and other infection control supplies, orders for these supplies will not be guaranteed for delivery and are subject to backorder status, allocation of amounts and/or cancellation. Orders placed cannot be canceled or returned.

**Please reach out to your private practice specialist or call customer service at 800.692.0804 with questions.**

# FOUR EASY WAYS TO ORDER

**GAINSpeed™**  
desktop ordering application



**GoetzeOnline.com**  
e-commerce ordering option



**800.692.0804**  
8am - 5pm CST



Fax:  
**800.552.7791**

