

PWR

Piezo Scaling & Air Polishing



Each patient provides unique challenges and only you have the power to choose how to apply your wealth of knowledge and skills. **Power ON with NEW PWR Piezo Scaling and Air Polishing products from HuFriedyGroup.** Choose Best in Practice calculus and stain removal, biofilm disruption, and implant cleaning with the **PWR Pair** and **PWR Piezo**. Built for your patient, designed for your power.

PWR Pair

Piezoelectric Scaler & Air Polisher



APPAIR

The power of all-in-one brings versatile air polishing to each clinician treating each patient. The **PWR Pair** brings you the power of supra- and subgingival calculus and stain removal, biofilm disruption, and implant cleaning.

PWRful Air Polishing Features

- Regulated, Fixed Air Pressure
- Continuous Flow of Air
- Slim Air Polishing Steel Handpieces
- Sub- and Supragingival Treatment
- Heated Water



PWRful Piezo Features

- Constant Performance Monitoring
- Soft Mode
- Three Different Power Zones
- Wide Array of PWR Piezo Tips
- Pulse Mode
- Slim LED Handpiece
- Irrigation System

PWR Piezo Piezoelectric Scaler



UPPIEZO

Put the PWR in power scaling with the **PWR Piezo**. Our Piezo technology allows for a comfortable experience for patient and clinician. The advanced feedback system gives the clinician perfect control over performance to improve results and increase comfort.

PWRful Piezo Features

- Constant Performance Monitoring
- Soft Mode
- Three Different Power Zones
- Wide Array of PWR Piezo Tips
- Pulse Mode
- Slim LED Handpiece
- Irrigation System



PWR Piezo Scaling

The PWR of piezo scaling is the power of choice. Power scale the way you choose with the expansive array of **PWR Piezo Tips**. From scaling to perio to endodontic and more, choose from an unrivaled variety of PWR Piezo Tips—for every indication, every treatment, and every patient you will face in your daily practice. Compatible with our PWR Pair and PWR Piezo devices.

PWR Piezo Tips | Scaling


S1 Piezo Tip	UPS1
S1-S Piezo Tip	UPS1-S
S2 Piezo Tip	UPS2
S3 Piezo Tip	UPS3
S4 Piezo Tip	UPS4
S5 Piezo Tip	UPS5

PWR Piezo Tips | Perio

P2 Piezo Tip	UPP2
P3 Piezo Tip	UPP3
P10 Piezo Tip	UPP10
P11 Piezo Tip	UPP11
P12 Piezo Tip	UPP12
P13 Piezo Tip	UPP13
P14 Piezo Tip	UPP14
P15 Piezo Tip	UPP15
P16L Piezo Tip	UPP16L
P16 R Piezo Tip	UPP16R

PWR Piezo Tips | Endodontic

ER1 Piezo Tip	UPER1
ER2 Piezo Tip	UPER2
ER3 Piezo Tip	UPER3
ER4 Piezo Tip	UPER4
ER5 Piezo Tip	UPER5
RD3 Piezo Tip - Diamond Coated	UPRD3
RD4 Piezo Tip - Diamond Coated	UPRD4
RD5 Piezo Tip - Diamond Coated	UPRD5



PWR Piezo K10 Tip Wrench
UPK10

Visit us online at [HuFriedyGroup.com/PWR-ON](https://www.HuFriedyGroup.com/PWR-ON)

PWR Air Polishing



PWR Air 90°
Handpiece
APHP90



PWR Air 120°
Handpiece
APHP120



PWR Air Perio
Handpiece
APHPPERIO

Achieve Best in Practice access with the **PWR Pair's 90° and 120° Angled Handpieces** for supra- and subgingival use in periodontal pockets up to 5 mm deep. All Air Polishing Handpieces easily attach with just a click and are safely secured, thanks to a dedicated security connection system.

The **PWR Pair's Perio Handpiece** equipped with the disposable sterile PWR Subgingival Perio Tip (included with the PWR Pair) provides optical access in pockets deeper than 5 mm. All Air Polishing Handpieces are resilient to cracking, being made from steel. Practice with the PWR of peace of mind: less interruptions and less time spent replacing handpieces.

PWR Prophylaxis Powders

remove supra- and subgingival biofilm, stains, and pigmentations along with treating peri-implantitis and cleaning of prosthesis and restorative materials. Choose from PWR Perio Powder (Glycine) and PWR Classic Powder (Sodium Bicarbonate), both coming in mint and neutral flavors. Our specially developed PWR Prophylaxis Powders are also hydrophobic, which helps lessen clogging in the PWR Pair.



Description/Flavor	Powder Type	SKU
Perio Powder <i>Mint</i>	Glycine	AP-PM
Perio Powder <i>Neutral</i>	Glycine	AP-PN
Perio Powder helps with removal of supra- and sub gingival biofilm, periodontal pocket detoxification, implant maintenance, peri-implantitis treatment, cleaning of fixed orthodontic appliances and cleaning of prosthesis and restorative materials.		
Classic Powder <i>Mint</i>	Sodium Bicarbonate	AP-CM
Classic Powder <i>Neutral</i>	Sodium Bicarbonate	AP-CN
Classic Powder helps with removal of supragingival biofilm, stains and pigmentations.		