

Planmeca PlanMill® 30 S

Efficient and cost-effective milling unit

High-accuracy
entry-level unit



High-accuracy entry-level milling unit for dental clinics

Enter the world of digital dental production with this cost-efficient single-spindle milling unit. **Planmeca PlanMill® 30 S** has been developed for dental businesses looking for efficient, accurate, and cost-effective milling solution. This mill guarantees a professional entry to the world of CAD/CAM technology.

- State-of-the-art design
- High accuracy
- High-speed single-spindle (max. 100k)
- Milling a crown takes approx. 16 to 18 minutes
- Rotary axis for milling blocks of your choice
- Fixed 5-station tool changer
- Block capacity 20 x 40 x 85 mm
- Dimensions when closed (w x h x d) 661 x 442 x 544 mm (26 x 17.4 x 21.4 in.)



Materials

Ceramic materials

- IPS e.max CAD
- IPS Empress CAD
- IPS Empress CAD Multi
- VITA SUPRINITY PC
- VITABLOCS Mark II
- VITABLOCS TriLuxe forte
- Straumann nIce

Hybrid ceramic materials

- VITA ENAMIC
- VITA ENAMIC multiColor
- GC CERASMART

Resin-based materials

- Telio CAD

Zirconium materials

- IPS e.max ZirCAD

PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, DentoVac, Digital perfection, Economat Plus, Elegant, Flexy, Mini-dent, Perio Fresh, PlanEasyMill, Planmeca 4D, Planmeca AINO, Planmeca ARA, Planmeca CAD/CAM, Planmeca CALM, Planmeca Chair, Planmeca Clarify, Planmeca Compact, Planmeca Creo, Planmeca Emerald, Planmeca FIT, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Lumo, Planmeca Maximity, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca mRomexis, Planmeca Noma, Planmeca Olo, Planmeca Online, Planmeca PlanCAD, Planmeca PlanCAM, Planmeca PlanClear, Planmeca PlanID, Planmeca PlanMill, Planmeca Planosil, Planmeca PlanPure, Planmeca PlanScan, Planmeca PlanView, Planmeca ProCeph, Planmeca ProFace, Planmeca ProID, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca Serenus, Planmeca SingLED, Planmeca SmartGUI, Planmeca Solanna, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Viso, Planmeca Verity, Planmeca Waterline Cleaning System, Planmeca Xtremity, Proline Dental Stool, ProTouch, Saddle Stool, SmartPan, SmartTouch, Trendy, and Ultra Relax are registered or non-registered trademarks of Planmeca in various countries.

