


Barracuda™ Water Ring Pumps

Single Water Ring Pump

Our Barracuda Single Water Ring Pump offers quiet operation and compact size designed for small to medium size practices.

- Large capacity, easy to clean filter, with clear bowl for easy inspection
- 1 & 2 HP Water Injected Pump models available to suit your specific requirements
- 10.55 CFM for 1 HP model and 23.5 CFM for 2 HP Model
- Quiet and compact
- Long life, chemical resistant brass impeller and housing
- Sealed bearings in motor
- Designed for easy servicing

NOW
with a **5 year**
limited warranty
from date of
installation.




Single Water Ring Pump
CV-101

Dual Water Ring Pump

For larger practices, the Barracuda Dual Water Ring Pump provides vacuum suction when operating in peak work periods. The dual pump offers back-up protection to minimize downtime and costly interruption to workflow.

- Same quality features found in RAMVAC® single pumps
- Two independent pumps and relay boxes to handle peak loads and back-up protection for 100% dual redundancy
- Optional water recirculator
- Optional top shelf with 300 lb. capacity for stacking compressors



Dual Water Ring Pump
MC-201

Specifications

	MC-201FS/FSW	MC-202FS/FSW	CV-101FS/FSW	CV-102FS/FSW
Horsepower	2 x 1 HP	2 x 2 HP	1 x 1 HP	1 x 2 HP
Volts	208-230	208-230	115 / 208-230	208-230
Amps <i>(each motor)</i>	7.5 each motor	15 each motor	13.2 / 7.2 - 6.6	11.9 - 10.7
Phase	Single	Single	Single	Single
RPM	3450	3450	3450	3450
Ambient Temperature	32° C / 104° F	32° C / 104° F	32° C / 104° F	32° C / 104° F
Cycle	60 Hz	60 Hz	60 Hz	60 Hz
Width	24"	24"	12"	14"
Depth	17"	17"	9"	11"
Height	24"	24"	16"	20"
Dimensions with Recirculator <i>(W x D x H)</i>	26.25" x 19.13" x 24"	26.25" x 19.13" x 24"	N/A	N/A
Shipping Weight	175 lbs.	215 lbs.	61 lbs.	78 lbs.
Max Sealed Vacuum	Approx. 20 - 25 inches Hg adjustable	Approx. 20 - 25 inches Hg adjustable	Approx. 20 - 25 inches Hg adjustable	Approx. 20 - 25 inches Hg adjustable
Users	4	6	2	3



For more information, please visit www.DentalEZ.com

Need a RAMVAC part? Visit www.DentalEZParts.com