



STATMATIC™

Handpiece Maintenance Unit

ULTIMATE CARE SYSTEM FOR HANDPIECES

and their chuck mechanisms.

STATMATIC™ is the automatic care and maintenance system that properly cleans, lubricates, and purges up to 3 handpieces prior to sterilization. The STATMATIC also incorporates a unique chuck cleaner to dislodge debris inside the chuck. You and your staff will profit from greater efficiency and fewer repairs of your handpieces.



Extends the Life of your Handpieces

It properly cleans, lubricates, and purges your handpieces to increase their working lifetime.



Optimizes Staff Productivity

Its automated process optimizes working procedures in your practice, allowing staff to focus on alternative tasks.



Time-Saving

Up to 3 handpieces can be cleaned, lubricated, and purged in less than a minute.



STATMATIC Technical Specifications

p/n S30000-NA

Dimensions (L x W x H)	28.5 x 19 x 40 cm / 11.4" x 7.6" x 16"	Sound Level	60 (65), dB (A)
Weight	8.6 kg / 18.9 lbs	Compressed Air	40 NL/min
Electrical Ratings	Rated voltage: 100-230 V Rated frequency: 50/60 Hz Rated power: max. 30 VA	Air Consumption	max. 85 %
Protection Class	1	Compressed-Air Quality	Residual oil content: max. 0,1 mg/m3 Residual dust: max 1µ; max. 1 mg/m3 Residual water: max. 0.1 g/m3 at -40°C
Ambient Temperature	5-40° C / 41-104 °F		



Dedicated chuck care port for maintaining the handpiece chuck



Easy spray can insertion with unique safety valve

p/n STATMATIC Spray Can: STM500-6NA



Absorption pad for excess lubricant and debris

Recommended Accessories



STATIS™

Air-driven high speeds

Made in Germany, SciCan's STATIS ML (medium head) and STATIS SL (small head) are powerful, efficient and ergonomic.



SANO™

E-Type handpieces

Electric and air handpieces designed for optimal control, power, and efficiency.



E-STATIS™

Motors

Fully customizable for high speed, low speed, and endodontic procedures and coupled with one of the smallest electric motors in the world.

Manufactured for:

SciCan Ltd.

1440 Don Mills Road,
Toronto, ON M3B 3P9

Tel: 416-445-1600

Fax: 416-445-2727

Toll free: 1-800-667-7733

customerservice@scican.com



Distributed by:

SciCan Inc.

701 Technology Dr.,

Canonsburg, PA 15317, USA

Tel: 724-820-1600

Toll Free: 1-800-572-1211

uscustomerservice@scican.com

EU Representative &

Distributed by:

SciCan GmbH

Wangener Strasse 78

88299 Leutkirch

Deutschland

Tel.: +49 (0)7561 98343 - 0

Fax: +49 (0)7561 98343 - 699

customerservice_eu@scican.com

Distributed by:

Office for Swiss Customers:

SciCan Medtech AG

Alpenstrasse 14

6300 Zug

Switzerland

Tel: +41 (0) 41-727 7027

Fax: +41 (0) 41-727 7029