



Power, Performance and Flexibility

The StarETorque Electric System combines the power and performance you've come to expect from StarDental. With customizable settings for multiple procedures, the StarETorque seamlessly combines power with performance. From caries removal to crown preparation, you can move from high speed to low speed settings with one touch.

The Power of Performance

- Compact Design for Stand-Alone or Detachable Display or Modular Mounting
- Textured Finish on Attachments for Comfortable Grip
- User Friendly Display with Customizable Operative Modes
- Multiple Attachments Cover a Broad Range of Operatory Needs



- Quiet, 60 W LED Brushless Motor
- Lightweight at 62 g to Minimize Fatigue
- 100–40,000 rpm Motor Speed
- Fully Autoclavable



StarETorque Attachments



1:5
Friction Grip
10,000 – 200,000 rpm
Cutting, crown and bridge



5:1
AutoLatch
non-Fiber-Optic
With No Light



1:1
AutoLatch
1,000 – 40,000 rpm
Caries removal



5:1
AutoLatch
Fiber-Optic
With Light



1:1
Straight
1,000 – 40,000 rpm
Chairside trimming and denture adjustment



StarDental[®] recommends lubricating your electric attachments after each use and prior to sterilization.

Personalize your handpiece with your group name or logo!
Ask your StarDental Customer Service Representative for more information.



For more information, please visit www.DentalEZ.com