

Hygiene Products Better Hygiene Made Easy

Hygiene handpieces from StarDental® feature flexibility and ergonomic designs to provide the very best in comfort and performance.

iStar™ Prophy Handpiece

Our iStar cordless prophy features a simple, compact, cordless design that dramatically simplifies the prophylaxis procedure. No more stress from handpiece tubing pulling or searching for a rheostat. The iStar is ideal for the modern dental operator and features thoughtful design elements that make work more efficient and comfortable.

- Cordless design eliminates potential tripping hazard from tubing
- Easy speed selection of five levels from 500-2,500 rpm
- Battery powered design has plenty of torque to get the job done
- Chargeable between patient use, iStar lasts a full day on a single charge
- On/Off finger tip control for intuitive operation

Select the disposable prophy angle of your choice!



Works with most prophy angles.



NOTES:

WARNING: Reproductive Harm - www.P65Warnings.ca.gov/product

Prophy Star® 3

Constructed with anodized aluminum and weighing just 62 grams, the lightweight Prophy Star® 3 helps minimize operator fatigue. The LubeFree design saves time and money. The 360-degree swivel provides better control for greater precision and comfort.



Titan® Blis-sonic™ Scaler

Lightweight Ergonomic
Sonic Scaler

Our scaler effectively removes hard calculus deposits and stains. It connects to existing air-driven tubings and foot control. Ergonomic silicone grips with comfortable positioning to reduce fatigue.



Scaler Tips

Scalers are packaged with five exclusive StarDental patented scaling tips: universal, sickle, perio, and perio right & left. Choose the tip that best suits your particular needs.



Please visit www.DentalEZ.com for more information.