

StarDental®

The Power of Performance




Why StarDental®?

For over 100 years, StarDental has manufactured easy to use, high quality handpieces that really perform. Our handpieces, hygiene products and electric systems are known the world over for dependability and performance. So when you think StarDental, think Power and Performance.



430 Series, Solara® & Star® Simplicity™ High-Speed Handpieces

StarDental's tradition of power and performance runs through all our high-speed handpieces, including the **430**, **Solara QT** and **StarSimplicity** Series of products. All StarDental products are designed with power and performance in mind.



430 Torque
The newest member of the 430 Family.



Solara QT
Lightweight titanium design.



StarSimplicity
Durable and economically designed.

The Perfect Fit

High Torque Design

- 27 Watts of power* for 430 Torque

Small Head Size for Better Visibility

- 10.5 mm head for Solara
- 11 mm for head 430 Series and Simplicity

Durable and Lightweight

- Titanium reduces hand and wrist fatigue for Solara

Ceramic Autochuck Turbines

- Smooth, concentric operation

Coaxial Water Spray

- Water spray technology keeps operatory site clear and cool

Bifurcated Fiber-Optics

- Dual beam glass design prevents shadowing in front of bur (*fiber-optic models only*)

The Power to Choose

- Available in LubeFree or Lubricated** for 430 and Solara

*Yapp, R., Cowen, M., & Powers, J. (2016). Handpiece Torque vs. Speed Performance (Rep.). Anne Arbor, MI: The Dental Advisor Microbiology Research Center. Data on file with StarDental. **In select models.

Star[®]ETorque[™] Electric System

The new **StarETorque** Electric System is designed to perform. The intuitive design and simple touch screen make every procedure easier, while delivering the flexibility you demand for high and low-speed procedures. The new shorter, brushless motor reduces overall weight, and is fully autoclavable. Performance, Power and Flexibility. StarETorque!



Performance and Flexibility

New Motor Design

- Short, lightweight, 62 g autoclavable brushless motor

Redesigned User Interface

- 16-Color touchpad display for intuitive setting changes

Attachments

- 5 attachments with a textured, satin finish for high and low-speed procedures

Work in any Mode

- 6 customizable presets for both ENDO and PREP modes

Flexibility

- From ENDO to PREP modes with just one touch

Titan® E, Titan-T™, Titan® 3 Low-Speed Handpieces

The StarDental Titan low-speed system provides all the flexibility and performance you'll need. The **Titan 3** offers legendary StarDental performance and durability. The **Titan-T** offers the same performance in a lightweight titanium body. The **Titan E** offers that same performance and e-coupling flexibility. So whichever you choose, you'll get the power of flexibility and performance.



Titan-T



Titan 3

Compliant with CDC Sterilization Guidelines

System Flexibility

- RPM range of 100-20,000 RPMs

Quick Connect Swivel

- 360° Quick Connect Swivel to reduce tubing drag and allow quick disconnect

Sterilization

- Fully autoclavable attachments and motors

Chairside Repairable Angles

- Minimize operator downtime

The Power to Choose

- Angles available in LubeFree or Lubricated

iStar™, Prophy Star® & Titan® Blis-sonic™ Scaler Hygiene Products

StarDental's hygiene line is focused on simplicity, performance and great design. From the **Titan Blis-sonic Scaler** and **Prophy Star 3** to the **iStar Cordless Handpiece**, StarDental has a simple solution that performs, no matter what your needs are. So when you think StarDental hygiene, think simplicity, performance and design.



Prophy Star 3



Titan Blis-sonic Scaler



iStar

Hygiene Made Simple

Prophy Star 3

- Lightweight, 62 g aluminum design and 360 degree swivel minimizes operator fatigue

Titan Blis-sonic Scaler

- Powerful Titan motor effectively removes hard calculus with this lightweight and ergonomic handpiece
- Available in swivel and fixed back end

iStar

- Simplifies procedures with an easy-to-use sleek, modern design
- Eliminate the fatigue of cord drag with this lightweight, fingertip controlled design
- Plenty of power capable of performing cleaning procedures



DentalEZ Integrated Solutions is committed to providing real solutions to everyday challenges in oral healthcare by uniquely combining innovation focused on simplification and efficiency in value based products and outstanding customer service and support.

For more information, please visit www.DentalEZ.com.

ColumbiaDentoform®

DentalEZ® Equipment

NevinLabs™

RAMVAC®

StarDental®

© 2018 DentalEZ, Inc. Dentoform, DentalEZ, Solara, Star, Simplicity, and Titan are registered trademarks and Blis-sonic, ETorque, iStar, NevinLabs, and Titan-T are trademarks of DentalEZ, Inc. RAMVAC is a registered trademarks of RAMVAC Dental Products, Inc.