

430 Series High-Speed Handpieces

The Power of Performance

The 430 series of handpieces from StarDental® delivers just the right combination of power and choice to drive your productivity. The 430 Torque has the highest power ratings tested at 27 watts*.

430 Series

The 430 Series handpieces are ergonomically designed with a stylish stainless steel or satin finish.

- Available in LubeFree, Lubricated, Fiber-Optic, and non-Fiber-Optic
- **Small head design** — Size and shape promote superior operatory field visibility
- **Double bucket rotor design** — Delivers powerful performance
- **Vortex washer** — Unique washer designed to minimize debris “suck back” into housing, preventing turbine damage
- **Bifurcated glass** — Dual beam glass design prevents shadowing in front of bur (*F.O. models only*)
- **Co-axial water spray** — Optimizes spray, cooling the bur and flushing the operatory site
- **100% stainless steel body** — Rugged durable construction



Choose the Right 430 for Your Practice



430 Torque Satin Finish

High Torque, high-speed handpiece with ceramic bearing autochuck turbine.

Available fiber-optic or non fiber-optic, LubeFree or Lubricated.

- 430SWL Torque LubeFree
- 430SWL Torque Lubricated
- 430SW Torque LubeFree
- 430SW Torque Lubricated



430SWL® Stainless or Satin Finish

High-speed fiber-optic handpiece with ceramic bearing autochuck turbine.

- 430SWL LubeFree
- 430SWL Lubricated
- 430SWL Satin LubeFree
- 430SWL Satin Lubricated



430SW™ Stainless or Satin Finish

High-speed non-fiber-optic handpiece with ceramic bearing autochuck turbine.

- 430SW LubeFree
- 430SW Lubricated
- 430SW Satin LubeFree
- 430SW Satin Lubricated



430SWL 45 Stainless Finish

High-speed handpiece with a 45 degree head. Available fiber-optic or non-fiber-optic, LubeFree.

- 430SWL 45 LubeFree
- 430SW 45 LubeFree



430K Stainless Finish

- 430K LubeFree
- 430K Lubricated

Personalize your handpiece with your group name or logo!

Ask your StarDental Customer Service Representative for more information.

LubeFree products include a 1-year limited warranty. Lubricated products include a 2-year limited warranty. Bifurcated solid cellular glass rod fiber optics include a 5-year limited warranty.



For more information, please visit www.DentalEZ.com